

General info

Protocol information

Protocol name	Tissue_RNA_Duo
Modified by	admin
Kit name	MagJET RNA Kit
Description	MagJET RNA Kit protocol for RNA purification from mammalian cultured cells, tissues, bacteria and yeast using KingFisher Duo Instrument.

Sample layout

1	
2	
3	
4	
5	
6	
7	
8	
9	
10	
11	
12	



Reagent info







A (DNase I)		RNA plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
1X Reaction Buffer with MgCl ₂ for DNase	200	-	Reagent	
DNase I (reconstituted)	5	-	Reagent	
B (Sample)		RNA plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Lysed sample	450	-	Reagent	
Magnetic Beads	40	-	Reagent	
Ethanol	400	-	Reagent	
C (Wash 1)		RNA plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Wash Buffer 1	700	-	Reagent	
D (Wash 2_1)		RNA plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Wash Buffer 2	700	-	Reagent	
E (Wash 2_2)		RNA plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Wash Buffer 2	700	-	Reagent	
F		RNA plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
-	-	-	-	
G (Tip Comb)		RNA plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
-	-	-	-	
H		RNA plate		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
-	-	-	-	
A (Elution Strip)		Elution strip		
Name	Well volume [µl]	Total reagent volume [µl]	Type	
Water, nuclease free	100	-	Reagent	

Dispensed reagents

A (DNase I)		RNA plate	
Name	Step	Well volume [µl]	Total reagent volume [µl]

Steps data

 Tip1	KingFisher Duo 12 tip comb	
 Pick-Up	RNA plate	(G) - Tip Comb
 Bind	RNA plate	(B) - Sample
Beginning of step	Precollect	No
	Release beads	No
Mixing / heating:	Mixing time, speed	00:05:00, Fast
	Heating during mixing	No
End of step	Postmix	No
	Collect count	3
	Collect time [s]	1
	Post-temperature	No
 Dry1	RNA plate	(B) - Sample
	Dry time	00:05:00
	Tip position	Outside well / tube
 DNase	RNA plate	(A) - DNase I
Beginning of step	Precollect	No
	Release time, speed	00:00:15, Bottom mix
Mixing / heating:	Mixing time, speed	00:15:00, Medium
	Heating temperature [°C]	37
End of step	Postmix	No
	Collect beads	No
	Post-temperature	No
 Dispense	RNA plate	(A) - DNase I
	Message	Add 200 ul Ethanol
	Dispensing volume [µl]	200
Reagent(s)	Name	Ethanol
	Volume [µl]	200
 Rebind	RNA plate	(A) - DNase I
Beginning of step	Precollect	No
	Release beads	No
Mixing / heating:	Mixing time, speed	00:05:00, Fast
	Heating during mixing	No
End of step	Postmix	No
	Collect count	3
	Collect time [s]	1
	Post-temperature	No

	Wash 1	RNA plate	(C) - Wash 1
	Beginning of step	Precollect	No
		Release time, speed	00:00:15, Bottom mix
	Mixing / heating:	Mixing time, speed	00:01:00, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	1
		Post-temperature	No
	Wash 2	RNA plate	(D) - Wash 2_1
	Beginning of step	Precollect	No
		Release time, speed	00:00:15, Bottom mix
	Mixing / heating:	Mixing time, speed	00:01:00, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	1
		Post-temperature	No
	Wash 3	RNA plate	(E) - Wash 2_2
	Beginning of step	Precollect	No
		Release time, speed	00:00:15, Bottom mix
	Mixing / heating:	Mixing time, speed	00:01:00, Fast
		Heating during mixing	No
	End of step	Postmix	No
		Collect count	3
		Collect time [s]	1
		Post-temperature	No
	Elution	Elution strip	(A) - Elution Strip
	Beginning of step	Precollect	No
		Release beads	Yes
	Mixing / heating:	Mixing time, speed	00:05:00, Fast
		Heating temperature [°C]	60
	End of step	Postmix	No
		Collect count	5
		Collect time [s]	30
		Post-temperature	No
	ReleaseBeads1	RNA plate	(C) - Wash 1
		Release time, speed	00:00:05, Fast
	Leave	RNA plate	(G) - Tip Comb

Lot info



No lot numbers have been defined.