


Corning Life Sciences Selection Guide

Issue 6



Trust Corning for Your Research

For superior quality and performance in life science research tools, you can count on Corning. Our comprehensive line of advanced products and technologies is designed to meet your evolving needs in cell-based and genomics research. To learn more, please visit www.corning.com/lifesciences, or call 1800.492.1110. Customers outside the U.S., please call 1978.442.2200 or contact your local support office (listed below).

CORNING

Corning Incorporated Life Sciences

Corning BV
Life Sciences
Koolhovenlaan 12
1119 NE Schiphol Rijk
The Netherlands

Corning Incorporated
Life Sciences
Tower 2, 4th Floor
900 Chelmsford St.
Lowell, MA 01851
t 800.492.1110
t 978.442.2200
f 978.442.2476

www.corning.com/lifesciences

Worldwide Support Offices

ASIA / PACIFIC

Australia
t 61 2-9416-0492
f 61 2-9416-0493

China
t 86 21-3222-4666
f 86 21-6288-1575

Hong Kong
t 852-2807-2723
f 852-2807-2152

India
t 91-124-235 7850
f 91-124-401 0207

Japan
t 81 (0) 3-3586 1996/1997
f 81 (0) 3-3586 1291/1292

Korea
t 82 2-796-9500
f 82 2-796-9300

Singapore
t 65 6733-6511
f 65 6861-2913

Taiwan
t 886 2-2716-0338
f 886 2-2716-0339

EUROPE

France
t 0800 916 882
f 0800 918 636

Germany
t 0800 101 1153
f 0800 101 2427

The Netherlands
t 31 20 655 79 28
f 31 20 659 76 73

United Kingdom
t 0800 376 8660
f 0800 279 1117



All Other European Countries
t 31 (0) 20 659 60 51
f 31 (0) 20 659 76 73

LATIN AMERICA

Brasil
t (55-11) 3089-7419
f (55-11) 3167-0700

Mexico
t (52-81) 8158-8400
f (52-81) 8313-8589

CORNING



New!

Online ordering
now available

CORNING

Possibilities Made Real

Corning is helping to make your research possibilities real with new and innovative products.

Introduction

Corning is pleased to present our Life Sciences Selection Guide. In this guide, you will find a selection of our newest, most innovative and most requested products.

For more than 150 years, Corning Incorporated has leveraged its materials science and process engineering expertise to collaborate closely with customers worldwide — turning what were once only possibilities into breakthrough realities.

One such reality is the Corning® Epic® System, a high-throughput label-free screening platform based on optical biosensor technology. The system performs both biochemical and cell-based drug discovery applications and offers drug developers the ability to evaluate promising new drug targets. It also allows for the observation of direct biological interactions not previously detectable in high-throughput applications.

For hard-to-attach cell lines, Corning offers a number of modified or synthetic surfaces including Corning CellBIND® Surface, Ultra-Low Attachment, and Ultra-Web™ Surfaces. If you are trying to prevent or reduce attachment, we offer plates, dishes, flasks and the CellSTACK® Culture Chambers with our Ultra-Low Attachment surface. We also offer two new vessel formats, the HYPERFlask™ Cell Culture Vessel and Low Profile flask, for conserving incubator space.


We have recently advanced our microplate line to include many enhancements and new products. Check out the new 384 well solid, low-volume and Poly-D-Lysine microplates, as well as 384 and 1536 well microplates with generic bar codes in the Microplates Section.

For up-to-date information on Corning Life Sciences' comprehensive range of products and services, go to www.corning.com/lifesciences where you can access:


- ▶ New Products Information
- ▶ Technical Information including:
 - Application Notes
 - Instruction Manuals
 - Product Bulletins
- ▶ Educational Opportunities
- ▶ Product Catalog Information
- ▶ Product Literature
- ▶ Complete Distributor Information

For additional product information, please visit www.corning.com/lifesciences, or contact your local support office (see back cover).


Try one of our newest innovations and see if it can help to make your research possibilities real.



Grow enough cells in one 1720 cm² HYPERFlask Cell Culture Vessel to seed 2,000 microplates.



Grow more *in vivo*-like cells on Ultra-Web™ synthetic surfaces.



Corning's Epic System, the new high-throughput, label-free detection technology allows researchers to bring drugs to market faster.



Improve cell attachment and yield with our Corning CellBIND surface – a nonbiological surface with high oxygen levels.

For up-to-date information on Corning Life Sciences' comprehensive range of products and services, go to www.corning.com/lifesciences.

Contents

WHAT'S NEW FROM CORNING	iii
--------------------------------------	-----

ONLINE CELL CULTURE AND ASSAY TRAINING	iv
---	----

PRODUCT ORDERING INFORMATION	v
---	---

DISTRIBUTORS	v
---------------------------	---

CELL CULTURE	1
---------------------------	---

Overview	2
----------------	---

Innovative Cell Culture Surfaces for the 21st Century	Foldout
--	---------

Corning® CellBIND® Surface	3
----------------------------------	---

Ultra-Low Attachment Surface	6
------------------------------------	---

Corning Labware with UltraWeb™ Surfaces	7
---	---

Cell Culture Flasks	9
----------------------------------	---

Cell Culture Dishes	13
----------------------------------	----

Multiple Well Plates	15
-----------------------------------	----

Transwell® Permeable Supports	20
--	----

Culture Tubes	25
----------------------------	----

Roller Bottles	26
-----------------------------	----

CellSTACK® Culture Chambers	28
--	----

CellCube® Systems	30
--------------------------------	----

E-Cube® Culture System	30
-------------------------------------	----

Spinner Flasks	31
-----------------------------	----

Erlenmeyer Flasks	36
--------------------------------	----

Cell Scrapers and Lifters	38
--	----

Pipets	99
---------------------	----

Spatulas and Microspatulas	112
---	-----

Centrifuge Tubes	113
-------------------------------	-----

Cryogenic Vials and Accessories	116
--	-----

Technical Appendix	39
---------------------------------	----

Corning Cell Culture Surfaces	39
-------------------------------------	----

Characteristics of Corning Plasticware	41
--	----

Chemical Compatibility of Corning Plasticware	42
---	----

Characteristics of Corning Centrifuge Tubes	42
---	----

MICROPLATES	45
--------------------------	----

Overview	46
----------------	----

Selecting the Best Microplate	47
--	----

96 Well Microplates	50
----------------------------------	----

Corning 96 Well Microplate Selection Guide	Foldout
--	---------

96 Well Assay Microplates	50
---------------------------------	----

96 Well Polystyrene Stripwell™ Microplates	54
--	----

96 Well Polypropylene Storage Microplates and Cluster Tubes	55
---	----

384 Well Microplates	57
-----------------------------------	----

Corning 384 Well Microplate Selection Guide	Foldout
---	---------

384 Well Assay Microplates	57
----------------------------------	----

384 Well Polystyrene Microplates	58
--	----

384 Well Polypropylene Storage Microplates	60
---	----

1536 Well Microplates	61
------------------------------------	----

Protein Crystallization Microplates	63
--	----

Microplate Accessories	65
-------------------------------------	----

Lids	65
------------	----

Storage Mats and Accessories	65
------------------------------------	----

Robolids	66
----------------	----

Sealing Tapes	66
---------------------	----

Technical Appendix	67
---------------------------------	----

Properties of Some Thermoplastics	67
---	----

Selected Corning Technical Literature	71
---	----

GENOMICS	73
-----------------------	----

Overview	74
----------------	----

Colony Picking, Bacterial Growth, and Storage	75
--	----

Purification	76
---------------------------	----

Quantitation and Detection	79
---	----

DNA Amplification	81
--------------------------------	----

Bar Code Customization	84
-------------------------------------	----

(continued next page)

Contents (CONTINUED)


MICROARRAY	85
Overview	86
Slide Selection Chart	87
Microarray Printing	88
Epoxide Coated Slides	88
UltraGAPS™ Coated Slides	89
GAPS™ II Coated Slides	90
Pronto!™ Universal Spotting Solution	91
Pronto! Epoxide Spotting Solution	92
Corning® Cover Glass	92
384 Well Microarray Printing Plates	93
Microarray Slide Mailers/Storage Boxes	93
Microarray Storage Pouches	94
Microarray Processing	94
Pronto! Universal Hybridization Kits	94
Pronto! Background Reduction Kit	96
Hybridization Chambers	96
LIQUID HANDLING	97
Overview	98
Pipets	99
Pipetting Aids	100
Pipettors	100
Pipet Tips	101
Reagent Reservoirs	104
Transtar-96® Well Liquid Transfer System	105
Aspirator	105
Vacuum Filters	105
Syringe Filters	108
Spin-X® Centrifuge Tube Filters	108
Storage Bottles	109
Containers	110
Cylinder	111
Erlenmeyer Flasks	111
Spatulas	112
Centrifuge Tubes	113
Microcentrifuge Tubes	115
Cryogenic Vials and Accessories	116
Technical Appendix	119
Selecting the Best Filter for Your Application	119
Characteristics of Corning Plasticware	122
Chemical Compatibility of Corning Plasticware	123
Characteristics of Corning Centrifuge Tubes	123
DISPOSABLE GLASS AND EQUIPMENT	127
CORNING® and PYREX® Disposable Glassware	128
Digital Stirring Hot Plates	135
INDEXES	137
Alphabetical Index	137
Catalog Number Index	139

What's New from Corning



CELL CULTURE & BIOPROCESS (Page 1)

- ▶ Corning® HYPERFlask™ Vessels
- ▶ Corning Low Profile 100 cm² Flask
- ▶ Ultra-Web™ Synthetic Surfaces
- ▶ Transwell® Permeable Supports Coated with Cultrex® Basement Membrane Extract
- ▶ Corning 1L and 3L Disposable Plastic Spinner Flasks
- ▶ Corning Aseptic Transfer Caps
- ▶ Corning Baffled Erlenmeyer flasks



HTS AND ASSAY MICROPLATES (Page 45)

- ▶ Corning Generic Bar Coded Polystyrene Microplates
- ▶ Corning Low Volume 384 Well Solid Black Flat Bottom Microplates
- ▶ Corning 1536 Well Echo™ Qualified Microplates
- ▶ Corning 384 Well Low Flange Solid Black and White Microplates



GENERAL LABWARE AND EQUIPMENT (Page 97)


- ▶ Corning Square Polycarbonate Storage Bottles
- ▶ PYREX® Wide Mouth Storage Bottles
- ▶ Clean Room Packaging for Pipets
- ▶ Corning CentriStar™ Centrifuge Tubes
- ▶ Corning Digital Hot Plates and Stirrers



PROTEOMICS AND GENOMICS (Page 73)

- ▶ Next Generation CrystalEX™ Protein Crystallization Microplates
- ▶ Corning Epoxide Coated Slides

Online Cell Culture and Assay Training



The Corning Scientific Seminar Series is a series of free online technical presentations that provide novel tips, best practices and proven techniques to help advance your research. Delivered by scientists, these one hour sessions have proven useful for technicians as well as for researchers who have been doing cell culture and assays for years.

Join us online for an upcoming session or download any of our previously recorded sessions. New topics are added monthly. Former topics include:

- ▶ Cell Culture Contamination – Every Researcher's Nightmare!
- ▶ Effects of Cell Culture Surfaces on Cellular Behavior
- ▶ Detecting, Removing and Managing Mycoplasma Contamination
- ▶ Grow More Cells! Scaling Up Cell Production
- ▶ Growing Happier Cells
- ▶ Growing Cells on Transwell Inserts – Tips and Techniques
- ▶ Life and Death *In Vitro* - Growth and Toxicity
- ▶ HeLa Cells - A Blessing or a Curse?
- ▶ More *In Vivo*-like Cell Cultures and Better Assays with Permeable Supports
- ▶ Optimizing Assay Performance through Microplate Attributes and Equipment Setting
- ▶ Primary Cell Culture – Tips and Techniques for Getting Started
- ▶ Solving Cell Culture Problems
- ▶ Using Frozen vs. Continuously Cultured Cells for HTS

Register at www.corning.com/lifesciences

What attendees had to say about past seminars:

"We are not getting information like this from anywhere. The seminar was amazing, very useful to my work. Thanks for organizing these sessions."

"Great tips! I'll pass along information I learned here to my colleagues to let them know how we should conduct cell culture properly. Many thanks again indeed."

"I use your seminars as training for new employees and estimate they save my company more than \$24,000 a year in training costs."



All attendees receive a certificate of completion.

Training is co-sponsored by:



CORNING



Product Ordering Information

Corning products are available as follows:

- ▶ For quick and easy ordering with up-to-the-minute product information, you can now order Corning products online. Please visit www.corning.com/lifesciences to register and begin placing your orders online.
- ▶ Directly from Corning's European Distribution Centre. To place your order, simply contact us as follows:

France

t 0800 916 882
f 0800 918 636
ServiceClients@
corning.com

United Kingdom

t 0800 376 8660
f 0800 279 1117
UKCustomerService@
corning.com

Germany

t 0800 101 1153
f 0800 101 2427
Kundenbetreuung@
corning.com

The Netherlands and All Other Countries

t +31 20 659 6051
f +31 20 659 7673
cceurnl@corning.com

- ▶ Through any Corning Distributor – please visit our website or see pages iii-vii for a full listing. To place an order, simply contact the Distributor of your choice. For each requested product, provide the Corning catalogue number, product description, and desired quantity.

Distributors

Visit www.corning.com/lifesciences for a complete listing of worldwide distributors.

Austria

SZABO-SCANDIC HandelsgmbH & Co KG

Quellenstr. 110
A-1100 Vienna
Austria
Phone: 43-1-489-3961-0
Fax: 43-1-489-3961-7
Email: mail@szabo-scandic.com
Web site: www.szabo-scandic.com
Products: Plastics

Sigma-Aldrich Handels GmbH

Vienna
Austria
Phone: 43-1-60-581-10
Fax: 43-1-60-581-20
Email: sigma@sigma.co.at
Web site: www.sigma-aldrich.com
Products: Plastics

Bahrain

Attieh Medico

P.O. Box 116105
21391 Jeddah
Saudi Arabia
Phone: 966-2-661-3613
Fax: 966-2-665-1358
Products: Plastics

Belgium

Elscolab NV

Industriepark Zwalubeek
B-9150 Kruibeke
Belgium
Phone: 32-3-250-1570
Fax: 32-3-252-8783
Email: elscolab@elscolab.com
Web site: www.elscolab.com
Products: Plastics

Fisher Scientific Bioblock Sprl

P.O. Box 567
B-7500 Tournai 1 / Doornik 1
Belgium
Phone: 056 260 260
Fax: 056 260 270
Email: belgium@bioblock.com
Web site: www.bioblock.be
Products: Plastics

Sigma-Aldrich NV/SA

Bornem
Belgium
Phone: 0800 14747
Fax: 0800 14745
Email: becusts@Europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

Czech Republic**BARIA s.r.o**

Jizni 393
252 44 Psary
Czech Republic
Phone: 800 555 056 / +420 244 911 228
Fax: +420 241 941 519
Email: info@baria.cz
Web site: www.baria.cz
Products: Plastics

GeneTica s.r.o.

Sluzeb 4
10852 Praha 10
Czech Republic
Phone: +420 272 701 055
Fax: +420 272 701 055
Email: kvapil@genetica.cz
Web site: www.genetica.cz
Products: Plastics

Sigma-Aldrich spol. S.r.o

Pobrezni 46
186 00 Praha 8
Prague
Czech Republic
Phone: +420-246 003 200
Fax: +420-246 003 291
Email: czecustsv@Europe.sial.com
Web site: www.sigma-aldrich.com/czech
Products: Plastics

Denmark**Fisher Biotech Line A/S**

Industrivej 3
DK-3550 Slangerup
Denmark
Phone: 45-70-27-99-20
Fax: 45-70-27-99-29
Email: biotechline@biotechline.dk
Web site: www.biotechline.dk
Products: Plastics

Sigma-Aldrich Denmark A/S

Kirkebjerg Alle 84
2605 Brondby, Copenhagen
Denmark
Phone: 45-43-56-59-10
Fax: 45-43-56-59-05
Email: denorder@europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

Egypt**Clinilab**

4,106 St., El-Etehad Square
21 El-Manial, Riham Tower, El-Maadi
11431 Cairo
Egypt
Phone: 20 2 525 7212
Fax: 20 2 525 7210
Email: clinilab@intouch.com
Products: Plastics

Estonia**Bio-Art**

Pepleri Str. 12-23
EE 51003+B16 Tartu
Estonia
Phone: 372-7-427085
Fax: 372-7-427085
Email: bioart@kodu.ee
Products: Plastics

Finland**LABNET**

Hermannin Rantatie 12B
00580 Helsinki
Finland
Phone: 358-9-4764-4300
Fax: 358-9-4764-4310
Email: ml.varis@rake.fi
Web site: www.labnet.fi
Products: Plastics

Sigma-Aldrich Finland

(YA-Kemia Oy)
Helsinki
Finland
Phone: 358 9 350 9250
Fax: 358 9 350 92555
Email: finorder@europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

France**Sigma-Aldrich Chimie S.a.r.l.**

Saint Quentin Fallavier (38)
France
Phone: 0800 211 408
Fax: 0800 031 052
Email: fradvsv@Europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

Fisher Scientific Bioblock

BP 50111
Parc d'innovattion
F-67403 ILLKIRCH CEDEX
France
Phone: 33-3-88-67-14-14
Fax: 33-3-88-67-85-11
Email: vente@bioblock.fr
Web site: www.bioblock.com
Products: Plastics

ATGC Biotechnologie

Rue Des Freres Montgolfiers
BP 5 Croissy
Marne Laval
F-77313
France
Phone: 33-1-60-95-51-00
Fax: 33-1-60-95-51-01
Email: contact@atgc.fr
Web site: www.ATGC.fr
Products: Plastics

Elvetec Services

ZAC Satolas Green
 Pusignan
 69881 MEYZIEUX
 France
 Phone: 33-4-72-22-43-00
 Fax: 33-4-72-22-43-00
 Email: elvetec@elvetec.fr
 Web site: www.elvetec.fr
 Products: Plastics

Germany**Fisher Scientific GmbH**

Im Heiligen Feld 17
 58239 Schwerte
 Germany
 Phone: 49 2304/932-5
 Fax: 49 2304/932950
 Web site: www.fishersci.com
 Products: Plastics

Sigma-Aldrich Chemie GmbH

Munich
 Germany
 Phone: 089-65 13-0
 Fax: 089-65 13-11 69
 Email: deorders@Europe.sial.com
 Web site: www.sigma-aldrich.com
 Products: Plastics

Omni Life Science GmbH & Co. KG

Schnackenburgallee 1160
 22525 Hamburg
 Germany
 Phone: 49 (0) 40 54774640
 Fax: 49 (0) 40 547746429
 Web site: www.ols-bio.de
 Products: Plastics

Schubert & Weiss OMNILAB GmbH & Co. KG

Fromundstrasse 34
 81547 Munchen
 Germany
 Phone: 49 89/6925718
 Fax: 49 89/6912413
 Products: Plastics

OMNILAB - LABORZENTRUM GmbH

Robert-Hooke-Str. 8
 28359 Bremen
 Germany
 Phone: 49 (0) 421 17 590
 Fax: 49 (0) 421 175 99 300
 Products: Plastics

Greece**Sigma-Aldrich (O.M.) LTD**

Athens
 Greece
 Phone: 0109948010
 Fax: 0109943831
 Email: grcustsv@Europe.sial.com
 Web site: www.sigma-aldrich.com
 Products: Plastics

ANTISEL A Selidis Bros S A

35 Anaximandrou
 6 Flem. Str.
 54250 Harilaou-Thessaloniki
 Greece
 Phone: 30-31-322525
 Fax: 30-31-321912
 Products: Plastics

Hungary**Sigma-Aldrich Kft**

Budapest
 Hungary
 Phone: (06-1) 269-6474
 Fax: (06-1) 235-9050
 Email: info@sigma.sial.hu
 Web site: www.sigma-aldrich.com
 Products: Plastics

Zenon Biotechnology Ltd.

Maros Str. 40
 H-6721 Szeged
 Hungary
 Phone: 36-62-424-290
 Fax: 36-62-476-558
 Email: zenonbio@deltav.hu
 Web site: www.zenonbio.hu
 Products: Plastics

Iceland**Lyfjadreifing ehf**

Icelandic Pharmaceuticals Ltd.
 Krokhsals 14
 Reykjavik 110
 Iceland
 Phone: 354-540-8000
 Fax: 354-540-8001
 Products: Plastics

Ireland**Fisher Scientific Ireland Ltd.**

Suite 3, Plaza 212
 Blanchards Town Corporate Park 2
 Ballycoolin
 Dublin 15
 Ireland
 Phone: +353 (0) 1885 5854
 Fax: +353 (0) 1899 1855

Sigma-Aldrich Ireland Ltd

Dublin
 Ireland
 Phone: 800 200 888
 Fax: 800 600 222
 Email: eicustsv@Europe.sial.com
 Web site: www.sigma-aldrich.com
 Products: Plastics

Israel

Getter (2000) Ltd.
P.O. Box 3500
52134 Ramat Gan
Israel
Phone: 972-3-5761520
Fax: 972-3-7523620
Products: Plastics

Italy

Celbio S.P.A.
Via Figino 20/22
20016 Pero (MI)
Italy
Phone: 39-02-3810-1465
Fax: 39-02-3819-557
Email: info.bio@celbio.it
Web site: www.celbio.it
Products: Plastics

Sigma-Aldrich S.r.l.
Via Gallarate, 154
20151 Milano
Italy
Phone: 02 33417310
Fax: 02 38010737
Email: itorder@Europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

Kenya

Jon.Achells & Sohne GmbH
P.O. Box 107727
D 28077 Bremen
Germany
Phone: 49-421-33670
Fax: 20-2-525-7210
Products: Plastics

Sciencescope LTD
Victor House
Kimathi street
PO Box 72963
0020 City Square
Nairobi, Kenya
Tel: 254 20 228 331
Fax: 254 20 318 221

Kuwait

Attieh Medico
P.O. Box 116105
21391 Jeddah
Saudi Arabia
Phone: 966-2-661-3613
Fax: 966-2-665-1358
Products: Plastics

Luxemburg

Fisher Bioblock Scientific Sprl
P.O. Box 567
B-7500 Tournai 1/Doornik 1
Belgium
Phone: 056 260 260
Fax: 056 260 270
Email: belgium@bioblock.com
Web site: www.bioblock.be
Products: Plastics

Netherlands

Fisher Emergo bv
A Fisher Scientific Company
P.O.Box 4
NL 1120 AA Landsmeer
The Netherlands
Phone: 020-4877000
Fax: 020-4877070
Email: info@emergofisher.com
Web site: www.fisheremergo.com
Products: Plastic

Omnilabo International BV
P.O.Box 3939
4800 DX Breda
The Netherlands
Phone: 31 76 5 795 795
Fax: 31 76 5 876 236
Email: info@omnilabo.nl
Web site: www.omnilabo.nl
Products: Plastics

Sigma-Aldrich Chemie BV
Zwijndrecht
The Netherlands
Phone: 0800 0229088
Fax: 0800 0229089
Email: nlcustsv@Europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

Norway

Sigma-Aldrich Norway AS
Tevlingveien 23
NO-1081 Oslo
Norway
Phone: 47-23-176000
Fax: 47-23-176010
Email: norsigma@europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

Oman

Attieh Medico
P.O. Box 116105
21391 Jeddah
Saudi Arabia
Phone: 966-2-661-3613
Fax: 966-2-665-1358
Products: Plastics

Poland**Sigma-Aldrich Sp.z.o.o.**

Poznan
Poland
Phone: 61 829 01 00
Fax: 61 829 01 20
Email: plcustsv@Europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

Biomibo

ul. Strazyglowska 15
04-872 Warszawa
Poland
Phone: 48-22-8720797 / 6126196
Fax: 48-22-8720797 / 6126196
Email: biomibo@supernet.pl
Products: Plastics

Portugal**Sigma-Aldrich Química, S.A.**

Sintra
Portugal
Free Phone: 800 20 21 80
Free Fax: 800 20 21 78
Email: poorders@Europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

Reagente 5

Rua Eng. Ezequil de Campos
330-B, Apartado 4710
4012 Porto Codex
Portugal
Phone: 351-2-6198090
Fax: 351-2-6198099
Products: Plastics

Qatar**Attieh Medico**

P.O. Box 116105
21391 Jeddah
Saudi Arabia
Phone: 966-2-661-3613
Fax: 966-2-665-1358
Products: Plastics

Russia**Corning Costar Moscow**

P.O. Box 33
Moscow 105484
Russia
Phone: 7-095-4612208
Fax: 7-095-4612208
Email: cosmos@orc.ru
Products: Plastics

Saudi Arabia**Attieh Medico**

P.O. Box 116105
21391 Jeddah
Saudi Arabia
Phone: 966-2-661-3613
Fax: 966-2-665-1358
Products: Plastics

Slovenia**Labormed d.o.o.**

Zg. Pirnice 96 c
Slo-1215 Medvode
Slovenia
Phone: 386-1-3621-414
Fax: 386-1-3621-415
Email: info@labormed.si
Web site: www.labormed.si
Products: Plastics

South Africa**Adcock Ingram Critical Care Pty Ltd**

P.O. Box 6888
2000 Johannesburg
South Africa
Phone: 27-11-494-8000
Fax: 27-11-494-8305
Email: adcock.scientific@tigerbrands.com
Web site: www.aiscientific.co.za
Products: Plastics

Spain**Cultek S.L.**

Avada. Cardenal Herrera
Oria 63
28034 Madrid
Spain
Phone: 34-91-729-0333
Fax: 34-91-358-1761
Email: cultek@cultek.com
Web site: www.cultek.com
Products: Plastics

Sweden**Fisher Scientific GTF**

S.Långebergsgatan 30
421 32 V Frölunda/Göteborg
Sweden
Phone: 46-3-168-94-00
Fax: 46-3-168-07-17
Email: GTF@gtf.se
Web site: www.gtf.se
Products: Plastics

DISTRIBUTORS

Sigma-Aldrich Sweden AB

Solkraftsvägen, 14C
SE-13570 Stockholm
Sweden
Phone: 46-8-742-42-00
Fax: 46-8-742-42-43
Email: sweorder@europe.sial.com
Web site: www.sigma-aldrich.com
Products: Plastics

Nordic Biolabs

Box 7293
Ritarslingan 16
18714 Taby
Sweden
Phone: 46-8-630-8500
Fax: 46-8-630-0905
Email: info@nordicbiolabs.se
Web site: www.nordicbiolabs.se
Products: Plastics

Switzerland

VITARIS AG

Blegistrasse 9
CH-6340 Baar
Switzerland
Phone: 41-41-769-00-00
Fax: 41-41-769-00-01
Email: info@vitaris.com
Web site: www.vitaris.com
Products: Plastics

Fisher Scientific AG

Wilstrasse 57
CH-5610 Wohlen
Switzerland
Phone: 41-56-618-41-11
Fax: 41-56-618-41-41
Email: info@ch.fishersci.com
Web site: www.ch.fishersci.com
Products: Plastics

United Arab Emirates

Attieh Medico

P.O. Box 116105
21391 Jeddah
Saudi Arabia
Phone: 966-2-661-3613
Fax: 966-2-665-1358
Products: Plastics

Emphor FZCO

P.O. Box 61232
Jebelau Free Zone
Dubai, UAE
Phone: 971 4 8830233
Fax: 971 50 457 3009
Email: ravi@emphor.biz
Web site: www.emphor.biz
Products: Plastics

United Kingdom

Fisher Scientific UK.

Bishop Meadow Road
Loughborough
Leicestershire
UK
Phone: 44-1509-231166
Fax: 44-1509-231893
Email: sales@fisher.co.uk
Web site: www.fisher.co.uk
Products: Plastics

Appleton Woods Limited

Lindon House
Heeley Road
Selly Oak
Birmingham
UK
Phone: 44-121-472-7353
Fax: 44-121-414-1075
Email: info@appletonwoods.co.uk
Web site: www.appletonwoods.co.uk
Products: Plastics


Sigma-Aldrich Company Ltd

Fancy Road
Poole
Dorset
BH12 4QH
UK
Phone: 44 1202 733114
Fax: 44 1202 715460
Email: ukcustsv@Europe.sial.com
Web site: www.sigma-aldrich.com/uk
Products: Plastics

Western Laboratory Consumables

Safeguard House
Hazleton Interchange
Horndean
Waterlooville
Hampshire
UK
Phone: 44 0844 8936501
Fax: 44 0870 7746540
Email: sales@consumables-uk.com
Web site: www.consumables-uk.com
Products: Plastics


Cell Culture



Low Profile Flasks, page 11



HYPERFlask™ Cell Culture Vessel,
page 12, 27



CellSTACK Culture Chambers, page 28

OVERVIEW	2
INNOVATIVE CELL CULTURE SURFACES FOR THE 21ST CENTURY	Foldout
CORNING® CELLBIND® SURFACE	3
ULTRA-LOW ATTACHMENT DISHES, PLATES, FLASKS, AND CELLSTACK® CULTURE CHAMBERS	6
CORNING LABWARE WITH ULTRA-WEB™ SURFACES	7
CELL CULTURE FLASKS	9
CELL CULTURE DISHES	13
MULTIPLE WELL PLATES	15
TRANSWELL® PERMEABLE SUPPORTS	20
CULTURE TUBES	25
ROLLER BOTTLES	26
CELLSTACK® CULTURE CHAMBERS	28
CELLCUBE® SYSTEMS	30
E-CUBE® CULTURE SYSTEM	30
SPINNER FLASKS	31
ERLENMEYER FLASKS	36
CELL SCRAPERS AND LIFTERS	38
PIPETS	99
SPATULAS AND MICROSPATULAS	112
CENTRIFUGE TUBES	113
CRYOGENIC VIALS AND ACCESSORIES	116
TECHNICAL APPENDIX	39
Corning Cell Culture Surfaces	39
Characteristics of Corning Plasticware	41
Chemical Compatibility of Corning Plasticware	42
Characteristics of Corning Centrifuge Tubes	42

Overview

DESIGNED FOR PERFORMANCE

Corning Life Sciences offers a full line of cell culture products that are manufactured under strict process controls guaranteeing consistent product performance. All Corning Life Sciences plastics manufacturing facilities are ISO 9001:9002 registered. ISO registration is recognized worldwide as a standard of excellence for quality systems.

In addition, customers can now obtain a Certificate of Compliance or product description for any Corning® or Costar® cell culture product from our website. This certificate details lot-specific information on component materials, sterility testing, pyrogen testing, cell attachment, and growth characteristics.

Also available are detailed drawings that highlight product dimensions. Drawings are available simply by calling your local Corning Life Sciences office.

ADDITIONAL QUALITY ASSURANCES

Nonpyrogenic Certification

Most Corning and Costar cell culture products are certified nonpyrogenic with a documented endotoxin level of equal to or less than 0.1 EU/mL. Endotoxins have been shown to cause variability in cell culture. Nonpyrogenic certification is just another way Corning helps ensure consistent cell culture results. Corning also offers a detailed technical bulletin on the effects of endotoxins in cell culture. This may be obtained by calling your local Corning Life Sciences office or by downloading the bulletin from the Corning web site www.corning.com/lifesciences.

Lot Number Traceability

To ensure accurate lot number traceability in biotechnology research and production facilities, most Corning and Costar cell culture flasks and most roller bottles feature a lot number individually printed on each product. Lot number traceability helps simplify quality assurance procedures for tracking and monitoring production and research processes.





Consistent Surface Chemistry

All Corning and Costar cell culture products are produced in ISO-certified facilities. Cell culture products are made from USP Class VI materials in accordance with documented manufacturing procedures. By carefully controlling both the materials we use and our manufacturing process, Corning is able to provide consistent surface chemistries across our entire line of cell culture products. This consistency increases the researcher's ability to produce reliable results.



Corning® CellBIND® Surface

A Novel Surface for Improved Cell Attachment, Serum Reduction, or the Elimination of Coatings



Corning CellBIND Surface is now available on flasks, CellSTACK® Culture Chambers, multiple well plates, 96 and 384 well plates, dishes, and roller bottles.

Increase Cell Growth and Yields with Corning CellBIND Surface

The Corning CellBIND surface enhances cell attachment under difficult conditions, such as reduced-serum or serum-free medium, resulting in higher cell yields. The first novel cell culture surface treatment in over 20 years.

Developed by Corning scientists, this patented technology (U.S. Patent No. 6,617,152) uses a microwave plasma process for treating the culture surface. This process improves cell attachment by incorporating significantly more oxygen into the cell culture surface, rendering it more hydrophilic (wetter) and increasing surface stability.

Benefits

- ▶ May eliminate the need for tedious, time-consuming, expensive and low stability biological coatings
- ▶ More quickly adapts cells to reduced-serum or serum-free conditions
- ▶ Increase cell survival following cryo-preservation
- ▶ Reduces premature cell detachment from confluent cultures especially in roller bottles
- ▶ Better cell attachment leads to increased cell growth and yields

- ▶ More consistent and even cell attachment
- ▶ Requires no refrigeration or special handling and is stable at room temperature

Same High Quality Standards as Other Corning Vessels

- ▶ Manufactured from optically clear polystyrene
- ▶ Rigorous QC testing for consistency and reproducibility
- ▶ Certified nonpyrogenic and sterile
- ▶ Lot numbers for quality assurance and tracking
- ▶ Corning CellBIND surface logo differentiates from standard treatment cell culture products and avoids mix-ups

Cell Dissociation Recommendations

Culture inoculating and harvesting should be performed in the same manner as methods currently being employed. Both enzymatic and nonenzymatic dissociating solutions have been successfully used to remove cells from Corning CellBIND surfaces. These include: Trypsin-EDTA, Accutase®, Versene®, Dispase®, and Citric Saline. Some dissociating agents, such as Dispase or Versene, should be removed by centrifugation prior to plating the cells.



HYPERFlask™ Cell Culture Vessel



CellSTACK Culture Chambers

Enhanced Attachment of LNCaP Cells to the Corning® CellBIND® Surface*




Figure 1. Left: Adherent cell recovery and growth of LNCaP cells 24 hours post-seeding. Data is average \pm standard error from 3 independent experiments. Right: Average \pm standard error from 3 independent experiments for 7 day growth after initial attachment.





Figure 2. Attachment of LNCaP cells. Cells were thawed and plated onto the Corning CellBIND Surface (right) or tissue culture treated (left) T25 flasks. 24 hours post seeding a random field was viewed by light microscopy (100X magnification).

*From *Enhanced Attachment of LNCaP Cells to the Corning CellBIND Surface*, Corning SnAPPShot publication CLS-AN-048.

Corning CellBIND Surface Product Ordering Information



Cat. No.	Description	Qty/Pk	Qty/Cs
<i>Roller Bottles</i>			
3907	Roller Bottle, 850 cm ² , Corning CellBIND Surface, Easy Grip Cap, Sterile	2	40
431134	Expanded Surface Roller Bottle, 1700 cm ² , Corning CellBIND Surface, PS, Easy Grip Cap, Sterile	20	20
431329	Roller Bottle, 850 cm ² , Corning CellBIND Surface, Vent Cap, Sterile	2	40
431344	Roller Bottle, 850 cm ² , Corning CellBIND Surface, Easy Grip Cap, Sterile	22	44
<i>Flasks</i>			
3289	Flask, 25 cm ² , Corning CellBIND Surface with Vent Cap, Sterile	20	200
3290	Flask, 75 cm ² , Corning CellBIND Surface with Vent Cap, Sterile	5	100
3073	Low Profile Flask, 100 cm ² , Corning CellBIND Surface with Vent Cap, Sterile	6	60
3291	Flask, 150 cm ² , Corning CellBIND Surface with Vent Cap, Sterile	5	50
3292	Flask, 175 cm ² , Corning CellBIND Surface with Vent Cap, Sterile	5	50
3293	Flask, 225 cm ² , Corning CellBIND Surface with Vent Cap, Sterile	5	25
3298	Flask, 175 cm ² , Corning CellBIND Surface with Phenolic Cap, Sterile	5	50
431328	Flask, 175 cm ² , Corning CellBIND Surface, Bar coded with Vent Cap, Sterile	7	84

Corning® CellBIND® Surface Product Ordering Information (Continued)

Cat. No.	Description	Qty/ Pk	Qty/ Cs
New 10010	HYPERFlask™ Vessel, 1720 cm ² , Corning CellBIND Surface, Bar Code, Sterile	4	4
New 10024	HYPERFlask Vessel, 1720 cm ² , Corning CellBIND Surface, Bar Code, Sterile	4	24
431346	Expanded Surface Flask, 235 cm ² , Corning CellBIND Surface, with Bar Code, Vent Cap, Sterile	7	42
3068	RoboFlask® Tissue Culture Vessel for Automation, 92.6 cm ² growth area, Corning CellBIND Surface, with Bar Code, Septum Cap, Sterile	10	50
3067	RoboFlask Tissue Culture Vessel for Automation, 92.6 cm ² growth area, Corning CellBIND Surface, with Bar Code, Septum Cap, Sterile	20	100

CellSTACK® Culture Chambers

3330	CellSTACK-1 Chamber, 636 cm ² growth area, Corning CellBIND Surface, Sterile	1	8
3310	CellSTACK-2 Chamber, 1,272 cm ² growth area, Corning CellBIND Surface, Sterile	1	5
3311	CellSTACK-5 Chamber, 3,180 cm ² growth area, Corning CellBIND Surface, Sterile	1	2
3312	CellSTACK-10 Chamber, 6,360 cm ² growth area, Corning CellBIND Surface, Sterile	1	2
3320	CellSTACK-10 Chamber, 6,360 cm ² growth area, Corning CellBIND Surface, Sterile	1	6
3321	CellSTACK-40 Chamber, 25,440 cm ² growth area, Corning CellBIND Surface, Sterile	1	2

Dishes

3294	Dish, 35 x 10 mm style, Corning CellBIND Surface, Sterile	10	210
3295	Dish, 60 x 15 mm style, Corning CellBIND Surface, Sterile	7	126
3296	Dish, 100 x 20 mm style, Corning CellBIND Surface, Sterile	5	40

Multiple Well Plates

3335	6 Well Plate, Corning CellBIND Surface, Clear, Sterile, with Lid	5	50
3336	12 Well Plate, Corning CellBIND Surface, Clear, Sterile, with Lid	5	50
3337	24 Well Plate, Corning CellBIND Surface, Clear, Sterile, with Lid	5	50
3338	48 Well Plate, Corning CellBIND Surface, Clear, Sterile, with Lid	5	50

Microplates

3300	96 Well Plate, Corning CellBIND Surface, Clear Bottom, Sterile, with Lid	5	50
3340	96 Well Plate, Corning CellBIND Surface, Black/Clear Bottom, Sterile, with Lid	5	50
3683	384 Well Plate, Corning CellBIND Surface, Black/Clear Bottom, Sterile, with Lid	10	50

Ultra-Low Attachment Dishes, Plates, Flasks, and CellSTACK® Culture Chambers

The Ultra-Low Attachment surface is a unique covalently bonded hydrogel surface that is hydrophilic and neutrally charged. It minimizes cell attachment, protein absorption and enzyme activation. The surface is noncytotoxic, biologically inert and nondegradable.



3261 and 3262 Ultra-Low Attachment Dishes



3814 T-75 Flask

Ultra-Low Attachment Tip

The Ultra-Low Attachment products may be useful for:

- ▶ Maintaining cells in a suspended, unattached state
- ▶ Preventing stem cells from attachment-mediated differentiation
- ▶ Preventing anchorage-dependent cells from dividing
- ▶ Reducing binding of attachment and serum proteins to the substrate

Ultra-Low Attachment Dishes Ordering Information

Cat. No.	Dish Style (mm)*	Height (mm)	Growth Area (cm ²)	Qty/Pk	Qty/Cs
3261	60	15	21	5	20
3262	100	20	55	5	20

*60 mm dish = 52.1 mm; 100 mm dish = 83.8 mm

Ultra-Low Attachment Plates Ordering Information


Cat. No.	Plate Type	Bottom Type	Diameter (mm)	Growth Area (cm ²)	Qty/Pk	Qty/Cs
3471	6 well plate	Flat	34.8	9.5	1	24
3473	24 well plate	Flat	15.6	1.9	1	24
3474	96 well plate	Flat	6.4	0.32	1	24
New 7007	96 well plate	Round	6.4	0.32	1	24

Ultra-Low Attachment Flasks Ordering Information

Cat. No.	Flask Type	Cap Style	Growth Area (cm ²)	Qty/Pk	Qty/Cs
3815	Rectangular	Vent	25	6	24
3814	Rectangular	Vent	75	4	24

Ultra-Low Attachment CellSTACK Ordering Information

Cat. No.	Flask Type	Cap Style	Growth Area (cm ²)	Qty/Pk	Qty/Cs
3303	CellSTACK Chamber, 1-Stack	Vent	636	1	8



Comparison of Cell Attachment in Ultra-Low vs. Standard Tissue Culture Treated Plates

Vero cells plated at 2.6×10^6 cells per well grown for 4 days at 37°C in a 5% CO₂ environment show a 99% reduction in cellular attachment vs. standard culture treated product.





Suggested working volumes for Ultra-Low Attachment products:

- ▶ 96 well plate: 0.1 to 0.2 mL/well
- ▶ 24 well plate: 0.4 to 0.8 mL/well
- ▶ 6 well plate: 1.9 to 3.8 mL/well
- ▶ 60 mm dish: 4.2 to 6.3 mL/dish
- ▶ 100 mm dish: 11.0 to 16.5 mL/dish
- ▶ 25 cm² flask: 5 to 7.5 mL/flask
- ▶ 75 cm² flask: 15 to 22.5 mL/flask
- ▶ CellSTACK Chamber, 1-STACK: 127 to 190 mL/stack

There are no special procedures that need to be followed in order to use this surface.

Corning® Labware with Ultra-Web™ Surfaces

Innovative synthetic surfaces which provide cells with 3-D *in vivo*-like nanofibrillar substrates for better performance and functionality



SEM cross sectional (top left) and top down (top right) images (10,000X) of the Ultra-Web Synthetic Surface. Pores between the nanofibers are so small that the cells grow on the surface rather than in it. Primary Rat Cortex Neurons Ultra-Web (bottom left): surfaces support neuronal attachment, neurite outgrowth, and promote earlier neuronal aggregation, and are compatible with fluorescent staining and visualization or photography of cells. Primary rat cortex neurons (20X) were stained with anti-Tuj-1 (Class III beta tubulin), Hoechst (nuclei) and anti-Synapsin I. Synapsin I binds synaptic vesicles containing neurotransmitters to components of the cytoskeleton and facilitate their release, and thus indicate preserved neuronal function. Synapsin staining is greatest in areas between bundle formation on Ultra-Web surface (red stain, indicated by arrows). Corning Labware with Ultra-Web synthetic surface (bottom right).

Corning® Labware with Ultra-Web™ Surfaces

Technical Information

Ultra-Web Synthetic Surfaces are composed of randomly orientated electrospun polyamide nanofibers with an average fiber diameter of ~180 nm. This creates a culturing substrate that mimics structural components within the basement membrane or extracellular matrix.

Ultra-Web synthetic surface is available with two surface chemistries:

- ▶ Untreated electrospun polyamide nanofibers with an uncharged slightly hydrophilic surface
- ▶ Polyamine treated electrospun polyamide nanofibers with a positively charged surface for enhanced cell attachment or binding and covalently linking biomolecules

Ultra-Web Synthetic Surface Applications

- ▶ Ideal for culturing liver, neuronal, kidney and stem cell lines or primary cultures where current surfaces do not provide the necessary culture performance or function
- ▶ Replacement for poly-lysine or animal-derived biological coatings
- ▶ Ideal substrate for binding cell attachment and growth factors to create more *in vivo*-like culture environments
- ▶ Compatible with cell-based luminescence-reporter gene and FLIPR calcium flux assays
- ▶ Promotes more *in vivo*-like morphology (spheroid and dome formation)

Ultra-Web Synthetic Surface Features

- ▶ Synthetic surfaces are more stable and consistent lot to lot than biological coatings
- ▶ Cells grow on the nanofiber surface, not in it, for easy harvesting
- ▶ Easy to view cells using phase contrast microscopy
- ▶ Ready to use and room temperature stable

- ▶ Animal component-free and irradiation sterilized (SAL 10⁻⁶) for extra security
- ▶ Compatible with most existing cell culture protocols and methods
- ▶ Easily coated with ECM molecules and growth factors

Imaging

Ultra-Web nanofibers provide a very thin three-dimensional surface on which the cells grow; consequently, focusing on and recognizing cells can be a little more difficult than traditional flat surfaces. Light microscopy, including phase contrast and differential interference contrast (DIC) can be used to view cells seeded on Ultra-Web surfaces.

Ultra-Web nanofibers should not interfere with imaging cells via fluorescence microscopy and has been tested successfully with Texas Red, Cy3, Cy5, FITC, and GFP filters. DAPI and Hoechst dyes demonstrate limited nonspecific binding to the Ultra-Web surface resulting in a slight increase in background staining dependent upon intensity of cell staining. Diluting DAPI or Hoechst staining solutions between 1:10 and 1:20 the recommended dose is suggested for routine staining to reduce this background staining. Testing with a cell-free control is advised.

Harvesting

Cells grown on Ultra-Web or Ultra-Web polyamine surfaces may be subcultured using standard cell dissociation techniques with trypsin, collagenase, or other enzymatic and nonenzymatic dissociation solutions or cell scraping (for recommended scrapers see Corning Cat. No. 3008 or 3010). Note: To aid cell detachment gentle pipetting or mechanical agitation by tapping the culture vessel may be used. Scraping can be used to detach any remaining cells. Gentle cell scraping results in minimal damage to the Ultra-Web surface. However, robust scraping will remove the Ultra-Web surface.



Corning Labware with Ultra-Web Surfaces Product Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
3870XX1	Corning 100 mm Dish with Ultra-Web Synthetic Surface, Sterile	5	20
3871XX1	Corning 100 mm Dish with Ultra-Web Synthetic Polyamine Surface, Sterile	5	20
3872XX1	Corning 96 Well Plate with Ultra-Web Synthetic Surface, Black/Clear Bottom, Sterile with Lid	5	10
3873XX1	Corning 96 Well Plate with Ultra-Web Synthetic Polyamine Surface, Black/Clear Bottom, Sterile with Lid	5	10


These products are covered by one or more of the following patents: U.S. Patent numbers 5002582, 5512329, 6121027, 6514734, 6924028, and 6955775. Additional U.S. and international patents pending.

Cell Culture Flasks

Corning® and Costar® flasks are available in a variety of sizes, designs and cap styles to meet your needs.

- ▶ **Corning CellBIND® Surface** is a novel cell culture treatment that increases surface wettability for more even and consistent cell attachment
- ▶ **Ultra-Low Attachment** flasks feature a covalently bound hydrogel layer that minimizes cell attachment, protein absorption and cellular activation
- ▶ Manufactured from optically clear virgin polystyrene
- ▶ Treated for optimal cell attachment
- ▶ Printed with lot numbers for ease in traceability
- ▶ 100% integrity tested
- ▶ Sterilized by gamma irradiation
- ▶ Certified nonpyrogenic

Flask Cap Styles



Plug seal caps feature one-piece linerless construction and are designed for use in closed systems, providing a liquid- and gas-tight seal. When loosened, this cap can also be used in open systems. This cap design was a Corning innovation that first appeared in 1974.



Phenolic style caps are designed (when loosened) for use in open systems requiring gas exchange. With the caps slightly loosened, gas is exchanged between the environments inside and outside of the flask.



Vent caps contain a 0.2 µm pore nonwetable membrane sealed to the cap, providing consistent, sterile gas exchange while minimizing the risk of contamination. These caps are highly recommended for use in all CO₂ incubators, especially for long-term use. The vent cap was a Corning innovation that first appeared in 1988.



Septum caps maintain a closed sterile environment within the RoboFlask™ vessel. The septum allows for adding or removing cells and solutions with a blunt tip cannula while reducing the opportunity for contamination. The septum is presplit to prevent coring of the septum by the cannula. The cap may also be removed to allow pipet access (up to 5 mL) or assist in harvesting of cells. This cap septum is validated for multiple entries.

Flask Neck Styles



Straight neck flasks are ideal for larger medium volumes since this design reduces medium sloshing into the cap.




Canted neck flasks allow easier pouring and improved access to the flask for pipetting or scraping. The canted neck design was a Corning innovation that first appeared in 1974.



Angled neck improves pipet access and reduces medium sloshing into the neck. This patented design was a Corning innovation that first appeared in 1988.

Flask Shapes


Choosing a flask shape is usually a matter of personal preference:



Low Profile flasks have reduced height for incubator space savings. The corner neck gives direct access to the flask corner.




Triangular and modified triangular flasks offer good pipet and cell scraper access to the corners. The wider base provides added stability.



Rectangular flasks have a ramp from the bottom to the canted neck for easier pouring and pipet access. Most canted neck flasks also have an antitip skirt to enhance stability.



Angled neck and traditional straight neck flasks utilize the entire bottom area for cell growth. Their design saves on space and reduces medium sloshing into the neck.



RoboFlask™ vessels are robotics compatible cell culture flasks offering 92.6 cm² cell growth surface area. The flasks are designed for use in automated cell culture systems utilizing a microplate-size format.



3056 25 cm² Triangular Flask with Vent Cap

Corning® Cell Culture Flask Ordering Information

25 cm² Growth Area Flasks

Cat. No.	Surface	Flask Style	Neck Style	Cap Style	Qty/Pk	Qty/Cs
430168	TC	Rectangular	Canted	Plug Seal	20	500
430372	TC	Rectangular	Canted	Phenolic-Style	20	500
430639	TC	Rectangular	Canted	Vent Cap	20	200
3055	TC	Triangular	Angled	Phenolic-Style	20	500
3056	TC	Triangular	Angled	Vent Cap	10	200
3289	Corning® CellBIND® Surface	Rectangular	Canted	Vent Cap	20	200
3815	Ultra-Low Attachment	Rectangular	Canted	Vent Cap	6	24

75 cm² Growth Area Flasks

Cat. No.	Surface	Flask Style	Neck Style	Cap Style	Qty/Pk	Qty/Cs
430641	TC	Rectangular	Canted	Vent Cap	5	100
430720	TC	Rectangular	Canted	Plug Seal	5	100
430725	TC	Rectangular	Canted	Phenolic-Style	5	100
3275	TC	Modified triangular	Straight	Phenolic-Style	5	100
3276	TC	Modified triangular	Straight	Vent Cap	5	100
3290	Corning CellBIND Surface	Rectangular	Canted	Vent Cap	5	100
3814	Ultra-Low Attachment	Rectangular	Canted	Vent Cap	4	24



430639 25 cm² Canted Neck Flask with Vent Cap

Cell Culture Tip

Check the Corning web site (www.corning.com/lifescience) for technical cell culture application bulletins.



430641 75 cm² Canted Neck Flask with Vent Cap




430725 75 cm² Canted Neck Flask with Phenolic Style Cap



3275 75 cm² Triangular Flask with Phenolic Style Cap

Cell Culture Flask Selection Tip


Corning's unique RoboFlask™ (U.S. Patent 7,078,228) Cell Culture Vessels are designed to an SBS standard microplate footprint for use in automated cell maintenance and assay systems.



3070 RoboFlask® Cell Culture Vessel with Septum Cap



430823 150 cm² Canted Neck Flask with Plug Seal Cap



431306 175 cm² Flask with Vent Cap and Bar Code

New!



3073 Low Profile Flask

Cell Culture Flask Selection Tip

The Low Profile 100 cm² flask:

- gives 33% more area in the footprint of a T75 flask
- has a 1/2 turn easy-opening cap
- saves 33% in incubator space
- uses 26% less plastic than a T75 flask

Cell Culture Flask Application Tip

Corning recommends 0.2 to 0.3 mL of medium per cm² of growth area.

92.6 cm² Growth Area RoboFlask™ Vessels

Cat. No.	Description	Qty/Pk	Qty/Cs
3070	RoboFlask Cell Culture Vessel for automation, tissue culture treated, with bar code, septum cap, sterile	20	100
3071	RoboFlask Cell Culture Vessel for manual use, tissue culture treated, with bar code, flat cap (without septum), sterile	20	100
3069	RoboFlask Cell Culture Vessel for automation, tissue culture treated, with bar code, septum cap, sterile	10	50
3059	RoboFlask Cell Culture Vessel for manual use, tissue culture treated, with bar code, flat cap (without septum), sterile	10	50
3067	RoboFlask Cell Culture Vessel for automation, Corning® CellBIND® surface treatment with bar code, septum cap, sterile	20	100
3068	RoboFlask Cell Culture Vessel for automation, Corning CellBIND surface treatment with bar code, septum cap, sterile	10	50

100 cm² Growth Area Low Profile Flask

Cat. No.	Description	Qty/Pk	Qty/Cs
3073	Low Profile Flask, 100 cm², Corning CellBIND Surface with Vent Cap, Sterile	6	60
3816	Low Profile Flask, 100 cm², tissue culture surface with Vent Cap, Sterile	6	60

150 cm² Growth Area Flasks

Cat. No.	Surface	Flask Style	Neck Style	Cap Style	Qty/Pk	Qty/Cs
430823	TC	Rectangular	Canted	Plug Seal	5	50
430824	TC	Rectangular	Canted	Phenolic Style	5	50
430825	TC	Rectangular	Canted	Vent Cap	5	50
3291	Corning CellBIND Surface	Rectangular	Canted	Vent Cap	5	50

162 cm² Growth Area Flasks

Cat. No.	Surface	Flask Style	Neck Style	Cap Style	Qty/Pk	Qty/Cs
3150	TC	Traditional	Straight	Phenolic Style	5	25
3151	TC	Traditional	Straight	Vent Cap	5	25



431082 225 cm² Angled Neck Flask with Vent Cap



10010 HYPERFlask Vessel

Cell Culture Flask Selection Tip

The novel HYPERFlask Vessel offers high yield and high performance with 10 growth surfaces and 1720 cm² growth area in the same footprint as the 175 cm² flask.



3001 225 cm² Canted Neck Flask with Vent Cap

175 cm² Growth Area Flasks

Cat. No.	Surface	Flask Style	Neck Style	Cap Style	Qty/Pk	Qty/Cs
431079	TC	Rectangular	Angled	Plug Seal	5	50
431080	TC	Rectangular	Angled	Vent Cap	5	50
431085	TC	Rectangular	Angled	Phenolic Style	5	50
431306*	TC	Rectangular	Angled	Vent Cap	7	84
431328*	Corning® CellBIND® Surface	Rectangular	Angled	Vent Cap	7	84
3292	Corning CellBIND Surface	Rectangular	Angled	Vent Cap	5	50
3298	Corning CellBIND Surface	Rectangular	Angled	Phenolic Style	5	50

*Flask prelabeled with bar code, validated for use with SelecT™ Robotic System.

225 cm² Growth Area Flasks

Cat. No.	Surface	Flask Style	Neck Style	Cap Style	Qty/Pk	Qty/Cs
431081	TC	Traditional	Angled	Plug Seal	5	25
431082	TC	Traditional	Angled	Vent Cap	5	25
3000	TC	Rectangular	Canted	Phenolic Style	4	24
3001	TC	Rectangular	Canted	Vent Cap	4	24
3293	Corning CellBIND Surface	Traditional	Angled	Vent Cap	5	25

235 cm² Expanded Growth Area Flask

Cat. No.	Surface	Flask Style	Neck Style	Cap Style	Qty/Pk	Qty/Cs
431346*	Corning CellBIND Surface	Rectangular	Angled	Vent Cap	7	42

*Flask prelabeled with bar code for use with SelecT Automation System with same footprint as the 175 cm² flask.

1720 cm² Growth Area HYPERFlask™ Vessel

Cat. No.	Description	Qty/Pk	Qty/Cs
10010*	HYPERFlask Vessel, 1720 cm ² , Corning CellBIND Surface, Bar Code, Sterile	4	4
10024*	HYPERFlask Vessel, 1720 cm ² , Corning CellBIND Surface, Bar Code, Sterile	4	24

*Flask prelabeled with bar code for use with SelecT Robotic System.


Cell Yields and Recommended Medium Volume

Corning and Costar® Flasks	Approximate Growth Area (cm ²)	Average Cell Yield*	Recommended Medium Volume (mL)	Maximum Working Volume (mL)†
25 cm ²	25	2.5 x 10 ⁶	5 - 7.5	10
75 cm ² Canted neck	75	7.5 x 10 ⁶	15 - 22.5	60
75 cm ² Straight neck	75	7.5 x 10 ⁶	15 - 22.5	90
RoboFlask™ Vessel	93	9.4 x 10 ⁶	20 - 30	70
100 cm ²	100	1.0 x 10 ⁷	20 - 30	40
150 cm ²	150	1.5 x 10 ⁷	30 - 45	210
162 cm ²	162	1.6 x 10 ⁷	32 - 48	175
175 cm ²	175	1.75 x 10 ⁷	35 - 52.5	250
225 cm ²	225	2.25 x 10 ⁷	45 - 67.5	370
235 cm ²	235	2.35 x 10 ⁷	47 - 70.5	250
1720 cm ²	1720	2.5 x 10 ⁸	565	565

*Assumes an average yield of 1 x 10⁵ cells/cm² from a 100% confluent culture. Yields from many cell types can be lower than this.

†Maximum working volume is the amount a flask can hold in the horizontal position when filled to the neck.

Cell Culture Dishes



3296 Corning CellBIND Surface 100 mm Dishes



430196 Gridded 60 mm Dish



431110 500 cm² Cell Culture Dish



3870XX1 100 mm Tissue Culture Dish with Ultra-Web Synthetic Surface

Cell Culture Tip

Check the Corning web site (www.corning.com/lifesciences) for technical cell culture application bulletins.

Corning Cell Culture Treated Dishes

- Corning® CellBIND® Surface is a novel cell culture treatment that increases surface wettability for more even and consistent cell attachment
- Ultra-Low Attachment dishes feature a covalently bound hydrogel layer that minimizes cell attachment, protein absorption and cellular activation
- Ultra-Web™ Synthetic Surface is composed of electrospun polyamide nanofibers, creating a culturing substrate that mimics structural components within the basement membrane or extracellular matrix.
- 6-pack carriers with only 10 dishes/bag are available for 100 mm dishes (Cat No. 430293)
- 245 mm square dishes offer 500cm² growth surface
- Manufactured from optically clear virgin polystyrene
- Sterilized by gamma radiation
- Certified nonpyrogenic
- Have stacking beads to aid in handling
- Supplied with vents to provide consistent gas exchange

Corning® Cell Culture Dish Ordering Information

Cat. No.	Surface	Dish Style* (mm)	Approx. Height (mm)	Growth Area (cm²)	Qty/Pk	Qty/Cs
3294	Corning CellBIND Surface	35	10	8	10	210
430165	TC	35	10	8	20	500
430166	TC	60	15	21	20	500
3295	Corning CellBIND Surface	60	15	21	7	126
3261	Ultra-Low Attachment	60	15	21	5	20
3262	Ultra-Low Attachment	100	20	55	5	20
430196	TC	60 with 2 mm grid	15	21	20	500
3296	Corning CellBIND Surface	100	20	55	5	40
3870XX1	Ultra-Web Synthetic Surface	100	20	55	5	20
3871XX1	Ultra-Web Synthetic Polyamine Surface	100	20	55	5	20
430167	TC	100	20	55	20	500
430293	TC	100	20	55	10	480
430599	TC	150	25	148	5	60
431110	TC	245	25	500	4	16

*Dish style (mm) = actual growth surface diameters: 35 mm dish = 34.4 mm; 60 mm dish = 52.1 mm; 100 mm dish = 83.8 mm; 150 mm dish = 139.1 mm. The square dishes have interior bottom dimensions of 224 mm x 224 mm.

Cell Culture Dish Application Tips

- ▶ The 150 and 245 mm culture dishes make excellent carriers and incubator trays for 35 and 60 mm dishes. This helps prevent spills and reduces opportunities for contamination.
- ▶ Corning recommends 0.2 to 0.3 mL of medium per cm² of growth area.

Corning Nontreated Cell Culture Dishes

- ▶ Manufactured from optically clear virgin polystyrene
- ▶ Not cell culture treated for applications where cell attachment is not desired
- ▶ Have stacking beads to aid in handling
- ▶ Supplied with vents to provide consistent gas exchange
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic

Corning Nontreated Cell Culture Dish Ordering Information

Cat. No.	Dish Style* (mm)	Height (mm)	Approx. Growth Area (cm ²)	Qty/Pk	Qty/Cs
430588	35	10	9	20	500
430589	60	15	21	20	500
430591	100	20	55	20	500
430597	150	25	152	5	60
431111†	245	25	500	4	16

*Note: Dish style (mm) = actual growth surface diameters: 35 mm dish = 34.4 mm; 60 mm dish = 52.1 mm; 100 mm dish = 83.8 mm; 150 mm dish = 139.1 mm.

†Cat. No. 431111 is a square dish with interior bottom plate dimensions of 224 mm x 224 mm.

Expected Cell Yields and Recommended Medium Volumes

Corning Dishes	Approximate Growth Area (cm ²)	Average Cell Yield*	Recommended Medium Volume (mL)†
35 mm	9	9.0×10^5	1.8 - 2.7
60 mm	21	2.1×10^6	4.2 - 6.3
100 mm	55	5.5×10^6	11 - 16.5
150 mm	152	1.52×10^7	30.4 - 45.6
245 mm (square)	500	5.0×10^7	100 - 150

*Assumes an average yield of 1×10^5 cells/cm² from a 100% confluent culture.


†Yields from many cell types can be lower than this.



3260 IVF Culture Dish

Costar® IVF Culture Dish


- ▶ 20 mm center well
- ▶ Inner well holds 3 mL of medium while the outer well holds 10 mL
- ▶ Treated for optimal cell attachment
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic
- ▶ For research use only



Costar IVF Culture Dish Ordering Information

Cat. No.	Size (mm)	Description (mm)	Center Well (mm)	Qty/Pk	Qty/Cs
3260	60	60 x 15	20	20	500


Multiple Well Plates




3516 6 Well Culture Plate



3513 12 Well Culture Plate



3524 24 Well Culture Plate



3548 48 Well Culture Plate

Costar® 6, 12, 24, and 48 Well Cell Culture Plates

- Corning® CellBIND® Surface is a novel cell culture treatment that increases surface wettability for more even and consistent cell attachment
- Ultra-Low Attachment plates feature a covalently bound hydrogel layer that minimizes cell attachment, protein absorption and cellular activation
- Nonreversible lids with condensation rings to reduce contamination, uniform footprint for ease in stacking
- Individual alphanumerical codes for well identification, flat bottoms
- Treated for optimal cell attachment (except where noted)
- Sterilized by gamma irradiation, certified nonpyrogenic

6, 12, 24, and 48 Well Plates Ordering Information


Cat. No.	Surface	Plate Type	Qty/Pk	Qty/Cs
6 Well Plates				
3335	Corning® CellBIND® Surface	Standard clear plate	5	50
3506	TC	Standard clear plate	5	100
3516	TC	Standard clear plate	1	50
3471	Ultra-Low Attachment	Standard clear plate with hydrogel*	1	24
12 Well Plates				
3336	Corning CellBIND Surface	Standard clear plate	5	50
3512	TC	Standard clear plate	5	100
3513	TC	Standard clear plate	1	50
24 Well Plates				
3337	Corning CellBIND Surface	Standard clear plate	5	50
3524	TC	Standard clear plate	1	100
3526	TC	Standard clear plate	1	50
3527	TC	Standard clear plate	5	100
3473	Ultra-Low Attachment	Standard plate with hydrogel*	1	24
48 Well Plates				
3338	Corning CellBIND Surface	Standard clear plate	5	50
3548	TC	Standard clear plate	1	100

*This covalently bonded hydrogel surface minimizes cell attachment, protein absorption, enzyme activation and cellular activation. The surface is noncytotoxic, biologically inert and nondegradable.


Well Dimensions, Expected Cell Yields, and Recommended Medium Volumes

Cell Culture Plates	Well Diameter (Bottom, mm)	Single Well Only				Entire Plate		
		Approx. Growth Area (cm ²)	Average Cell Yield*	Total Well Volume (mL)	Working Volume (mL)	Approx. Growth Area (cm ²)	Average Cell Yield*	Working Volume (mL)
6 well	34.8	9.5	9.5×10^5	16.8	1.9 - 2.9	57	5.7×10^6	11.4 - 17.1
12 well	22.1	3.8	3.8×10^5	6.9	0.760 - 1.14	45.6	4.56×10^6	9.1 - 13.7
24 well	15.6	1.9	1.9×10^5	3.4	0.380 - 0.570	45.6	4.56×10^6	9.1 - 13.7
48 well	11	0.95	9.5×10^4	1.6	0.19 - 0.285	45.6	38.4×10^6	9.1 - 13.7

*Assumes an average yield of 1×10^5 cells/cm² from a 100% confluent culture. Yields from many cell types can be lower than this.



3596 96 Well Culture Plate

3610 and 3603
96 Well Clear Bottom Plates3917 and 3916
96 Well Solid Plates3872XX1
96 Well Plate with Ultra-Web
Synthetic Surface

Corning® and Costar® 96 Well Cell Culture Plates


- ▶ **Corning CellBIND® Surface** is a novel cell culture treatment that increases surface wettability for more even and consistent cell attachment
- ▶ **Ultra-Low Attachment** plates feature a covalently bound hydrogel layer that minimizes cell attachment, protein absorption and cellular activation
- ▶ **Ultra-Web™ Synthetic Surface** is composed of electrospun polyamide nanofibers, creating a culturing substrate that mimics structural components within the basement membrane or extracellular matrix.
- ▶ **Corning Poly-D-Lysine (PDL)** microplates are coated with PDL (molecular weight range of 70 to 150 kDa) giving the surface a net positive charge for better cell attachment.
- ▶ Nonreversible lids with condensation rings to reduce contamination (except where noted)
- ▶ Treated for optimal cell attachment (except where noted)
- ▶ Sterilized by gamma radiation, certified nonpyrogenic
- ▶ Individual alphanumeric codes for well identification, flat bottoms (except where noted)

Black plates are designed to lower background in fluorescent assays and reduce crosstalk. White plates are designed for luminescent assays. Some plates have the Corning CellBIND surface or a poly-D-lysine coating to enhance cell attachment. Corning offers many other 96 well plate types for applications other than cell culture; for a complete listing, check the catalog at www.corning.com/lifesciences.


Corning® Assay Surface Properties and Applications

Corning Surface	Applications	Binding Interaction	Sample Properties
Standard Tissue Culture Surface	<ul style="list-style-type: none"> ▶ Assays using standard attachment dependent cell lines 	Hydrophilic and ionic interactions (negatively charged)	Allows cell attachment and binding
Corning CellBIND Surface	<ul style="list-style-type: none"> ▶ Assays for difficult to attach cells ▶ Help cells stay attached during washing steps 	Hydrophilic and ionic interactions (negatively charged)	Enhances cell attachment uniformity and binding to polystyrene
Poly-D-Lysine-Coated Surface	<ul style="list-style-type: none"> ▶ Assays for difficult to attach cells ▶ Help cells stay attached during washing steps 	Hydrophilic and ionic interactions (positively charged)	Enhances cell attachment and binding
Ultra-Low Attachment Surface	<ul style="list-style-type: none"> ▶ Assays where preventing cell attachment is required ▶ Hybridoma production and clonal isolation by limiting dilution 	Nonionic hydrogel layer reduces or eliminates ionic and hydrophobic binding	Prevents or reduces cell attachment and binding
Ultra-Web Surface	<ul style="list-style-type: none"> ▶ Assays where cell attachment or performance is enhanced by using a 3D surface 	Hydrophilic interactions on a 3D surface	Enhances cell attachment and performance
Ultra-Web Polyamine Surface	<ul style="list-style-type: none"> ▶ Assays where cell attachment or performance is enhanced by using a 3D surface 	Hydrophilic and ionic interactions (positively charged) on a 3D surface	Enhances cell attachment and performance


Well Geometry




Flat Bottom



Round Bottom



V Bottom



Half Area

Cell Culture Tip

Check the Corning web site (www.corning.com/lifesciences) for technical cell culture application bulletins.

96 Well Plate Ordering Information

Cat. No.	Surface	Description	Qty/ Pk	Qty/ Cs
Clear Plates				
3300	Corning® CellBIND® Surface	Standard clear plate	5	50
3596	TC	Standard clear plate	1	50
3997	TC	Standard clear plate	10	50
3598	TC	Standard clear plate	5	100
3599	TC	Standard clear plate	1	100
3585	TC	Standard clear plate with special low evaporation lid	5	50
3595	TC	Standard clear plate with special low evaporation lid	1	50
3696	TC	96 well half area, flat bottom clear plate	1	50
3697	TC	96 well half area clear plate	20	100
3790	Not Treated	96 well round bottom, polypropylene plate with polystyrene lid	1	50
3799	TC	96 well round bottom clear plate	1	50
3894	TC	96 well V-bottom clear plate	1	50
3665	Poly-D-lysine	Standard clear plate, coated	20	100
9102	TC	8-well strip plate, assembled 12 strips per plate	1	50
3474	Ultra-Low Attachment	Standard clear plate with hydrogel*	1	24
7007	Ultra-Low Attachment	96 well round bottom plate with hydrogel*	1	24
White Plates				
3917	TC	Solid white plate	20	100
3362	TC	Solid white plate without lid	25	100
3688	TC	96 well half area solid white plate	20	100
3885	TC	96 well half area white plate with clear bottom	20	100
3610	TC	White plate with clear bottom	1	48
3903	TC	White plate with clear bottom	20	100
3666	Poly-D-lysine	White plate with clear bottom	20	100
Black Plates				
3340	Corning CellBIND Surface	Black plate with clear bottom	5	50
3872XX1	Ultra-Web™ Synthetic Surface	Black plate with clear bottom, with lid	5	10
3873XX1	Ultra-Web Synthetic Polyamine Surface	Black plate with clear bottom, with lid	5	10
3916	TC	Solid black plate	20	100
3875	TC	96 well half area solid black plate	20	100
3882	TC	96 well half area black plate with clear bottom	20	100
3603	TC	Black plate with clear bottom	1	48
3904	TC	Black plate with clear bottom	20	100
3667	Poly-D-lysine	Black plate with clear bottom	25	100
3614	TC	Black plate with special optics, ultrathin, clear bottom, without lid	20	100
Lids and Tape				
3099	—	Universal lid, sterile	25	50
3345	—	Breathable Sealing tape, sterile	50	500
3930	—	Rigid styrene lid with condensation rings, sterile	1	100
3931	—	Rigid styrene lid with condensation rings, sterile	25	50

*This covalently bonded hydrogel surface minimizes cell attachment, protein absorption, enzyme activation and cellular activation. The surface is noncytotoxic, biologically inert and nondegradable.

96 Well Cell Culture Plates

Well Dimensions, Expected Cell Yields, and Recommended Medium Volume

Cell Culture Plates	Well Diameter (Bottom, mm)	Single Well Only				Entire Plate		
		Approx. Growth Area (cm ²)	Average Cell Yield*	Total Well Volume (mL)	Working Volume (mL)	Approx. Growth Area (cm ²)	Average Cell Yield*	Working Volume (mL)
96 well flat bottom	6.4	0.32	3.2 x 10 ⁴	0.36	0.100 - 0.200	30.7	3.07 x 10 ⁶	9.6 - 19.2
96 well round bottom	6.4	NA [†]	NA [†]	0.33	0.100 - 0.200	NA [†]	NA [†]	9.6 - 19.2
96 well V bottom	6.4	0.38	3.8 x 10 ⁴	0.29	0.100 - 0.200	36.5	3.65 x 10 ⁶	9.6 - 19.2
96 half area	4.5	0.16	1.6 x 10 ⁴	0.19	0.050 - 0.100	15.4	1.54 x 10 ⁶	4.8 - 9.6

*Assumes an average yield of 1 x 10⁵ cells/cm² from a 100% confluent culture. Yields from many cell types can be lower than this.

[†]Because these wells are round, the surface area available for cell attachment is dependent on the medium volume used.


384 Well Cell Culture Plates

- Flat bottoms and lids
- New low volume plates have only a 50 µL total well volume, with recommended working volume of 5 to 40 µL
- Treated for optimal cell attachment
- Sterilized by gamma radiation
- Certified nonpyrogenic

Black plates are designed to lower background in fluorescent assays and reduce crosstalk. White plates are designed for luminescent assays. Some plates have the Corning® CellBIND® Surface or a poly-D-lysine coating to enhance cell attachment. Corning offers many other 384 well plate types for applications other than cell culture; for a complete listing, check the catalog at www.corning.com/lifesciences.

384 Well Cell Culture Plate Ordering Information

Cat. No.	Description	Surface	Qty/Pk	Qty/Cs
Clear Plates				
3701	Standard clear plate, low flange	TC	20	100
3662	Standard clear plate	Poly-D-lysine	25	100
White Plates				
3570	Solid white plate	TC	10	50
3707	White plate with clear bottom	TC	20	100
3663	White plate with clear bottom	Poly-D-lysine	25	100
3826	Solid white plate, low volume	TC	10	50
Black Plates				
3571	Solid black plate, low flange	TC	10	50
3712	Black plate with clear bottom	TC	20	100
3664	Black plate with clear bottom	Poly-D-lysine	25	100
3683	Black plate with clear bottom	Corning CellBIND Surface	10	50
3542	Low volume, black plate with clear bottom	TC	10	50
3822	Low volume, solid black plate	TC	10	50
3985	Black optical imaging plate with clear bottom	TC	20	100



3707 and 3712
384 Well Clear Bottom Plates

Well Dimensions, Expected Cell Yields, and Recommended Medium Volumes

Cell Culture Plates	Well Diameter (Bottom, mm)	Single Well Only				Entire Plate		
		Approx. Growth Area (cm ²)	Average Cell Yield*	Total Well Volume (mL)	Working Volume (mL)	Approx. Growth Area (cm ²)	Average Cell Yield*	Working Volume (mL)
Standard 384 Well	2.7 x 2.7 [†]	0.056	5.6 x 10 ³	0.125	0.025 - 0.050	21.5	2.15 x 10 ⁶	9.6 - 19.2
Low Volume 384 Well	2.0	0.031	3.1 x 10 ³	0.050	0.005 - 0.040	12.0	1.2 x 10 ⁶	1.9 - 15.3

*Assumes an average yield of 1 x 10⁵ cells/cm² from a 100% confluent culture. Yields from many cell types can be lower than this.

[†]These wells are square.


1536 Well Cell Culture Plates

- Superior performance compared to competitor plates: lower CVs, higher signal-to-noise ratios, and lower background fluorescence
- Compatible with bar coding, standard readers and automation
- Recommended working volume of up to 8 µL
- Treated for optimal cell attachment
- Flat bottoms and lids
- Sterilized by gamma radiation
- Certified nonpyrogenic

Black plates are designed to lower background in fluorescent assays and reduce crosstalk. White plates are designed for luminescent assays. Corning offers other 1536 well plate types for applications other than cell culture; for a complete listing, check the catalog at www.corning.com/lifesciences.

1536 Well Cell Culture Plate Ordering Information

Cat. No.	Description	Surface	Qty/Pk	Qty/Cs
<i>Clear Plates</i>				
3853	Standard clear plate	TC	20	100
<i>White Plates</i>				
3727	Solid white plate	TC	10	50
3855	Solid white plate, low volume	TC	20	100
<i>Black Plates</i>				
3726	Solid black plate	TC	10	50
3893	Black clear bottom plate	TC	10	50
3854	Solid black plate, low volume	TC	20	100




3893 1536 Well Culture Plates

Well dimensions, Expected Cell Yields, and Recommended Medium Volumes

Cell Culture Plates	Well Diameter (Bottom, mm)	Single Well Only				Entire Plate		
		Approx. Growth Area (cm ²)	Average Cell Yield*	Total Well Volume (µL)	Working Volume (µL)	Approx. Growth Area (cm ²)	Average Cell Yield*	Working Volume (mL)
1536 Well Clear Flat Bottom	1.63 x 1.63	0.025	2.5 x 10 ³	12.5	5 - 8	38.3	3.8 x 10 ⁶	7.7 - 15.4
1536 Well Solid Flat Bottom	1.53 x 1.53	0.023	2.3 x 10 ³	12.5	5 - 8	35.3	3.5 x 10 ⁶	7.7 - 15.4

*Assumes an average yield of 1 x 10⁵ cells/cm² from a 100% confluent culture. Yields from many cell types can be lower than this.

Transwell® Permeable Supports



Transwell cell culture inserts are convenient, easy-to-use permeable support devices for the study of both anchorage-dependent and anchorage-independent cell lines


- ▶ Designed to produce a cell culture environment that closely resembles the *in vivo* state
- ▶ Allows polarized cells to feed basolaterally and thereby carry out metabolic activities in a more natural fashion
- ▶ Unique patented self-centered hanging design prevents medium wicking between the insert and outer well
- ▶ Permits access to the lower compartment through windows in the insert wall
- ▶ Suspended design allows for undamaged co-culturing of cells in the lower compartment
- ▶ Available in a range of pore sizes and different membranes to satisfy diverse experimental requirements

Characteristics of Transwell Membranes

Characteristics	Polyester (PET)	Polycarbonate	PTFE/Collagen
Optical properties	Clear	Translucent	Clear when wet
Cell visibility	Good	Poor	Cell outlines
Tissue culture treated	Yes	Yes	No
Membrane thickness	10 µm	10 µm	50 µm
Matrix/ECM coatable	Yes	Yes	Yes
Collagen treated	No	No	Yes
Available Pore Sizes (µm)	0.4, 1.0, 3.0, 8.0	0.4, 3.0, 5.0, 8.0	0.4, 3.0

Chemical Compatibility

All of the Transwell membranes are compatible with histological fixatives including methanol and formaldehyde. The polyester Transwell membranes have the best overall chemical resistance. These membranes (but not the polystyrene housings) are compatible with many alcohols, amines, esters, ethers, ketones, oils and some solvents, including many halogenated hydrocarbons and DMSO but are not recommended for use with strong acids and bases.



Transwell Permeable Supports Tip

Check the Corning web site (www.corning.com/lifesciences) for an extensive list of references, listed by application, citing the use of Transwell permeable supports in cell culture research.

Pore Density

Of the three types of Transwell membranes, only the PTFE does not have a defined pore density because it is a tortuous path membrane. The two membranes with a nominally defined pore density are polycarbonate and polyester. The polyester Transwell membranes do not have as high a pore density as the polycarbonate Transwell but have better optical clarity as a result. The nominal pore densities for Corning® Polycarbonate and Polyester (PET) membranes are given in the following table.

Packaging of Transwell Polycarbonate and Polyester Membrane Inserts

The cell culture inserts come prepackaged in the appropriate multiple well plate as follows:

- 24 mm diameter inserts are packaged 6 inserts in a 6 well plate, four 6 well plates per case, for a total of 24 inserts.
- 12 mm diameter inserts are packaged 12 inserts in a 12 well plate, four 12 well plates per case, for a total of 48 inserts.
- 16.5 mm diameter inserts are packaged 12 inserts in a 24 well plate, four 24 well plates per case, for a total of 48 inserts.

All Transwell-COL collagen coated inserts are individually packaged and each case includes the appropriate multiple well plate.



3401 12 mm Polycarbonate Transwell Insert



3419 75mm Polycarbonate Transwell Insert

Nominal Pore Densities for Transwell Polyester and Polycarbonate Membranes

Pore Size	Nominal Pore Density	
	Polycarbonate Membrane Transwell (pores/cm ²)	Transwell-Clear Polyester Membrane (pores/cm ²)
0.4 µm	1 x 10 ⁸	4 x 10 ⁶
1.0 µm	n/a	1.6 x 10 ⁶
3.0 µm	2 x 10 ⁶	2 x 10 ⁶
5.0 µm	4 x 10 ⁵	n/a
8.0 µm	1 x 10 ⁵	1 x 10 ⁵

Growth Areas and Recommended Medium Volumes for Transwell Permeable Supports

Multiple Well Plate or Dish Type	Transwell Insert Diameter (mm)	Insert Membrane Growth Area (cm ²)	Volume Added per Plate Well	Volume Added to Inside of Transwell Insert (mL)
HTS 96	4.26	0.143	0.235	0.075
HTS 24	6.5	0.33	0.6	0.1
24 well	6.5	0.33	0.6	0.1
12 well	12	1.12	1.5	0.5
6 well	24	4.67	2.6	1.5
100 mm dish	75	44	13	9


Transwell® Polycarbonate Membrane Insert

- 10 µm thick translucent membrane
- Pore sizes ranging from 0.4 µm to 8 µm diameters
- Treated for optimal cell attachment
- Supplied in multiple well plates
- Membrane must be stained for cell visibility
- Sterilized by gamma radiation

Transwell Polycarbonate Membrane Permeable Support Ordering Information

Cat. No.	Membrane Diameter (mm)	Growth Surface Area (cm ²)	Membrane Pore Size (µm)	Tissue Culture Treated	Inner Packaging*	Inserts/Cs
3413	6.5	0.33	0.4	Yes	12/plate*	48
3415	6.5	0.33	3.0	Yes	12/plate*	48
3421	6.5	0.33	5.0	Yes	12/plate*	48
3422	6.5	0.33	8.0	Yes	12/plate*	48
3401	12	1.12	0.4	Yes	12/plate	48
3402	12	1.12	3.0	Yes	12/plate	48
3412	24	4.67	0.4	Yes	6/plate	24
3414	24	4.67	3.0	Yes	6/plate	24
3428	24	4.67	8.0	Yes	6/plate	24
3419	75	44	0.4	Yes	1/dish	12
3420	75	44	3.0	Yes	1/dish	12

*6.5 mm membrane diameter are packaged 12 inserts in a 24 well plate, 4 plates per case.



3458 6.5 mm Polycarbonate Transwell Cultrex Insert

Corning® Transwell® Invasion Inserts

Transwell Polycarbonate Membrane Insert Coated with Cultrex® Basement Membrane Extract

- ▶ 8 µm pore size membrane, 10 µm thick
- ▶ Coated with Cultrex basement membrane extract
- ▶ For cell invasion assays
- ▶ Supplied in 24-well plates, 12 inserts per plate

Transwell Cultrex Polycarbonate Membrane Permeable Support Ordering Information

Cat. No.	Membrane Diameter (mm)	Growth Surface Area (cm ²)	Membrane Pore Size (µm)	Tissue Culture Treated	Inner Packaging	Inserts/Cs
3458	6.5	0.33	0.8	Yes	12/plate	24

Transwell-Clear Polyester Membrane Insert

- ▶ 10 µm transparent membrane
- ▶ Treated for optimal cell attachment
- ▶ Excellent visibility under phase contrast microscopy
- ▶ Supplied in multiple well plates
- ▶ Sterilized by gamma radiation

Transwell-Clear Insert Ordering Information

Cat. No.	Membrane Diameter (mm)	Growth Surface Area (cm ²)	Membrane Pore Size (µm)	Inner Packaging*	Inserts/Cs
3450	24	4.67	0.4	6/plate	24
3452	24	4.67	3.0	6/plate	24
3460	12	1.12	0.4	12/plate	48
3462	12	1.12	3.0	12/plate	48
3470	6.5	0.33	0.4	12/plate*	48
3472	6.5	0.33	3.0	12/plate*	48

*6.5 mm membrane diameter are packaged 12 inserts in a 24 well plate, 4 plates per case.


Transwell-COL Collagen-Coated Membrane Insert

- ▶ Transparent collagen treated PTFE membrane
- ▶ Promotes cell attachment and spreading
- ▶ Equimolar mixture of types I and III collagen
- ▶ Individually packaged
- ▶ Multiple well plates included in each case
- ▶ Supplied sterile

Transwell-COL Insert Ordering Information

Cat. No.	Membrane Diameter (mm)	Growth Surface Area (cm ²)	Membrane Pore Size (µm)	Inner Packaging	Multiple Well Plate	Inserts/Cs
3491	24	4.67	0.4	Individual	6 well	24
3492	24	4.67	3.0	Individual	6 well	24
3493	12	1.12	0.4	Individual	12 well	24
3494	12	1.12	3.0	Individual	12 well	24
3495*	6.5	0.33	0.4	Individual	24 well	24
3496*	6.5	0.33	3.0	Individual	24 well	24

*Includes twenty-four 6.5 mm inserts packaged separately with two 24 well plates.



3450 24 mm Transwell-Clear Insert



3491 24 mm Transwell-COL Collagen-Coated Insert



3407 12 mm Snapwell Inserts

Snapwell™ Inserts


- ▶ A modified Transwell permeable support containing a 12 mm diameter membrane supported by a detachable ring
- ▶ Once cells are grown to confluence on the Snapwell insert, the ring can be placed in a vertical or horizontal diffusion chamber*
- ▶ Sterilized by gamma radiation
- ▶ Packaged in 6 well plates

Snapwell Insert Ordering Information

Cat. No.	Membrane Pore Size (µm)	Membrane	Inner Packaging	Inserts/Cs
3407	0.4	Polycarbonate	6/plate	24
3802	3.0	Polycarbonate	6/plate	24
3801	0.4	Clear Polyester	6/plate	24

*Diffusion Chambers are available through Harvard Apparatus (www.harvardapparatus.com)

Corning® HTS Transwell®-24 Well Permeable Supports




HTS Transwell-24 Well Permeable Support

- ▶ The HTS Transwell-24 Well Permeable Support has an array of 24 wells with membrane inserts connected by a rigid, robotics-friendly tray that enables all 24 Transwell supports to be handled as a single unit
- ▶ Cell growth area is 0.33 cm²/well
- ▶ Choice of either polyester (PET) membrane (0.4 µm pore size) or polycarbonate (PC) membrane (0.4 µm, 3.0 µm pore sizes)
- ▶ Treated for optimal cell attachment
- ▶ Individual pack has 2 HTS Transwell-24 units loaded into two open reservoir trays and two individually wrapped 24 well plates
- ▶ Bulk pack has 12 HTS Transwell-24 units loaded into 24 well plates only. Reservoirs may be purchased separately
- ▶ Sterilized by gamma radiation

HTS Transwell-24 Well Permeable Supports Ordering Information

Cat. No.	Description	Membrane Pore Size (µm)	Membrane	Qty/Pk	Plates/Cs
3396	HTS Transwell-24, individual	0.4	PC	1	2
3397	HTS Transwell-24, bulk	0.4	PC	12	12
3398	HTS Transwell-24, individual	3.0	PC	1	2
3399	HTS Transwell-24, bulk	3.0	PC	12	12
3395	HTS Transwell nontreated reservoir	—	—	12	48
3378	HTS Transwell-24, bulk	0.4	PET	12	12
3379	HTS Transwell-24, individual	0.4	PET	1	2



HTS Transwell-96 System

Corning® HTS Transwell®-96 Well Permeable Support Systems and Plates

- ▶ The HTS Transwell-96 Well Permeable Support has an array of 96 wells with membrane inserts connected by a rigid, robotics-friendly tray that enables all 96 inserts to be handled as a single unit
- ▶ Choice of either polyester (PET) membrane (1.0 μm , 8.0 μm pore sizes) or polycarbonate (PC) membrane (0.4 μm , 3.0 μm , 5.0 μm pore sizes)
- ▶ 0.143 cm^2 membrane area per well, providing 20 to 50% more surface area for cell growth than other commercially available systems
- ▶ Large apical and basolateral access ports allow efficient media sampling and facilitate automated or manual access
- ▶ Optimized for automation, with multichannel feeder ports, improved gripping surface, and standard bar codes
- ▶ The reservoir plate allows for simultaneous feeding of 96 wells and comes with a removable media stabilizer to reduce the risk of spills during handling
- ▶ The receiver plate isolates each well to enable 96 individual assays
- ▶ Sterilized by gamma radiation
- ▶ The HTS Transwell-96 Systems (0.4 μm PC and 1.0 μm PET) are packaged with the 96 well insert plate in a reservoir plate and includes the 96 well receiver plate with lid.
- ▶ The HTS Transwell-96 Well Plates (3.0 and 5.0 μm PC, 8.0 μm PET) are packaged with the 96 well insert plate in the 96 well receiver plate with lid. Reservoir plates may be purchased separately.

HTS Transwell®-96 Well Permeable Supports Ordering Information

Cat. No.	Description	Membrane Pore Size (μm)		Qty/ Pk	Qty/ Cs
			Membrane		
3381	HTS Transwell-96 System, reservoir and receiver plates with 2 lids	0.4	PC	1	1
3391	HTS Transwell-96 System, reservoir and receiver plates with 2 lids	0.4	PC	1	5
3380	HTS Transwell-96 System, reservoir and receiver plates with 2 lids	1.0	PET	1	1
3392	HTS Transwell-96 System, reservoir and receiver plates with 2 lids	1.0	PET	1	5
3385	HTS-Transwell-96 Well Plate, receiver plate and lid, individual	3.0	PC	1	2
3386	HTS-Transwell-96 Well Plate, receiver plate and lid, bulk	3.0	PC	4	8
3387	HTS-Transwell-96 Well Plate, receiver plate and lid, bulk	5.0	PC	4	8
3388	HTS-Transwell-96 Well Plate, receiver plate and lid, individual	5.0	PC	1	2
3374	HTS-Transwell-96 Well Plate, receiver plate and lid, individual	8.0	PET	1	2
3384	HTS-Transwell-96 Well Plate, receiver plate and lid, bulk	8.0	PET	4	8
3382	HTS Transwell-96 receiver plate with lid, tissue culture treated	n/a	n/a	10	10
3383	HTS Transwell-96 reservoir plate with removable media stabilizer and lid, not treated	n/a	n/a	10	10
3583	HTS Transwell-96 black receiver plate with lid, tissue culture treated	n/a	n/a	10	10



Netwell Inserts

Netwell™ Inserts

- ▶ Costar® Netwell inserts have polyester mesh bottoms attached to polystyrene rings or housing
- ▶ They are used as tissue carriers, supports and strainers for culture of small organs, tissue slices or explants at the air-media interface
- ▶ Handy carrier for immunocytochemical staining of tissue slices (see accessories below)
- ▶ Provides coarse filtration of tissue homogenates, cell suspensions and microcarriers
- ▶ Available in two mesh sizes and diameters
- ▶ Supplied sterile and preloaded in 6- or 12-well plates
- ▶ 24 mm Netwell inserts fit in Corning® 50 mL plastic centrifuge tubes

Netwell Inserts Ordering Information

Cat. No.	Membrane Dia. (mm)	Polyester Membrane Mesh Size (µm)	Sterile	Inner Packaging	Inserts/Cs
3477	15	74	Yes	12/plate	48
3478	15	500	Yes	12/plate	48
3479	24	74	Yes	6/plate	48
3480	24	500	Yes	6/plate	48

Netwell Accessories

- ▶ Specially designed Netwell carriers and handles allow simultaneous processing of up to 12 samples per carrier
- ▶ Polystyrene reagent trays are available in white for colorimetric reaction contrast, or black for better visibility of tissue sections
- ▶ Each carrier kit contains eight carriers and eight handles

Netwell Accessories Ordering Information

Cat. No.	Description	Qty/Cs
3517	Netwell Reagent Tray, black	200
3519	Netwell Reagent Tray, white	200
3520	Netwell Carrier Kit, 15 mm	8
3521	Netwell Carrier Kit, 24 mm	8



Netwell Accessories

Culture Tubes



430172 Culture Tube

Culture Tubes

- ▶ Manufactured from optically clear polystyrene
- ▶ Threaded plug seal caps prevent leakage
- ▶ Cell culture treated tubes supplied racked
- ▶ Untreated tubes provided bulk packed
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic

Culture Tube Ordering Information

Cat. No.	Treated	Size (mm)	Cap Style	Qty/Pk	Qty/Cs
430157	No	16 x 125	Screw Top	25	500
430172	Yes	16 x 125	Screw Top	50	500

Roller Bottles



430849 850 cm² Roller Bottle

Roller Bottle Application Tips

- Corning recommends 0.2 to 0.3 mL of medium per cm² of growth area.
- Corning recommends setting roller rack speeds to provide 0.5 to 1.0 RPM.



430852 Expanded Surface Roller Bottle

Roller Bottles

- Manufactured from virgin polystyrene
- Treated for optimal cell attachment
- One piece seamless construction
- Most bottles have graduations.
- All bottles have printed lot numbers to aid in product traceability.
- Sterilized by gamma radiation
- Certified nonpyrogenic

Roller Bottle Ordering Information

Cat. No.	Surface	Surface Area (cm ²)	Cap Style	Graduations	Qty/Pk	Qty/Cs
430195	TC	490	Plug Seal	No	2	40
430699	TC	1,750	Easy Grip	Yes	10	20
430849	TC	850	Easy Grip	Yes	2	40
431133	TC	850	Easy Grip	Yes	20	20
431198	TC	850	Easy Grip Vent	Yes	2	40
430851	TC	850	Easy Grip	Yes	5	40
431318	TC	850	Easy Grip	No	20	80
431321	TC	850	Easy Grip	Yes	22	44
3907	Corning® CellBIND® Surface	850	Easy Grip	Yes	2	40
431329	Corning CellBIND Surface	850	Easy Grip Vent	Yes	2	40
431344	Corning CellBIND Surface	850	Easy Grip	Yes	22	44

Expanded Surface Roller Bottles

- Same features as standard roller bottles
- Ribbed design provides twice the surface area with the same exterior dimensions

Expanded Surface Roller Bottle Ordering Information

Cat. No.	Surface	Surface Area (cm ²)	Cap Style	Graduations	Qty/Pk	Qty/Cs
430852	TC	1,700	Easy Grip	Yes	2	40
430853	TC	1,700	Easy Grip	Yes	5	40
431134	Corning CellBIND Surface	1,700	Easy Grip	Yes	20	20
431135	TC	1,700	Easy Grip	Yes	20	20
431191	TC	1,700	Easy Grip Vent	Yes	20	20

Expected Cell Yields and Recommended Medium Volumes

Corning Roller Bottles	Approximate Growth Area (cm ²)	Average Cell Yield*	Recommended Medium Volume (mL)
490 cm ² roller bottle	490	4.9 x 10 ⁷	100 - 150
850 cm ² roller bottle	850	8.5 x 10 ⁷	170 - 255
1700 cm ² roller bottle	1,700	1.7 x 10 ⁸	340 - 510
1750 cm ² roller bottle	1,750	1.75 x 10 ⁸	350 - 525

*Assumes an average yield of 1 x 10⁵ cells/cm² from a 100% confluent culture. Yields from many cell types can be lower than this.



Easy Grip Cap features large knurls designed for ergonomic handling.



Easy Grip Vent Cap is designed for applications requiring consistent gas exchange.



Plug Seal Cap, designed for use in closed systems, provides a liquid- and gas-tight seal. When loosened, this cap can be used in open systems.

Polyethylene Roller Bottle Caps

Caps are sold separately and are available individually wrapped in either Easy Grip or Easy Grip Vent Cap designs.

Cat. No.	Cap Style	Qty/Pk	Qty/Cs
430698	Easy Grip	1	100
431132	Easy Grip Vent	1	300

Corning is committed to partnering with you, our customer, to provide solutions that increase your efficiency and productivity. We offer the ability to customize packaging and cap design to meet your specific requirements. Minimum order quantities apply. Please call us or contact your local Corning Office for more details. See back cover for contact information.

New!

HYPERFlask™ Cell Culture Vessel

The new Corning HYPERFlask Vessel offers 1720 cm² growth area in the footprint of a traditional 175 cm² flask. This high yield, high performance flask utilizes a multilayered gas permeable growing surface for efficient gas exchange.

Cat. No.	Description	Qty/Pk	Qty/Cs
10010	HYPERFlask Vessel, 1720 cm ² , Corning CellBIND Surface, Bar Code, Sterile	4	4
10024	HYPERFlask Vessel, 1720 cm ² , Corning CellBIND Surface, Bar Code, Sterile	4	24

Corning® CellSTACK® Culture Chambers



Corning CellSTACK Culture Chambers

- ▶ Available in Five Sizes
 - 1-Stack with 636 cm² cell growth area
 - 2-Stack with 1,272 cm² cell growth area
 - 5-Stack with 3,180 cm² cell growth area
 - 10-Stack with 6,360 cm² cell growth area
 - 40-Stack with 25,440 cm² cell growth area
- ▶ Choice of traditional surface treatment, new **Corning CellBIND® Surface** for enhanced cell attachment, or **Ultra-Low Attachment Surface** for reduced cell attachment
 - Great for reducing serum levels
 - Better attachment increases cell yields
 - May eliminate need for expensive coatings
- ▶ Greater Chamber Durability
 - Superior mechanical strength and structural integrity
 - Self venting caps prevent pressure build-up during transport
 - 100% leak tested prior to shipping
- ▶ Greater Cleanliness
 - Improved assembly procedures reduce particulates
 - Certified nonpyrogenic and sterilized by gamma irradiation
- ▶ Continuous Supply Reliability
 - Manufactured in USA under GMP conditions
- ▶ Easier to Use
 - Larger openings with threaded closures and vented caps
 - Footprint identical to competitor's product



CellSTACK Chamber, 40-Stack

Corning CellSTACK Culture Chambers Ordering Information

Cat. No.	Surface	Growth Area (cm ²)	Description	Qty/ Pk	Pk/ Cs
3330	Corning CellBIND Surface	636	CellSTACK-1 Chamber	1	8
3268	TC	636	CellSTACK-1 Chamber	1	8
3310	Corning CellBIND Surface	1,272	CellSTACK-2 Chamber	1	5
3269	TC	1,272	CellSTACK-2 Chamber	1	5
3311	Corning CellBIND Surface	3,180	CellSTACK-5 Chamber	1	2
3319	TC	3,180	CellSTACK-5 Chamber	1	2
3313	TC	3,180	CellSTACK-5 Chamber	1	8
3320	Corning CellBIND Surface	6,360	CellSTACK-10 Chamber	1	6
3312	Corning CellBIND Surface	6,360	CellSTACK-10 Chamber	1	2
3270	TC	6,360	CellSTACK-10 Chamber	1	2
3271	TC	6,360	CellSTACK-10 Chamber	1	6
3321	Corning CellBIND Surface	25,440	CellSTACK-40 Chamber	1	2
3272	TC	25,440	CellSTACK-40 Chamber	1	2
3303	Ultra-Low Attachment	636	CellSTACK-1 Chamber	1	8

Corning® CellSTACK® Accessories are Simply Better!

Corning offers a variety of accessories to simplify handling and reduce contamination risks when processing CellSTACK Chambers.

For Better Filling

A variety of optional filling caps are available to allow direct aseptic transfer of media and cells via pumping or gravity feed. Several coupling devices are available on these filling caps with or without integrally sealed USP Class VI certified C-Flex® tubing. Optional filling caps with attached filters with hydrophobic membranes provide for gas exchange and faster aseptic venting during liquid transfers. Extra sterile vented or unvented 33 mm replacement caps are also available.

For Better Stacking

Reusable stacking devices fit between CellSTACK Chambers to keep them level and optimize incubator space while providing clearance for gas exchange.

For Better Options

Sometimes, currently available accessories just don't fit a customer's needs. This is why Corning will work with you to design a CellSTACK Chamber accessory that will make your work flow process more efficient and reliable.

For large scale production using CellSTACK-40 Chambers there are automated systems that can save on labor while increasing reliability and efficiency.

Call us to discuss your specific requirements.

Corning CellSTACK Accessories Ordering Information

Cat. No.	Description	Qty/ Pk	Qty/ Cs
3331	Stacking device, ABS, nonsterile	1	5
3332	Universal cap*, with vented overcap, sterile	1	4
3969	Solid cap, sterile	1	6
3968	Vent cap, 0.2 mm membrane, sterile	1	6
3281	Vent cap, 3/8" (9.5 mm) ID tubing, 7 cm length, Pall® Acro 50, PVDF filter, sterile	1	5
3282	Fill cap, 1/8" (3.2 mm) ID tubing, female luer lock with male luer plug, sterile	1	5
3283	Fill cap, 3/8" (9.5 mm) ID tubing and 5/16" (7.94 mm) barbed fitting, sterile	1	5
3284	Vent cap, 3/8" (9.5 mm) ID tubing, 7 cm length, Pall Bacterial Air Vent, sterile	1	4
3324	Two vented over caps and one solid over cap for the Universal Cap, sterile	5	100
3333	Fill cap*, 1/4" (6.4 mm) ID tubing, 70 cm length, male MPC coupling with female end cap, sterile	1	4
3328	Fill cap, female MPC coupling, 1/4" (6.4 mm) ID barbed fitting with male end cap, sterile	1	4
3329	Fill cap, female MPC coupling, 3/8" (9.5 mm) ID barbed fitting with male end cap, sterile	1	4
3334	Fill cap, male MPC coupling, 1/4" (6.4 mm) ID barbed fitting with female end cap, sterile	1	4
3339	Fill cap, male MPC coupling with male end cap, 3/8" (9.5 mm) ID barbed fitting with female end cap, sterile	1	4

*All caps are 33 mm thread caps.



CellSTACK Accessories



3328 Fill Cap, Female MPC Coupling



3281 0.2 µm Vent Cap



3284 Bacterial Air Vent



3333 Fill Cap, Male MPC Coupling



3331 Stacking Device



3332 Universal Cap

CellCube® Systems



The CellCube System provides a fast, simple, and compact method for the mass culture of attachment-dependent cells. It uses a tissue culture treated growth surface for cell attachment, and continually perfuses the cells with fresh medium for increased cell productivity. The CellCube System is comprised of four pieces of capital equipment: the system controller, oxygenator, pump tower, and circulation pump, and is designed to use disposable CellCube Modules. Performance data from the CellCube System can be easily scaled to the production system. Please inquire about CellCube System pricing. Corning provides on-site technical support for the CellCube System.

The CellCube Modules provide a traditional tissue culture treated surface or new Corning® CellBIND® Surface for the growth of attachment dependent cells. The CellCube System provides an environment which more closely simulates *in vivo* conditions and reliably distributes nutrients and oxygen with low differential gradients across all cells within the modules.

CellCube System Ordering Information

Cat. No.	Description	Qty/Cs
3040	CellCube 6 Liter Oxygenator, Process Scale	1
3041	CellCube 6 Liter Oxygenator, Production Scale	1
3020	CellCube Set Up Kit, Bioprene® Tubing, Process Scale	1
3021	CellCube Set Up Kit, Bioprene Tubing, Production Scale	1
3022	CellCube Set Up Kit, Sta-Pure Tubing, Process Scale	1
3023	CellCube Set Up Kit, Sta-Pure Tubing, Production Scale	1
3101	CellCube Single Module System 6 Liter Oxygenator, Complete	1
3139	CellCube Single Module System Secondary Oxygen Probe (25 x 70 mm)	1
3138	CellCube Single Module System Secondary Oxygen Probe Holder	1
3144	CellCube Single Module System Oxygen Probe Cable	1
3165	CellCube Single Module System 12 mm Dissolved Oxygen Probe Membrane Kit	1
3166	CellCube System 25 mm Dissolved Oxygen Probe Membrane Kit	1
3136	CellCube Single Module System Stainless Steel Stand	1
3135	CellCube Single Module System Setup Kit	1
3200	CellCube 10-Stack Module (8,500 cm ²), Tissue Culture Treated	2
3201	CellCube 25-Stack Module (21,250 cm ²), Tissue Culture Treated	1
3304	CellCube 25-Stack Module (21,250 cm ²), Corning CellBIND Surface	1
3264	CellCube 100-Stack Module (85,000 cm ²), Tissue Culture Treated	1
3302	CellCube 100-Stack Module (85,000 cm ²) Corning CellBIND Surface	1

Corning E-Cube™ Culture System



Corning E-Cube Culture System

The E-Cube system provides a simple method to determine if your cells will grow in the CellCube module prior to providing in the resources and funding that would be necessary for the larger, automated CellCube system.

Corning E-Cube Culture System Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
3286	E-Cube System Kit (without CellCube module)	1	1
3200	CellCube Module 10-Stack	1	2

Corning E-Cube Culture System Accessories Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
430518	1 L Storage Bottle with cap	2	24
401654	45 mm Cap with 2 stainless steel ports	1	1
3287	E-Cube Fittings	1	1

Spinner Flasks



3152 and 3153
Disposable Spinner Flasks



3561 and 3563
Disposable Spinner Flasks

Corning® Disposable Spinner Flasks

- ▶ The Corning disposable spinner flask system comes ready-to-use with paddle and integrated magnet, eliminating the need for time-consuming assembly or cleaning and reassembly
- ▶ Molded from virgin polystyrene and gamma-irradiated, each spinner flask system assures a clean sterile unit. No more concerns with detergent residues or contamination
- ▶ Made of ISO 10993 compliant polystyrene, the vessel is comparable to conventional glass spinner flasks for growth of suspension cell lines and any attachment-dependent cultures using microcarrier beads. The 1L and 3L impellers are made of ISO 10993 compliant polypropylene.
- ▶ The paddle size and height is optimized for each vessel size. A unique integrated magnet provides smooth, even rotation at required speeds on slow-speed stirrers. Heat build-up in the vessel is reduced by means of a specially designed flange that raises the vessel off the stir-plate surface for the 125 mL and 500 mL flasks only.

Corning Disposable Spinner Flasks Ordering Information


Cat. No.	Description	Capacity (mL)	Center Neck (mm)	Sidearm Neck (mm)	Qty/Cs
3152	Disposable Spinner Flask	125	70	25	12
3153	Disposable Spinner Flask	500	100	45	12
3561	Disposable Spinner Flask	1000	-	45	6
3563	Disposable Spinner Flask	3000	-	45	4
3569	Disposable Spinner Flask with accessory attached	1000	-	45	6

Replacement Caps & Aseptic Transfer Caps

Cat. No.	Description	Capacity (mL)	Sidearm Neck (mm)	Qty/Cs
3567	Vent Cap, 0.2 µm Vent	500, 1000, 3000	45	4
3565	Side Arm Aseptic Transfer Cap, DipTube w/ 0.2 µm Vent, MLL	500	45	2
3562	Side Arm Aseptic Transfer Cap, DipTube w/ 0.2 µm Vent, MLL	1000	45	2
3564	Side Arm Aseptic Transfer Cap, DipTube w/ 0.2 µm Vent, MLL	3000	45	2



1L and 3L Disposable Spinner Flasks with Accessories



3565, 3562 and 3564
Aseptic Transfer Cap



3567 Vent Cap



4500-1L and 4500-250
Spinner Flasks



ProCulture Spinner Flasks

ProCulture® Glass Spinner Flask with Angled Sidearms

- ▶ Baffles enhance aeration and agitation of contents of the flask.
- ▶ Unique impeller design ensures optimal stirring.
- ▶ Sidearm designs permit easy access of 25 mL and 50 mL pipettes
- ▶ Visit www.corning.com/lifesciences to view additional Corning spinner flask accessories


ProCulture Spinner Flasks with Angled Sidearms Ordering Information

Cat. No.	Description	Capacity	Center Neck (mm)	Sidearm Neck (mm)	Qty/Cs
4500-125	Spinner	125 mL	70	32	1
4500-250	Spinner	250 mL	70	32	1
4500-500	Spinner	500 mL	100	45	1
4500-1L	Spinner	1L	100	45	1
4500-3L	Spinner	3L	100	45	1
4500-6L	Spinner	6L	100	45	1
4500-8L	Spinner	8L	100	45	1
4500-15L	Spinner	15L	100	45	1
4500-36L	Spinner	36L	100	45	1
4502-3L	Spinner	3L	120	45	1
4502-6L	Spinner	6L	120	45	1
4502-8L	Spinner	8L	120	45	1
4502-15L	Spinner	15L	120	45	1
4502-36L	Spinner	36L	120	45	1
4504-3L	Spinner	3L	140	45	1
4504-6L	Spinner	6L	140	45	1
4504-8L	Spinner	8L	140	45	1
4504-15L	Spinner	15L	140	45	1
4504-36L	Spinner	36L	140	45	1


Retrofit Kits are available for converting older Corning® ProCulture Spinner Flasks to fit newer dual-bearing impellers.

ProCulture Spinner Flasks with Vertical Sidearms Ordering Information


Cat. No.	Capacity	Center Neck (mm)	Number of Vertical Sidearms	Sidearm Neck (mm)	Qty/Cs
4510-8L	8L	100	4	45	1
4510-15L	15L	100	4	45	1
4510-36L	36L	100	6	45	1
4512-8L	8L	120	4	45	1
4512-15L	15L	120	4	45	1
4512-36L	36L	120	6	45	1
4514-15L	15L	140	4	45	1
4514-36L	36L	140	6	45	1



Vertical Sidearm Fittings,
Gas Delivery and Venting



Vertical Sidearm Fittings,
Media Handling



Dual Angled Sidearm Fittings

Gas Handling Fittings, Vertical Sidearm Flasks

- Used to provide gases into larger spinner flasks with vertical sidearms
- Fittings are comprised of a PET insert with a Viton® O-Ring and a polypropylene sealing cap
- Gas filters are PTFE, 0.2 µm porosity
- The 316 stainless steel tubes are held in place by Noryl® nuts with integrated ferrules
- The fittings are completely autoclavable

Cat. No.	Description	Dimension	Qty/Cs
4519-100	Sidearm fitting, gas delivery	1/8" Inlet	1
4519-102	Sidearm fitting, gas delivery	1/4" Inlet	1
4519-104	Sidearm fitting, delivery and vent	1/8" and 1/4"	1
4519-106	Sidearm fitting, vent cap, 0.2 µm	50 mm filter	1
4519-177	Sidearm fitting, vent cap, 0.2 µm, Sanitary	50 mm filter	1

Media Handling Fittings, Vertical Sidearm Flasks


- Used to introduce medium aseptically into large spinner flasks with vertical sidearms
- Fittings are comprised of a PET insert with a Viton O-Ring and a polypropylene sealing cap
- Gas filters are PTFE, 0.2 µm porosity
- The 316 stainless steel tubes are held in place by Noryl nuts with integrated ferrules
- The fittings are completely autoclavable

Cat. No.	Description	Fits Flask Size	Tubing O.D. (inches)	Qty/ Case
4519-112	Sidearm fitting, dual, media handling	8L, 15L	1/8"	1
4519-114	Sidearm fitting, dual, media handling	36L	1/8"	1
4519-116	Sidearm fitting, dual, media handling	8L, 15L	1/4"	1
4519-118	Sidearm fitting, dual, media handling	36L	1/4"	1
4519-120	Sidearm fitting, combo, media handling	8L, 15L	1/8", 1/4"	1
4519-122	Sidearm fitting, combo, media handling	36L	1/8", 1/4"	1
4519-124	Sidearm fitting, single, media handling	8L, 15L	1/2"	1
4519-126	Sidearm fitting, single, media handling	36L	1/2"	1
4519-176	Sidearm fitting, dual, media handling, EPDM	8L, 15L	1/4"	1


Gas or Media Handling Fittings, Angled Sidearm Flasks, Dual Style

- Dual angled sidearm fittings can be used for aseptically transferring medium into or out of angled sidearm spinner flasks or for sparging the medium with gases
- Fittings are comprised of a PET insert with a Viton O-ring and a polypropylene sealing cap
- Two 316 stainless steel tubes which extend to the bottom of the flask, are held in place by Noryl nuts with integrated ferrules
- The fittings are completely autoclavable

Cat. No.	Description	Flask Size	Tubing O.D. (inches)	Qty/ Case
4519-150	SA fitting, Dual	1L	1/8"	1
4519-151	SA fitting, Dual	3L	1/8"	1
4519-152	SA fitting, Dual	6L	1/8"	1
4519-153	SA fitting, Dual	8L	1/8"	1
4519-173	SA fitting, Dual	1L	1/8", 1/4"	1
4519-121	SA fitting, Dual	8L	1/8", 1/4"	1
4519-174	Sidearm fitting, Dual	500 mL	1/8" angled to 125 mL level, 1/4"	1
4519-154	Sidearm fitting, Dual	1L	1/4"	1
4519-155	Sidearm fitting, Dual	3L	1/4"	1
4519-156	Sidearm fitting, Dual	6L	1/4"	1
4519-157	Sidearm fitting, Dual	8L	1/4"	1
4519-170	Sidearm fitting, Dual	15L	1/4"	1



Combination and Triple Angled Sidearm Fittings



Sidearm Fittings for Sensors



Impeller Assembly

Gas or Media Handling Fittings, Angled Sidearm Flasks, Combination Style

- Used to aseptically transfer medium, sparge the cell culture medium directly or add gases to the head space above the cell culture medium
- Fittings are comprised of a PET insert with a Viton® O-ring and a polypropylene sealing cap
- One or two 316 stainless steel tubes extend to the bottom of the flask; the other is a shorter 6" length
- Both tubes are held in place by Noryl® nuts with integrated ferrules
- The fittings are completely autoclavable

Cat. No.	Description	Flask Size	Tubing O.D. (inches)	Qty/Case
4519-158	Sidearm fitting, combination	1L	1/8"	1
4519-159	Sidearm fitting, combination	3L	1/8"	1
4519-160	Sidearm fitting, combination	6L	1/8"	1
4519-161	Sidearm fitting, combination	8L	1/8"	1
4519-162	Sidearm fitting, combination	1L	1/4"	1
4519-163	Sidearm fitting, combination	3L	1/4"	1
4519-164	Sidearm fitting, combination	6L	1/4"	1
4519-165	Sidearm fitting, combination	8L	1/4"	1
4519-171	Sidearm fitting, combination	15L	1/4"	1
4519-166	Sidearm fitting, combination, triple	1L	1/8"	1
4519-167	Sidearm fitting, combination, triple	3L	1/8"	1
4519-168	Sidearm fitting, combination, triple	6L	1/8"	1
4519-169	Sidearm fitting, combination, triple	8L	1/8"	1

Fittings for Insertion Probes, Vertical Sidearm Flasks

- Used to secure pH, O₂, or temperature sensors in large spinner flasks with vertical sidearms
- Fittings are comprised of a PET insert with a Viton O-ring and a polypropylene sealing cap
- The 316 sensors are held in place by Noryl nuts with integrated ferrules
- The fittings are completely autoclavable

Cat. No.	Description	Sensor O.D. (mm)	Qty/Cs
4519-108	Sidearm fitting, sensor, O ₂ probes	12	1
4519-128	Sidearm fitting, sensor, temperature probes	12	1
4519-110	Sidearm fitting, sensor, pH probes	12	1
4519-172	Sidearm fitting, sensor, pH or O ₂	18	1

Impeller Assembly for Magnetically-Driven Bioreactor


Stainless steel impeller shaft with modified impeller blade for use with probes to create a small bioreactor.

Cat. No.	Description	Qty/Cs
402648	Impeller assembly, stainless steel, dual bearing, modified for probes, 3L	1
402649	Impeller assembly, stainless steel, dual bearing, modified for probes, 6L	1
401392	Impeller assembly, stainless steel, dual bearing, modified for probes, 8L	1
401661	Impeller assembly, stainless steel, dual bearing, modified for probes, 15L	1
402650	Impeller assembly, stainless steel, dual bearing, modified for probes, 36L	1

Cap Assembly for Magnetically-Driven Bioreactor

Cap assembly for small bioreactor with various fitting arrangements.


Cat. No.	Description	Qty/Cs
402579	Cap Assembly, 120 mm, Glass Filled PBT, 3 (3/8"), 1 (1/4") fittings	1
402576	Cap Assembly, 120 mm, Glass Filled PBT, 2 (12 mm), 2 (1/4") fittings	1
402577	Cap Assembly, 120 mm, Glass Filled PBT, 2 (12 mm), 2 (1/4"), 1 (3/8") fittings	1



Spare Parts for Sidearm Fittings

Securing Caps

Corning Cat. No.	Description	Qty/Cs
1395-32LTC	Cap, solid, 32 mm, orange	1
1395-45LTC	Cap, solid, 45 mm, orange	1
1395-45LTR	Drip ring, 45 mm, clear	1
1395-45LTMC	Cap, vented, securing, 45 mm, .22 PTFE, grey	10



Direct Drive Motor

Direct Drive Motors

- High torque, low rpm stirrer designed to maintain constant low speed
- Gearhead stirrer delivers 14.5 in-lbs of torque
- Maximum speed is 350 rpm
- Weight of motor is 9 lbs (4.1kg)
- Available with 120VAC 60Hz or 230VAC 50Hz

Cat. No.	Description	Qty/Cs
400640	120VAC, 60 Hz Motor	1
402645	230VAC, 50 Hz Motor	1



Direct Drive Shaft/Cap Assembly

Direct Drive Shaft/Cap Assemblies

- For 8L, 15L, or 36L paddle assemblies
- Used on all series 4510 and 4512 Spinner flasks

Cat. No.	Description	Qty/Cs
402614	For 100 mm Neck Flasks	1
400649	For 120 mm Neck Flasks	1

Direct Drive Paddle Assemblies

- For series 4510, 4512, and 4514 Spinner flasks when coupled to a direct drive motor
- Paddle assemblies will couple to 100 mm and 120 mm cap assemblies

Cat. No.	Description	Qty/Cs
4515-8L	Paddle assembly only for 8L flask	1
4515-15L	Paddle assembly only for 15L flask	1
4515-36L	Paddle assembly only for 36L flask	1



Direct Drive Paddle Assembly

Erlenmeyer Flasks



Sizes range from 125 mL to 3L for plain and baffled Erlenmeyer flasks



Unique baffled design with a molded-in "1/3 Fill" line for convenience on all baffled Erlenmeyer flasks.

Shaker Flask Application Tip

Corning recommends starting with a shaking rate of 75-125 RPM (orbital shaker) and a medium volume of 30-40% of the nominal flask capacity.

Corning® Erlenmeyer Flasks

Corning baffled and plain Erlenmeyer and Fernbach culture flasks are ideal for shaker culture applications and storage. Like all Corning flasks, the Erlenmeyer flasks are certified nonpyrogenic and sterile.

- Polycarbonate construction: USP Class VI material provides excellent optical clarity and mechanical strength
- Sizes range from 125 mL to 3L
- Baffled or plain bottom options in every size
- Molded-in graduations for accuracy
- Vent cap option for continuous gas exchange while ensuring sterility and preventing leakage
- Individually packaged and radiation sterilized for ease of use
- All flasks have the highest Sterility Assurance Level (SAL) of 10^{-6}
- Certified nonpyrogenic

Corning Polycarbonate Erlenmeyer Flasks and Caps Ordering Information

Baffled Bottom Erlenmeyer Flasks

Cat. No.	Description	Sterile	Qty/Cs
431405	Erlenmeyer Flask, Baffled, 125 mL, Vent Cap	Yes	50
431404	Erlenmeyer Flask, Baffled, 125 mL, Plug Seal Cap	Yes	50
431407	Erlenmeyer Flask, Baffled, 250 mL, Vent Cap	Yes	50
431406	Erlenmeyer Flask, Baffled, 250 mL, Plug Seal Cap	Yes	50
431401	Erlenmeyer Flask, Baffled, 500 mL, Vent Cap	Yes	25
431408	Erlenmeyer Flask, Baffled, 500 mL, Plug Seal Cap	Yes	25
431403	Erlenmeyer Flask, Baffled, 1L, Vent Cap	Yes	25
431402	Erlenmeyer Flask, Baffled, 1L, Plug Seal Cap	Yes	25
431256	Erlenmeyer Flask, Baffled, 2L, Vent Cap	Yes	6
431253	Fernbach Culture Flask, Baffled, 3L, Vent Cap	Yes	4

Plain Bottom Erlenmeyer Flasks

Cat. No.	Description	Sterile	Qty/Cs
431143	Erlenmeyer Flask, 125 mL, Vent Cap	Yes	50
430421	Erlenmeyer Flask, 125 mL, Plug Seal Cap	Yes	50
431144	Erlenmeyer Flask, 250 mL, Vent Cap	Yes	50
430183	Erlenmeyer Flask, 250 mL, Plug Seal Cap	Yes	50
431145	Erlenmeyer Flask, 500 mL, Vent Cap	Yes	25
430422	Erlenmeyer Flask, 500 mL, Plug Seal Cap	Yes	25
431147	Erlenmeyer Flask, 1L, Vent Cap	Yes	25
431146	Erlenmeyer Flask, 1L, Plug Seal Cap	Yes	25
431255	Erlenmeyer Flask, 2L, Vent Cap	Yes	6
431252	Fernbach Culture Flask, 3L, Vent Cap	Yes	4

Replacement Erlenmeyer Flask Caps

Corning® Polypropylene Erlenmeyer Flask Caps are also available separately. They are sterile, individually packaged and available for the 500 mL and 1L*, 2L and 3L flask sizes.

Cat. No.	Description	Sterile	Qty/Cs
431372*	43 mm Vent Cap, 500 mL and 1L Erlenmeyer Flask*	Yes	50
431339	48 mm Vent Cap, 2L Erlenmeyer Flask	Yes	24
431364	48 mm Flat Cap, 2L Erlenmeyer Flask	Yes	24
431340	70 mm Vent Cap, 3L Erlenmeyer Flask	Yes	24
431363	70 mm Flat Cap, 3L Erlenmeyer Flask	Yes	24

*The 43 mm cap for the 500 mL and 1L sizes are available Made to Order only with a 5 case minimum.

Aseptic Transfer Caps

Corning Erlenmeyer Flask Aseptic Transfer Caps are available separately. They are sterile and arrive individually doubled bagged. They are available for the 1L, 2L and 3L flask sizes. The transfer caps have two ports. One port ends in a 0.2 µm Acro 50 mm disk and the other port is C-Flex Tubing ending in either a male luer lock or a male polycarbonate quick connect. The Diptube reaches all the way to the bottom of the flask for easy aseptic transfer of your liquid handling processes.

Cat. No.	Description	Sterile	Qty/Cs
431444	43 mm Cap, 1L, Diptube with 0.2 µm Vent, MLL	Yes	5
431445	43 mm Cap, 1L, Diptube with 0.2 µm Vent, MPC	Yes	5
431446	48 mm Cap, 2L, Diptube with 0.2 µm Vent, MLL	Yes	6
431447	48 mm Cap, 2L, Diptube with 0.2 µm Vent, MPC	Yes	6
431448	70 mm Cap, 3L, Diptube with 0.2 µm Vent, MLL	Yes	4
431449	70 mm Cap, 3L, Diptube with 0.2 µm Vent, MPC	Yes	4



431363 Flat Cap



431340 Vent Cap



Aseptic Transfer Caps, MLL



Aseptic Transfer Caps, MPC

Cell Scrapers and Lifters



3008 Cell Lifter



3010 Small Cell Scraper

Cell Scrapers and Cell Lifters

- › Useful for the manual harvesting of cells
- › Blade design minimizes cell damage and ensures even contact with the growth surface
- › Cell lifter is useful for harvesting cells (especially stem cells) in dishes
- › Scrapers designed for use in flasks
- › Individually wrapped
- › Sterilized by gamma radiation
- › Certified nonpyrogenic

Cell Scraper and Lifter Ordering Information

Cat. No.	Description	Blade Length (cm)	Handle Length (cm)	Qty/Pk	Qty/Cs
3008	Cell lifter	1.9	18	1	100
3010	Small scraper	1.8	25	1	100
3011	Large scraper	3.0	39	1	100

Technical Appendix

CORNING® CELL CULTURE SURFACES

Introduction

For over eighty years Corning has been developing products and surfaces for cell culture. Corning currently offers six polystyrene-based surfaces (Table 1) for growing cells including the most recent technology revolution, the patented Corning CellBIND® surface (U.S. Patent 6,617,152):

Most of these early plastic vessels were made from polystyrene, a long carbon chain polymer with benzene rings attached to every other carbon. Polystyrene was chosen because it has excellent optical clarity, is easy to mold and is relatively inexpensive. However, it also has one significant drawback: it is a very hydrophobic (nonwetable) polymer to which cells have difficulty attaching. Fortunately, the surface of polystyrene can be easily modified by a variety of chemical (sulfuric acid) and physical (corona discharge, gas-plasma or irradiation) methods. Using these methods, hydroxyl, ketone, aldehyde, carboxyl and amine groups can readily be grafted onto the polymer (Figure 1). These groups modify the surface characteristics changing the uncharged hydrophobic surface into a more ionic hydrophilic surface. Polystyrene can also be modified through chemical reactions to allow the covalent attachment of a variety of reactive groups that can be used for the subsequent covalent immobilization of biomolecules. For additional information, please check the References.

Corning CellBIND® Surface

The Corning CellBIND culture surface, the first novel cell culture surface treatment in over 20 years, is designed to improve cell attachment under difficult conditions, such as reduced-serum or serum-free medium, resulting in higher cell yields. It is also useful for growing “difficult” cells such as primary cultures or transfected cells over expressing proteins. Developed by Corning scientists, this patented technology (U.S. Patent 6,617,152) uses a novel microwave plasma process for treating the culture surface. This process improves cell attachment by incorporating significantly more oxygen into the cell culture surface than traditional plasma or corona discharge treatments, rendering it more hydrophilic (wetable) and increasing the stability of the surface.




Figure 1. Polystyrene can be surface modified by the addition of a variety of different chemical groups, by breaking the carbon chain backbone, or by opening the benzene ring (not shown).

Unlike biological coatings, the Corning CellBIND surface is a nonbiological surface that requires no special handling or storage. Because the polymer is treated, rather than coated, the surface is more consistent and stable. This enhanced cell performance has already led to a major biotechnology company choosing Corning roller bottles with the Corning CellBIND surface for producing a new FDA approved protein therapeutic.

Corning CellBIND surface benefits:

- ▶ Quickly adapts cells to reduced serum or serum-free conditions
- ▶ May eliminate the need for tedious, time-consuming, expensive and low stability biological coatings
- ▶ Stable at room temperature, requires no refrigeration or special handling
- ▶ Gives more consistent and even cell attachment for difficult to attach cell lines, especially transfected cells
- ▶ Reduces premature cell detachment from confluent cultures especially in roller bottles and during cell-based assays

The Corning CellBIND® surface is available on flasks, multiple well plates, CellSTACK® Culture Chambers, roller bottles, 96 well plates, 384 well plates, and dishes.

Table 1. Corning Cell Culture Surfaces

Corning Surface	Binding Interaction	Sample Properties
Corning CellBIND Surface modified polystyrene surface	Hydrophilic and ionic (negatively charged)	Improves cell attachment and binding to polystyrene
Standard Tissue culture treated polystyrene	Hydrophilic and ionic (negatively charged)	Allows cell attachment and binding to polystyrene
Untreated polystyrene	Hydrophobic	Significantly reduces the attachment of most cells
Ultra Low Attachment coated polystyrene	Hydrophilic and nonionic	Hydrogel layer prevents the attachment of almost all cells
Poly-D-lysine coated polystyrene	Hydrophilic and ionic (positively charged)	Improves cell attachment and binding to polystyrene
Ultra-Web™ Synthetic Surface	Hydrophilic and ionic interactions (positively charged) on a 3D surface	Promotes more <i>in vivo</i> -like morphology

Ultra-Web™ Synthetic Surfaces

Ultra-Web Synthetic Surfaces are composed of randomly orientated electrospun polyamide nanofibers with an average fiber diameter of ~180 nm. This creates a culturing substrate that mimics structural components within the basement membrane or extracellular matrix.

Ultra-Web synthetic surface is available with two surface chemistries:

- ▶ Untreated electrospun polyamide nanofibers with an uncharged slightly hydrophilic surface
- ▶ Polyamine treated electrospun polyamide nanofibers with a positively charged surface for enhanced cell attachment or binding and covalently linking biomolecules

Ultra-Web Synthetic Surface Applications:

- ▶ Ideal for culturing liver, neuronal, kidney and stem cell lines or primary cultures where current surfaces do not provide the necessary culture performance or function
- ▶ Replacement for poly-lysine or animal-derived biological coatings
- ▶ Ideal substrate for binding cell attachment and growth factors to create more *in vivo*-like culture environments
- ▶ Compatible with cell-based luminescence-reporter gene and FLIPR calcium flux assays
- ▶ Promotes more *in vivo*-like morphology (spheroid and dome formation)

Ultra-Web Synthetic Surface Benefits:

- ▶ Synthetic surfaces are more stable and consistent lot to lot than biological coatings
- ▶ Cells grow on the nanofiber surface, not in it, for easy harvesting
- ▶ Easy to view cells using phase contrast microscopy
- ▶ Ready to use and room temperature stable

Ultra-Low Attachment Coated Polystyrene Surface

The Corning Ultra-Low Attachment surface is a covalently bound hydrogel layer that is hydrophilic and neutrally charged. Since proteins and other biomolecules passively adsorb to polystyrene surfaces through either hydrophobic or ionic interactions, this hydrogel surface naturally inhibits nonspecific immobilization via these forces, thus inhibiting subsequent cell attachment. This surface is very stable, non-cytotoxic, biologically inert and nondegradable. Corning offers the Ultra-Low Attachment surfaces on dishes, plates, flasks, and CellSTACK® Culture Chamber 1-Stack.

This Ultra-Low Attachment surface has been shown to successfully inhibit attachment of anchorage dependent MDCK, VERO, and C6 cells grown for a period of time equal to that necessary to obtain confluent cell growth on the control surface (standard tissue culture treated polystyrene; Figure 2). This surface has also been shown to inhibit the attachment and activation of macrophages and neutrophils.

Ultra-Low Attachment culture vessels are useful for:

- ▶ Studying tissue-specific functions of certain cancer cells (i.e., MCF-7 breast cancer cells)




Figure 2. Single cell derived colonies of C6 glioma cells grow as flattened attached colonies in standard tissue culture treated surface (left panel) but form unattached spherical colonies on the ultra low attachment surface (right panel).

- ▶ Preventing stem cells from attachment-mediated differentiation
- ▶ Selectively culturing tumor or virally transformed cells as unattached colonies (substitute for soft agar assays)

Poly-D-lysine Coated Surface

Some assays and procedures require enhanced binding of cells to polystyrene. Corning poly-D-lysine (PDL) microplates are coated with PDL (molecular weight range of 70 to 150 kDa) by a proprietary method. This synthetic polymeric coating creates a uniform net positive charge on the plastic surface which, for some cell types, can enhance cell attachment, growth and differentiation, especially in serum-free and low serum conditions. PDL surfaces often improve attachment and growth of primary neurons, glial cells, neuroblastomas, and a variety of transfected cell lines, including HEK-293. Corning offers poly-D-lysine coated 96 and 384 well microplates.

Standard Tissue Culture Treated Polystyrene Surface

Standard Corning® polystyrene cell culture vessels are surface modified using either corona discharge (flasks, dishes and microplates) or gas-plasma (roller bottles and culture tubes). These processes generate highly energetic oxygen ions which graft onto the surface polystyrene chains (Figure 1) so that the surface becomes hydrophilic and negatively charged when placed in medium. Corning offers the standard tissue culture treated surface on flasks, dishes, multiple well plates, CellSTACK® Culture Chambers, roller bottles and culture tubes.

Untreated Polystyrene Surface

Natural, unmodified polystyrene surfaces are hydrophobic and only bind cells and biomolecules through passive hydrophobic interactions. Corning offers untreated polystyrene culture dishes and microplates for growing cells in stationary suspension or other applications where reduced cell attachment is desired. However, these untreated vessels are sterilized by low dose gamma irradiation, which slightly increases the wettability of the surface. Since some transformed cell lines (CHO-k1, for example) and macrophages will attach and grow on these hydrophobic surfaces, Corning also offers an Ultra-Low Attachment Surface (see below) for use in situations where cell attachment must be kept to an absolute minimum.

CHARACTERISTICS OF CORNING® PLASTICWARE

	Polystyrene	Polyethylene (High Density)	Polypropylene	Polycarbonate	Nylon	P.T.F.E. (Teflon®)	Polyethylene Terephthalate (PET)
PHYSICAL CHARACTERISTICS							
Basic Properties	Biologically inert, hard, excellent optical qualities	Biologically inert, high chemical resistance toughness	Biologically inert, high chemical resistance, exceptional resistance	Clear, very tough, inert, high temperature transmission	Tough, heat resistant, machinable, high moisture vapor slippery surface	Biologically and chemically inert, high resistant	Biologically inert, hard, tough, excellent optical qualities
Clarity	Clear	Opaque	Translucent	Clear	Opaque	Opaque	Clear
Autoclave Results	Melts	May	Withstands distort	Withstands several cycles	OK one cycle	OK	Melts
Heat Distortion Point	147-175°F 64-80°C	250°F 121°C	275°F 135°C	280-290°F 138-143°C	300-356°F 150-180°C	250°F 121°C	250°F 70°F
Burning Rate	Slow	Slow	Slow	Self-extinguishing	Self-extinguishing	None	–
EFFECTS OF LABORATORY REAGENTS							
Weak Acids	None	None	None	None	None	None	None
Strong Acids	Oxidizing acids attack	Oxidizing acids attack	Oxidizing acids attack	May be attacked	Attacked	None	Oxidizing acids attack
Weak Alkalies	None	None	None	None	None	None	None
Strong Alkalies	None	None	None	Slowly attacked	None	None	Attacked
Organic Solvents	Soluble in aromatic chlorinated hydrocarbons	Resistant below 80°C	Resistant below 80°C	Soluble in chlorinated hydrocarbons; partly soluble in aromatics	Resistant	Resistant	Soluble in aromatic or chlorinated hydrocarbons
GAS PERMEABILITY OF THIN WALL PRODUCTS*							
O ₂	Low	High	High	Very low	Very low	–	Very low
N ₂	Very low	Low	Low	Very low	Very low	–	Very low
CO ₂	High	Very high	Very high	Low	–	–	Low

Portions of this table courtesy of Modern Plastics Encyclopedia. Most data are from tests by A.S.T.M. methods. Tables show averages or ranges. Many properties vary with manufacturer, formulation, testing laboratory, and the specific operating conditions.

*Obtained from a table which lists gas permeability in CC/100 sq. inches per 24 hrs./mil.

CHEMICAL COMPATIBILITY OF CORNING® PLASTICWARE

	PS	PP	PVC	CA	PC	CN	NY	MCE	PTFE	PET
<i>Acids</i>										
Hydrochloric acid (25%)	G	G	G	N	R	R	N	O	R	R
Hydrochloric acid (concentrated)	F	G	F	N	R	N	N	N	R	O
Nitric acid (concentrated)	P	P	P	N	R	N	N	N	O	N
Nitric acid (25%)	P	G	F	N	R	L	N	O	R	R
<i>Alcohols</i>										
Butanol	G	G	G	R	R	R	R	R	R	R
Ethanol	G	G	G	R	R	N	R	O	R	R
Methanol	G	G	G	R	R	N	R	O	R	R
<i>Amines</i>										
Aniline	G	G	P	N	N	R	R	N	R	O
Dimethylformamide	P	G	F	N	N	N	R	N	R	N
<i>Bases</i>										
Ammonium hydroxide (25%)	F	G	G	R	N	R	R	O	N	O
Ammonium hydroxide (1N)	F	G	G	N	N	R	R	O	N	N
Sodium hydroxide	G	G	G	N	N	N	R	N	R	N
<i>Hydrocarbons</i>										
Hexane	P	G	F	R	R	R	R	R	R	R
Toluene	P	G	P	R	O	R	R	R	R	N
Xylene	P	F	P	R	R	R	R	R	R	N
Dioxane	P	G	P	N	N	N	R	N	R	R
Dimethylsulfoxide (DMSO)	P	G	P	N	N	N	R	N	R	O*
<i>Halogenated Hydrocarbons</i>										
Chloroform	P	N	P	N	N	R	R	N	R	R
Methylene chloride	P	F	P	N	N	R	R	N	R	N
<i>Ketones</i>										
Acetone	P	G	P	N	O	N	R	N	R	R
Methyl ethyl diketone	P	G	P	N	O	N	R	O	R	R

*Can be used with aqueous solutions containing up to 20% DMSO.

R = Recommended, L = Limited Resistance, N = Not Recommended, O = Testing Advised, F = Fair, G = Good, P = Poor, PP = Polypropylene, PVC = Polyvinyl Chloride, CA = Cellulose Acetate, PC = Polycarbonate, PTFE = Polytetrafluoroethylene PS = Polystyrene, CN = Cellulose Nitrate, NY = Nylon, MCE = Mixed Cellulose Esters, PET = Polyethylene Terephthalate.

CHARACTERISTICS OF CORNING CENTRIFUGE TUBES

The following information is provided to serve as a general guideline for determining suitability of Corning centrifuge tubes for your applications. In addition, Corning recommends following the procedures outlined by the centrifuge manufacturer, as well as conducting a trial run to determine proper conditions before beginning any critical applications.

Corning centrifuge tubes are tested for leakage. They should not break or leak if used in a properly balanced rotor with suitable carriers, holders, and adapters that fully support the tubes when run in accordance with the guidelines in this section. These tubes are intended for one-time use only; reuse is not recommended as breakage or leakage may occur.

The recommended working temperature range for Corning centrifuge tubes is 0 to 40°C. The suitability of these tubes for storage below 0°C depends on both the solution and the

storage conditions. In general, the polypropylene and PET tubes are more resistant to stress at low temperatures than polystyrene. It is strongly recommended that a trial run be performed under actual conditions to test the suitability of the tubes for frozen storage.

Suggestions for Safe Centrifugation

- ▀ **Caution:** When centrifuging pathogenic organisms, clinical specimens known or suspected of being infectious, or any other potentially biohazardous materials, approved safety containment systems should be used. Contact your centrifuge manufacturer for appropriate accessories or recommendations.
- ▀ Read protocols and instruction manuals carefully. Do not confuse speed or revolutions per minute (RPM) with relative centrifugal force (RCF). Instructions for centrifuging a sample at a given RPM and time are incomplete unless the rotor or radius is specified. Protocols should always state the time and RCF value for centrifuging a sample.

- Proper balancing and distribution of the load in a centrifuge is critical for optimum performance and to prevent damage to the tubes or centrifuge. Opposing buckets or loads should always be balanced within the range specified by the manufacturer. Tubes should always be distributed in the buckets with respect to the center of rotation as well as the pivotal axis of the bucket. Failure to do this may prevent the bucket from achieving a horizontal position during the centrifugation run. Uneven separations or tube failure may result.

These centrifuge tubes are intended for use by persons knowledgeable in safe laboratory practices. Failure can result from surface damage, exceeding the specified RCF values, using unsuitable support systems, improper temperatures, or incompatible chemicals.

The RCF ratings for Corning® disposable centrifuge tubes have been established at room temperature using tubes filled to nominal capacity with water and spun in a horizontal rotor

centrifuge for 5 minutes. The centrifuge must be equipped with the recommended carriers, adapters, and cushions that fully support the tubes. If an angle head rotor is used or proper support is not provided, RCF values will be lower. Use of liquid other than water may also lower RCF values. Please consult your centrifuge specifications and the nomogram table (page 44) to determine speeds at which maximum RCF is achieved.

Chemical Compatibility of Disposable Plastic Centrifuge Tubes

The mechanical strength, flexibility, color, weight and dimensional stability of all plastic centrifuge tubes are affected to varying degrees by the chemicals with which they come in contact. Specific operating conditions, especially temperature, RCF, rotor type, carrier design, and run length will also affect tube performance.

Physical Properties of Disposable Plastic Centrifuge Tubes

	Clear Polypropylene	New Polyethylene Terephthalate
Recommended Working Temp*	0-40°	0-40°
Heat Distortion Point	121°	70°
Flexibility	Moderate	Rigid
Transparency	Clear	Clear
Maximum RCF:		
15 mL Tube	12,000 x g	3,600 x g
50 mL Tube	15,500 x g	3,600 x g
250 mL Tube	—	—
500 mL Tube	—	—

*At room temperature for 24 hours.

Chemical Resistance of Disposable Plastic Centrifuge Tubes*


Chemical Class	Polyethylene Terephthalate	Polypropylene	Polyethylene Caps
Acids (weak)	1	1	1
Acids	3	1	1
Alcohols	1	1	1
Aldehydes	3 ^a	2 ^a	1
Bases	3	1	1
Esters	2	2	2
Hydrocarbons:			
Aliphatic	1	2	3
Aromatic	3	3 ^b	3
Halogenated	3	3	3
Ketones	2	2 ^c	2

*At room temperature for 24 hours.

1 = Recommended; 2 = Suitable for most applications. However, a trial run under specific operating conditions is recommended; 3 = Not recommended.

Note: a = Formaldehyde, rated 1; b = Phenol, rated 1; c = Acetone, rated 1.

Nomogram for Computing Relative Centrifugal Force



To calculate the RCF value at any point along the tube or bottle, measure the radius, in mm, from the center of the centrifuge spindle to the particular point. Draw a line from the radius value on the right hand column to the appropriate centrifuge speed on the left-hand column. The RCF value is the point where the line crosses the center column. The nomogram is based on the formula:

$$RCF = (11.17 \times 10^{-7}) RN^2$$

where:

R = Radius in mm from centrifuge spindle to point in tube bottom

N = Speed of spindle in RPM

Fold out for

Innovative Cell Culture Surfaces for the 21st Century



Innovative Cell Culture Surfaces for the 21st Century

Corning® Surfaces

For over 30 years, Corning culture vessels have been modified using corona discharge and vacuum plasma to generate better surfaces for growing attached cells.

Today's new culture technologies, such as stem cells and tissue engineering, require new surfaces with new capabilities. Corning's investments in developing surface technologies are paving the way for these new cell culture applications. See for yourself why Corning is the first and only name to trust for surfaces that are backed with a performance guarantee.

Surfaces for Enhancing Cell Attachment

Corning CellBIND® Surface

The unique Corning CellBIND surface uses a patented microwave process for incorporating significantly more oxygen into the cell culture surface, rendering it better for cell attachment especially under difficult conditions.

- Quickly adapts cells to reduced serum or serum-free conditions
- Improves attachment and yield
- No special handling or storage required

Corning Labware with Ultra-Web™ Synthetic Surfaces

Innovative Ultra-Web synthetic nanofiber surfaces offer cells

a more *in vivo*-like 3-D fibrillar topography for cells where current surfaces do not allow the desired cell culture performance or function.

- Growing cells, such as stem cells, liver, neuronal and primary cultures
- Helps promote a more *in vivo*-like cell morphology not obtainable with plastic 2-D surface
- Easy to view and harvest cells using standards methods
- Compatible with cell-based luminescence-reporter gene and FLIPR calcium flux assays



Corning's research center at Sullivan Park, Corning, New York

Corning Microplates with Poly-D-Lysine Coated Surface

Corning Poly-D-Lysine (PDL) microplates are coated with PDL (molecular weight range of 70 to 150 kDa) giving the surface a net positive charge for better cell attachment.

- Improves differentiation of primary neurons, glial cells, neuroblastomas
- Enhances attachment of transfected cell lines, including HEK-293
- Helps cells stay attached during assay processing

Surfaces for Reducing or Preventing Cell Attachment

Corning Ultra-Low Attachment Coated Polystyrene Surface


The Corning Ultra-Low Attachment surface uses a covalently bound hydrogel layer to inhibit cell attachment.

- Growing primary cultures of tumor or adult stem cells as unattached spheroids
- Preventing anchorage-dependent cells, such as fibroblasts, from attaching and dividing
- Promoting embryoid body formation from ES cells


Corning Cell Culture Surfaces	Cell Culture Formats							
	Flasks	Dishes	Multiple Well Plates	Microplates	Roller Bottles	CellSTACK® Chambers	CellCube® Chambers	Culture Tubes
For enhancing cell attachment:								
Original Tissue Culture Surface	■	■	■	■	■	■	■	■
Corning CellBIND Surface	■	■	■	■	■	■	■	
Ultra-Web Surfaces		■		■				
Poly-D-Lysine Coated Surface				■				
For reducing or preventing cell attachment:								
Ultra-Low Attachment Surface	■	■	■	■		■		
Untreated Surface		■		■				■

For more information or product numbers, reference the format categories within the Cell Culture section of this catalog.


Microplates



1536 Well Echo™ Qualified Microplates,
page 62



1536 Well Black Clear Bottom Microplates,
page 61



384 Well Low Flange Microplate, page 58

OVERVIEW	46
SELECTING THE BEST MICROPLATE	47
96 WELL MICROPLATES	50
96 Well Assay Microplates	50
96 Well Polystyrene Stripwell™ Microplates	54
96 Well Polypropylene Storage Microplates and Cluster Tubes	55
384 WELL MICROPLATES	57
384 Well Assay Microplates	57
384 Well Polystyrene Microplates	58
384 WELL POLYPROPYLENE STORAGE MICROPLATES	60
1536 WELL MICROPLATES	61
PROTEIN CRYSTALLIZATION MICROPLATES	63
MICROPLATE ACCESSORIES	65
Lids	65
Storage Mats and Accessories	65
Robolids	66
Sealing Tapes	66
TECHNICAL APPENDIX	67
Corning Assay Surface Properties and Applications	67
Selected Corning Technical Literature	71

Overview

DESIGNED FOR PERFORMANCE

Corning has been setting the standard for excellence in life science labware for over 85 years. With our comprehensive line of plasticware, including assay products, we continue to be an industry leader. Corning strives for the highest standards in product design and plastics molding.

Corning Life Sciences microplates and accessories are manufactured under strict process controls guaranteeing consistent product performance. Our manufacturing facility is located in Kennebunk, Maine, registered to the ISO 9001 2000 standards. ISO registration is recognized worldwide as a standard of excellence for quality systems.

Customers can request a Certificate of Compliance for any Corning® microplate. Also available are detailed product descriptions and drawings that highlight product dimensions and testing procedures. All are available by contacting your local Corning Life Sciences office. See the back cover of this guide for a listing.

CORNING MICROPLATE EQUIPMENT COMPATIBILITY PROGRAM

The increasing use of automated laboratory equipment demands consumables that have been qualified for fit and function. Corning microplates are designed with automation compatibility in mind and meet industry standards. In addition, Corning has a comprehensive equipment compatibility program in which leading equipment manufacturers certify the compatibility of our microplates with their instruments.

For the most up-to-date information on equipment compatibility, Corning maintains a *Microplate Equipment Compatibility Guide* on our web site at www.corning.com/lifesciences. This on-line guide is searchable by instrument type, plate type, and by manufacturer name.

LIFE SCIENCES EARLY ACCESS TO DEVELOPMENT – CORNING'S L.E.A.D. PROGRAM

Corning is committed to meeting the rapidly evolving needs of the life sciences laboratory. We are continually developing new and innovative products that are compatible with the latest advances in technology and instrumentation. Our L.E.A.D. program gives researchers access to these products and special pricing prior to their full market release. Contact your local Corning Life Sciences office or representative for more information about the products currently available through this program.

SELECTING THE BEST CORNING® MICROPLATE FOR YOUR APPLICATION

Corning offers a range of microplates in a variety of well designs and sizes, polymer materials and colors, and surface treatments. This guide includes 96, 384, and 1536 well microplates. Information on Corning plates in lower density formats (e.g., 24 and 48 well plates) can be found in our on-line product catalog at www.corning.com/lifesciences.

There are three simple steps for selecting the best Corning microplate for your application:

- 1 Choose the Corning microplate format and well design
- 2 Choose the Corning microplate material and color
- 3 Choose the Corning surface treatment

1 Choose the Corning Microplate Format and Well Design

Corning microplate dimensions meet industry standards, ensuring compatibility with all microplate equipment and automation. Our microplates feature an A-1 corner notch design. The A-1 corner notch allows for quick visual orientation of plates when setting up automation runs, thereby reducing chances for robotics problems and lost productivity.

Corning microplates are available in several well shapes, optimized to meet different application requirements.

- **Flat bottom** – for bottom reading plate readers and cell culture applications
- **Round bottom** – for improved mixing and washing
- **V-bottom** – for easier removal of total well contents
- **Easy Wash™ bottom** (round to narrowed flat well bottom) – for improved washing in immunoassays

In addition, Corning offers Half Area microplates for the 96 well format and Low Volume microplates for the 384 well format. These microplates are ideal for assays using reduced working volumes and can provide savings in reagent and compound use.

Well Shape Selection Chart

Well Shape	Microplate Format						
	96 Well	96 Well Stripwell™	Half Area 96 Well	384 Well	Low Volume 384 Well	1536 Well	2 µL 1536 Well
Flat bottom	■	■	■	■	■	■	
Round bottom	■				■	■	■
V-bottom	■						
Easy Wash bottom	■						

Detailed information about well volume, working volumes, and plate dimensions for Corning 96, 384, and 1536 well microplates are provided throughout this guide.

2 Choose the Corning Microplate Material and Color

Corning uses different polymers for microplates to support various application requirements. Selection of the appropriate polymer material and color can improve assay performance. Additional technical information on key polymers can be found in the appendix at the end of this guide.

Material Selection Chart

Plate Material	Microplate Format						
	96 Well	96 Well Stripwell	Half Area 96 Well	384 Well	Low Volume 384 Well	1536 Well	2 µL 1536 Well
Clear polystyrene	■	■	■	■			■
Solid black or white polystyrene	■	■	■	■	■	■	■
Clear bottom black or white polystyrene	■		■	■	■	■	
Polypropylene	■			■			
Solid black or white polypropylene	■			■*			
Flexible vinyl (PVC)	■						
UV	■		■	■			

*Only available in black polypropylene

Corning® microplates are available in different materials:

- ▶ **Clear polystyrene microplates** are used for cell culture and colorimetric (absorbance) assays.
- ▶ **Black and white polystyrene microplates** can be used for fluorescent and luminescent assays. Solid black polystyrene plates are designed to reduce well-to-well crosstalk and background for fluorescent assays. Solid white polystyrene plates are designed to reduce well-to-well crosstalk, enhance luminescent signals and reduce background for luminescent assays. Both black and white plates are available with clear bottoms for use in cell-based assays and microscopy applications, and allow top or bottom reading capabilities.
- ▶ **Polypropylene microplates** are ideal for compound storage or assays that require high resistance to solvents including DMSO and ethanol. The Corning ClearPro™ 96 well microplate is also available and has greater clarity than standard polypropylene for easier visual inspection of samples.
- ▶ **Black and white polypropylene microplates** can be used for fluorescent and luminescent assays and reduce nonspecific binding problems observed with polystyrene plates. The polypropylene material is also highly resistant to many commonly used solvents.
- ▶ **Flexible vinyl (PVC) microplates** are economical, nonsterile general assay 96 well plates. Due to their flexible nature, these microplates are not compatible with automation.
- ▶ **UV microplates** allow UV absorbance readings with low background especially at 260 to 280 nm, and are ideal for determining protein or nucleic acid concentration.

3 Choose the Corning Surface Treatment

Corning offers polystyrene microplates with a variety of modified surfaces. These surfaces can support binding or covalent immobilization of cells, proteins, nucleic acids, and other biomolecules. Additional information on these surfaces can be found in the Technical Appendix at the end of this guide.

Surface Treatment Selection Chart

Surface Treatment	Microplate Format						
	96 Well	96 Well Stripwell™	Half Area 96 Well	384 Well	Low Volume 384 Well	1536 Well	2 µL 1536 Well
<i>For General Assay</i>							
Not Treated (medium binding)	■	■	■	■	■	■	■
High Binding	■	■	■	■	■		■
Nonbinding (NBS™)	■		■	■	■	■	
Sulfhydryl (Sulfhydryl-BIND™) Binding	■	■					
Carbohydrate (Carbo-BIND™) Binding	■	■					
Photo-reactive (Universal-BIND™) Binding	■	■					
Amine Binding		■					
<i>For Cell Culture</i>							
Tissue Culture (TC) Treated	■	■	■	■	■	■	■
Ultra-Low Attachment	■						
Corning® CellBIND® Surface	■			■			
Poly-D-Lysine	■			■	■		

Corning offers various surface treatments for microplates:

- ▶ **Not treated (or medium binding) polystyrene surface** is hydrophobic in nature and binds biomolecules through passive interactions. It is suitable primarily for the immobilization of large molecules, such as antibodies, that have large hydrophobic regions that can interact with the surface.
- ▶ **High binding surface** is capable of binding medium (>10 kD) and large biomolecules that possess ionic groups and/or hydrophobic regions.
- ▶ **Nonbinding surface (NBS)** is a Corning proprietary treatment technology used on polystyrene microplates to create a nonionic hydrophilic surface (polyethylene oxide-like) that minimizes molecular interactions. Ideal for reducing protein and nucleic acid binding at low concentrations, and increasing assay signal to noise.

- ▶ **Corning® CellBIND® Surface** is a Corning proprietary treatment which provides improved consistency and even cell attachment.
- ▶ **Tissue culture treated (TC-Treated) surface** is used for the attachment and growth of anchorage-dependent cells.
- ▶ **Ultra-Low Attachment surface** has a covalently bonded hydrogel designed to minimize cell attachment, protein absorption, enzyme activation and cellular activation. This surface is noncytotoxic, biologically inert and nondegradable.
- ▶ **Poly-D-lysine coated surface** can improve attachment of difficult-to-attach cells.
- ▶ **Sulfhydryl (Sulfhydryl-BIND™) binding surface** has covalently-linked maleimide groups that covalently couple to sulfhydryl groups via SH moieties. Ideal for assays requiring site-directed orientation of a biomolecule, especially antibodies.
- ▶ **Carbohydrate (Carbo-BIND™) binding surface** has hydrazide groups covalently coupled to carbohydrate groups. Ideal for assays requiring site-directed orientation of a biomolecule (oxidized antibodies, carbohydrates, and glycosylated proteins) while maintaining enzymatic or immunological activity.
- ▶ **Photo-reactive (Universal-BIND™) surface** covalently immobilizes biomolecules via abstractable hydrogens using UV illumination, resulting in a carbon-carbon bond. Although linkage is nonspecific and does not allow for site-directed orientation of a biomolecule, this surface may be useful for immobilization of double stranded DNA, antigens of unknown structure, and mixtures of biomolecules (e.g., cell lysates).
- ▶ **Amine surface** has positively charged amine groups (2×10^{13} reactive sites/cm²) that can be used for covalent immobilization via bifunctional crosslinkers.


BAR CODE CUSTOMIZATION

New

Generic Bar Codes

Corning now offers a line of generic bar coded plates to better meet the demands of your screening needs (see list of available plates on back).

- ▶ No lead time: Plates are in stock and ready to ship
- ▶ Surface identification: The surface treatment of the microplate is identified in the human readable portion of the bar code:
 - NT = Non Treated
 - TC = Tissue Culture Treated
 - NB = Non Binding Surface
- ▶ Labels applied to all 4 sides of the microplate to ensure usability regardless of scanner location
- ▶ Each microplate is specially treated to reduce the impact of static build-up
- ▶ Code 128 bar code format ensures compatibility with most bar code scanning and software systems



Generic Bar Code Microplate

Custom Designed Bar Codes

Corning will assist in designing and implementing a bar code label to meet your exact specifications. We can provide bar code label test samples at the front end of a project, to confirm decodability and ensure flawless performance in your end-use process. Our other customization features include:

- ▶ Superior print quality and resolution
- ▶ Flexible bar code label positioning
- ▶ Resistant to most commonly used organic solvents

Dependable Durability

Bar codes have been quality tested for optimal readability, chemical resistance, and temperature variation.

Expert Advice

Most Corning® microplates are suitable for bar code customization. Contact Corning Life Sciences or your local representative for more information.

96 Well Microplates

Corning offers a complete line of 96 well microplates for laboratory miniaturization and automation. These microplates are available for different applications:

- ▶ 96 well assay microplates
 - General assays – Not treated, NBS™, covalent binding, high binding, flexible vinyl (PVC), and UV microplates
 - Cell-based assays – Tissue culture treated, Corning® CellBIND® Surface, poly-D-lysine, and Ultra-Low Attachment polystyrene microplates
 - Immunoassays – EIA/RIA polystyrene plates (medium and high binding)
- ▶ 96 well polystyrene Stripwell™ microplates
- ▶ 96 well polypropylene storage microplates and cluster tubes

This selection guide does not include 96 well microplates for PCR and genomics. Please refer to the Corning Genomics Selection Guide for information on these products (page 73).

For additional microplate information, refer to *Selecting the Best Corning Microplate for Your Application* in the Overview section of this guide (page 47).

96 WELL ASSAY MICROPLATES


Corning offers a wide variety of assay microplates. They are organized into five groups:

- ▶ 96 Well Clear Polystyrene Microplates
- ▶ 96 Well Solid Black and White Polystyrene Microplates
- ▶ 96 Well Clear Bottom Black and White Polystyrene Microplates
- ▶ 96 Well UV Microplates
- ▶ 96 Well Clear Flexible Vinyl (PVC) Microplates

Corning 96 well polystyrene plates are offered in standard volume formats and in lower volume format (called Corning half area plates). Corning 96 well polystyrene microplates have plate dimensions (length x width x height) of 127.76 x 85.48 x 14.22 mm that meet proposed industry standards.

96 Well Plate Types	Well Bottom Shape	Total Well Volume (µL)	Recommended Working Volume (µL)
Standard 96 Well	Flat	360	75 to 200
Standard 96 Well	Round	330	75 to 200
Standard 96 Well	V	320	75 to 200
Standard 96 Well	Easy Wash™	360	75 to 200
Half area 96 Well, Solid	Flat	190	25 to 125
Half area 96 Well, Clear Bottom	Flat	205	25 to 125


96 Well Geometry and Dimensions



Corning® tissue culture treated microplates have the same surface treatment used on other Corning culture vessels. In addition to this traditional surface, Corning offers three additional surfaces: Corning® CellBIND® surface treatment for improving consistency and even cell attachment, a poly-D-lysine coating for enhancing attachment of difficult-to-attach cell lines, and an Ultra-Low Attachment surface for minimizing cell attachment.



96 Well Clear Microplates



96 Well EIA/RIA Microplates

Corning® CellBIND® Surface for Optimizing Cell-Based Assay Performance

- Available in 96 and 384 well black clear bottom microplates and 96 well clear solid microplates
- Surface treatment improves consistency and more even cell attachment, and may improve attachment of difficult-to-attach cell lines
- Not a coating, requires no special handling, and is stable at room temperature
- Sterilized by gamma radiation and certified nonpyrogenic

96 Well Clear Polystyrene Microplates

- Cell culture plates are sterilized by gamma radiation and certified nonpyrogenic
- Lids available where indicated (Information on lids and other microplate accessories can be found beginning on page 65.)


96 Well Clear Polystyrene Microplate Ordering Information

Cat. No.	Plate Format	Well Bottom	Surface Treatment	Sterile	Qty/Pk	Qty/Cs
3360	Standard Plate, no Lid	Round	TC-Treated	Yes	25	100
3366	Standard Plate	Round	High Bind	No	25	100
3367	Standard Plate	Round	Not Treated	Yes	1	50
3788	Standard Plate, with Lid	Round	Not Treated	Yes	20	100
3795	Standard Plate	Round	Not Treated	Yes	25	100
3798	Standard Plate	Round	Not Treated*	No	25	100
3797	Standard Plate	Round	Not Treated	No	25	100
3799	Standard Plate, with Lid	Round	TC-Treated	Yes	1	50
3894	Standard Plate, with Lid	V	TC-Treated	Yes	1	50
3896	Standard Plate	V	Not Treated	Yes	1	48
3897	Standard Plate	V	Not Treated	No	25	100
3898	Standard Plate	V	Not Treated*	No	25	100
2503	Standard Plate	Flat	Universal-BIND™	No	1	50
2507	Standard Plate	Flat	Carbo-BIND™	No	1	50
2509	Standard Plate	Flat	Sulfhydryl-BIND™	No	1	50
New 3300	Standard Plate, with Lid	Flat	Corning® CellBIND® Surface	Yes	5	50
3361	Standard Plate, with Lid	Flat	High Bind	Yes	20	100
3370	Standard Plate, with Lid	Flat	Not Treated	Yes	20	100
3474	Standard Plate, with Lid	Flat	Ultra-Low Attachment	Yes	1	24
3585	Standard Plate, with Lid**	Flat	TC-Treated	Yes	5	50
3590	Standard Plate	Flat	High Bind	No	1	100
3591	Standard Plate	Flat	Not Treated	No	1	50
3595	Standard Plate, with Lid**	Flat	TC-Treated	Yes	1	50
3596	Standard Plate, with Lid	Flat	TC-Treated	Yes	1	50
3598	Standard Plate, with Lid	Flat	TC-Treated	Yes	5	100
3599	Standard Plate, with Lid	Flat	TC-Treated	Yes	1	100
3628	Standard Plate, with Lid	Flat	TC-Treated	Yes	20	100
3641	Standard Plate	Flat	NBS™	No	25	100
3665	Standard Plate, with Lid	Flat	Poly-D-Lysine	Yes***	20	100
3997	Standard Plate, with Lid	Flat	TC-Treated	Yes	10	50
9017	Standard Plate	Flat	Not Treated	No	25	100
9018	Standard Plate	Flat	High Bind	No	25	100
3690	Half Area Plate	Flat	High Bind	No	25	100
3695	Half Area Plate	Flat	Not Treated	No	25	100
3696	Half Area Plate, with Lid	Flat	TC-Treated	Yes	1	50
3697	Half Area Plate, with Lid	Flat	TC-Treated	Yes	20	100
3368	Standard Plate	Easy Wash™	Not Treated	No	25	100
3369	Standard Plate	Easy Wash	High Bind	No	25	100

*Processed to improve hydrophilicity for hemagglutination and similar assays.

**Special low evaporation lid

***Aseptically manufactured



96 Well Black and White Polystyrene Microplates

96 Well Solid Black and White Polystyrene Microplates

- Designed to reduce well-to-well crosstalk
- White plates enhance luminescent signals and have low background luminescence and fluorescence
- Black plates have low background fluorescence and minimize light scattering

96 Well Solid Black and White Polystyrene Microplate Ordering Information


Cat. No.	Plate Format	Plate Color	Well Bottom	Surface Treatment	Sterile	Qty/Pk	Qty/Cs
3605	Standard Plate	White	Round	NBS	No	25	100
3789	Standard Plate	White	Round	Not Treated	No	25	100
3792	Standard Plate	Black	Round	Not Treated	No	25	100
3362	Standard Plate, no Lid	White	Flat	TC-Treated	Yes	25	100
3600	Standard Plate	White	Flat	NBS™	No	25	100
3650	Standard Plate	Black	Flat	NBS	No	25	100
3912	Standard Plate	White	Flat	Not Treated	No	25	100
3915	Standard Plate	Black	Flat	Not Treated	No	25	100
3916	Standard Plate, with Lid	Black	Flat	TC-Treated	Yes	20	100
3917	Standard Plate, with Lid	White	Flat	TC-Treated	Yes	20	100
3922	Standard Plate	White	Flat	High Bind	No	25	100
3925	Standard Plate	Black	Flat	High Bind	No	25	100
3990	Standard Plate	White	Flat	NBS	No	5	25
3991	Standard Plate	Black	Flat	NBS	No	5	25
3642	Half Area Plate	White	Flat	NBS	No	25	100
3686	Half Area Plate	Black	Flat	NBS	No	25	100
3688	Half Area Plate, with Lid	White	Flat	TC-Treated	Yes	20	100
3693	Half Area Plate	White	Flat	Not Treated	No	25	100
3694	Half Area Plate	Black	Flat	Not Treated	No	25	100
3875	Half Area Plate, with Lid	Black	Flat	TC-Treated	Yes	20	100
3992	Half Area Plate	White	Flat	NBS	No	5	25
3993	Half Area Plate	Black	Flat	NBS	No	5	25

96 Well Clear Bottom Black and White Polystyrene Microplates

- Bottoms are 60% thinner than conventional polystyrene plates, resulting in lower background fluorescence and enabling readings down to 340 nm
- Opaque walls prevent well-to-well crosstalk
- Optically clear flat bottom permits direct microscopic viewing

96 Well Clear Bottom Black and White Polystyrene Microplate Ordering Information

Cat. No.	Plate Format	Plate Color	Well Bottom	Surface Treatment	Sterile	Qty/Pk	Qty/Cs
3340	Standard Plate, with Lid	Black	Flat	Corning® CellBIND® Surface	Yes	5	50
3372	Standard Plate, with Lid	Black	Flat	Poly-D-Lysine	Yes	10	50
3601	Standard Plate	Black	Flat	High Bind	No	25	100
3603	Standard Plate, with Lid	Black	Flat	TC-Treated	Yes	1	48
3604	Standard Plate	White	Flat	NBS	No	25	100
3610	Standard Plate, with Lid	White	Flat	TC-Treated	Yes	1	48
3614	Special Optics Plate, no Lid	Black	Flat	TC-Treated	Yes	25	100
3615	Special Optics Plate, with Lid	Black	Flat	Not Treated	No	25	100
3631	Standard Plate	Black	Flat	Not Treated	No	25	100
3632	Standard Plate	White	Flat	Not Treated	No	25	100
3651	Standard Plate	Black	Flat	NBS™	No	25	100
3666	Standard Plate, with Lid	White	Flat	Poly-D-Lysine	Yes*	20	100
3667	Standard Plate, with Lid	Black	Flat	Poly-D-Lysine	Yes*	20	100



96 Well Clear Bottom Black and White Microplates

96 Well Clear Bottom Black and White Polystyrene Microplate Ordering Information (Continued)

Cat. No.	Plate Format	Plate Color	Well Bottom	Surface Treatment	Sterile	Qty/Pk	Qty/Cs
3903	Standard Plate, with Lid	White	Flat	TC-Treated	Yes	20	100
3904	Standard Plate, with Lid	Black	Flat	TC-Treated	Yes	20	100
3995	Standard Plate	White	Flat	NBS™	No	5	25
3998	Standard Plate, with Lid	Black	Flat	Poly-D-Lysine	Yes	5	25
3682	Half Area Plate, with Lid	Black	Flat	Poly-D-Lysine	Yes	10	50
3721	Half Area Plate	Black	Flat	TC-Treated	Yes	5	25
3880	Half Area Plate	Black	Flat	Not Treated	No	25	100
3881	Half Area Plate	Black	Flat	NBS	No	25	100
3882	Half Area Plate, with Lid	Black	Flat	TC-Treated	Yes	20	100
3883	Half Area Plate	White	Flat	Not Treated	No	25	100
3884	Half Area Plate	White	Flat	NBS	No	25	100
3885	Half Area Plate, with Lid	White	Flat	TC-Treated	Yes	20	100
3886	Half Area Plate, no Lid	White	Flat	TC-Treated	Yes	25	100
3887	Half Area Plate, no Lid	Black	Flat	TC-Treated	Yes	25	100
3994	Half Area Plate	White	Flat	NBS	No	5	25

*Aseptically manufactured

Tip for Improving Optical Performance in Fluorescent Assays

Corning® Special Optics 96 Well Microplates have black walls with ultra thin, clear bottoms for sharp, clear images and minimal background in fluorescent assays.



96 Well UV Microplate – Certified DNase- and RNase-free

96 Well UV Microplates

The Corning® 96 well UV microplate has a UV-transparent well bottom and is ideal for determining protein and/or nucleic acid concentrations.

- Certified DNase- and RNase-free
- UV-transparent bottom is molded directly to an acrylic base for greater strength and maximum leak resistance
- Total well volume: flat bottom – 360 µL; recommended working volume of 75 to 200 µL
- UV half area microplate has well volume of 205 µL; working volume of 25 to 125 µL
- Allows UV absorbance readings with low background, especially at 260 to 280 nm
- Lids are available separately. (Information on lids and other microplate accessories can be found beginning on page 65.)

96 Well UV Microplate Ordering Information

Cat. No.	Plate Format	Well Bottom	Sterile	Qty/Pk	Qty/Cs
3635	Standard Plate	Flat	No	25	50
3679	Half Area Plate	Flat	No	25	50

96 Well Clear Flexible Vinyl (PVC) Microplates

- Untreated PVC microplates are economical plates for solution-based assays, serial dilutions, and general storage applications.
- Well volume of 250 µL (260 µL for V-bottom); working well volume of 50 to 150 µL
- Lids are not available.

96 Well Clear Flexible PVC Microplate Ordering Information

Cat. No.	Plate Format	Well Bottom	Sterile	Qty/Pk	Qty/Cs
2797	Standard Plate	Round	No	25	100
2897	Standard Plate	V	No	25	100
2595	Standard Plate	Flat	No	25	100

Tip for Reducing Reagent Use

Corning 96 Well Half Area Microplates can save on valuable reagents by reducing the amount of reagent needed per well, while still retaining the ability to be read in standard plate readers. These microplates have a recommended working volume of 25 µL to 125 µL and are available untreated or with tissue culture, high bind, or NBS treatment.



Stripwell Microplates



Standard vs. Low Volume

Stripwell Low Volume Microplates

Big Cost Savings!

- Save 70% or more on antibody costs
- Save 50% or more on reagent costs

Features

- Total well volume: 190 μ L
- Recommended working volume: 75 to 125 μ L
- Same height/path length as a standard strip
- Standard 96 well center-to-center spacing

Custom Stripwell Microplate Colors



96 WELL POLYSTYRENE STRIPWELL® MICROPLATES

Corning® Stripwell plates are designed for *in vitro* diagnostic assays. The flat bottom strips are designed to easily break apart and are pre-assembled in an “egg-crate” style strip holder that allows each individual well to be positioned back into the plate once broken.

- Stripwell plates have 96 well flat bottom polystyrene format
- Low volume and standard Stripwell microplates have well volumes of 190 μ L and 360 μ L, respectively
- 1 x 8 strips are designed to fit only one way into the strip holder, eliminating the chance of misorientation
- Accessories can be found beginning on page 65.

Stripwell Microplates Ordering Information

Stripwell Low Volume Microplates

Cat. No.	Color	Binding Property	Qty/Pk	Qty/Cs
2480	Clear	Medium	25	100
2481	Clear	High	25	100
2482	Black	Medium	25	100
2483	Black	High	25	100
2484	White	Medium	25	100
2485	White	High	25	100

Standard Stripwell Microplates

Cat. No.	Color	Binding Property	Qty/Pk	Qty/Cs
2592*	Clear	High	25	100
2593*	Clear	Medium	25	100
2580**	Clear	High	200	800
9102***	Clear	TC-Treated, Sterile	1	50
3913	White	Medium	25	100
3923	White	High	25	100
3914	Black	Medium	25	100
3924	Black	High	25	100

*Product has a certified surface chemistry

**Individual 1 x 8 Strips without frame, bulk packed

***Microplates individually packaged with lid

Surface Modified Stripwell Microplates, Clear

Cat. No.	Description	Surface Chemistry	Well Volume	Qty/Pk	Qty/Cs
2388	Amine	Amine	360 μ L	1	50
2504	Universal-BIND™ Surface	Universal	360 μ L	1	50
2506	DNA-BIND™ Surface	N-oxysuccinimide	360 μ L	1	50
2508	Carbo-BIND™ Surface	Hydrazide	360 μ L	1	50
2510	Sulfhydryl-BIND™ Surface	Maleimide	360 μ L	1	50

Strip Accessories

Cat. No.	Description	Sterile	Qty/Pk	Qty/Cs
2572	Strip Holder “egg crate”	No	5	20
2578	96 Well Strip Ejector	No	5	5

Color Coding

Corning offers customers the ability to color code their Stripwell microplates. Currently there are 14 colors available from which to choose on both our certified high and medium binding Stripwell plates. In addition to the clear strip, two other colors can be applied to the same plate. Color-coded Stripwell microplates are made to order and minimum order requirements do apply. If interested, please contact your local Corning representative.

96 WELL POLYPROPYLENE STORAGE MICROPLATES AND CLUSTER TUBES

96 Well Polypropylene Microplates and Storage Blocks

Corning polypropylene microplates offer both small volume and large volume (blocks) well formats to meet assay and storage requirements.


- Flat, round or V-shaped well bottom
- Feature uniform skirt heights for greater robotic gripping surface
- Solvent resistant polypropylene provides compatibility with many common organic solvents (e.g., DMSO, ethanol, methanol)
- Certified DNase- and RNase-free
- Available sterile or nonsterile
- Refer to the Microplate Accessories section for information about microplate accessory products including sealing tapes and mats.

96 Well Polypropylene Microplate Dimensions and Well Volumes

Well Shape	Total Well Volume (μL)	Well Depth (mm)	Well Diameter (mm)	Plate Dimensions (L x W x H) (mm)
96 Well Flat Bottom	360	10.67	6.86	127.76 x 85.48 x 14.22
96 Well Round Bottom	360	11.3	6.86	127.76 x 85.48 x 14.22
96 Well V-bottom	320	11.13	6.86	127.76 x 85.48 x 14.22
96 Well V-bottom, Expanded Volume	450	12.43	8.50	127.76 x 85.48 x 14.35
96 Well 0.5 mL Block	500	25.3	6.86	127.76 x 85.48 x 27.18
96 Well 1 mL Block	1000	39.9	6.86	127.76 x 85.09 x 41.66
96 Well 2 mL Block	2000	42.04	8.13 (width)	128.27 x 85.85 x 43.94



Corning ClearPro Microplate (Cat. No. 3371) has higher clarity than standard polypropylene plates and allows users to visually inspect their samples in each well.



96 Well Polypropylene Storage Blocks with Storage Mat

96 Well Polypropylene Microplate Ordering Information

Cat. No.	Plate Format	Color	Well Bottom	Sterile	Qty/Pk	Qty/Cs
3355	Standard Plate	White	Round	No	25	100
3356	Standard Plate	Black	Round	No	25	100
3359	Standard Plate*	Clear	Round	Yes	25	100
3365	Standard Plate*	Clear	Round	No	25	100
3371	Corning® ClearPro™ Plate	Clear	Round	No	25	100
3364	Standard Plate	Clear	Flat	No	25	100
3343	Expanded Volume Plate	Clear	V	No	10	50
3344	Expanded Volume Plate	Clear	V	Yes	10	50
3357	Standard Plate	Clear	V	Yes	25	100
3363	Standard Plate	Clear	V	No	25	100

*Upgraded features include: superior clear polypropylene, lowered perimeter ridge for improved sealing, and added rigidity and dimensional stability for improved automated handling.

96 Well Polypropylene Storage Block Ordering Information

Cat. No.	Plate Format	Well Volume	Well Bottom	Sterile	Qty/Pk	Qty/Cs
3958	1 mL Round Well Block	1 mL	Round	Yes	5	25
3959	1 mL Round Well Block	1 mL	Round	No	5	100
3956	0.5 mL Round Well Block	0.5 mL	V	Yes	10	50
3957	0.5 mL Round Well Block	0.5 mL	V	No	100	100
3960	2 mL Square Well Block	2 mL	V	Yes	5	25
3961	2 mL Square Well Block	2 mL	V	No	5	100



Cluster Tube Systems

96 Well Cluster Tubes

- ▶ Composed of 96 polypropylene tubes in a standard microplate format
- ▶ Feature 1.2 mL tubes that are available individually or in strips of eight tubes
- ▶ Polyethylene tube caps are available in 8-cap strips

96 Well Cluster Tube Ordering Information

Cat. No.	Format	Sterile	Rack	Qty/Pk	Qty/Cs
4401	Individual	No	No	960/Bag	960
4408	8-Tube Strip	No	No	120/Bag	120
4410	Individual	No	Yes	96/Rack	960
4411	Individual	Yes	Yes	96/Rack	960
4412	8-Tube Strip	No	Yes	12/Rack	120
4413	8-Tube Strip	Yes	Yes	12/Rack	120
4418	8-Cap Strip	Yes	No	12/Bag	120

384 Well Microplates

Corning offers a variety of 384 well microplates for high throughput assays and storage. Microplates are grouped by application:

- ▶ 384 well assay microplates
 - General assays – Not treated, NBS™, high binding, and UV microplates
 - Cell-based assays – Tissue culture treated, Corning® CellBIND® Surface, and poly-D-lysine coated polystyrene microplates
- ▶ 384 well polypropylene storage microplates

This selection guide does not include 384 well microplates for PCR and genomics. Please refer to the Corning Genomics Selection Guide (page 73) or web site (www.corning.com/lifesciences) for further information on these products.


For additional microplate information, refer to *Selecting the Best Corning Microplate for Your Application* in the Overview section of this guide (page 47).

384 WELL ASSAY MICROPLATES

Corning offers a wide variety of assay microplates. They are organized into five groups:

- ▶ 384 Well Clear Polystyrene Microplates
- ▶ 384 Well Solid Black and White Polystyrene Microplates
- ▶ 384 Well Clear Bottom Black and White Polystyrene Microplates
- ▶ 384 Well UV Microplates

For assays performed in reduced volumes, Corning 384 well low volume polystyrene plates are available in solid round bottom and in black clear bottom formats.



Low Volume 384 Well Solid Round Bottom Microplates

Unique well design for optimal assay performance

- ▶ Raised well bottom for higher sensitivity
- ▶ Raised rim for decreased wicking and contamination
- ▶ Round bottom for better Z factor and minimized trapped air
- ▶ Conical well molded in the shape of a light cone for efficiency


384 well plate types	Well Bottom Shape	Total Well Volume (μL)	Recommended Working Volume (μL)
Standard 384 well	Flat	112	20 to 80
Low volume 384 well, solid	Round	35	5 to 20
Low volume 384 well, clear bottom	Flat	50	5 to 40

- ▶ Corning 384 well polystyrene microplates have plate dimensions (length x width x height) of 127.76 x 85.48 x 14.22 mm that meet proposed industry standards



384 Well Geometry and Dimensions

Corning 384 well microplates for cell culture include tissue culture treated, Corning CellBIND Surface, and poly-D-lysine coated microplates. The tissue culture treated microplates have the same surface treatment used on other Corning cell culture vessels while the poly-D-lysine treatment improves attachment of anchorage-dependent cells. The new Corning CellBIND Surface treatment can provide improved consistency and even cell attachment.



384 Well Clear Microplates

384 Well Clear Polystyrene Microplates

- ▶ Total well volume of 112 µL; working well volume of 20 to 80 µL
- ▶ Cell culture plates are sterilized by gamma radiation and certified nonpyrogenic
- ▶ The 384 well universal optics NBS™ plate is manufactured using an advanced polymer with high clarity and improved chemical resistant properties.
- ▶ Lids available as indicated. (Information on lids and other microplate accessories can be found beginning on page 65.)

384 Well Clear Polystyrene Microplate Ordering Information

Cat. No.	Plate Format	Well Bottom	Surface Treatment	Sterile	Qty/ Pk	Qty/ Qty/Cs
3640	Standard Plate	Flat	NBS	No	25	100
New 3640BC	Standard Plate with Bar Code Labels	Flat	NBS	No	25	100
3662	Clear Plate, with Lid	Flat	Poly-D-Lysine	Yes*	20	100
3680	Standard Plate, with Lid	Flat	Not Treated	Yes	20	100
3700	Standard Plate	Flat	High Bind	No	25	100
3701	Clear Plate, with Lid	Flat	TC-Treated	Yes	20	100
3702	Standard Plate	Flat	Not Treated	No	25	100
New 3702BC	Standard Plate with Bar Code Labels	Flat	Not Treated	No	25	100
3723	Universal Optics Plate (Standard)	Flat	NBS	No	25	100


*Aseptically manufactured.

384 Well Solid Black and White Polystyrene Microplates


- ▶ Designed to reduce well-to-well crosstalk during fluorescent and luminescent assays

384 Well Solid Black and White Polystyrene Microplate Ordering Information


Cat. No.	Plate Format	Plate Color	Well Bottom	Surface Treatment	Sterile	Qty/ Pk	Qty/ Cs
3673	Low Volume Plate	White	Round	NBS	No	25	100
3674	Low Volume Plate	White	Round	Not Treated	No	25	100
3676	Low Volume Plate	Black	Round	NBS	No	25	100
3677	Low Volume Plate	Black	Round	Not Treated	No	25	100
3678	Low Volume Plate	Black	Round	High Bind	No	25	100
New 3570	Solid White Plate, with Lid	White	Flat	TC-Treated	Yes	10	50
New 3571	Solid Black Plate, with Lid	Black	Flat	TC-Treated	Yes	10	50
New 3572	Standard Plate, Low Flange	White	Flat	Not Treated	No	10	50
New 3573	Standard Plate, Low Flange	Black	Flat	Not Treated	No	10	50
New 3574	Standard Plate, Low Flange	White	Flat	NBS	No	10	50
New 3574BC	Standard Plate, Low Flange with Bar Code Labels	White	Flat	NBS	No	10	50
New 3575	Standard Plate, Low Flange	Black	Flat	NBS	No	10	50
New 3575BC	Standard Plate, Low Flange with Bar Code Labels	Black	Flat	NBS	No	10	50
New 3820	Low Volume Plate	Black	Flat	NBS	No	10	50
New 3821	Low Volume Plate	Black	Flat	Not Treated	No	10	50
New 3821BC	Low Volume Plate with Bar Code Labels	Black	Flat	Not Treated	No	10	50
New 3822	Low Volume Plate, with Lid	Black	Flat	TC-Treated	Yes	10	50
New 3824	Low Volume Plate	White	Flat	NBS	No	10	50
New 3824BC	Low Volume Plate with Bar Code Labels	White	Flat	NBS	No	10	50
New 3826	Low Volume Plate with Lid	White	Flat	TC-Treated	Yes	10	50
New 3826BC	Low Volume Plate with Lid and Bar Code Labels	White	Flat	TC-Treated	Yes	10	50




384 Well Solid Low Flange Microplates



384 Well Low Volume Solid Microplates



384 Well Clear Bottom Black and White Microplates



384 Well Low Volume Black Clear Bottom Microplates

384 Well Clear Bottom Black and White Polystyrene Microplates

- Suited for fluorescent and luminescent assays using either top or bottom detection plate readers

384 Well Clear Bottom Black and White Microplate Ordering Information

Cat. No.	Plate Format	Plate Color	Well Bottom	Surface Treatment	Sterile	Qty/Pk	Qty/Cs
3540	Low Volume	Black	Flat	Not Treated	No	10	50
3542	Low Volume, Clear Bottom Plate with Lid	Black	Flat	TC-Treated	Yes	10	50
3544	Low Volume	Black	Flat	NBS	No	10	50
3643	Low Volume	Black	Flat	Poly-D-Lysine	Yes	10	50
3653	Standard Plate	White	Flat	NBS	No	25	100
3663	Clear Bottom Plate with Lid	White	Flat	Poly-D-Lysine	Yes*	20	100
3664	Clear Bottom Plate with Lid	Black,	Flat	Poly-D-Lysine	Yes*	20	100
3655	Standard Plate	Black	Flat	NBS	No	25	100
3683	Clear Bottom Plate with Lid	Black	Flat	CellBIND® Corning® Surface	Yes	10	50
3706	Standard Plate	White	Flat	Not Treated	No	25	100
3707	Clear Bottom Plate with Lid	White	Flat	TC-Treated	Yes	20	100
3711	Standard Plate	Black	Flat	Not Treated	No	25	100
3712	Clear Bottom Plate with Lid	Black	Flat	TC-Treated	Yes	20	100
3985	Optical Imaging Plate with Clear Bottom and Lid	Black	Flat	TC-Treated	Yes	20	100
3895BC	Optical Imaging Plate with Clear Bottom, Lid and Bar Code Labels	Black	Flat	TC-Treated	Yes	20	100

*Aseptically manufactured

384 Well UV Microplate

- Offers certified performance at 260 to 280 nm
- Provides consistently low background and well to well uniformity
- Performance approaches that of quartz cuvettes. Certified DNase- and RNase-free

384 Well UV Microplate Ordering Information

Cat. No.	Plate Format	Well Bottom	Sterile	Qty/Pk	Qty/Cs
3675	Standard Plate	Flat	No	5	25

384 WELL POLYPROPYLENE STORAGE MICROPLATES

384 Well Polypropylene Storage Microplates

Corning polypropylene microplates offer both small volume and large volume (blocks) well formats to meet assay and storage requirements.

384 Well Polypropylene Microplate Dimensions and Well Volumes

Well Shape	Total Well Volume (µL)	Well Depth (mm)	Well Diameter (mm)	Plate Dimensions (L x W x H) (mm)
384 Well Low Volume Low Profile Plate	20	6.30	3.30	127.76 x 85.48 x 10.00
384 Well Round Bottom Plate	95	11.56	3.63	127.76 x 85.48 x 14.22
384 Well Round Bottom Block	180	25.11	3.63	127.76 x 85.48 x 27.81
384 Well V-Bottom Block	240	22.31	3.30*	127.76 x 85.48 x 24.73

*Width of square well.

- ▶ Resistant to many common organic solvents (e.g., DMSO, ethanol, methanol)
- ▶ Black polypropylene microplate (Cat. No. 3658) is ideal for fluorescent assays requiring solvent resistance
- ▶ Certified DNase- and RNase-free
- ▶ Refer to the Microplate Accessories section for information about microplate accessory products including sealing tapes and mats.

384 Well Polypropylene Microplate Ordering Information

Cat. No.	Plate Format	Well Bottom	Well Volume (µL)	Sterile	Qty/ Pk	Qty/ Cs
3656	Standard Plate, Clear	Round	95	Yes	25	100
3657	Standard Plate, Clear	Round	95	No	25	100
3658	Standard Plate, Black	Round	95	No	25	100
3672	Low Volume, Low Profile, Clear	Conical	20	No	10	50

384 Well Polypropylene Storage Block Ordering Information

Cat. No.	Plate Format	Well Bottom	Well Volume (µL)	Sterile	Qty/ Pk	Qty/ Cs
3964	384 Well Block, Clear	Round	180	Yes	5	25
3965	384 Well Block, Clear	Round	180	No	5	100
3342	384 Well Block, Clear	V	240	Yes	5	50
3347	384 Well Block, Clear	V	240	No	5	50



384 Well Polypropylene Storage Microplates


1536 Well Microplates

Corning® 1536 well microplates are our highest density microplates available for high throughput screening. The microplates conform to standard microplate footprint and dimensions. These microplates are offered in solid black and white polystyrene plates, with round or flat bottoms, and in black clear bottom formats.


Corning also offers a ultra-thin 1536 well microplate with a total well volume of 2 µL. This uniquely designed plate represents leading edge technology in assay miniaturization, with the length and width dimensions and microplate footprint meeting industry standards.

1536 Well Standard Polystyrene Microplates

- ▶ Total well volume of 10 µL for round well plates and 12.8 µL for flat bottom plates
- ▶ Recommended working volume of up to 8 µL
- ▶ Round well bottoms for reduced air entrapment and improved CVs and Z factor
- ▶ Raised well bottoms for higher sensitivity
- ▶ Flood reservoir on four sides to reduce instrument contamination
- ▶ Lids are available separately. Corning lid Cat. No. 3098 is compatible with these plates.
(Information on lids and other microplate accessories can be found beginning on page 65.)




1536 Well Solid Round Bottom Microplates



1536 Well Black Clear Bottom Microplates

1536 Well Polystyrene Microplate Ordering Information

Cat. No.	Plate Format	Plate Color	Well Bottom	Surface Treatment	Sterile	Qty/ Pk	Qty/ /Cs
3936	Standard Plate	Black	Round	Not Treated	No	10	50
3937	Standard Plate	White	Round	Not Treated	No	10	50
3724	Standard Plate	Black	Flat	Not Treated	No	10	50
New 3724BC	Standard Plate with Bar Code Labels	Black	Flat	Not Treated	No	10	50
3725	Standard Plate	White	Flat	Not Treated	No	10	50
New 3725BC	Standard Plate with Bar Code Labels	White	Flat	Not Treated	No	10	50
3726	Standard Plate with Lid	Black	Flat	TC-Treated	Yes	10	50
3727	Standard Plate with Lid	White	Flat	TC-Treated	Yes	10	50
New 3727BC	Standard Plate, with Lid and Bar Code Labels	White	Flat	TC-Treated	Yes	10	50
3728	Standard Plate	Black	Flat	NBS™ Surface	No	10	50
New 3728BC	Standard Plate with Lid and Bar Code Labels	Black	Flat	NBS Surface	No	10	50
3729	Standard Plate	White	Flat	NBS Surface	No	10	50
New 3729BC	Standard Plate with Bar Code Labels	White	Flat	NBS Surface	No	10	50
3891	Clear Bottom	Black	Flat	Not Treated	No	10	50
New 3891BC	Clear Bottom with Bar Code Labels	Black	Flat	Not Treated	No	10	50
3893	Clear Bottom with Lid	Black	Flat	TC-Treated	Yes	10	50
New 3893BC	Clear Bottom with Lid and Bar Code Labels	Black	Flat	TC-Treated	Yes	10	50
3895	Clear Bottom	Black	Flat	NBS Surface	No	10	50



1536 Well 2 µL Polystyrene Microplates

1536 Well 2 µL Polystyrene Microplates

- ▶ A variety of assays, including enzyme assays, receptor-ligand assays, and cell-based assays have been effectively performed in these plates.
- ▶ Recommended working volume of up to 1.5 µL
- ▶ The plates are demarcated in a 8 x 12 array with each square containing 16 wells
- ▶ Eight extra wells on both the left and right sides of the plate that can be used to run controls
- ▶ Series of notches that allow stacked plates to be easily separated from one another
- ▶ Lids are available separately, Cat. No. 3849. (Information on lids and other microplate accessories can be found beginning on page 65.)

1536 Well 2 µL Polystyrene Microplate Ordering Information


Cat. No.	Plate Format	Color	Well Bottom	Surface Treatment	Sterile	Qty/Pk	Qty/Cs
3850	Low Volume Plate	Clear	Flat	Not Treated	No	20	100
3851	Low Volume Plate	Black	Flat	Not Treated	No	20	100
3852	Low Volume Plate	White	Flat	Not Treated	No	20	100
3853	Low Volume Plate	Clear	Flat	TC-Treated	Yes	20	100
3854	Low Volume Plate	Black	Flat	TC-Treated	Yes	20	100
3855	Low Volume Plate	White	Flat	TC-Treated	Yes	20	100
3857	Low Volume Plate	White	Flat	High Bind	No	20	100
3858	Low Volume Plate	Clear	Flat	High Bind	No	20	100

1536 Well Echo™ Qualified Microplate

- ▶ Corning-Labcyte joint development delivers optimal acoustic performance on the Labcyte Echo 550 Compound Reformatter
- ▶ Plates lot tested and certified to meet performance specifications
- ▶ Enhanced flatness provides low intra- and inter-plate CVs
- ▶ Low flange base is designed for bar code customization and robotic handling


Corning® 1536 Well Echo Qualified COC Microplate Ordering Information

Cat No.	Description	Well Bottom	Surface	Sterile	Qty/Pk	Qty/Cs
3730	1536 Well Clear COC Plate	Flat	Not Treated	No	10	50



1536 Well Echo Microplate

Protein Crystallization Microplates



96 Well CrystalEX Microplates


- Corning® 96 and 384 well crystallization microplates are optimized for high throughput protein crystal growth and screening
- Designed for sitting drop applications
- Meet 96 and 384 well microplate standards for automation

Next Generation CrystalEX™ Microplates for 96 Well High Throughput Sitting Drop Protein Crystallization

- Conforms to SBS specifications for full compatibility in automated crystal screening
- Multiple formats and versatility for custom options to maximize crystal formation, identification and analysis, and harvesting
 - Choose from five unique protein well shapes
 - Available in two materials, including a special zero polarization polymer (PZero)
 - Options include 1, 3, or 5 protein wells per reservoir well
- PZero polymer is superior for zero background polarization and nonbirefringence
- COC polymer offers strong chemical compatibility and good optical clarity
- Reservoir numbers are embossed on each individual well for easy identification

Next Generation CrystalEX Microplate Designs


One reservoir well is flanked by either one, three, or five protein wells, with SBS-standard spacing between the centers of adjacent well clusters.



Alphanumeric markers in each well cluster for easy identification under the microscope.



Five different protein well shapes are available:




*The crystal cup facilitates collection and centering of the protein crystals after incubation.

Corning Next Generation CrystalEX Microplates Ordering Information

Cat. No.	Protein Well Size	Protein Well Shape	No. of Protein Wells	Material	Treated	Qty/ Pk	Qty/ Cs
3556	4 µL	Round	1	COC	No	10	50
3551	4 µL	Conical flat	1	COC	Yes*	10	50
3840	2 µL	Conical flat	3	COC	No	10	50
3552	2 µL	Round	3	PZero	No	10	50
3553	2 µL	Conical flat	3	PZero	No	10	50
3554	2 µL	Flat	3	PZero	No	10	50
3555	2 µL	Conical flat – crystal cup	3	PZero	No	10	50
3550	1 µL	Conical flat – crystal cup	3	PZero	No	10	50
3557	1 µL	Conical flat – crystal cup	5	PZero	No	10	50

*Surface processed for hydrophilicity.



96 and 384 Well Protein Crystallization Microplates

96 Well CrystalEX™ Crystallization Microplates

- Features 96 large reservoir (reagent) wells and 96 corresponding protein wells
- Conical bottom protein wells allow for improved centering of the protein drop
- Compatible with manual pipettors and automation
- Novel merged well design provides efficient vapor space for protein crystallization


384 Well CrystalEX Crystallization Microplates

- Meets industry standards for 384 well microplate footprint and well locations
 - Ideal for fully automated crystal screening
- Features 192 reservoir wells and 192 corresponding protein wells
- Flat bottom protein wells are optimized for imaging of crystals
- Reservoir and protein wells are positioned to be compatible with multi-head dispensing equipment (up to 96 well heads)

96 and 384 Well CrystalEX Crystallization Microplate Ordering Information

Cat. No.	Plate Format	Reservoir Well Volume (μL)	Protein Well Volume (μL)	Sterile	Qty/Pk	Qty/Cs
3773	96 Well Plate, Conical Bottom	210	10	No	10	50
3785*	96 Well Plate, Conical Flat Bottom, Treated	210	7	No	10	50
3775	384 Well Plate, Flat Bottom	100	3.4	No	10	50

*Surface processed for hydrophilicity



96 Well Crystallization Microplate with Universal Optical Sealing Tape


Universal Optical Sealing Tape for Next Generation CrystalEX and CrystalEX Microplates

- High optical quality, pressure-sensitive tape ensures tight sealing to minimize evaporation
- Ideal for microscopic observation of crystals
- Suitable for use between -70°C and 100°C
- Compatible with commonly used aqueous solutions and organic solvents

Accessory for Next Generation CrystalEX and CrystalEX Microplates

Cat. No.	Description	Qty/Pk	Qty/Cs
6575	Universal Optical Sealing Tape	100	100

Microplate Accessories



Lids

Optimizing Sealing Conditions on Corning Polypropylene Microplates

Corning offers an application note (Corning Literature No. ALSP-AN-011) describing effective sealing with the ABgene® ALPS-100 automated plate sealer.



Corning Storage Mat Applicator

Lids

- ▶ All lids are made of rigid polystyrene except where indicated
- ▶ All lids have a corner notch on the A1 corner (except where indicated) to correspond to the corner notches found on all Corning® microplates
- ▶ The Universal Lid without a corner notch (Cat. No. 3098) does not need to be oriented in any particular direction to be placed on Corning plates. The lid also has a shorter skirt than standard lids
- ▶ The black Universal Lid (Cat. No. 3935) is suitable for fluorescent and other light-sensitive assays
- ▶ The DMSO-resistant cyclic-olefin lid (Cat. No. 3085) is tinted amber in color for light-sensitive assays and is 100% DMSO-resistant

Microplate Lid Ordering Information

Cat. No.	Description	Plate Compatibility	Sterile	Qty/ Pk	Qty/ Cs
3930	Low Evaporation Lid with Corner Notch and Condensation Rings	96 well microplates only (not 2 mL block)	Yes	1	100
3931	Low Evaporation Lid with Corner Notch and Condensation Rings	96 well microplates only (not 2 mL block)	Yes	25	50
3098	Universal Lid without Corner Notch	All microplates	Yes	25	100
3099	Universal Lid with Corner Notch	All microplates	Yes	25	50
3935	Black Universal Lid with Corner Notch	All microplates	Yes	25	50
3085	DMSO-resistant Cyclic-olefin Lid without Corner Notch	All microplates	No	25	50
3849	1536 Well 2 μ L Lid	2 μ L 1536 Well Microplates only	Yes	20	100

Storage Mats and Accessories

- ▶ Multiple formats are offered for specific and precise fit on 96 and 384 well plates and blocks
- ▶ Storage Mats Cat. Nos. 3080 and 3083 are manufactured from DMSO-resistant EVA (ethyl vinyl acetate) polymer
- ▶ Certified DNase- and RNase-free
- ▶ Can be applied manually or with Storage Mat Applicator

Storage Mats and Accessories Ordering Information

Cat. No.	Description	Sterile	Qty/Pk	Qty/Cs
3080	Round Well Storage Mat for 96 Well Plates and Blocks	No	25	100
3083	Square Well Storage Mat for Corning 2 mL Square Blocks	No	1	50
3346	Storage Mat for Expanded Volume 96 Well Microplates	No	10	50
3341	Storage Mat for 384 Well V-Bottom Blocks	No	10	50
3081	Storage Mat Applicator	N/A	1	1



96 and 384 Well Robolids


Robolids

- Combines the sealing ability of a storage mat with the rigidity of a plastic lid
- Designed for repeated application and removal by automation and to prevent short-term evaporation
- Silicone sealing plugs for organic solvent resistance and low extractables
- Can be used manually or with automation


Robolid Ordering Information

Cat. No.	Description	Sterile	Qty/Pk	Qty/Cs
3090	96 Well Robolid with Corner Notch	No	25	50
3089	384 Well Robolid with Corner Notch	No	25	50

Moisture Exchange with Corning® Robolids



Sealing Tapes



Sealing Mats and Tapes

- Easy application and removal for short- and long-term storage
- Provide tight seal to minimize evaporation and condensation
- Acetate Sealing Tape (Cat. No. 3095) is suitable for use between -16°C and 38°C, is transparent, and is not pierceable
- Aluminum Sealing Tape (Cat. No. 6569, 6570) is suitable for use between -80°C and 150°C, is not transparent, and is pierceable
- Breathable Sealing Tape (Cat. No. 3345) allows gas exchange across the surface
- Universal Optical Sealing Tape (Cat. No. 6575) is suitable for use between -70°C and 100°C, and is transparent

Sealing Tape Ordering Information

Cat. No.	Description	Sterile	Qty/Pk	Qty/Cs
3095	Acetate Sealing Tape for all Microplates	No	100	100
6524	Polyethylene Sealing Tape	No	100	100
6570	Aluminum Sealing Tape for 96 Well Microplates	No	100	100
6569	Aluminum Sealing Tape for 384 Well Microplates	No	100	100
3345	Breathable Sealing Tape	Yes	50	500
6575	Universal Optical Sealing Tape	No	100	100

Technical Appendix


Corning® Assay Surface Properties and Applications

Corning Surface	Applications	Binding Interaction	Sample Properties	Performance Criteria
FOR BIOCHEMICAL ASSAYS				
NBS™ coated polystyrene surface	<ul style="list-style-type: none"> SPA assays Homogeneous assays 	None – Inhibits hydrophobic and ionic interactions	Significantly reduces (<2 ng/cm ²) protein and nucleic acid binding	95% reduction of nonspecific binding of protein compared to untreated polystyrene
Medium Binding (Untreated) modified polystyrene surface	<ul style="list-style-type: none"> Homogeneous and heterogeneous assays 	Hydrophobic	Large biomolecules >20kD with large or abundant hydrophobic regions	96w clear: Well to well CV ≤ 5% 96w black: Well to well CV ≤ 15% (HT) Well to well CV ≤ 3% (HO) 96w white: Well to well CV ≤ 8% (HT) Well to well CV ≤ 5% (HO) 384w clear: Well to well CV ≤ 10% (HT) 384w black and white: Well to well CV ≤ 15% (HT) Well to well CV ≤ 5% (HO)
High Binding modified polystyrene surface	<ul style="list-style-type: none"> ELISA and other heterogeneous assays 	Hydrophobic and ionic (negatively charged)	Improves binding of medium to large biomolecules (>10kD) that are positively charged with or without hydrophobic regions.	96w clear: Well to well CV ≤ 3% 96w black: Well to well CV ≤ 8% 96w white: Well to well CV ≤ 10% 384w clear: Well to well CV ≤ 10% 384w black and white: Well to well CV ≤ 15%
Aminated-modified polystyrene surface	<ul style="list-style-type: none"> Used with bifunctional crosslinkers (i.e., glutaraldehyde, carbodiimide) to covalently couple to functional groups (primary amines, thiols, and carboxyls) on biomolecules. 	Hydrophilic and ionic (positively charged); allows covalent immobilization	Small negatively charged biomolecules OR biomolecules possessing an appropriate functional amine, carboxyl or thiol group.	CV ≤ 5% Percent Covalent Coupling ≥ 95%
DNA-BIND® modified polystyrene surface	<ul style="list-style-type: none"> Immobilization of aminated DNA for use in nucleic acid hybridization assays and solid-phase PCR Immobilization of peptides and other small primary amine containing molecules 	Allows covalent immobilization to amine groups via binding to NOS groups	Small to medium biomolecules, especially DNA, possessing an available amine group.	CV ≤ 15%
Sulphydryl-BIND™ modified polystyrene surface	<ul style="list-style-type: none"> Assays requiring site-directed orientation of a particular biomolecule, especially antibodies 	Allows covalent immobilization via SH moieties on maleimide groups	Biomolecules possessing an accessible sulphydryl group or reducible disulfide bond.	CV ≤ 15% Activated/NonActivated ≥ 2.0 Activated = reduced disulfide bonds
Carbo-BIND™ modified polystyrene surface	<ul style="list-style-type: none"> Assays requiring site-directed orientation of a particular biomolecule (oxidized antibodies, carbohydrates and glycosylated proteins) while maintaining enzymatic or immunological activity 	Allows covalent immobilization via binding to hydroxide groups	Biomolecules possessing carbohydate moieties available for periodate activation.	CV ≤ 15% Activated/Non-activated ≥ 3.0 Activated = periodate activation
Universal-BIND™ modified polystyrene surface	<ul style="list-style-type: none"> Immobilization of double-stranded DNA Immobilization of antigens of unknown structure (available functional groups unidentified) Immobilization of samples containing a mixture of biomolecules, such as cell lysate samples Immobilization of other nonproteinaceous molecules, such as glycolipids 	Allows covalent immobilization via UV cross-linking to abstractable hydrogen	Biomolecules with abstractable hydrogen.	CV ≤ 15% Activated/Non-Activated ≥ 2.0 Activated – by UV
FOR CELL-BASED ASSAYS				
Standard Tissue Culture Surface	<ul style="list-style-type: none"> Assays using standard attachment dependent cell lines 	Hydrophilic and ionic interactions (negatively charged)	Allows cell attachment and binding	≥95% confluency (attachment dependent cell line)
Corning® CellBIND® Surface	<ul style="list-style-type: none"> Assays for difficult to attach cells Help cells stay attached during washing steps 	Hydrophilic and ionic interactions (negatively charged)	Enhances cell attachment uniformity and binding to polystyrene	96 Well Plates: CV ≤10%; Wells with cells/wells without cells – 2X signal from MTS assay 384 Well Plates: CV ≤20%; Wells with cells/wells without cells – 2X signal from MTS assay
Poly-D-Lysine-Coated Surface	<ul style="list-style-type: none"> Assays for difficult to attach cells Help cells stay attached during washing steps 	Hydrophilic and ionic interactions (positively charged)	Enhances cell attachment and binding	96 Well: CV ≤15%; PDL/TCT ≤2.0 serum free HEK cells 384 Well: CV ≤ 20%; PDL/TCT ≥1.5 serum free HEK cells
Ultra-Low Attachment Surface	<ul style="list-style-type: none"> Assays where preventing cell attachment is required Hybridoma production and clonal isolation by limiting dilution 	Nonionic hydrogel layer reduces or eliminates ionic and hydrophobic binding	Prevents or reduces cell attachment and binding	≥95% cell attachment inhibition
Ultra-Web™ Surface	<ul style="list-style-type: none"> Assays where cell attachment or performance is enhanced by using a 3D surface 	Hydrophilic interactions on a 3D surface	Enhances cell attachment and performance	
Ultra-Web Polyamine Surface	<ul style="list-style-type: none"> Assays where cell attachment or performance is enhanced by using a 3D surface 	Hydrophilic and ionic interactions (positively charged) on a 3D surface	Enhances cell attachment and performance	

Technical Appendix (CONTINUED)

MICROPLATES

Corning® Ultra-Low Attachment Microplate (Cat. No. 3474) has a covalently bonded hydrogel layer to minimize cell attachment, protein absorption, enzyme activation and cellular activation. The surface is noncytotoxic, biologically inert, and nondegradable.



Comparison of Cell Attachment in Ultra-Low vs. Standard Tissue Culture Treated Plates
Vero cells plated at 2.6 x 10⁶ cells per well grown for 4 days at 37°C in a 5% CO₂ environment show a 99% reduction in cellular attachment vs. standard culture treated product.

High Binding Plate Certification of Corning EIA/RIA Microplates

Corning offers 96 well EIA/RIA plates and Stripwell™ microplates manufactured from a special medical grade polystyrene for uniform binding, high optical clarity, and low background absorption.

Certification Standards	High Binding	Medium Binding (Not Treated)
Well-to-well coefficient of variation (CV)	≤3%	≤5%
Average high and low wells from the mean	≤8%	≤15%
Background absorbance units from the mean	±0.005	±0.005

Corning high binding plates have a binding capacity of approximately 500 ng of mouse IgG/cm². The nontreated plates have a binding capacity of approximately 250 ng of Mouse IgG/cm². Corning tests its EIA/RIA plates on a lot-to-lot basis and the certification results for each lot are made available upon request by contacting your local Corning Life Sciences office. In addition, five ELISA Technical Bulletins are available at www.corning.com/lifesciences.


NBS Binding Performance

NBS microplates have a nonionic hydrophilic well surface, and are ideal for minimizing protein binding in homogeneous assays.

Binding in ng/cm ²	¹²⁵ I-IgG	¹²⁵ I-BSA	¹²⁵ I-Insulin	³² P-oligo DNA	³² P-λ phage DNA
Polystyrene	400	450	310	22	6
Polypropylene	380	440	370	3	<2
NBS on Polystyrene	<2.5	<2.5	5	<2	<2

Benefits of NBS™ on Homogeneous Assays


Fluorescence-based Assay Performance with Corning® NBS™ Low Volume Microplates



Higher Sensitivity for Fluorescence Polarization Assays with 384 Well Corning NBS Low Volume Microplates (Cat. No. 3676)

Data demonstrates *Streptomyces griseus* protease activity on BODIPY fluorescent labeled (FL) casein substrate. Protease activity is measured as a reduction in millipolarization (mP) units. A significant reduction in fluorescence polarization was observed at the lowest concentration of enzyme in a 10 µL volume.


Reduced Nonspecific Protein Binding with Corning NBS Microplates



NBS Surface Significantly Reduces Nonspecific Binding of a BODIPY FL Casein Substrate to Corning Microplates

Dilutions of BODIPY FL casein in digestion buffer were incubated for 30 min at room temperature in black Corning untreated and NBS microplates (Cat. No. 3654). Control wells contained digestion buffer only. Microplates were washed 3 times with PBS, pH 7.4, and 200 µL/well of digestion buffer alone was added to the wells. Fluorescence intensity was measured.


Miniaturization of Calcium Mobilization Assay in Corning 384 Well Low Volume Black Clear Bottom Tissue Culture Treated Microplate (Cat. No. 3542)



The chromatograms shown here are the rapid increase of calcium signals in Transfected CHO-K1 cells upon the addition of carbachol (n=5 wells). Transfected CHO-K1 cells of M1WT2 (ATCC, CRL-1984) were seeded at 5,000 cells per well in 10 µL medium and then grown in standard CO₂ incubator overnight (37°C). After the addition of 10 µL calcium dye solution per well, the plates were incubated in 37°C for 30 min. After equilibrating to RT for 30 min, plates were loaded to Flexstation reader (Molecular Devices, Inc.). Five µL of 50 µM carbachol solution (final concentration 10 µM) was transferred to induce the response (or 5 µL of plain buffer for the negative controls). The calcium signal was monitored for 60 sec. Assay was performed with Calcium 3 kit (Molecular Devices, Inc.).

Technical Appendix (CONTINUED)

Performance of Corning® 384 Well Poly-D-Lysine Microplate (Cat. No. 3664)




Corning 384 Well Poly-D-Lysine (PDL) Microplates

have over 60% more cell attachment capacity than those of a leading competitor. Comparison of cell attachment capacity with Corning PDL coated plates to competitor's PDL coated plates and uncoated TC-treated plates. BHK-21 cells (1×10^4 cells/well) were incubated in 25 μ L of DMEM F-12 media in 8 replicate wells for 1 hour (37°C, 5% CO₂) on 384 well black/clear bottom microplates.


Data provided by Sigma-Aldrich Corporation. Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.

Improved Kinase Performance with Corning 1536 Well Solid Black NBS Microplate (Cat. No. 3728)



The fluorescence of the fluorogenic substrate is quenched in this assay. Upon phosphorylation, the quenching mechanism is released, resulting in a significant increase in fluorescence intensity (FI), and therefore, kinase activity can be monitored continuously. The total reaction volume was 8 μ L and contained 20 mM Tris-HCl (pH 7.6), 5 mM MgCl₂, 5 mM DTT, 10% Lipid Activator, 6 μ M fluorogenic substrate, 10 μ M ATP and 50 pg/mL PKC β -II. Signals were measured by Acquest™ reader (Molecular Devices, Inc.). The PKC assay was developed by Applied Biosystems, Inc.

Performance of Corning 1536 Well 10 μ L Round Well Microplate (Cat. No. 3936)



Fluorescent Polarization Assay on Corning 1536 10 μ L Assay Microplate

10 ng/ μ L, 5 ng/ μ L and 1 ng/ μ L of *Streptomyces griseus* protease were incubated with 2.0 ng/ μ L of BODIPY FL casein substrate in 5 μ L volumes for 10 minutes at room temperature. (Corning 1536 Well 10 μ L black microplate, untreated, Cat. No. 3936).

Selected Corning Technical Literature

All literature is available in PDF file format at www.corning.com/lifesciences.

Assay Microplates

Binding Comparison of Polymer Surfaces: Introducing Non-Binding Surface Microplates

Corning® 96-well NBS™ microplates are ideal for homogeneous assays in high throughput screening. Studies of protein and nucleic acid binding to the NBS, when compared to polystyrene and polypropylene surfaces, demonstrate significant reduction in nonspecific binding.

Chemiluminescent HRP-Based Assay Using Corning White Microplate

A comparison of the performance of white microplates from several microplate manufacturers to that of Corning 96 well white microplate using a model HRP based luminescent assay system.

Corning Non-Binding Surface Microplates for Fluorescent HTS Assays

This 4-page technical note evaluates the efficacy of the Corning NBS microplate for use in a homogeneous fluorescence polarization protease assay.

Corning Non-Binding Surface Treatment to Reduce Non-Specific Binding To Microplates

This 2-page technical note evaluates Corning NBS microplates for Scintillation Proximity Assays.

Corning 384 Well Low Volume Microplate Performance in Miniaturized Assays

This technical note describes the performance of Low Volume microplates using a homogeneous fluorescence polarization assay at low volumes.

Design and Performance of the Corning 2 µL 1536 Well Plate

This 2-page technical note describes the design features and performance criteria for Corning 2 µL 1536 well microplates.

Fluorescent Polarization Kinase Assay Miniaturization in Corning 96 Well Half Area and 384 Well Microplates

This 4-page technical note examines assay miniaturization in Corning 96 well, 96 well Half Area, and 384 well microplates using fluorescence polarization tyrosine kinase assays.

Cell Culture Microplates

Helpful Hints to Manage Edge Effects of Cultured Cells for High Throughput Screening

This technical note is a compendium of techniques, collected from Corning Cell Culture facilities and customers, to reduce the occurrence of irregular patterns of cell adhesion or “edge effect” in microplates.

Poly-D-Lysine Coated Microplates

This 2-page application report describes binding and performance characteristics, and provides operating protocols for Corning’s poly-D-lysine microplates.

Immunoassay Microplates

Corning offers five ELISA Technical Bulletins:

- ▶ Immobilization Principles — Selecting the Surface
- ▶ Optimizing the Immobilization of Protein and other Biomolecules
- ▶ Effective Blocking Procedures
- ▶ Optimizing the Separation Step on 96 Well Microplates
- ▶ Selecting the Detection System – Colorimetric, Fluorescent, Luminescent

Storage Applications

Corning ClearPro™ 96 Well Polypropylene Microplates

This 4-page technical note describes the heat sealing and storage performance characteristics for Corning ClearPro microplates.

New Storage Mat Applicator System Meets Customers’ Strict Storage Requirements

This 2-page application note describes the performance characteristics of the Corning Storage Mat Applicator and the Corning products with which it is compatible.

Recommendations for Heat Sealing Corning Polypropylene Storage Products Using the ABgene® Automated Laboratory Plate Sealer

This 3-page application note describes the critical parameters for sealing Corning microplates with the ABgene Automated Laboratory Plate Sealer.

Fold out for

96 Well Microplate Selection Guide



CORNING 96 WELL MICROPLATE SELECTION GUIDE


[illegible]

Fold out for

384 Well Microplate Selection Guide



Genomics



96 and 384 Well UV Microplates, page 79

- OVERVIEW 74
- COLONY PICKING, BACTERIAL GROWTH, AND STORAGE 75
- PURIFICATION 76
- QUANTITATION AND DETECTION 79
- DNA AMPLIFICATION 81
- BAR CODE CUSTOMIZATION 84

Overview

FROM START TO FINISH – MEETING THE NEEDS OF THE GENOMICS LABORATORY

Corning's dedication to quality and technology has produced this comprehensive line of products for the genomics laboratory. Whatever aspect of research you are involved in – from culturing cells or microorganisms to printing and hybridizing DNA arrays, Corning's quality and breadth of line delivers reliable results. All of Corning's products are manufactured under stringent quality guidelines as an assurance of consistent performance from device to device and lot to lot. Featured in this brochure are our newest products for the high throughput genomics laboratory:

- ▶ Thermowell® Gold PCR reaction vessels for conventional and real-time PCR, and cycle sequencing
- ▶ 96 well half area UV plates for nucleic acid quantitation
- ▶ Low profile BioAssay dishes that are robotic friendly and maximize incubator and storage space

THE EQUIPMENT COMPATIBILITY PROGRAM

The increasing use of automated laboratory equipment demands laboratory disposables whose fit and function have been qualified. Our microplates are designed with automation compatibility in mind and they meet industry standards for plate dimensions. In addition, Corning Life Sciences maintains a comprehensive equipment compatibility program in which leading equipment manufacturers certify the compatibility of our products with their instruments. This information is continually updated with our new products as well as new instruments. For the most current information, visit our website: www.corning.com/lifesciences.

LIFE SCIENCES EARLY ACCESS TO DEVELOPMENT – THE L.E.A.D. PROGRAM

Corning is committed to meeting the rapidly evolving needs of the life sciences laboratory. We are continually developing innovative new products that are compatible with the latest advances in technology and instrumentation. Our L.E.A.D. program gives researchers access to these products and special pricing prior to their full market release. Contact your local Corning Life Sciences office or representative for more information about the products currently available through this program.

EXPERT ASSISTANCE IS JUST A TELEPHONE CALL OR E-MAIL AWAY

Customer service and technical representatives are available to answer any question – from pricing and product availability to protocols and applications advice. Our offices around the world are able to respond promptly to your inquiry regardless of your location. Contact us at your local office (listed on the back cover).

Colony Picking, Bacterial Growth, and Storage

245 mm Square BioAssay Dishes



245 mm Square BioAssay Dish

Square bioassay dishes are made from polystyrene and are certified nonpyrogenic. They are packed with lids and are designed with a stacking bead so that they will stack securely without slipping. The dishes are compatible with automated colony picking instruments.

Cat. No.	Description	Automation Compatibility	Qty/Pk	Qty/Cs
431111	245 mm x 245 mm, Square, 18 mm Deep Nontreated Dish, Sterile	PBA Flexys™ and the Genetix “Q” Bot® automated colony picking and gridding robots	4	16
431272	245 mm x 245 mm, Square, 18 mm Deep Nontreated Dish, Sterile	AutoGen AutoGenesys, BioRobotics BioPick, BioGrid, TAS and MicroGrid II high volume automated colony picking systems	4	16
New 431301	245 mm x 245 mm, Low Profile, Sterile, Nontreated Dish	PBA Flexys, Genetix “Q” Bot, BioRobotics, BioPick	5	20

96 and 384 Well Polypropylene Blocks for Growth and Storage



384 Well Polypropylene Blocks

96 and 384 well deep well blocks feature well designs for optimal liquid handling and are certified DNase- and RNase-free.

Cat. No.	Description	Well Shape	Sterile	Well Volume	Qty/Pk	Qty/Cs
3956	96 Well	Round V	Yes	0.5 mL	10	50
3957	96 Well	Round V	No	0.5 mL	10	100
3958	96 Well	Round	Yes	1 mL	5	25
3959	96 Well	Round	No	1 mL	5	100
3960	96 Well	V-Bottom	Yes	2 mL	5	25
3961	96 Well	V-Bottom	No	2 mL	5	100
3964	384 Well	Square-Round	Yes	180 µL	5	25
3965	384 Well	Square-Round	No	180 µL	5	100
3342	384 Well	Square V	Yes	240 µL	5	50
3347	384 Well	Square V	No	240 µL	5	50

Disposable Culture Flasks




Disposable Plastic Erlenmeyer Flasks


Disposable plastic Erlenmeyer flasks are made from optically clear polycarbonate and feature a wide, easy access mouth. The polycarbonate construction also delivers mechanical strength for shaker culture applications. Each flask is individually packaged and radiation sterilized. The polypropylene plug seal caps offer two positions: open to allow gas exchange or closed for a liquid-tight seal. The vent caps allow free gas exchange while offering a liquid-tight, contamination-free seal.

Cat. No.	Capacity (mL)	Graduations (mL)	Neck Diameter (mm)	Cap Style	Sterile	Qty/Pk	Qty/Cs
430421	125	25	26	Plug Seal	Y	1	50
431143	125	25	26	Vent Cap	Y	1	50
430183	250	25	31	Plug Seal	Y	1	50
431144	250	25	31	Vent Cap	Y	1	50
430422	500	50	43	Plug Seal	Y	1	25
431145	500	50	43	Vent Cap	Y	1	25
431146	1000	50	43	Plug Seal	Y	1	25
431147	1000	50	43	Vent Cap	Y	1	25

Purification



384 FiltrEX Filter Plates



96 FiltrEX Filter Plates

FiltrEX™ 96 and 384 Well Filter Plates

Corning® FiltrEX filter plates meet the industry standards for plate dimensions. The rigid side walls make the plate ideal for automation and the wide skirt accepts bar codes. Individual filter disks are encapsulated in the plate by a patented* process that ensures 100% integrity of each well. The design of the nozzle prevents sample cross-contamination and wicking. Glass fiber filter plates can be used for a variety of applications, such as plasmid isolation, DNA purification, PCR[†] clean-up or receptor/ligand binding assays. They are a cost-saving alternative to expensive DNA prep kits. Use the low-binding hydrophilic PVDF membrane for lysate clarification, protein kinase assays, or bead- or resin-based separation assays. Visit the Technical Information Center at our web site for additional application information.

384 Well FiltrEX Filter Plates

Cat. No.	Membrane	Pigment	Sterile	Well Volume (µL)	Qty/Pk	Qty/Cs
3531	0.45 µm PVDF	White	No	180	5	25
3533	0.66 mm Glass Fiber	White	No	180	5	25

96 Well FiltrEX Filter Plates

Cat. No.	Description	Sterile	Well Volume (µL)	Qty/Pk	Qty/Cs
3504	0.2 µm PVDF Membrane, Hydrophilic	No	350	10	50
3505	0.2 µm PVDF Membrane, Hydrophilic	Yes	350	10	50
3510	0.25 mm Glass Fiber Filter	No	350	10	50
3511	0.66 mm Glass Fiber Filter	No	350	10	50
3514	Fluid Guard for FiltrEX 96 Well Filter Plates	No	—	100	100

Please contact us for customized membranes.

*U.S. Patent No. 6,391,241

[†]PCR is covered by patents owned by Hoffman-LaRoche Inc., Nutley, NJ. Use of the PCR process requires a license.

Volume Adapter and Applicator

A volume adapter allows larger volumes (up to 1 mL) to be applied to the 96 well filter plates. The applicator easily assembles and disassembles the filter plate and adapter, and ensures a perfect, leak-free fit.

Cat. No.	Description	Qty/Pk	Qty/Cs
3584	Volume Adapter, Nonsterile	10	50
3507	Applicator	1	1

Collection Microplates


FiltrEX 96 and 384 well filter plates meet industry standards for plate dimensions and can be used with a broad range of collection plates. Polystyrene and polypropylene plates are available with a variety of well geometries. Commonly used collection plates are listed below. For information about other compatible collection plates, please contact us.

Cat. No.	Description	Well Volume (µL)	Qty/Pk	Qty/Cs
3371	96 Well, Round Bottom Polypropylene ClearPro™ Plate	360	25	100
3795	96 Well, Round Bottom Polystyrene Plate	360	25	100
3897	96 Well, V-Bottom Polystyrene Plate	320	25	100
3657	384 Well, Square Well, Round Bottom Polypropylene Plate	95	25	100
3965	384 Well, Square Well, Round Bottom Polypropylene Block	180	5	100
3702	384 Well, Flat Bottom Polystyrene Plate	125	25	100




3584 Volume Adapter

FiltrEX™ 96 and 384 Well Filter Plate Construction




The proprietary nozzle design and individual, integrally-sealed filter disks prevent filtrate cross contamination and wicking. The rigid construction and wide skirt allow for robotic handling and bar coding.

FiltrEX Filter Plate Performance




1 2 3 4 5 6 7 8

Agarose gel electrophoresis of Plasmid DNA prepared using Corning® FiltrEX 384 well filter plates.
Plasmid DNA samples isolated with glass fiber filter plates were separated in a 1% agarose gel in 1X TAE buffer. Lane 1 contains 10 µL of Hi-Lo™ markers (total DNA = 1 µg). Lanes 2-8 contain plasmid preparations isolated using 7 different wells of the filter plate. Plasmid DNA was recovered in approximately 55 µL total volume. Two microliters (2 µL) of recovered plasmid were loaded in each lane of the gel.

Marker	+	-	+	-	+	-	+	-	+	-	+	-	+
													

Integrally sealed wells eliminate sample cross contamination.
Alternating wells of negative control (-) or plasmid DNA (+) filtered through FiltrEX 384 well Filter Plates were analyzed for cross-contamination by PCR. PCR products were not detected in the negative control wells, indicating the absence of contaminating DNA.

Unpurified	Purified
	
	-1 kb
	-250 bp

PCR Clean-Up
PCR products were purified using a FiltrEX 384 well glass fiber filter plate. Primer dimers were efficiently removed with good recovery of the PCR products.



Spin-X Centrifuge Tube Filters

Spin-X® Centrifuge Tube Filters

Spin-X centrifuge tube filters consist of a membrane-containing filter unit within a centrifuge tube. They filter by centrifugation for bacteria removal, particle removal, HPLC sample preparation, removal of cells from media and DNA removal from agarose or acrylamide gels. Maximum RCF** is 16,000 x g.

Cat. No.	Membrane	Well Volume (working μL)	Pore Size (μm)	Tube Size (mL)	Qty/Cs
8160	CA	500	0.22	2.0	96
8161*	CA	500	0.22	2.0	100
8162	CA	500	0.45	2.0	96
8163*	CA	500	0.45	2.0	100
8169*	NY	500	0.22	2.0	200
8170*	NY	500	0.45	2.0	200

CA = Cellulose Acetate, NY = Nylon

*Indicates that the product is nonsterile and certified nonpyrogenic.

**RCF = Relative Centrifugal Force.

Spin-X Tube Purification of DNA from Agarose Gels

Introduction

Purification of DNA from an agarose gel with the Spin-X tube is quick and efficient, unlike the electroelution, dialysis, and “freeze-squeeze” methods. The Spin-X method consists of two simple steps: excision of the band from the gel and centrifugation in the Spin-X tube. Yields range from 30 to 80%.

Protocol*

1. Electrophorese DNA in an agarose gel containing ethidium bromide.
2. After electrophoresis, illuminate the gel under long wavelength UV light, then, using a sharp instrument, carefully excise the band of interest (30-15,000 bp).
3. Place the gel slice into the filter cup of the Spin-X tube (Cat. No's. 8160, 8161, 8162, 8163) and mix with 100 to 200 μL of distilled water or Tris-EDTA.
4. Spin the tube at about 13,000 x g for 5 to 20 minutes at room temperature.
5. Collect the DNA from the microcentrifuge tube; the agarose gel will be retained on the Spin-X membrane. If needed, ethanol precipitate the DNA to remove any EDTA present.

Note: DNA yield may increase with the incorporation of one or all of the following steps:


1. Macerate the gel slice prior to placement in the Spin-X tube.
2. Prior to centrifugation in step #4, freeze the gel slice at -70°C in a separate tube, then allow to thaw.
3. After the initial centrifugation, add an additional 200 μL of buffer to the Spin-X tube and centrifuge again.
4. Spin for a longer period of time.

*Schwarz, Herbert and Whitton, J. Lindsay, 1992. A Rapid, Inexpensive Method for Eluting DNA from Agarose or Acrylamide Gel Slices Without Using Toxic or Chaotropic Materials. Vol. 13, No. 2, Biotechniques.

Quantitation and Detection

96 and 384 Well UV Microplates



These plates have a unique UV-transparent bottom; ideal for determining protein and/or nucleic acid concentrations. The UV-transparent bottom is molded to the top without adhesives for greater strength and maximum leak resistance. Plates are certified for low background and consistent performance at 260 and 280 nm. Their broad linear detection range allows reliable detection of high and low concentrations of biomolecules.



96 Well Half Area, 96 and 384 Well UV Microplates


Cat. No.	Format	Bottom	Well Volume (μL)	Qty/Pk	Qty/Cs
3635	96 Well	Flat	370	25	50
3675	384 Well	Flat	125	5	25
3679	96 Well Half Area	Flat	205	25	50

384 Well UV Microplate Performance




Background absorbance of three samples of the Corning® 384 well UV plate bottom material.

This material features consistently low absorbance over a broad wavelength range, including well into the UV. The three samples showed negligible background absorbance across the entire visible spectrum (left panel) and very low background in the UV range (right panel).



Comparison of DNA detection with the Corning 384 well UV plate to individual samples read in a quartz cuvette.

For each indicated DNA concentration, triplicate 100 μL samples were read in a quartz cuvette with a Beckman DU® spectrophotometer. Six samples (90 μL) were read in a Corning 384 well UV plate at each concentration with a Tecan ULTRA™ reader. These sample volumes were chosen in order to maintain a 1 cm path length (smaller volume samples can be read in the UV plate). The Corning UV plate demonstrates a broad linear range enabling the reliable detection of high and low concentrations as well as good sample to sample consistency (CV's of <2% at 50 μg/mL DNA).



DNA-BIND Assay Microplates


DNA-BIND™ Assay Microplates

DNA-BIND surface covalently couples to amine groups, providing a convenient method to immobilize aminated single-stranded DNA by either the 5' or 3' end for hybridization, amplification, or other DNA-based assays. 96 well plates and 1 x 8 Stripwell™ plates come without lids. Protocols and application information are available on our web site:


www.corning.com/lifesciences.

Cat. No.	Format	Well Shape	Qty/Pk	Qty/Cs
2497	1 x 8 Stripwell Plate, White	Flat	1	50
2505	96 Well Plate, Clear	Flat	1	50
2506	1 x 8 Stripwell Plate, Clear	Flat	1	50
2498	96 Well Plate, Black	Flat	1	50
2499	96 Well Plate, White	Flat	1	50
2525	96 Well Plate, Clear	Flat	1	10


DNA-BIND Surface Performance




Reaction of N-oxy succinimide with an aminated biomolecule



Kinetics of oligonucleotide coupling



Detection of hybridization at increasing oligonucleotide concentrations



Effect of post-synthetic purification method on signal strength and background


DNA Amplification

Thermowell® GOLD PCR[†] Reaction Vessels from Corning

Thermowell GOLD PCR reaction vessels exemplify Corning's commitment to innovation: to develop superior quality, reliable, and versatile products to complement today's dynamic changes in technology. The wide variety of options offered by Thermowell GOLD provides researchers the choices necessary for complete compatibility with laboratory equipment. Look to Thermowell GOLD for PCR, sequencing, and real-time PCR.

Thermowell GOLD 384 Well Polypropylene PCR Microplates and Accessories

Thermowell GOLD 384 well PCR microplates feature exceptional dimensional stability following thermocycling, and are fully compatible with automation, commonly used thermal cyclers, and Applied Biosystems® sequencing adapters (see compatibility table).




Thermowell GOLD 384 Well PCR Microplates

Cat. No.	Description	Qty/Pk	Qty/Cs
3757	384 Well Polypropylene PCR Microplate, clear	10	50
3756	384 Well Polypropylene PCR Microplate, black	10	50
3699 ^a	Silicone Rubber Sealing Mat – 384 Well Microplates	1	25
6569	Aluminum Sealing Tape-384 Well Microplates	100	100
6575	Universal Optical Sealing Tape	100	100

^aThermowell Sealing Mats, available for 384 Well PCR Plates, are easy to apply and remove, fully autoclavable and reusable (at least five times). These silicone rubber mats offer a cost effective alternative to other sealing methods and provide 100% sealing when used in conjunction with clamp or screw-down heated lid thermal cyclers.

Thermowell GOLD and Thermowell 96 Well Polypropylene PCR Microplates and Accessories

Thermowell GOLD 96 well PCR microplates are offered in five formats to ensure maximum flexibility and a perfect match for your applications. The original Thermowell 96 well PCR microplates are universal fit and can be cut into 3 x 8 well segments.



Thermowell GOLD 96 Well Polypropylene PCR Microplates

Cat. No.	Description	Qty/Pk	Qty/Cs
6551	96 Well Microplate, Clear – Thermowell	25	25
3752	96 Well Microplate, Full Skirt, Clear – Thermowell GOLD	10	50
3751	96 Well Microplate, Full Skirt, Black – Thermowell GOLD	10	50
3753	96 Well Microplate, Half Skirt, Clear – Thermowell GOLD	10	50
3755	96 Well Microplate, Half Skirt, Black – Thermowell GOLD	10	50
3754 ^b	96 Well Microplate, Elevated Skirt, Clear – Thermowell GOLD	10	50

^bFully compatible with ABI 3700 and 3730.

Polycarbonate PCR Microplates

Cat. No.	Format	Model Name	Well Volume (µL)	Qty/Pk	Qty/Cs
6509	96 Well	Model P	200	1	25
6511	96 Well	Model M	200	1	25

Thermal Cycler Compatibility Guide for Polycarbonate PCR Microplates

Cat. No.	Name	Compatible Thermal Cyclers
6509	Model P	Applied Biosystems GeneAmp® PCR System 9600 ^c , Barnstead Thermolyne Amplitron II®, Techne® Cyclogene, and Gene E with 96 x 0.2 mL block
6511	Model M	MJ Research PTC-100-96V, PTC-200 DNA Engine™, Biometra Uno - Thermoblocker™, Coy Corporation Temp Cycler II, Corbett Research FTS-960, Hybaid OmniGene with Microblock, Quatro BioSystems T-C-40

^cRequires the use of the Spacer Block and Frame (Cat. No. 6527).

[†]PCR is covered by patents owned by Hoffman-LaRoche Inc., Nutley, NJ. Use of the PCR process requires a license.





Plate Well Dimensions



Sealing Tape, Sealing Mats, and Cap Strips for PCR



Thermowell GOLD PCR Tubes



Thermowell GOLD 8 Well PCR Tube Strips

PCR Sealing Tape and Sealing Mats

Sealing tapes prevent evaporation and enable oil-free operation when used with thermal cyclers with heated lids. The universal optical sealing tape can be used in detection coupled with PCR systems (real-time PCR).

Cat. No.	Description	Qty/Pk	Qty/Cs
6569	Aluminum Sealing Tape-384 Well Microplates	100	100
6570	Aluminum Sealing Tape-96 Well Microplates	100	100
6575	Universal Optical Sealing Tape for Real Time PCR	100	100
6555	Thermowell™ Sealing Mat-96 Well Microplates	1	25
3699	Thermowell Sealing Mat-384 Well Microplates	1	25
3087	Silicone Rubber Septa Mat	10	50

Thermowell® GOLD and Thermowell PCR Tubes

Individual PCR tubes are made of thin wall polypropylene and designed for precise fit in heat blocks to optimize heat transfer. Tubes are tested and certified to be free of DNase and RNase, are autoclavable at 121°C and withstand centrifugation to 10,000 x g.

Cat. No.	Volume (mL)	Cap Style	Color	Qty/Pk	Qty/Cs
<i>Individual Polypropylene PCR Tubes – Thermowell GOLD</i>					
3745	0.2	Flat	Clear	500	1,000
3744	0.2	Flat	Assorted	500	1,000
3747	0.2	Dome	Clear	500	1,000
3746	0.2	Dome	Assorted	500	1,000
3750	0.5	Flat	Clear	500	1,000
3749	0.5	Flat	Assorted	500	1,000

Individual Polypropylene PCR Tubes – Thermowell

Cat. No.	Volume (mL)	Cap Style	Color	Qty/Pk	Qty/Cs
6530	0.5	Flat	Clear	250	1,000
6531	0.2	Dome	Clear	96	960
6571	0.2	Flat	Clear	96	960

Thermowell GOLD and Thermowell 8 Well PCR Tube Strips

Tube strips consist of eight 0.2 mL thin wall polypropylene tubes connected together. Dual connectors between adjacent tubes eliminate inadvertent breakage during sample handling. Tube strips are designed for precise fit in thermal cyclers to optimize heat transfer. Thermowell GOLD cap strips are sold separately from Thermowell GOLD tube strips. Original Thermowell tube strips and cap strips are packaged together. Tube strips are tested and certified to be free of DNase/RNase contamination and are autoclavable at 121°C.




Cat. No.	Description	Qty/Pk	Qty/Cs
3741	0.2 mL 1 x 8 Tube Strips, Clear – Thermowell GOLD	125	1,250
3740	0.2 mL 1 x 8 Tube Strips, Assorted Colors - Thermowell GOLD	125	1,250
6542	0.2 mL 1 x 8 Tube Strips, Clear – Thermowell	60	300
6547*	0.2 mL 1 x 8 Tube Strips, Assorted – Thermowell	60*	300
3743	1 x 8 Cap Strips, Domed, Clear – Thermowell GOLD	125	1,250
3748	1 x 8 Cap Strips, Domed, Assorted Colors – Thermowell GOLD	125	1,250
3742	1 x 8 Optically Clear Flat Cap Strips, for RT-PCR** – Thermowell GOLD	125	1,250

*60 of each color per bag; 1 bag of each color per case.

**Optically Clear Flat Cap Strips are designed for real-time PCR. Suitable for use with Thermowell GOLD 0.2 mL 1 x 8 PCR tube strips and 96 well microplates.

See page 83 for Compatibility Guide and Volume Reference table.


Compatibility Guide for Thermocyclers, Sequencers, and Real Time PCR

		Thermowell® GOLD Microplates		
				
		96 Well Half Skirt	96 Well Full Skirt	384 Well
Thermal Cyclers				
Applied BioSystems®	GeneAmp® 9600	■		
	GeneAmp 9700	■		■
Biometra®	Uno®	■	■	
	Uno II®	■		■
	T1 Thermocycler®	■	■	■
	Tgradient®	■	■	
	Trobot®	■	■	■
Bio-Rad®	iCycler™	■		■
Eppendorf	MasterCycler®	■	■	
Ericomp	SingleBlock®	■		
	TwinBlock®	■		
	Deltacycler I®	■		
Flexi	Gene		■	■
	Genius		■	
ThermoHybaid	PCR Sprint®		■	
	PCR Express®	■	■	■
	MultiBlock System	■	■	■
	Touchdown®	■	■	■
	Omnigene®	■	■	
	Omn-E®	■	■	
MJ Research™	PTC 200 DNA Engine™	■	■	■
	PTC 225 DNA Tetrad®	■	■	■
	PTC 100®	■	■	
MWG™	Primus 96®	■	■	
	Primus 384®			■
Stratagene®	Robocycler®		■	■
TaKaRa	TP 240®		■	
	TP 3000®	■	■	
Techne®	Touchgene X®	■	■	■
RT-PCR Thermal Cyclers				
Applied BioSystems	ABI PRISM® 7000	■		■
	ABI PRISM 7700	■	■	■
	ABI PRISM 7900 HT	■		■
Bio-Rad	iCycler®		■	
Stratagene®	Mx 4000®		■	
Sequencers				
Applied BioSystems	ABI PRISM 3100	■		■
	ABI PRISM 3700	■		■
	ABI PRISM 3730	■		■
Amersham Biosciences	MegaBACE™ 500		■	
	MegaBACE 1000 Mark II		■	
	MegaBACE 4000			■
MJ Research™	BaseStation®		■	
	Transgenomic		■	

Thermowell® GOLD PCR Microplates Volume Reference Table

Format	Total Volume	Working Volume
384 Well PCR Microplates	55 µL	50 µL
96 Well PCR Microplates, Full Skirt	240 µL	200 µL
96 Well PCR Microplates, Half Skirt	340 µL	300 µL
96 Well PCR Microplates, Elevated Skirt	340 µL	300 µL

Bar Code Customization



Dependable Durability

Bar codes have been quality tested for optimal readability, chemical resistance and temperature variation.

What is a Bar Code*?

The same kind of bar codes you see in stores and supermarkets can be very useful to your lab. Consisting of a series of black bars and light spaces representing letters and/or numerals, a bar code is an easy-to-use vehicle for data collection. The specific arrangement of these bars and spaces follows strict rules known as a "symbology."

How Does a Bar Code Work?

Bar codes reflect spots of light into a scanner in varying amounts. These differences in reflection are translated into electrical signals by a light detector inside the scanner. The signals are converted into binary ones and zeros, which are used in various combinations to stand for specific numbers and letters.

Common Characteristics of a Bar Code

The Quiet Zones

The areas immediately adjacent to the beginning and the end of the bar code symbol. These zones define the parameters of the code. As a rule of thumb, zones should be 0.25" or larger to prevent misreads.

Start and Stop Characters

Found at the beginning and end of the bar code symbol. They tell the scanner from which direction information is being received.

Interpretation Line

Appears above or beneath a bar code where human readable information appears.

Corning, Beyond the Common Bar Code

- ▶ 2.75" x 0.25" label size
- ▶ Linear (1-D) bar codes: Code 128, Code 3 of 9, Interleaved 2 of 5
- ▶ 10 Mil Narrow Bar Element (X-dimension = 0.010")
- ▶ Multiple bar code labels on a single plate
- ▶ Label placement on any side of a Corning microplate
- ▶ Customer sequence is electronically stored and can be maintained even if plates or projects change.

Custom Designed Bar Codes

Corning will assist in designing and implementing a bar code label to meet your exact specifications. We will provide bar code label test samples at the front end of a project, to confirm decodability and ensure flawless performance in your end-use process. Our other customization features include:

- ▶ Flexible bar code and corresponding human readable layout/orientation on the bar code label, for compatibility with the internal bar code scanner inside your automated instruments
- ▶ Color coding
- ▶ Superior print quality and resolution
- ▶ Flexible bar code label positioning
- ▶ Resistant to most commonly used organic solvents

Expert Advice

Most Corning genomics products are suitable for bar code customization. Contact Corning Life Sciences or your local representative for more information.

*Information provided by Computype, Inc.

Microarray



UltraGAPS Coated Slides, page 89



Corning Cover Glass, page 92

OVERVIEW	86
SLIDE SELECTION CHART	87
MICROARRAY PRINTING	88
Epoxide Coated Slides	88
UltraGAPS™ Coated Slides	89
GAPS™ II Coated Slides	90
Pronto!™ Universal Spotting Solution	91
Pronto! Epoxide Spotting Solution	92
Corning® Cover Glass	92
384 Well Microarray Printing Plates	93
Microarray Slide Mailers/Storage Boxes	93
Microarray Storage Pouches	94
MICROARRAY PROCESSING	94
Pronto! Universal Hybridization Kits	94
Pronto! Background Reduction Kit	96
Hybridization Chambers	96

Overview

The quality and reliability of microarray results largely depend on the quality and consistency of both the glass substrate and the reagents used to manufacture and process the arrays. Corning has a history rich in science and technology, with expertise in glass and surface modification, optics, biochemistry and molecular biology, which has led to many innovations for life science research. Using this broad-based knowledge, Corning provides complete solutions to customers' complex problems and enables the achievement of breakthrough discoveries.

TOOLS FOR EVERY STEP IN THE PROCESS

- ▶ Premium glass substrates for printing microarrays
- ▶ Optimized reagents for the highest possible performance and control throughout the microarray process
- ▶ Storage products to facilitate the process and preserve sample content

SUPERIOR TECHNICAL AND APPLICATIONS INFORMATION

- ▶ Protocols included with every case of product, optimized from Corning's vast research experience
- ▶ Expert assistance that is just an e-mail or phone call away
- ▶ Field Applications support with a direct link to our Applications Group

UNPARALLELED MANUFACTURING PROCESSES

Corning® glass slides are manufactured using a proprietary coating process in a Class 100 cleanroom and undergo numerous quality control tests. Every slide is meticulously inspected for the presence of contaminating particulates, scratches and other defects before and after coating, ensuring a substrate of unmatched cleanliness, consistency, reliability and integrity.


The reagents in the Pronto!TM Microarray Kits are quality controlled to deliver consistency at every step in the process. They are optimized for use with Corning microarray slides, allowing the highest possible level of performance, standardization, and control.

Slide Selection Chart

Slide	Attachment Chemistry	Probe Types	Recommended Spotting Media	Applications
Epoxide	Covalent, Epoxysilane	Oligonucleotides	150 mM sodium phosphate, pH 8.5, 0.005% SDS	<ul style="list-style-type: none"> Transcriptional profiling SNP analysis
UltraGAPS™	Ionic, aminosilane	Double-stranded DNA	Pronto!™ Universal Spotting Solution 30 to 50% DMSO 3xSSC 150 mM sodium phosphate, pH 7.5	<ul style="list-style-type: none"> Transcriptional profiling Array CGH ChIP on Chip
GAPS™ II	Ionic, aminosilane	Proteins	20% Glycerol in PBS (ligand dependent)	<ul style="list-style-type: none"> Antibody screening Functional assays



Microarray Printing




Epoxide Coated Slides

Corning® Epoxide Coated Slides provide the optimal, uniform surface chemistry for covalent attachment of **unmodified or amino-modified short oligonucleotides** (~30-mer), as well as long oligonucleotides (>50-mer) and cDNA. Use Corning Epoxide Coated Slides with the Pronto! Universal Hybridization Kit (see page 94) to achieve the highest possible level of overall microarray performance.

- Versatility**
 - ▶ Ideal for short oligonucleotides, long oligonucleotides, and cDNA
 - ▶ Print with unmodified or amino-modified oligonucleotides
 - ▶ No UV crosslinking or baking step required for DNA coupling
- Reproducibility**
 - ▶ Minimal contribution to interarray variability (less than 3% CV)
- Sensitivity**
 - ▶ Detect 1 pg RNA spiked into 4 µg of total RNA sample
- Specificity**
 - ▶ Differentiate between 90% homologous oligos (3 mismatches in 30-mer oligonucleotides)

Epoxide Coated Slides Ordering Information

Cat. No.	Description	Slides/Pk	Slides/Cs
40040	Epoxide Slide Starter Kit (10 Epoxide Coated Slides, 0.8 mL Short Oligo and 0.8 mL cDNA/Long Oligo Hybridization Solution)	5	10
40041	Epoxide Coated Slides with Bar Code	5	25
40042	Epoxide Coated Slides without Bar Code	5	25
40043	Epoxide Coated Slides with Bar Code, Bulk Pack	25	25
40044	Epoxide Coated Slides without Bar Code, Bulk Pack	25	25



UltraGAPS™ Coated Slides




The Gamma Amino Propyl Silane surface on UltraGAPS Coated Slides is ideal for printing long (>50-mer) oligonucleotides, as well as cDNA. UltraGAPS Coated Slides have a more hydrophobic surface than competitors' slides, resulting in smaller, more consistent spot size. Each lot is tested for consistent spot morphology, signal intensity, and low background in a hybridization assay. Some of the applications for which UltraGAPS Coated Slides are ideally suited include: gene expression analysis, genotyping, and CGH (comparative genomic hybridization).

The Pronto!™ Universal Spotting Solution (see page 91) has been optimized for use with the UltraGAPS Coated Slides and provides excellent spot morphology for microarray printing. Use the Pronto! Universal Hybridization Kit (see page 94) in conjunction with these slides to achieve the highest level of microarray performance.

- Reproducibility
- Dynamic Range
- Manufacturing Excellence
- Minimal contribution to interarray variability (less than 5% CV)
 - Low background autofluorescence
 - Consistent spot morphology
 - Uniform surface treatment
 - Higher hydrophobicity

UltraGAPS Coated Slides Ordering Information

Cat. No.	Description	Slides/Pk	Slides/Cs
40015	UltraGAPS Coated Slides with Bar Code	5	25
40016	UltraGAPS Coated Slides without Bar Code	5	25
40017	UltraGAPS Coated Slides with Bar Code, Bulk Pack	25	25
40018	UltraGAPS Coated Slides without Bar Code, Bulk Pack	25	25
40019	UltraGAPS Slide Starter Kit (Includes 10 UltraGAPS Coated Slides, 5 mL Universal Spotting Solution)	5	10
40024	Pronto! Universal Validation Kit (Includes 10 UltraGAPS Coated Slides, 15 mL Universal Spotting Solution, Pronto! Hybridization Kit for 10 arrays)	5	10
40025	Pronto! Universal Printing Kit (Includes 25 UltraGAPS Coated Slides, 50 mL Universal Spotting Solution)	25	25



27,000 Feature Array on UltraGAPS Coated Slides (three magnifications)
A 27,000 feature array was printed on UltraGAPS slides, processed and hybridized using the reagents in the Pronto! Universal Hybridization Kit. The inset shows a magnification to highlight the low background, uniform spot morphology, and signal intensity of a the array.

Data courtesy of A. Borg, Ph.D., Lund University, Sweden.




GAPS™ II Coated Slides

GAPS II Coated Slides are manufactured from a proprietary ultraflat glass that enhances microarray performance, enabling more accurate reading of microarrays by confocal laser scanners. GAPS II Coated Slides are manufactured using the same coating process and attachment chemistry as the original GAPS amino-silane coated slides, enabling researchers to use the same protocols that they optimized for GAPS slides. Use GAPS II Coated Slides with the Pronto!™ Universal Hybridization Kit (see page 94) to achieve maximum microarray performance.


- Flexibility
- Binding Capacity
- Dynamic Range
- Recommended for both DNA and protein arrays
- High DNA retention for maximum signal strength
- Low background autofluorescence

GAPS II Coated Slides Ordering Information

Cat. No.	Description	Slides/Pk	Slides/Cs
40003	GAPS II Coated Slides with Bar Code	5	25
40004	GAPS II Coated Slides without Bar Code	5	25
40005	GAPS II Coated Slides with Bar Code, Bulk Pack	25	25
40006	GAPS II Coated Slides without Bar Code, Bulk Pack	25	25




GAPS II Coated Slide



Ordinary Silane Coated Slide

Spot Morphology on GAPS II Coated Slides

Note uniform spot morphology, high signal strength, and ultra-low background with the GAPS II Coated Slide. Images courtesy of Dr. John Quackenbush of the Institute for Genomic Research (TIGR), Rockville, MD.



Functional Peptide Array on GAPS II Coated Slides

The Caspase-3 substrate NH₂-DEVDA-Biotin was suspended in Corning® Epoxide Spotting Solution and printed in quadruplicate onto anhydride-derivitized GAPS II Coated Slides. Peptide arrays were incubated with avidin-Cy[®]3 in the absence or presence of Caspase-8 or Caspase-3 (as indicated), and scanned at 532 nm. The printed DEVDA peptide retained function on the array, as indicated by the reduced fluorescence seen in the spots treated with Caspase-3, but not Caspase-8. Note: GAPS II Coated Slides have also been used successfully for protein arraying without derivatization.

Data generated by Corning R&D.



Pronto!™ Universal Spotting Solution


Pronto! Universal Spotting Solution is optimized for both long oligonucleotides (>50-mer) and cDNA printed on UltraGAPS™ Coated Slides. The proprietary formulation provides excellent spot morphology and has an extremely low evaporation rate. Pronto! Universal Spotting Solution is available in bulk as well as part of both the UltraGAPS Slide Starter Kit and Pronto! Universal Printing Kit.

- ▶ Low background autofluorescence
- ▶ Low evaporation rate
- ▶ Ensures consistent DNA printing concentration
- ▶ Eliminates need for volume adjustments
- ▶ Provides for even distribution of spotted DNA across entire array

Pronto! Universal Spotting Solution Ordering Information


Cat. No.	Description	Qty/Pk	Qty/Cs
40027	Pronto! Universal Spotting Solution, 250 mL	1	1

Cat. No.	Description	Slides/Pk	Slides/Cs
40019	UltraGAPS Slide Starter Kit (Includes 10 UltraGAPS Coated Slides, 5 mL Universal Spotting Solution)	5	10
40025	Pronto! Universal Printing Kit (Includes 25 UltraGAPS Coated Slides and 50 mL Universal Spotting Solution)	25	25



Pronto! Universal Spotting Solution – Spot Uniformity

Quality control testing for Pronto! Universal Spotting Solution requires consistent spots when using 12 pins printed 25 times.



Pronto! Universal Spotting Solution – Low Evaporation

Pronto! Universal Spotting Solution evaporative losses are <5% over 4 hours, as compared to evaporative losses of >25% with other commercial spotting solutions.




Pronto!™ Epoxide Spotting Solution

Pronto! Epoxide Spotting Solution should be used for printing all types of DNA, including short oligonucleotides (~30-mer), long oligonucleotides (>50-mer), and cDNA printed on Corning® Epoxide Coated Slides. When used with Corning Epoxide Coated Slides, this spotting solution provides spot size control for printing high density arrays without contributing to background fluorescence. Pronto! Epoxide Spotting Solution is available in bulk (250 mL) or as part of the Corning Epoxide Slide Starter Kit.

- Provides controlled spot size for high density arrays
- No significant contribution to background fluorescence of arrays
- Low evaporation rate
- Enhanced spot morphology

Pronto! Epoxide Spotting Solution Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
40047	Pronto! Epoxide Spotting Solution, 250 mL	1	1
Cat. No.	Description	Slides/Pk	Slides/Cs
40040	Epoxide Slide Starter Kit (10 Corning Epoxide Coated Slides, 5 mL Pronto! Epoxide Spotting Solution, 0.8 mL Short Oligo Hybridization Solution)	5	10



Varying Spotting Solution Formulations to Adjust Spot Size

Pronto! Epoxide Spotting Solution formulation can be adjusted to alter spot size at will. DNAs were dissolved in Pronto! Epoxide Spotting Solution to which varying amounts of sodium dodecyl sulfate (SDS) had been added, and were printed in quadruplicate onto Epoxide Coated slides using 120 μm solid pins. The top row (80 μm feature diameter) had no addition, whereas adding increasing amounts of SDS resulted in correspondingly larger feature diameters.


Corning® Cover Glass

Corning Cover Glass is manufactured from special, optically clear glass. The cover glass is resistant to surface attack or weathering and will remain clear for extended periods of time. The flatness is controlled by a machine process resulting in a trouble-free fit to slides for a wettable and bubble-free mount.

The thickness of No. 1½ cover glass is 0.16 to 0.19 mm. Cover glass is packaged in plastic boxes for protection and convenience. Cover glasses in sizes and thicknesses other than those listed are available.

Cover Glass Ordering Information

Cat. No.	Description	Approx. Pcs/Oz	Qty/Cs
2870-22	Corning Cover Glass, Square, 22 x 22 mm, No. 1½	135	10 oz
2940-223	Corning Cover Glass, Rectangular, 22 x 30 mm, No. 1½	97	10 oz
2940-224	Corning Cover Glass, Rectangular, 22 x 40 mm, No. 1½	73	10 oz
2940-225	Corning Cover Glass, Rectangular, 22 x 50 mm, No. 1½	58	10 oz
2940-243	Corning Cover Glass, Rectangular, 24 x 30 mm, No. 1½	89	10 oz
2940-244	Corning Cover Glass, Rectangular, 24 x 40 mm, No. 1½	67	10 oz
2940-245	Corning Cover Glass, Rectangular, 24 x 50 mm, No. 1½	54	10 oz
2940-246	Corning Cover Glass, Rectangular, 24 x 60 mm, No. 1½	45	10 oz



384 Well Microarray Printing Plates

Corning® 384 well polypropylene microplates are available in both low and full volume well formats to meet source plate requirements for printing DNA content onto microarray slides. The plates are manufactured from solvent resistant, virgin polypropylene that is compatible with many organic solvents including DMSO. The plates feature rigid, full length skirts for full compatibility with automation.

The 384 Well Low Volume Microarray Printing Plate (Cat. No. 3672), with a working volume of 2 to 20 μL , has a conical V-bottom, square well geometry that provides for maximum sample recovery. The 384 Well Full Volume Storage Plate (Cat. No. 3656) has a total well volume of 95 μL .

- ▶ Well design provides for maximum sample recovery
- ▶ Resistant to many organic solvents including DMSO
- ▶ Certified DNase- and RNase-free
- ▶ Automation compatible

384 Well Microarray Printing Plates Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
3672	384 Well Microarray Printing Plate, Polypropylene, Low Volume	10	50
3656	384 Well Storage Plate, Polypropylene, Full Volume	25	100
6569	Aluminum Sealing Tape for 384 Well Microplates	100	100
3099	Universal Lid for 384 Well Microplates	25	50
3085	DMSO Resistant Lid for 384 Well Microplates	25	50


Microarray Slide Mailers/Storage Boxes

The plastic containers in which UltraGAPS™ Coated Slides are shipped also function as storage boxes for printed arrays. These containers are available as either 5 slide mailers or 25 slide storage boxes. The trans-orange plastic has low transmittance in the 500 to 700 nm wavelength range which helps protect Cy®3 and Cy®5 dyes from photobleaching. These rigid plastic containers do not shed particles or outgas volatile chemicals that may contaminate microarray slides.

The Corning® 25 Slide Storage Box has a lift off lid which is easy to open and close. The 5 Slide Mailer has a hinged lid that snaps closed tightly to prevent slides from accidentally falling out.

Microarray Slide Mailers/Storage Boxes Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
40082	5 Slide Mailer for Microarrays	50	50
40081	25 Slide Storage Box for Microarrays	10	20



Transmittance Through Corning Trans-Orange Slide Mailers

Low transmittance (500-700 nm) helps protect Cy3 and Cy5 from photobleaching.



Microarray Storage Pouches

Corning® Microarray Storage Pouches for 5- and 25-slide holders are the same pouches in which Corning UltraGAPS™ and Epoxide Coated Slides are shipped. These tear-resistant, foil-laminated pouches can be used to store and ship microarrays.

When heat-sealed, the pouches protect microarrays from light, humidity, and environmental contaminants. Each pouch comes affixed with a 3" x 4" white marking label.

Microarray Storage Pouches Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
40085	5 Slide Storage Pouch	50	50
40086	25 Slide Storage Pouch	50	50

Microarray Processing



Pronto!™ Universal Hybridization Kits

Pronto! Universal Hybridization Kits (Cat. Nos. 40026 and 40028) provide all of the reagents necessary to perform hybridizations of fluorescently labeled cDNA to microarrays printed on Corning Epoxide, UltraGAPS, or GAPS™ II Coated Slides. The Pronto! Universal Validation Kit (Cat. No. 40024) contains all of the reagents from above as well as 10 UltraGAPS Coated Slides and 15 mL of Universal Spotting Solution.




- | | |
|-----------------------------------|---|
| Pre-Soak Solution | ▶ Remove existing autofluorescence from printed microarrays |
| Pre-Hybridization Solution | ▶ Block background fluorescence during array hybridization |
| Hybridization Solutions | ▶ Solutions compatible with cDNA, long oligonucleotide, and short oligonucleotide content |
| | ▶ Ready to use (no dilution required) |
| | ▶ Contain blockers to increase specificity |
| Wash Solutions | ▶ Quality tested to ensure manufacturing consistency |

Pronto! Universal Hybridization Kit Ordering Information

Cat No.	Product	Reactions
40024	Pronto! Universal Validation Kit (Includes 10 UltraGAPS Coated Slides, 15 mL Universal Spotting Solution)	10
40028	Pronto! Universal Hybridization Kit	10
40026	Pronto! Universal Hybridization Kit	25
40030*	Pronto! Hybridization Kit without Pre-soak	25
40090	Pronto! cDNA Long Oligo Hybridization Solution, 20 mL	
40048	Pronto! Short Oligo Hybridization Solution, 4 mL	

*Larger volumes of reagents specially designed for use in automated hybridization stations.

A



Superior Reproducibility



Four separate 4K arrays were processed using the Pronto!™ Universal Hybridization Kit to demonstrate reproducibility (A). Differential gene expression patterns as represented by ratios of normalized Cy³/Cy⁵ in the cluster diagram were found to be very consistent between the four arrays (B). Interslide CVs were shown to be <10% for each of 3 separate experiments performed (C).

Pronto!™ Background Reduction Kit

The Pronto! Background Reduction Kit is designed to eliminate background autofluorescence and prepare printed arrays for hybridization. It also can be used as the final step in the array fabrication process. The strong reducing effect of this treatment leads to increased sensitivity and specificity by removing autofluorescent background due to oxidation. The kit includes liquid sodium borohydrate and 1 L of Pre-Soak Solution which provides enough reagents for the treatment of at least 50 arrays.

Pronto! Background Reduction Kit Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
40029	Pronto! Background Reduction Kit	1	1



With Presoak

Without Presoak

	With Presoak	Without Presoak
Background cutoff (RFU)	100.6	183.3
Number of features ≥ 2X background	1221	891

Detect Low Expressing Genes

Use of the Pronto! Background Reduction Kit results in the increased detection of low expressing genes (see table). Reduction of background autofluorescence is evident when 4K human arrays were processed using the presoak reagents in the Pronto! Background Reduction Kit. Arrays that were processed with the presoak reagents (left image) had a lower background detection cutoff than those processed without presoak (right image).

Hybridization Chambers


Corning® Hybridization Chambers are designed to hold microarray slides (25 x 75 mm) at constant humidity during hybridization incubations. The O-ring and retaining clips ensure that the reusable chambers remain watertight when submerged in waterbaths and airtight in hybridization ovens. Wells in the base hold 10 to 15 µL of water to maintain optimal interior humidity.

The original Corning Hybridization Chamber (Cat. No. 2551) provides the ideal interior height and volume for use with one slide of the standard 1 mm thickness and a standard coverglass. The Corning Hybridization Chamber II (Cat. No. 40080) has an increased interior depth which not only allows for single slide hybridizations, but also allows the user to place two arrays face-to-face and hybridize using a single labeled target. This chamber can also be used with raised-edge cover-slips (Erie Scientific M-Series Lifter Slips™) that are thicker or taller than standard thin coverglass.

Hybridization Chambers Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
2551	Hybridization Chamber	1	5
40080	Hybridization Chamber II with Increased Depth	1	5
40001	Replacement O-rings (fit both chambers)	5	5

Liquid Handling



Stripettes, page 99



Vacuum Filters, page 105



Cryogenic Vials, page 116

OVERVIEW	98
PIPETS	99
PIPETTING AIDS	100
PIPETTORS	100
PIPET TIPS	101
REAGENT RESERVOIRS	104
TRANSTAR-96® WELL LIQUID TRANSFER SYSTEM	105
ASPIRATOR	105
VACUUM FILTERS	105
SYRINGE FILTERS	108
SPIN-X® CENTRIFUGE TUBE FILTERS	108
STORAGE BOTTLES	109
CONTAINERS	110
CYLINDER	111
ERLENMEYER FLASKS	111
SPATULAS	112
CENTRIFUGE TUBES	113
MICROCENTRIFUGE TUBES	115
CRYOGENIC VIALS AND ACCESSORIES	116
TECHNICAL APPENDIX	119
Selecting the Best Filter for Your Application	119
Characteristics of Corning® Plasticware	122
Chemical Compatibility of Corning Plasticware	123
Characteristics of Corning Centrifuge Tubes	123

Overview

DESIGNED FOR PERFORMANCE

Corning Life Sciences offers a full line of liquid handling products that are manufactured under strict process controls guaranteeing consistent product performance. All Corning Life Sciences plastics manufacturing facilities are ISO 9001:2000 registered. ISO registration is recognized worldwide as a standard of excellence for quality systems.

In addition, customers can now request a Certificate of Quality for any Corning® or Costar® liquid handling product. Certificates are available at www.corning.com/lifesciences. This certificate details lot-specific information on component materials, sterility testing and pyrogen testing. Also available are detailed product descriptions and drawings that highlight product dimensions and testing procedures. All are available simply by calling your local Corning Life Sciences office.

NONPYROGENIC CERTIFICATION

Most Corning and Costar liquid handling products are certified non-pyrogenic with a documented endotoxin level of equal to or less than 0.1 EU/mL. Endotoxins have been shown to cause variability in cell culture. Nonpyrogenic certification is another way Corning helps ensure consistent cell culture results. Corning also offers a detailed technical bulletin on the effects of endotoxins in cell culture. This may be obtained by calling your local Corning Life Sciences office or by downloading the bulletin from the Corning web site www.corning.com/lifesciences.



Pipets




Stripette Serological Pipets



Clear Plastic Wrap



Paper/Plastic Wrap



Bulk Pack



Clean Room Packaging

Stripette® Serological Pipets

- ▶ Stripette pipets are sterile, nonpyrogenic, and RNase-/DNase-free.
- ▶ Exclusive antidrip tip assures accurate delivery.
 - Available in 25, 50 and 100 mL sizes
- ▶ Color-coded magnifier stripes make volume reading easier.
- ▶ Bidirectional graduations provide choice of ascending and descending scales
- ▶ Negative graduations allow additional working volume.
- ▶ Four packaging options:
 - Individually wrapped, clear plastic
 - Individually wrapped, paper/plastic
 - Bulk packed for large-scale sterile and nonsterile liquid handling applications
 - Clean room packed; individually wrapped, paper/plastic, triple bagged, SAL 10⁻⁶

Stripette Pipets Ordering Information

Cat. No.	Capacity (mL)	Graduations (mL)	Negative Grads. (mL)	Color Coded Stripe	Qty/Pk	Qty/Cs
<i>Individually Wrapped, Clear Plastic Wrap</i>						
4011	1	1/100	0.2	Yellow	100/bag	1,000
4012	1	1/100	0.2	Yellow	100/bag	200
4021	2	1/100	0.2	Green	100/bag	1,000
4051	5	1/10	2.5	Blue	50/bag	200
4101	10	1/10	3.0	Orange	50/bag	200
4492*	10	1/10	3.0	Orange	50/bag	200
4251	25	2/10	10.0	Red	50/bag	200
4501	50	1/2	10.0	Purple	25/bag	100
4484	100	1	N/A	Aqua	10/bag	100
<i>Individually Wrapped, Paper/Plastic Wrap</i>						
4485	1	1/100	0.2	Yellow	50/bag	1,000
4486	2	1/100	0.2	Green	50/bag	1,000
4487	5	1/10	2.5	Blue	50/bag	200
4488	10	1/10	3.0	Orange	50/bag	200
4489	25	2/10	10.0	Red	25/bag	200
4490	50	1/2	10.0	Purple	25/bag	100
4491	100	1	N/A	Aqua	10/bag	100
<i>Bulk Packed in Bags</i>						
4010	1	1/100	0.2	Yellow	50/bag	1,000
4020	2	1/100	0.2	Green	50/bag	1,000
4050	5	1/10	2.5	Blue	50/bag	500
4100	10	1/10	3.0	Orange	50/bag	500
4250	25	2/10	10.0	Red	25/bag	200
4500	50	1/2	10.0	Purple	25/bag	100
<i>Clean Room Pack, Individually Wrapped, Paper/Plastic, Triple Bagged</i>						
7041	1	1/100	0.2	Yellow	50/bag	1,000
7042	2	1/100	0.2	Green	50/bag	1,000
7045	5	1/10	2.5	Blue	50/bag	200
7015	10	1/10	3.0	Orange	50/bag	200
7016	25	2/10	10.0	Red	25/bag	200
7017	50	1/2	10.0	Purple	25/bag	100
7000	100	1	N/A	Aqua	10/bag	100

*Cat. No. 4492 features a wide tip for handling viscous fluids.



Aspirating Pipets

Aspirating Pipets

Aspirating pipets are sterile, ungraduated and unplugged polystyrene pipets for aspirating liquid using vacuum suction.

Aspirating Pipets Ordering Information

Cat. No.	Volume (mL)	Packaging	Qty/Pk	Qty/Cs
4975	1	Individually wrapped, bulk packed	50	1,000
9186	2	Individually wrapped, clear plastic wrap	100	1,000
9099	5	Individually wrapped, clear plastic wrap	50	200

Pipetting Aids



Stripettor Pipetting Aid

Stripettor™ Pipetting Aids

- ▶ Lightweight, adjustable speed control, and designed for use with all serological pipets
- ▶ Nose cones are autoclavable and have a replaceable 0.2 µm hydrophobic sterilizing filter
- ▶ Operates on a rechargeable 9V nickel hydride battery and features an LED light on the handle that lets the user know when to recharge
- ▶ Unit is fully operational while recharging

Stripettor Pipetting Aid Ordering Information

Cat. No.	Product Description	Qty/Cs
4910	Stripettor with sterile filter, rechargeable battery and recharger/adaptor	1
4911	Grommet replacement (silicone pipet holder)	1
4922	0.2 µm hydrophobic replacement filter	4
4923	0.2 µm hydrophobic replacement filter	25
4914	Recharger/adaptor for 4910	1

Pipettors



Lambda Single Channel Pipettor

Lambda® Single Channel Pipettor

- ▶ Corning® Lambda pipettors have a contoured handgrip and hook-style hand rest for greater comfort and less fatigue during prolonged use
- ▶ Quick-turn volume adjustment knob and easy-to-read digital volume display makes volume selection easier
- ▶ Volume ranges include 0.1 to 2 µL, 0.5 to 10 µL, 2 to 20 µL, 10 to 100 µL, 20 to 200 µL, and 100 to 1000 µL
- ▶ Bottom part of unit is autoclavable
- ▶ Backed by a three-year warranty

Lambda Single Channel Pipettor Ordering Information

Cat. No.	Volume Range (µL)	Qty/Cs
4959	0.1-2.0	1
4960	0.5-10	1
4961	2-20	1
4962	10-100	1
4963	20-200	1
4964	100-1,000	1
4958	Pipettor Stand	1



8 and 12-Pette Multichannel Pipettors

8-Pette® and 12-Pette® Multichannel Pipettors

- ▶ Costar® 8-Pette and 12-Pette multichannel pipettors feature a unique, ergonomic trigger-style aspiration and dispense control mechanism designed to reduce thumb fatigue during repetitive pipetting
- ▶ Volume range is 20 to 200 µL
- ▶ Volume is adjusted with a vernier-scale spindle
- ▶ Pipettors are entirely autoclavable

8-Pette and 12-Pette Multichannel Pipettors Ordering Information

Cat. No.	Volume Range (µL)	Channels	Qty/Cs
4880	20-200	12	1
4888	20-200	8	1

Pipet Tips



Universal Fit Pipet Tips

Universal Fit 200 and 1000 µL Pipet Tips

- ▶ Corning® universal fit tips are designed to provide a reliable fit with all major brand pipettors. (A Pipet Tip Compatibility Guide can be requested or downloaded from the Corning website.)
- ▶ Beveled orifice helps ensure accurate fluid delivery
- ▶ 1-200 µL universal fit tips are graduated at the 10, 50, and 100 µL volumes
- ▶ Select from three packaging options:
 - Racked tips are certified RNase-/DNase-free and nonpyrogenic
 - Stack rack tips feature a stack of five racks, each containing 96 tips, for a total of 480 tips in a space-saving design
 - Bulk packed tips are nonsterile and very economical

Universal Fit 200 and 1,000 µL Pipet Tips Ordering Information

Cat. No.	Volume Range (µL)	Format	Color	Sterile	Racks/Cs	Tips/Cs
<i>Racked Tips</i>						
4860	1-200	96 Tips/Rack	Yellow	Yes	10	960
4863	1-200	96 Tips/Rack	Natural	No	10	960
4864	1-200	96 Tips/Rack	Natural	Yes	10	960
4865	1-200	96 Tips/Rack	Yellow	No	10	960
4867	100-1,000	100 Tips/Rack	Blue	No	10	1,000
9032	100-1,000	100 Tips/Rack	Blue	Yes	10	1,000



Universal Fit Hinged Rack
Pipet Tips

Universal Fit 200 and 1000 μ L Pipet Tips (Continued)

Universal Fit 200 and 1,000 μ L Pipet Tips Ordering Information

Cat. No.	Volume Range (μ L)	Format	Color	Sterile	Racks/Cs	Tips/Cs
<i>Stack Rack Tips</i>						
4803	1-200	480 Tips/Stack Rack	Natural	No	10	4,800
4804	1-200	480 Tips/Stack Rack	Natural	Yes	10	4,800
4806	1-200	480 Tips/Stack Rack	Natural	No	2	960
<i>Bulk Packed Tips</i>						
4844	1-200	Bulk Pack	Natural	No	1,000	10,000
4845	1-200	Bulk Pack	Yellow	No	1,000	10,000
4862	1-200	Bulk Pack	Natural	No	1,000	1,000
4866	1-200	Bulk Pack	Yellow	No	1,000	1,000
4846	100-1,000	Bulk Pack	Blue	No	1,000	10,000
4868	100-1,000	Bulk Pack	Blue	No	1,000	1,000
<i>Universal Fit Hinged Rack Pipet Tips</i>						
4711	1-200 μ L	96 Tip Hinged Rack	Yellow	Yes	10	960
4712	1-200 μ L	96 Tip Hinged Rack	Yellow	No	10	960
4710	1-200 μ L	96 Tip Insert for Hinged Rack	Yellow	No	10 Inserts	960
4714	100-1000 μ L	100 Tip Hinged Rack	Blue	Yes	10	1,000
4713	100-1000 μ L	100 Tip Hinged Rack	Blue	No	10	1,000
4715	100-1000 μ L	100 Tip Insert for Hinged Rack	Blue	No	10 Inserts	1,000

Smart Rack Pipet Tip Refill System

- ▶ Corning® Smart Rack makes refilling pipet tip racks easier than ever
- ▶ Tips are contained on an autoclavable plastic reload card and transferred to a rack with a disposable reloading device (included)
- ▶ Compatible with many popular brand 200 μ L 96-tip racks
- ▶ Two configurations are available – 94-tip and 96-tip. In the 94-tip configuration, each reload card contains 94 tips and two corner anchoring pins that secure the card to the rack
- ▶ The 96-tip configuration does not include the corner anchoring pins.
- ▶ Smart Rack tips are nonsterile, autoclavable, RNase-/DNase-free, and nonpyrogenic

Smart Rack Pipet Tip Refill System Ordering Information

Cat. No.	Tip Volume (μ L)	Color	Tips/Pack	Packs/Cs	Tips/Cs
4786	200	Natural	940	5	4,700
4787	200	Natural	960	5	4,800

Pipet Tip Loading System


- ▶ The Corning pipet tip loading system makes reloading racks economical and effortless
- ▶ The system consists of a one-piece base that attaches to a “magazine” consisting of 10 layers of 96 tips
- ▶ Simply place the loader and magazine over an empty rack, lower the magazine, and “click” – the rack is loaded and ready for use
- ▶ Tip loading system tips are nonsterile, autoclavable, DNase/RNase-free, and nonpyrogenic

Pipet Tip Loading System Ordering Information

Cat. No.	Description
4780	Starter Kit, natural 1-200 μ L tips; includes 1 tip loader, 1 magazine with 960 tips and 10 empty racks
4781	Starter Kit, yellow 1-200 μ L tips; includes 1 tip loader, 1 magazine with 960 tips and 10 empty racks



Smart Rack Pipet Tip
Refill System



Pipet Tip Loading System

Pipet Tip Loading System Refill Magazines

- ▶ Tip loading system refills consist of magazines containing 10 layers of 96 tips
- ▶ Each magazine refills 10 racks

Pipet Tip Loading System Refill Magazines Ordering Information

Cat. No.	Tip Volume (µL)	Color	Tips/Magazine	Magazines/Cs	Tips/Cs
4783	1-200	Natural	960	5	4,800
4785	1-200	Yellow	960	5	4,800

IsoTip Filtered Pipet Tips

- ▶ IsoTip filtered pipet tips feature an inert, hydrophobic barrier that prevents aerosolized contaminants from coming in contact with pipettor shafts
- ▶ Ideal for applications where avoiding cross contamination is critical, such as DNA amplification and radioisotope handling
- ▶ Packaged sterile
- ▶ Certified RNase-/DNase-free and nonpyrogenic
- ▶ A Pipet Tip Compatibility Guide can be requested or downloaded from the Corning website.

IsoTip Filtered Pipet Tips Ordering Information

Cat. No.	Volume Range (µL)	Precise Fit	Tips/Rack	Racks/Cs	Tips/Cs
4801	0.1-2.0	Gilson® and other popular ultra-micropipettors	96	10	960
4807	0.2-10	Gilson and other popular ultra-micropipettors	96	10	960
4808	0.5-10	Eppendorf® and other popular ultra-micropipettors	96	10	960
4821	1-30	All popular research-grade pipettors	96	10	960
4823	1-200	All popular research-grade pipettors	96	10	960
4810*	1-200	All popular research-grade pipettors	96	10	960
4809	100-1,000	All popular research-grade pipettors	100	10	1,000

*Extended length tip designed for use with 2-20 µL, 10-100 µL, and 20-200 µL pipettors.

1 to 200 µL Gel-Loading Pipet Tips

- ▶ Corning® gel-loading pipet tips feature a capillary end that allows easy access into vertical and horizontal electrophoresis gels
- ▶ Total capacity of 200 µL
- ▶ Certified RNase-/DNase-free and nonpyrogenic
- ▶ Tips are 83 mm in length

1 to 200 µL Gel-Loading Pipet Tips Ordering Information

Cat. No.	Tip Shape	End Thickness (mm)	Sterile	Tips/Rack	Racks/Cs	Tips/Cs
4853	Round	0.5	No	200	2	400
4854	Flat	0.4	No	200	2	400
4884	Flat	0.2	No	200	2	400



IsoTip Filtered Pipet Tips



Gel-Loading Pipet Tips and Microvolume Gel-Loading Pipet Tips

Microvolume Gel-Loading Pipet Tips

- ▶ Corning microvolume gel-loading tips feature a capillary end for gel-loading and are designed for use with Gilson and other popular ultra-micropipettors
- ▶ Working volume of 0.2 to 10 μ L
- ▶ Certified RNase-/DNase-free and nonpyrogenic

Microvolume Gel-Loading Pipet Tips Ordering Information

Cat. No.	Tip Shape	End Thickness (mm)	Sterile	Tips/Rack	Racks/Cs	Tips/Cs
4815	Flat	0.2	No	200	2	400

Microvolume Pipet Tips

- ▶ Microvolume tips provide accurate, reliable performance in the 0.1-10 μ L range for major brand ultra-micropipettors
- ▶ All racked tips are certified RNase-/DNase-free and nonpyrogenic

Microvolume Pipet Tips Ordering Information

Cat. No.	Volume Range (μ L)	Fit	Sterile	Qty/Pk	Tips/Cs
4826	0.1-10	Gilson® and other popular ultra-micropipettors	No	96/rack	960
4894	0.1-10	Gilson and other popular ultra-micropipettors	Yes	96/rack	960
4840	0.1-10	Gilson and other popular ultra-micropipettors	No	1,000/bag	10,000
4830	0.5-10	Eppendorf® and other popular ultra-micropipettors	Yes	96/rack	960
4834	0.5-10	Eppendorf and other popular ultra-micropipettors	No	96/rack	960
4901	0.5-10	Eppendorf and other popular ultra-micropipettors	No	1,000/bag	10,000



Microvolume Pipet Tips


Reagent Reservoirs

Costar® Reagent Reservoirs are ideal for repetitively filling multichannel pipettors

- ▶ Manufactured from modified polystyrene
- ▶ Sterile
- ▶ Disposable

Reagent Reservoirs Ordering Information

Cat No.	Volume (mL)	Color	Qty/Pk	Qty/Cs
4870	50	Natural	5/bag	200
4871	50	Natural	1/bag	100
4872	100	White	5/bag	200
4873	100	White	1/bag	100



4870 and 4871
50 mL Reagent Reservoir



4872 and 4873
100 mL Reagent Reservoir

Transtar-96® Well Liquid Transfer System



Transtar-96 System

- ▶ The Costar® Transtar-96 System is a portable, autoclavable liquid handling device for use with 96 well plates
- ▶ A sterile 96 tip disposable cartridge, which loads into the Transtar system, enables liquids to be aspirated, transferred and dispensed over a volume range of 25 to 200 µL in 5 µL increments
- ▶ The Transtar-96 System is ideal for changing cell culture media and screening monoclonal antibodies
- ▶ Transtar-96 System accuracy is rated at ±5% at all volume levels

Transtar-96 Well Liquid Transfer System Ordering Information

Cat. No. Qty/Cs	Description	Sterile	Qty/Pk	
7605	Transtar-96, adjustable-volume pipettor	N/A	1	1
7606	Transtar elevator	N/A	1	1
7610	Transtar disposable cartridges	Yes	1	24
4876	Transtar disposable reservoir liner, open	Yes	1	100
4877	Transtar disposable reservoir liner, 12-channel	Yes	1	100
4878	Transtar disposable reservoir liner, 8-channel	Yes	1	100

The Costar aspirator is an aspirating device for safe liquid removal/disposal from a variety of laboratory vessels using standard disposable pipet tips.



Transtar Disposable Cartridge

Aspirator



Aspirator

Aspirator Ordering Information

Cat. No.	Description
4930	Aspirator device (includes hand piece, grommet for accessory attachment, and single-channel adapter for use with disposable pipet tips)
4931	8-channel adapter for use with disposable pipet tips

Vacuum Filters



Vacuum Filter

Corning offers a variety of filter systems, membranes, pore sizes, and materials. For help in selecting the best filter combination for your research, please refer to the Technical Appendix for *Selecting the Best Filter for Your Application* on page 119.

115 mL Vacuum Filters

- ▶ 60 mm diameter membrane
- ▶ Low center of gravity and wide base for stability
- ▶ Separate pour spout to remove filtered sample which minimizes contamination
- ▶ Individually packaged, sterile, certified nonpyrogenic

115 mL Vacuum Filters Ordering Information

Cat. No.	Membrane	Volume (mL)	Pore Size (µm)	Qty/Cs
430944	CA	115	0.22	24
430945	CA	115	0.45	24

CA = Cellulose Acetate



Tube Top Vacuum Filter

150 mL Tube Top Vacuum Filters

- ▶ 50 mm diameter membrane
- ▶ Minimizes unnecessary transfers by filtering directly into 50 mL centrifuge tube
- ▶ Includes two centrifuge tube stands with each case
- ▶ Each polypropylene centrifuge tube is supplied with an individually wrapped cap for storage
- ▶ Individually packaged, sterile, certified nonpyrogenic

150 mL Tube Top Vacuum Filters Ordering Information

Cat. No.	Membrane	Funnel Size/ Tube Size (mL)	Pore Size (µm)	Qty/Cs
430314	CA	150/50	0.45	12
430320	CA	150/50	0.22	12

CA = Cellulose Acetate

Vacuum Filter Systems

- ▶ Four sizes: 150 mL; 250 mL, 500 mL, and 1L
- ▶ Adapters are color coded by membrane type for easy product identification
- ▶ Angled hose connector simplifies vacuum line attachment
- ▶ Receiver bottles feature easy grip sides for improved handling
- ▶ Individually packaged, sterile, certified nonpyrogenic
- ▶ Caps for receiver bottles are sterile and individually packaged
- ▶ Extra plastic storage bottles are available, see page 109



Vacuum Filter Systems

Vacuum Filter Systems Ordering Information

Cat. No.	Membrane	Funnel/Bottle Volume (mL)	Pore Size (µm)	Color-Coded Adapter	Qty/Cs
<i>150 mL Capacity, 50 mm Diameter Membrane</i>					
431153	PES	150/150	0.22	Yellow	12
431154	CA	150/150	0.22	Orange	12
431155	CA	150/150	0.45	Orange	12
<i>250 mL Capacity, 50 mm Diameter Membrane</i>					
430756	CN	250/250	0.22	Blue	12
430767	CA	250/250	0.22	Orange	12
430768	CA	250/250	0.45	Orange	12
430771	NY	250/250	0.2	Red	12
431096	PES	250/250	0.22	Yellow	12
<i>500 mL Capacity, 70 mm Diameter Membrane</i>					
430758	CN	500/500	0.22	Blue	12
430769	CA	500/500	0.22	Orange	12
430770	CA	500/500	0.45	Orange	12
430773	NY	500/500	0.2	Red	12
431097	PES	500/500	0.22	Yellow	12
<i>1,000 mL Capacity, 90 mm Diameter Membrane</i>					
430186	CN	1,000/1,000	0.22	Blue	12
430515	NY	1,000/1,000	0.2	Red	12
430516	CA	1,000/1,000	0.45	Orange	12
430517	CA	1,000/1,000	0.22	Orange	12
431098	PES	1,000/1,000	0.22	Yellow	12
431205*	CA	500*/1,000	0.22	Orange	12
431206*	CA	500*/1,000	0.45	Orange	12

*500 mL Funnel with 70 mm membrane.

PES = Polyethersulfone, CA = Cellulose Acetate, CN = Cellulose Nitrate, NY = Nylon.



Bottle Top Vacuum Filters

Bottle Top Vacuum Filters

- Individually packaged, sterile and certified nonpyrogenic
- Adaptors are color coded by membrane type
- Available in 33 mm and 45 mm neck sizes to fit most glass and plastic media storage bottles
- 45 mm neck sizes fit on Corning® plastic storage bottles, see page 109

Bottle Top Vacuum Filters Ordering Information

Cat. No.	Membrane	Volume (mL)	Neck Size (mm)	Pore Size (µm)	Color-Coded Adapter	Qty/Cs
<i>150 mL Capacity, 50 mm Diameter Membrane</i>						
430624	CA	150	33	0.22	Orange	48
430625	CA	150	33	0.45	Orange	48
430626	CA	150	45	0.22	Orange	48
430627	CA	150	45	0.45	Orange	48
431160	PES	150	33	0.22	Yellow	48
431161	PES	150	45	0.22	Yellow	48
<i>500 mL Capacity, 70 mm Diameter Membrane</i>						
430049	NY	500	45	0.2	Red	12
430512	CA	500	33	0.45	Orange	12
430513	CA	500	45	0.22	Orange	12
430514	CA	500	45	0.45	Orange	12
430521	CA	500	33	0.22	Orange	12
431117	PES	500	33	0.22	Yellow	12
431118	PES	500	45	0.22	Yellow	12
<i>1,000 mL Capacity, 90 mm Diameter Membrane</i>						
430015	CA	1,000	45	0.22	Orange	12
431174	PES	1,000	45	0.22	Yellow	12

PES = Polyethersulfone, CA = Cellulose Acetate, CN = Cellulose Nitrate, NY = Nylon.

Syringe Filters



Syringe Filters

- ▶ A variety of membranes are available to meet your needs: Polyethersulfone (PES) – low protein binding and faster flow rates; surfactant-free cellulose acetate (SFCA) – lowest protein binding; Teflon (PTFE) – chemical resistance; regenerated cellulose (RC) – best choice for DMSO compatibility; Nylon (NY) – hydrophilic, surfactant-free and lowest extractable.
- ▶ 100% integrity tested, certified nonpyrogenic and noncytotoxic, manufactured in accordance with ISO 9002 standards

Syringe Filters Ordering Information

Cat. No.	Diameter (mm)	Pore Size (µm)	Membrane Material	Housing Material	Sterile	Inlet/Outlet	Packaging	Qty/Cs
431212	4	0.2	RC	PP	Yes	LL/LS	Ind	50
431215	15	0.2	RC	PP	Yes	LL/LS	Ind	50
431218	26	0.2	SFCA-PF	AC	Yes	LL/LS	Ind	50
431219	26	0.2	SFCA	AC	Yes	LL/LS	Ind	50
431220	26	0.45	SFCA	AC	Yes	LL/LS	Ind	50
431221	26	0.8	SFCA	AC	Yes	LL/LS	Ind	50
431222	25	0.2	RC	PP	Yes	LL/LS	Ind	50
431224	25	0.2	NY	PP	Yes	LL/LS	Ind	50
431225	25	0.45	NY	PP	Yes	LL/LS	Ind	50
431227*	50	0.2	PTFE	PP	Yes	HB/HB	Ind	12
431229	26	0.2	PES	AC	Yes	LL/LS	Ind	50
431231	25	0.45	PTFE	PP	No	LL/LS	Bulk	50

PP = Polypropylene, AC = Acrylic Copolymer, LL = Luer Lock/Female, LS = Luer Slip/Male, HB = Hose Barb, NY = Nylon, PES = Polyethersulfone, PTFE = Teflon, RC = Regenerated Cellulose, SFCA = Surfactant Free Cellulose Acetate, SFCA-PF = Surfactant Free Cellulose Acetate with Prefilter.

*Recommended as in-line air filter.

Spin-X® Centrifuge Tube Filters



Spin-X Centrifuge Tube Filters

- ▶ Costar® Spin-X centrifuge tube filters consist of a membrane-containing filter unit within a microcentrifuge tube.
- ▶ Uses:
 - Removing bacteria, cells and particles from liquids
 - HPLC sample preparation
 - DNA removal from agarose or acrylamide gels. Maximum RCF (Relative Centrifugal Force [x g]) is 16,000

Spin-X Centrifuge Tube Filters Ordering Information

Cat. No.	Membrane Material	Working Volume (µL)	Pore Size (µm)	Sterile	Tube Size (mL)	Qty/Cs
8160	CA	500	0.22	Yes	2.0	96
8161	CA	500	0.22	No	2.0	100
8162	CA	500	0.45	Yes	2.0	96
8163	CA	500	0.45	No	2.0	100
8169	NY	500	0.22	No	2.0	200
8170	NY	500	0.45	No	2.0	200

CA = Cellulose Acetate, NY = Nylon.

Storage Bottles – Polystyrene



- ▶ Disposable polystyrene bottles for storage of media, buffers and other aqueous solutions
- ▶ Two styles:
 - Low profile, easy grip style has sides that facilitate handling
 - Traditional style has smooth sides
- ▶ Plug seal caps (45 mm) provide an airtight seal and help minimize the risk of contamination.
- ▶ Bottles can be used with Corning® Vacuum Filter Systems, see page 105
- ▶ Sterile, certified nonpyrogenic



Corning Easy Grip Style Storage Bottles Ordering Information

Cat. No.	Volume (mL)	Neck Size (mm)	Qty/Pk	Qty/Cs
431175	150	45	2	24
430281	250	45	2	24
430282	500	45	2	24
430518	1,000	45	2	24

Costar® Traditional Style Storage Bottles Ordering Information

Cat. No.	Volume (mL)	Neck Size (mm)	Qty/Pk	Qty/Cs
8388	125	45	1	24
8390	250	45	1	12
8393	500	45	1	12
8396	1,000	45	1	12

Storage Bottles – Square, Polycarbonate



Corning square polycarbonate storage bottles are easier to handle, require less space (13-20%) on the shelf or in the autoclave and are ideal for mixing, sampling and storage.

- ▶ Strong polycarbonate bottles are more break-resistant than other glass or plastic bottles.
- ▶ Screened white enamel graduations are easier to see than molded graduations.
- ▶ Sterilized by gamma radiation
- ▶ Large white marking spot for easier identification
- ▶ Bottles can be autoclaved once at 121°C and 15 psi. Repeated autoclaving weakens polycarbonate and is not recommended.
- ▶ Optional reusable caps with silicone septa (Corning Cat. Nos. 1395-45HTSC, 1395-45SS, 1395-45TS) for syringe sampling or introduction of reagents are available for these caps.
- ▶ Store up to -80°C

Corning Square Polycarbonate Storage Bottles Ordering Information

Cat. No.	Capacity (mL)	Shape	Bottle Material	Neck Dia. (mm)	Qty/Bag	Qty/Cs
431430	150	Square	Polycarbonate	45	1	24
431431	250	Square	Polycarbonate	45	1	24
431432	500	Square	Polycarbonate	45	1	24
431433	1000	Square	Polycarbonate	45	1	24

Containers



- Flexible polypropylene bottom with snap-on polyethylene lid serves as a beaker or storage container
- Graduated in both milliliters and ounces
- Certified nonpyrogenic

Containers Ordering Information

Cat. No.	Description	Sterile	Capacity (mL)	Qty/Pk	Qty/Cs
430179	Container and Lid	Yes	250	1	100
430180	Container Only	Yes	250	20	500
430181	Lid Only	Yes	n/a	20	500

DISPOSABLE SAMPLE CONTAINERS

1700 Corning® Coliform Water Test Sample Container, Sterile* with Sodium Thiosulfate Tablet



Sterile container used in testing for the presence of coliform, a microbiological contaminant in drinking water. Manufactured from pure polypropylene in a sterile environment. The one-piece container has attached lid to reduce chance of contamination. Locking arrow assures sterility has not been compromised. The EPA fill line of 100 mL \pm 2.5% makes it easy to use. A sodium thiosulfate tablet has been added to each container thus saving lab prep time and expense. Leak tight when sealed properly. An added benefit is the tie-down to protect from accidental opening and also serves as a custody seal. Sample label and instructions for use are supplied with each. A low cost, convenient product which meets EPA requirements.

Cat. No.	Description	Capacity (mL)	Approx. Diam. x Height (mm)	Qty/Cs
1700-100	Container w/tablet	100-120	65 x 120	100

*Sterile-by-process.

1705 Corning Water Test Sample Container, Sterile* without Sodium Thiosulfate Tablet



Sterile container used in the testing of non-chlorinated drinking water. Manufactured from pure, recyclable polypropylene. The one-piece container has attached lid to reduce chance of contamination. Locking arrow assures sterility has not been compromised. Leak tight when sealed properly. An added benefit is the tie-down which protects against accidental opening.

Cat. No.	Description	Capacity (mL)	Approx. Diam. x Height (mm)	Qty/Cs
1705-100	Container w/o Tablet	100-120	65 x 120	100

*Sterile-by-process.

1730 Corning Snap-Seal Plastic Sample Containers



Designed for a wide variety of applications, these containers provide a reliable leak-tight seal when closed properly. The Snap-Seal locking device keeps the cap closed and secure. The specially designed hinged cap stays in place in use, reducing the chance of sample contamination. The containers are made of recyclable polypropylene, in a translucent style for normal usage. The containers are graduated in both milliliters and ounces, and the cap has a rough surface for marking.

Cat. No.	Capacity	Color	Approx. Diam. x Height (mm)	Qty/Cs
1730-5X	0.45 oz. (13 mL)	Natural	16 x 94	500
1730-2C	1.5 oz. (45mL)	Natural	30 x 84	400
1730-4H	4 oz. (120 mL)	Natural	45 x 91	200
1730-4L	4 oz. (120 mL)	Natural	68 x 52	200
1730-8	8 oz. (240 mL)	Natural	80 x 75	100
1730-10	10 oz. (300 mL)	Natural	63 x 112	100

Cylinder



- ▶ Optically clear polystyrene
- ▶ Sterile
- ▶ Graduated for accurate dispensing
- ▶ A polyethylene dust cover is included

Cylinder Ordering Information

Cat. No.	Capacity (mL)	Graduation (mL)	Sterile	Qty/Pk	Qty/Cs
430182	100	1	Yes	1	50


Erlenmeyer Flasks




1L Erlenmeyer Flask

Polycarbonate Erlenmeyer Flasks

- ▶ Made from optically clear polycarbonate
- ▶ Ideal for shaker culture applications
- ▶ Two-position polypropylene plug seal caps can be open for gas exchange or closed for liquid-tight seal
- ▶ Vent cap option for continuous gas exchange while ensuring sterility and preventing leakage
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic



Breathable two-position plug seal caps feature one-piece linerless construction with a flexible plug for a gas- and liquid-tight seal. In addition, the unique breathable cap design allows use in either an open or closed mode.



Vent caps contain a 0.2 µm nonwetttable membrane sealed to the cap, providing consistent, sterile gas exchange while minimizing the risk of contamination.

Shaker Flask Application Tip

Corning recommends starting with a shaking rate of 75-125 RPM (orbital shaker) and a medium volume of 30-40% of the nominal flask capacity.

Polycarbonate Erlenmeyer Flasks Ordering Information

Cat. No.	Capacity (mL)	Graduation (mL)	Neck Diameter (mm)	Cap Style	Qty/Pk	Qty/Cs
430421	125	25	26	Plug Seal	1	50
431143	125	25	26	Vent Cap	1	50
430183	250	25	31	Plug Seal	1	50
431144	250	25	31	Vent Cap	1	50
430422	500	50	43	Plug Seal	1	25
431145	500	50	43	Vent Cap	1	25
431146	1,000	50	43	Plug Seal	1	25
431147	1,000	50	43	Vent Cap	1	25



431255 2L Erlenmeyer Flask



431252 3L Fernbach Culture Flask

Polycarbonate 2L and 3L Flasks

- ▶ Made from optically clear polycarbonate
- ▶ Ideal for shaker and suspension culture applications
- ▶ Available with or without baffled bottoms
- ▶ Features a hydrophobic membrane in cap for applications requiring sterile gas exchange
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic
- ▶ Solid and vent caps available separately

Polycarbonate 2L and 3L Flasks Ordering Information

Cat. No.	Description	Sterile	Qty/Cs
431255	Erlenmeyer Flask, 2L, Polycarbonate	Yes	6
431256	Erlenmeyer Flask, 2L, Polycarbonate, Baffled Bottom	Yes	6
431252	Fernbach Culture Flask, 3L, Polycarbonate	Yes	4
431253	Fernbach Culture Flask, 3L, Polycarbonate, Baffled Bottom	Yes	4
431339	Cap, Vented, 48 mm for 2L Flask	Yes	24
431340	Cap, Vented, 70 mm for 3L Flask	Yes	24
431364	Cap, Solid, 48 mm for 2L	Yes	24
431363	Cap, Solid, 70 mm for 3L	Yes	24

Spatulas



Spatulas



Microspatulas

- ▶ Corning® spatulas are designed to save time and to provide contamination-free samples
- ▶ Individually packaged, certified RNase-/DNase-free, nonpyrogenic, antistatic and sterile
- ▶ Eliminates the recycling and resterilizing necessary with reusable spatulas
- ▶ Available in five different configurations
- ▶ Microspatulas are available in two configurations

Spatulas Ordering Information

Cat. No.	Description	Qty/Cs
3003	Spatula, Tapered Blade/Spoon	100
3004	Spatula, Small Spoon/Spoon	100
3005	Spatula, Round End/Spoon	100
3006	Spatula, V-Scoop/Spoon	100
3007	Spatula, Flat End/Spoon	100
3012	Microspatula, Tapered End/Scoop	50
3013	Microspatula, Rounded End/Scoop	50

Centrifuge Tubes



15 mL Centrifuge Tube with CentriStar Cap

15 mL Centrifuge Tubes

- ▶ Corning® 15 mL centrifuge tubes feature black printed graduations and a large white marking spot
- ▶ Available with your choice of cap styles; the advanced CentriStar™ cap or the original plug seal cap
- ▶ Available in racks or bulk packed in ziplock, resealable sleeves
- ▶ Sterile, certified nonpyrogenic, and RNase-/DNase-free
- ▶ 95 k PA (14 psi) pressure tested – plug seal cap only
- ▶ Foam racks also available separately

15 mL Centrifuge Tubes Ordering Information

Cat. No.	Material	Cap Style	Max. RCF	Qty/Pk	Qty/Cs
430053	PET	Plug Seal Cap	3,600	25/Sleeve	500
430055	PET	Plug Seal Cap	3,600	50/Rack	500
430052	PP	Plug Seal Cap	12,000	50/Rack	500
430766	PP	Plug Seal Cap	12,000	25/Sleeve	500
430790	PP	CentriStar Cap	12,000	50/Rack	500
430791	PP	CentriStar Cap	12,000	25/Sleeve	500
431355	Foam Centrifuge Tube Rack, 15 mL				20

PP = Polypropylene, PET = Polyethylene Terephthalate, RCF = Relative Centrifugal Force (x g).

50 mL Centrifuge Tubes

- ▶ Corning 50 mL centrifuge tubes feature black printed graduations and a large white marking spot
- ▶ Available with your choice of cap styles: the advanced CentriStar cap or the original plug seal cap
- ▶ Available in racks or bulk packed in ziplock, resealable sleeves
- ▶ Sterile, certified nonpyrogenic, and RNase-/DNase-free
- ▶ 95 k PA (14 psi) pressure tested
- ▶ Foam racks also available separately

50 mL Centrifuge Tubes Ordering Information

Cat. No.	Material	Cap Style	Max. RCF	Qty/Pk	Qty/Cs
430290	PP	Plug Seal Cap	15,500	25/Rack	500
430291	PP	Plug Seal Cap	15,500	25/Sleeve	500
430304	PET	Plug Seal Cap	3,600	25/Rack	500
430828	PP	CentriStar Cap	15,500	25/Rack	500
430829	PP	CentriStar Cap	15,500	25/Sleeve	500
4558	PP	CentriStar Cap	15,500	25/Universal Rack*	300
4365	Foam Centrifuge Tube Rack, 50 mL			–	20
4366	Universal Foam Centrifuge Tube Rack, 15 mL and 50 mL			–	20

PP = Polypropylene, PET = Polyethylene Terephthalate, RCF = Relative Centrifugal Force (x g).

*New innovative universal rack can hold 50 mL and 15 mL tubes securely, allowing researchers to work with and store both size tubes in the same rack, saving bench and storage space.



50 mL Centrifuge Tube with CentriStar Cap



Bulk Pack – Ziplock Bag



Universal Rack



Foam Centrifuge Tube Racks



Self Standing 50 mL Centrifuge Tube with Flat Cap

Self-Standing 50 mL Centrifuge Tubes

- ▶ Corning® 50 mL centrifuge tubes feature black printed graduations and a large white marking spot
- ▶ Available with your choice of flat or the original plug seal cap
- ▶ Tubes are bulk packed in ziplock, resealable sleeves
- ▶ 95 k PA (14 psi) pressure tested – plug seal cap only
- ▶ Sterile, certified nonpyrogenic, and RNase-/DNase-free

Self-Standing 50 mL Centrifuge Tubes Ordering Information


Cat. No.	Material	Cap Style	Max. RCF	Qty/Sleeve	Qty/Cs
430897	PP	Plug Seal Cap	3,000	25	500
430921	PP	Flat Cap	3,000	25	500

PP = Polypropylene, RCF = Relative Centrifugal Force (x g).

CentriStar™ Cap

Corning 15 mL and 50 mL centrifuge tubes are now available with the advanced CentriStar cap. The CentriStar cap has an easy-on/easy-off flat top and offers advanced ergonomics with its wider knurls and roll-over edge design for easier gripping. This design includes a revolutionary plug feature that minimizes the risk of seepage when used under recommended conditions.

Dimensions of Corning 15 mL and 50 mL Centrifuge Tube



250 mL and 500 mL Centrifuge Tubes and Support Cushions

- ▶ Corning 250 mL and 500 mL polypropylene tubes are ideal for applications requiring large-volume centrifugation
- ▶ Each case of tubes contains a rack to facilitate handling
- ▶ Support cushions must be used with this product unless the rotor has appropriately shaped V-bottom holders
- ▶ Tubes are sterile and certified nonpyrogenic

250 mL and 500 mL Centrifuge Tubes Ordering Information

Cat. No.	Description	Material	Cap Style	Max RCF	Qty/Pk	Qty/Cs
430776	250 mL Tube	PP	Plug	6000	6	102
430236	250 mL Support Cushion	PEI	n/a	n/a	n/a	6
431123	500 mL Tube	PP	Plug	6000	6	36
431124	500 mL Support Cushion	PEI	n/a	n/a	n/a	6

PP = Polypropylene, PEI = Polyetherimide, RCF = Relative Centrifugal Force (x g).



500 and 250 mL Centrifuge Tubes

Microcentrifuge Tubes

Corning offers two styles of microcentrifuge tubes: traditional snap cap tubes for quick access or screw cap tubes for greater sealing security.



Microcentrifuge Tubes

Snap Cap Polypropylene Microcentrifuge Tubes

- ▶ Costar® microcentrifuge tubes are certified RNase-/DNase-free
- ▶ Supplied nonsterile and are autoclavable
- ▶ External graduations and frosted writing spot for easy sample identification
- ▶ Positive seal design allows for repeated opening and closing
- ▶ Flat cap surface for convenient labeling
- ▶ Withstands a maximum RCF of 17,000 x g
- ▶ Costar low binding microcentrifuge tubes feature a bonded polymer technology that reduces protein and nucleic acid binding, resulting in better sample recovery

Snap Cap Polypropylene Microcentrifuge Tubes Ordering Information

Cat. No.	Volume (mL)	Color	Qty/Pk	Qty/Cs
<i>Snap Cap Microcentrifuge Tubes</i>				
3208	0.65	Natural	500	1,000
3209*	0.65	Rainbow*	200	1,000
3620	1.7	Natural	500	500
3621	1.7	Natural	500	5,000
3622*	1.7	Rainbow*	100	500
3213	2.0	Natural	500	1,000
<i>Low Binding Snap Cap Microcentrifuge Tubes Ordering Information</i>				
3206	0.65	Natural	500	500
3207	1.7	Natural	250	250


*Rainbow pack includes one bag each of blue, green, yellow, red, and orange tubes.




Microcentrifuge Tubes

Screw Cap Polypropylene Microcentrifuge Tubes

- ▶ Corning® polypropylene microcentrifuge tubes feature screw caps that provide a tight secure seal
- ▶ Choice of attached cap with silicone O-ring or unattached rim seal cap
- ▶ All tubes have a large white marking spot.
- ▶ Withstands a maximum RCF of 13,000 x g
- ▶ Sterile



▶ Attached loop cap allows for optimum one-handed convenience. Silicone O-ring gasket provides a snug seal, safeguarding samples against leakage.




▶ Easy-to-use unattached rim seal cap design twists on or off in a single turn.

Screw Cap Polypropylene Microcentrifuge Tubes Ordering Information

Cat. No.	Volume (mL)	Cap Style	O-ring	Self Standing	Qty/Cs
430909	1.5	Attached	Yes	No	500
430915	2.0	Attached	Yes	Yes	500
430917	2.0	Unattached	No	Yes	500


Cryogenic Vials and Accessories

Corning offers three styles of cryogenic vials as well as storage racks and boxes.




External Thread Cryogenic Vial

- Color-coded polypropylene cap inserts simplify vial identification. Available in variety packs of white, blue, green, red, and yellow.
- Silicone washer provides a secure seal.
- Easy-to-read black graduations for partial volumes
- Self-standing base, self-locking skirt



Internal Thread Cryogenic Vial

- Color-coded polypropylene cap inserts simplify vial identification. Available in variety packs of white, blue, green, red, and yellow.
- Silicone washers or rubber O-rings provide a secure seal.
- Easy-to-read black graduations for partial volumes
- Self-standing base, self-locking skirt



External Thread Plug Seal Cap

- Sure-grip plug seal screw cap
- Inner cap ring assures a tight seal.



External Thread Cryogenic Vials

External Thread Cryogenic Vials

- Manufactured from polypropylene to withstand temperatures down to -196°C
- Larger marking spot
- Black graduations
- Certified RNase-/DNase-free
- Vials have a silicone washer for a secure seal.
- Vials may be color coded with inserts, see page 118
- Self-standing vials have a special base design allowing them to be locked into cryogenic rack and tray (Cat. No. 430525 or 431131) for single-handed manipulation
- Sterilized by gamma radiation
- Certified nonpyrogenic
- Free foam rack with each case

External Thread Cryogenic Vials Ordering Information

Cat. No.	Capacity (mL)	Style	Self-Standing	Qty/Pk	Qty/Cs
430658	1.2	Conical Bottom	Yes	50	500
430659	2.0	Round Bottom	Yes	50	500
430661	2.0	Round Bottom	No	50	500
430662	4.0	Round Bottom	Yes	50	500
430663	5.0	Round Bottom	Yes	50	500

Warning! Do not use cryogenic vials for storage in the liquid phase of liquid nitrogen. Only store vials in the vapor phase above the liquified gas. Always use appropriate safety equipment when removing vials from cryogenic storage.

Cryogenic Vial Safety Tip

Appropriate safety equipment (gloves, face shields, biological safety cabinets, hoods, etc.) should always be used to protect personnel when removing vials or ampules from cryogenic storage systems.



Internal Thread
Cryogenic Vials

Internal Thread Cryogenic Vials

- ▶ Manufactured from polypropylene to withstand temperatures down to -196°C
- ▶ Larger marking spot
- ▶ Black graduations
- ▶ Certified RNase-/DNase-free
- ▶ Vials have a silicone washer or rubber O-ring for a secure seal
- ▶ Vials may be color coded with inserts, see page 118
- ▶ Self-standing vials have a special base design allowing them to be locked into cryogenic rack and tray (Cat. No. 430525 or 431131) for single-handed manipulation
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic
- ▶ Free foam rack with each case

Internal Thread Cryogenic Vials Ordering Information

Capacity Cat. No.	(mL)	Style	Self- Standing	Seal Type	Qty/Pk	Qty/Cs
430487	1.2	Conical Bottom	Yes	Washer	50	500
430488	2.0	Round Bottom	Yes	Washer	50	500
430489	2.0	Round Bottom	No	Washer	50	500
430490	4.0	Round Bottom	No	Washer	50	500
430491	4.0	Round Bottom	Yes	Washer	50	500
430492	5.0	Round Bottom	No	Washer	50	500
430656	5.0	Round Bottom	Yes	Washer	50	500
431386	2.0	Round Bottom	Yes	Washer	50	250

Warning! Do not use cryogenic vials for storage in the liquid phase of liquid nitrogen. Only store vials in the vapor phase above the liquified gas. Always use appropriate safety equipment when removing vials from cryogenic storage.

External Thread Cryogenic Vials with Plug Seal Cap

- ▶ Manufactured from polypropylene to withstand temperatures down to -196°C
- ▶ Vials feature an external thread with a traditional plug seal cap design for a secure seal
- ▶ Cap does not accept color-coded inserts
- ▶ Sterilized by gamma radiation
- ▶ Certified nonpyrogenic

External Thread Cryogenic Vials with Plug Seal Cap Ordering Information

Cat. No.	Capacity (mL)	Style	Self-Standing	Qty/Pk	Qty/Cs
430289	2.0	Round Bottom	No	50	500

Warning! Do not use cryogenic vials for storage in the liquid phase of liquid nitrogen. Only store vials in the vapor phase above the liquified gas. Always use appropriate safety equipment when removing vials from cryogenic storage.



External Thread Cryogenic
Vials with Plug Seal Cap



Cap Inserts

Cap Inserts for Cryogenic Vials

- ▶ Cap inserts provide color coding for easy sample identification
- ▶ Inserts are packaged in resealable bags
- ▶ Nonsterile
- ▶ Cap inserts fit all Corning® cryogenic vials except Cat. No. 430289

Cryogenic Vials Cap Inserts Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
430499	Assorted colors, polypropylene cap inserts: 100 each of white, blue, red, green, and yellow	50	500
2015	White polypropylene cap inserts	50	500
2016	Blue polypropylene cap inserts	50	500
2017	Red polypropylene cap inserts	50	500
2018	Green polypropylene cap inserts	50	500
2019	Yellow polypropylene cap inserts	50	500

Cryogenic Vial Racks and Storage Boxes

- ▶ Reusable racks are designed for use with most cryogenic vials
- ▶ Cat. No. 430525 has a locking feature for use with all Corning self-standing vials


Cryogenic Vial Racks and Storage Boxes Ordering Information

Cat. No.	Description	Qty/Pk	Qty/Cs
430525	Polycarbonate rack and tray, holds 30 vials; self-locking design in ice/water bath	1	1
430526	Polycarbonate rack only, holds 30 vials; self-locking design	1	1
431131	Reusable orange polypropylene vial rack, holds 50 vials; self-locking design	2	2
431119	81 count (9 x 9 array) Cryogenic Box, for 1-2 mL vials	5	10
431120	81 count (9 x 9 array) Cryogenic Box, for 4-5 mL vials	5	10
431121*	100 count (10 x 10 array) Cryogenic Box, for 1-2 mL vials	5	10

*431121 accepts internally threaded cryogenic vials only.




Cryogenic Vial Racks



431119 Cryogenic Storage Box



431120 Cryogenic Storage Box



431121 Cryogenic Storage Box

Technical Appendix

SELECTING THE BEST FILTER FOR YOUR APPLICATION

Choosing a filter does not have to be complicated - Corning has simplified the process. Just follow these four easy steps:

Step 1: Match your application with the best pore size.

Step 2: Select the best membrane and housing material for your application.

Step 3: Select the correct membrane diameter to optimize flow rate and throughput.

Step 4: Choose the best filter design for your application.

Step 1: Match Your Application with the Best Pore Size

The pore size is usually determined by your application or objective.

- ▶ Routine laboratory sterilization of most media, buffers, biological fluids and gases is usually done with 0.2 or 0.22 μm pore filter membranes
- ▶ Clarification and prefiltration of solutions and solvents is best accomplished with 0.45 μm or larger filter membranes
- ▶ Prefiltration to improve filter performance can also be accomplished by the use of glass fiber prefilters sold separately.

Use Table 1 to match your applications with a recommended membrane and pore size.

Step 2: Select the Best Membrane and Housing Material for Your Application

Your filter unit must be fully compatible with the chemical characteristics of your sample.

- ▶ Some filter membranes contain nontoxic wetting agents that may interfere with some applications
- ▶ Other membranes may bind proteins or other macromolecules leading to premature filter clogging or loss of valuable samples

Therefore, it is very important to understand their characteristics and the potential effects filter membranes can have on the solutions they contact. The following four graphs (Figure 1)

Table 1. Selecting the Pore Size

Application	Pore Size (μm)	Membrane Availability
Sterilization and Ultracleaning of Aqueous Solutions	0.20 to 0.22	All Membranes except Teflon™
Ultracleaning of Solvents (HPLC)	0.20 to 0.22	RC*, Teflon, Nylon
Clarification of Aqueous Solutions	0.45	All Membranes except Teflon
Clarification of Solvents (HPLC)	0.45	RC, Teflon, Nylon
Coarse Particle Removal	0.8	SFCA*, Glass Fiber Prefilters

*RC = Regenerated Cellulose, SFCA = Surfactant-Free Cellulose Acetate.

Figure 1. Important Performance Characteristics of Corning® Filter Membranes



PES = Polyethersulfone, CA = Cellulose Acetate, CN = Cellulose Nitrate, NY = Nylon

compare the flow rates, levels of extractable materials, and relative amounts of protein binding of four of the most popular membranes used in Corning® filters. Combining this with the information from Tables 2 and 3 (page 120) will help you choose the best Corning membranes for your applications.

Corning Filter Membrane Materials

Polyethersulfone (PES) membranes are the best for filtering cell culture media. PES has very low protein binding and extractables. PES also demonstrates faster flow rates than cellulose or nylon membranes.

Cellulose acetate (CA) membranes have a very low binding affinity for most macromolecules and are especially recommended for applications requiring low protein binding, such as filtering culture media containing sera. However, both cellulose acetate and cellulose nitrate membranes are naturally hydrophobic and have small amounts (less than 1%) of non-toxic wetting agents added during manufacture to ensure proper wetting of the membrane. If desired, these agents can be easily removed prior to use by filtering a small amount of warm purified water through the membrane or filter unit. Surfactant free cellulose acetate membranes, with very low levels of extractables, are available on some Corning® syringe filters.

Cellulose nitrate (CN) membranes are recommended for filtering solutions where protein binding is not a concern. They are recommended for use in general laboratory applications such as buffer filtration. Corning's cellulose nitrate membranes are Triton X-100®-free and noncytotoxic.

Nylon membranes are naturally hydrophilic and are recommended for applications requiring very low extractables since they do not contain any wetting agents, detergents or surfactants. Their greater chemical resistance makes them better for filtering more aggressive solutions, such as alcohols and

DMSO. However, like cellulose nitrate membranes, they may bind greater amounts of proteins and other macromolecules than do the cellulose acetate or PES membranes. They are recommended for filtering protein-free culture media.

Regenerated cellulose (RC) membranes are hydrophilic and have very good chemical resistance to solvents, including DMSO. They are used to ultraclean and de-gas solvents and mobile phases used in HPLC applications.

Teflon™ (PTFE; polytetrafluoroethylene) membranes are naturally and permanently hydrophobic. They are ideal for filtering gases, including humidified air. The extreme chemical resistance of Teflon membranes makes them very useful for filtering solvents or other aggressive chemicals for which other membranes are unsuitable. Because of their hydrophobicity, Teflon membranes must be prewetted with a solvent, such as ethanol, before aqueous solutions can be filtered.

Glass fiber filters are used as depth filters for prefiltering solutions. They have very high particle loading capacity and are ideal for prefiltering dirty solutions and difficult to filter biological fluids such as sera.

Table 2. Characteristics of Corning Filter Membranes

	Cellulose Nitrate	Cellulose Acetate	Nylon	Polyether-Sulfone	Regenerated Cellulose	Teflon (PTFE)
Flow rates for medium with 10% serum	Good	Very Good	Poor	Best	NA	NA
Wetting Agents	Yes	Yes	No, naturally hydrophilic	No	Yes	Does not wet
Protein Binding	Very high	Very low	Low to moderate	Very low	Low	NA
DNA Binding	High	Very low	Very high	Very low	Low	NA
Chemical Resistance	Low	Low	Moderate to high	Low	Very high	Very high

Table 3. Chemical Resistance Guide for Corning Filters

This information has been developed from a combination of laboratory tests, technical publications, or material suppliers. It is believed to be reliable. Due to conditions outside of Corning's control, such as variability in temperatures, concentrations, duration of exposure and storage conditions, no warranty is given or is to be implied with respect to this information.

Chemical Class	Filter Membranes							Housing Materials				
	CN	CA	PC	NY	PES	RC	PTFE	PS	PP	AC	PYR	PVC
Weak Acids	2	2	1	2	3	1	1	1	1	2	1	1
Strong Acids	3	2	3	3	3	3	1	2	1	3	2	1
Alcohols	3	1	1	1	1	1	1	2	1	3	1	1
Aldehydes	2	3	2	2	3	2	1	3	1	3	1	3
Aliphatic Amines	3	3	3	1	1	1	1	3	1	3	1	2
Aromatic Amines	3	3	3	2	3	1	1	3	1	3	1	3
Bases	3	3	3	2	3	2	1	1	1	2	2	1
Esters	3	3	2	1	3	1	1	3	2	2	1	3
Hydrocarbons	2	2	2	2	3	1	1	3	2	2	1	2
Ketones	3	3	2	2	3	1	1	3	2	3	1	3

Key: 1 = Recommended, 2 = May be suitable for some applications, a trial run is recommended, 3 = Not recommended, CN = Cellulose Nitrate, CA = Cellulose Acetate, NY = Nylon, PYR = PYREX Glass, PC = Polycarbonate, PES = Polyethersulfone, RC = Regenerated Cellulose, PS = Polystyrene, PTFE = Polytetrafluoroethylene (Teflon), PP = Polypropylene, PVC = Polyvinylchloride, AC = Acrylic Copolymer.

Corning® Filter Housing Materials

The filter housing materials also must be compatible with the solutions being filtered.

Polystyrene (PS) is used in the filter funnels and storage bottles for all of the Corning plastic vacuum filters. This plastic polymer should only be used in filtering and storing nonaggressive aqueous solutions and biological fluids. Refer to Table 3 (page 120) for more chemical compatibility information.

Acrylic copolymer (AC) and **Polyvinyl chloride (PVC)** are used in some of the Corning syringe filter housings. These plastics should only be used in filtering less aggressive aqueous solutions and biological fluids. Refer to Table 3 for more chemical compatibility information.

Polypropylene (PP) is used in the Spin-X® centrifuge filters and some of the syringe and disc filter housings. This plastic polymer has very good resistance to many solvents. Refer to Table 3 for more chemical compatibility information.

Chemical Compatibility

The mechanical strength, color, appearance, and dimensional stability of Corning filters are affected to varying degrees by the chemicals with which they come into contact. Specific operating conditions, especially temperature and length of exposure, will also affect their chemical resistance. Table 3 provides basic information on the chemical resistance of Corning filter membranes and housings.

Step 3: Select the Correct Membrane Diameter to Optimize Flow Rate and Throughput

The third step is selecting a filter that will have enough volume capacity or throughput to process your entire sample quickly and efficiently. This is primarily determined by the effective surface area of the membrane. Table 4 shows the relationship between filter diameter, effective filtration surface area and expected throughput volumes. The lower values are typical of viscous or particle-laden solutions; the higher values are typical of buffers or serum-free medium.

Step 4: Choose the Best Filter Design for Your Application

Disposable Plastic Vacuum Filters

These sterile filters are available in four styles: complete filter/storage systems, bottle top filters, centrifuge tube top filters, or one-piece filter systems. Four membranes are available to meet all of your filtration needs: cellulose acetate, cellulose nitrate, nylon, or polyethersulfone.

Table 4. Typical Expected Throughput Volumes

Filter Diameter and Description	Effective Filtration Surface (cm ²)	Expected Throughput (mL)*
4 mm syringe/disc	0.07	0.05-3
15 mm syringe/disc	1.7	3-15
25 mm syringe/disc	4.8	15-100
26 mm syringe/disc	5.3	15-100
50 mm disc	19.6	100-750
50 mm vacuum system	16.6	100-750
60 mm vacuum system	24.6	200-1,000
70 mm vacuum system	38.5	300-1,500
90 mm vacuum system	58.1	500-2,000

*These values assume an aqueous solution and a 0.2 µm membrane. Solutions containing sera or other proteinaceous materials will be at the lower end of the range. Use of prefilters may extend the throughput 50 to 100% above the values shown.

Disposable Syringe/Disc Filters

The smaller conventional Corning syringe disc-type filters (4, 15, 25, and 26 mm diameter) are used with syringes which serve as both the fluid reservoir and the pressure source. The HPLC certified nonsterile syringe filters are available with nylon, regenerated cellulose or Teflon® (PTFE) membranes in polypropylene housing for extra chemical resistance. The sterile tissue culture tested syringe filters are available in PES, regenerated cellulose (ideal for use with DMSO-containing solutions) or surfactant-free cellulose acetate membranes in either polypropylene or acrylic copolymer housings.

The larger 50 mm diameter disc filter has a Teflon (PTFE) membrane and polypropylene housing with hose barb connectors. This product is ideal for filtering aggressive solvents or gases and applications requiring sterile venting of gases. Because they have a hydrophobic (will not pass aqueous solutions) membrane, they are also ideal for protecting vacuum lines and pumps.

Spin-X® Disposable Centrifuge Tube Filters

Costar® Spin-X centrifuge tube filters consist of a membrane-containing (either cellulose acetate or nylon) filter unit within a polypropylene microcentrifuge tube. They filter small sample volumes by centrifugation for bacteria removal, particle removal, HPLC sample preparation, removal of cells from media, and purification of DNA from agarose and polyacrylamide gels. (See Corning Technical Bulletin: *Spin-X Purification of DNA from agarose gels* at www.corning.com/lifesciences.)

Corning FiltrEX™ 96 and 384 Well Filter Plates

Information on Corning FiltrEX 96 and 384 well filter plates can be found in the Corning Genomics Selection Guide or on the Corning Life Sciences web site www.corning.com/lifesciences.

Table 5. Corning® Filter Designs

Design	Sterile	Filter Diameters (mm)	Available Membrane Materials	Pore Sizes (µm)	Special Features
Syringe Filters	Some	4, 15, 25, and 26	RC, PES, SFCA, NY, and PTFE	0.2, 0.45, and 0.8	Ideal for small volume pressure filtration
Disc Filters	Yes	50	PTFE	0.2	Ideal for filtering solvents and gases
Vacuum Filter Storage Systems	Yes	50, 70 and 90	PES, CA, CN, and Nylon	0.2, 0.22, and 0.45	Easy grip bottles for storing filtrate
Bottle Top Vacuum Filters	Yes	50, 70 and 90	PES, CA, CN, and Nylon	0.2, 0.22, and 0.45	2 neck widths to fit most glass bottles
Tube Top Vacuum Filters	Yes	50	CA	0.22 and 0.45	Minimizes unnecessary transfers by filtering into a 50mL centrifuge tube
115 One Piece Vacuum Filters	Yes	60	CA and CN	0.2 and 0.45	Very economical with separate pour spout
Spin-X® Centrifuge Filters	Some	7.7	CA and Nylon	0.22 and 0.45	Ideal for purifying DNA from agarose gels
FiltrEX™ 96 and 384 Well Filter Plates	Some	6.4, 3.2	PVDF, GlassFiber, PES, NC, and UF	0.2, 0.45, 1.2 and others	Clear, opaque, or solvent resistant*

*Call for specific details; several custom-made products available.

CN = Cellulose Nitrate, CA = Cellulose Acetate, PES = Polyethersulfone, RC = Regenerated Cellulose, PTFE = Polytetrafluoroethylene (Teflon), SFCA = Surfactant-Free Cellulose Acetate.

CHARACTERISTICS OF CORNING PLASTICWARE

		Polystyrene	Polyethylene (High Density)	Polypropylene	Polycarbonate	Nylon	P.T.F.E. (Teflon®)
Physical Characteristics	Basic Properties	Biologically inert, hard, excellent optical qualities	Biologically inert, high chemical resistance	Biologically inert, high chemical resistance, exceptional toughness	Clear, very tough, inert, high temperature resistance	Tough, heat resistant, machinable, high moisture vapor transmission	Biologically and chemically inert, high resistant slippery surface
	Clarity	Clear	Opaque	Translucent	Clear	Opaque	Opaque
	Autoclave Results	Melts	May distort	Withstands several cycles	Withstands one cycle	OK	OK
	Heat Distortion Point	147-175°F 64-80°C	250°F 121°C	275°F 135°C	280-290°F 138-143°C	300-356°F 150-180°C	250°F 121°C
	Burning Rate	Slow	Slow	Slow	Self-extinguishing	Self-extinguishing	None
Effects of Laboratory Reagents	Weak Acids	None	None	None	None	None	None
	Strong Acids	Oxidizing acids attack	Oxidizing acids attack	Oxidizing acids attack	May be attacked	Attacked	None
	Weak Alkalies	None	None	None	None	None	None
	Strong Alkalies	None	None	None	Slowly attacked	None	None
	Organic Solvents	Soluble in aromatic chlorinated hydrocarbons	Resistant below 80°C	Resistant below 80°C	Soluble in chlorinated hydrocarbons; partly soluble in aromatics	Resistant	Resistant
Gas Permeability of Thin Wall Products*	O ₂	Low	High	High	Very low	Very low	—
	N ₂	Very low	Low	Low	Very low	Very low	—
	CO ₂	High	Very high	Very high	Low	—	—

Portions of this table courtesy of Modern Plastics Encyclopedia. Most data are from tests by A.S.T.M. methods. Tables show averages or ranges. Many properties vary with manufacturer, formulation, testing laboratory, and the specific operating conditions.

*Obtained from a table which lists gas permeability in CC/100 sq. inches per 24 hrs./mil.

CHEMICAL COMPATIBILITY OF CORNING PLASTICWARE

	PS	PP	PVC	CA	PC	CN	NY	MCE	PTFE	PET
<i>Acids</i>										
Hydrochloric acid (25%)	G	G	G	N	R	R	N	O	R	R
Hydrochloric acid (concentrated)	F	G	F	N	R	N	N	N	R	O
Nitric acid (concentrated)	P	P	P	N	R	N	N	N	O	N
Nitric acid (25%)	P	G	F	N	R	L	N	O	R	R
<i>Alcohols</i>										
Butanol	G	G	G	R	R	R	R	R	R	R
Ethanol	G	G	G	R	R	N	R	O	R	R
Methanol	G	G	G	R	R	N	R	O	R	R
<i>Amines</i>										
Aniline	G	G	P	N	N	R	R	N	R	O
Dimethylformamide	P	G	F	N	N	N	R	N	R	N
<i>Bases</i>										
Ammonium hydroxide (25%)	F	G	G	R	N	R	R	O	N	O
Ammonium hydroxide (1N)	F	G	G	N	N	R	R	O	N	N
Sodium hydroxide	G	G	G	N	N	N	R	N	R	N
<i>Hydrocarbons</i>										
Hexane	P	G	F	R	R	R	R	R	R	R
Toluene	P	G	P	R	O	R	R	R	R	N
Xylene	P	F	P	R	R	R	R	R	R	N
Dioxane	P	G	P	N	N	N	R	N	R	R
Dimethylsulfoxide (DMSO)	P	G	P	N	N	N	R	N	R	O*
<i>Halogenated Hydrocarbons</i>										
Chloroform	P	G	P	N	N	R	R	N	R	N
Methylene chloride	P	F	P	N	N	R	R	N	R	N
<i>Ketones</i>										
Acetone	P	G	P	N	O	N	R	N	R	R
Methyl ethyl diketone	P	G	P	N	O	N	R	O	R	R

*Can be used with aqueous solutions containing up to 20% DMSO.

R = Recommended, L = Limited Resistance, N = Not Recommended, O = Testing Advised, F = Fair, G = Good, P = Poor, PP = Polypropylene, PVC = Polyvinyl Chloride, CA = Cellulose Acetate, PC = Polycarbonate, PTFE = Polytetrafluoroethylene PS = Polystyrene, CN = Cellulose Nitrate, NY = Nylon, MCE = Mixed Cellulose Esters, PET = Polyethylene Terephthalate.

CHARACTERISTICS OF CORNING® CENTRIFUGE TUBES

The following information is provided to serve as a general guideline for determining suitability of Corning centrifuge tubes for your applications. In addition, Corning recommends following the procedures outlined by the centrifuge manufacturer, as well as conducting a trial run to determine proper conditions before beginning any critical applications.

Corning centrifuge tubes are tested for leakage. They should not break or leak if used in a properly balanced rotor with suitable carriers, holders, and adapters that fully support the tubes when run in accordance with the guidelines in this section. These tubes are intended for one-time use only; reuse is not recommended as breakage or leakage may occur.

The recommended working temperature range for Corning centrifuge tubes is 0 to 40°C. The suitability of these tubes for storage below 0°C depends on both the solution and the

storage conditions. In general, the polypropylene and PET tubes are more resistant to stress at low temperatures than polystyrene. It is strongly recommended that a trial run be performed under actual conditions to test the suitability of the tubes for frozen storage.

Suggestions for Safe Centrifugation

- Caution:** When centrifuging pathogenic organisms, clinical specimens known or suspected of being infectious, or any other potentially biohazardous materials, approved safety containment systems should be used. Contact your centrifuge manufacturer for appropriate accessories or recommendations.
- Read protocols and instruction manuals carefully. Do not confuse speed or revolutions per minute (RPM) with relative centrifugal force (RCF). Instructions for centrifuging a sample at a given RPM and time are incomplete unless the rotor or radius is specified. Protocols should always state the time and RCF value for centrifuging a sample.

- Proper balancing and distribution of the load in a centrifuge is critical for optimum performance and to prevent damage to the tubes or centrifuge. Opposing buckets or loads should always be balanced within the range specified by the manufacturer. Tubes should always be distributed in the buckets with respect to the center of rotation as well as the pivotal axis of the bucket. Failure to do this may prevent the bucket from achieving a horizontal position during the centrifugation run. Uneven separations or tube failure may result.

These centrifuge tubes are intended for use by persons knowledgeable in safe laboratory practices. Failure can result from surface damage, exceeding the specified RCF values, using unsuitable support systems, improper temperatures, or incompatible chemicals.

The RCF ratings for Corning® disposable centrifuge tubes have been established at room temperature using tubes filled to nominal capacity with water and spun in a horizontal rotor

centrifuge for 5 minutes. The centrifuge must be equipped with the recommended carriers, adapters, and cushions that fully support the tubes. If an angle head rotor is used or proper support is not provided, RCF values will be lower. Use of liquid other than water may also lower RCF values. Please consult your centrifuge specifications and the nomogram table (page 125) to determine speeds at which maximum RCF is achieved.

Chemical Compatibility of Disposable Plastic Centrifuge Tubes

The mechanical strength, flexibility, color, weight and dimensional stability of all plastic centrifuge tubes are affected to varying degrees by the chemicals with which they come in contact. Specific operating conditions, especially temperature, RCF, rotor type, carrier design, and run length will also affect tube performance.

Physical Properties of Disposable Plastic Centrifuge Tubes

	Clear Polypropylene	New Polyethylene Terephthalate
Recommended Working Temp*	0-40°	0-40°
Heat Distortion Point	121°	70°
Flexibility	Moderate	Rigid
Transparency	Clear	Clear
Maximum RCF:		
15 mL Tube	8,400 x g	3,600 x g
50 mL Tube	9,400 x g	3,600 x g
250 mL Tube	—	—
500 mL Tube	—	—

Chemical Resistance of Disposable Plastic Centrifuge Tubes*


Chemical Class	Polyethylene Terephthalate	Polypropylene	Polyethylene Caps
Acids (weak)	1	1	1
Acids	3	1	1
Alcohols	1	1	1
Aldehydes	3 ^a	2 ^a	1
Bases	3	1	1
Esters	2	2	2
Hydrocarbons:			
Aliphatic	1	2	3
Aromatic	3	3 ^b	3
Halogenated	3	3	3
Ketones	2	2 ^c	2

*At room temperature for 24 hours.

1 = Recommended; 2 = Suitable for most applications. However, a trial run under specific operating conditions is recommended; 3 = Not recommended.

Note: a = Formaldehyde, rated 1; b = Phenol, rated 1; c = Acetone, rated 1.

Nomogram for Computing Relative Centrifugal Force



To calculate the RCF value at any point along the tube or bottle, measure the radius, in mm, from the center of the centrifuge spindle to the particular point. Draw a line from the radius value on the right hand column to the appropriate centrifuge speed on the left-hand column. The RCF value is the point where the line crosses the center column. The nomogram is based on the formula:

$$RCF = (11.17 \times 10^{-7}) RN^2$$

where:

R = Radius in mm from centrifuge spindle to point in tube bottom

N = Speed of spindle in RPM

Disposable Glassware and Equipment



Digital Stirring Hot Plates, page 135

CORNING® AND PYREX® DISPOSABLE GLASSWARE 128

DIGITAL STIRRING HOT PLATES 135

Corning® and PYREX® Disposable Glassware

Corning Life Sciences is pleased to present our Disposable Glass Selection Guide. In this guide, you will find a selection of Corning's newest and most requested products. Corning's PYREX and Corning Brand glassware is made to the same high quality standards as reusable glassware but are disposable to make them more convenient and economical to use.

Corning Quality

- ▶ PYREX low expansion borosilicate glass has been the standard in chemistry labs for over 90 years
- ▶ More heat and chemical resistant than plastics
- ▶ Choice of a wide selection of products, sizes and styles designed to meet your needs

Corning Convenience

- ▶ Wide variety of reusable glassware, disposable glassware and plastics to help you get your research done
- ▶ An economical alternative to cleaning and recycling reusable glassware
- ▶ Disposal is easy and safer
- ▶ Serological and Bacteriological pipets are sterile and ready to use
- ▶ All glassware products can be sterilized by autoclaving

For up-to-date information on Corning Life Sciences' comprehensive range of products and services, go to **www.corning.com/lifesciences** where you can access:

- ▶ New Products Information
- ▶ Technical Information including:
 - Instruction Manuals
 - Product Bulletins
- ▶ Educational Opportunities
- ▶ Product Catalog Information
- ▶ Product Literature
- ▶ Complete Distributor Information

For additional product information, please visit **www.corning.com/lifesciences**, or call +31 (0) 20 659 60 51 or contact your local support office listed on the back cover.

Ordering Information

Corning products are available directly from Corning (order online, by phone, fax or email).

Or, to place an order through a distribution partner, simply contact the distributor of your choice. For each requested product, provide the Corning catalog number, product description, and desired quantity.

PYREX® Disposable Glass Serological Pipets

These disposable PYREX serological pipets are calibrated “to deliver” (TD) with blow-out and have negative graduations.

- ▶ Available sterile and nonsterile; nonsterile are unplugged and bulk packed
- ▶ Choose shorty pipets for working in hoods or tight spaces
- ▶ Choose wide orifice tips for use with cells or viscous samples
- ▶ Long, slender, tapered tips to make pipetting go faster and easier

PYREX Disposable Glass Serological Pipets, TD, Individually Wrapped, Sterile, Plugged Ordering Information

Cat. No.	Capacity	Negative Grad.	Grad. Interval	Qty/ Pk	Approx. Qty/Cs
7077-1N	1.0 mL	0.2 mL	.01 mL	200	800
7077-2N	2.0 mL	0.2 mL	.0 mL	180	720
7077-5N	5.0 mL	1.0 mL	.1 mL	120	720
7077-10N	10.0 mL	2.0 mL	.1 mL	100	600

PYREX Disposable Glass Serological Pipets, TD, Flip-Top Canister Pack, Sterile, Plugged Ordering Information

Cat. No.	Capacity	Negative Grad.	Grad. Interval	Qty/ Pk	Approx. Qty/Cs
7078D-1	1.0 mL	0.2 mL	.01 mL	50	500
7078D-5	5.0 mL	1.0 mL	. mL	25	400
7078D-10	10.0 mL	2.0 mL	.1 mL	25	400

PYREX Disposable Glass Serological Pipets, TD, Multi-Pack, Sterile, Plugged Ordering Information

Cat. No.	Capacity	Negative Grad.	Grad. Interval	Qty/ Pk	Approx. Qty/Cs
7078-5X	0.5 mL	0.2 mL	.01 mL	25	500
7078-1N	1.0 mL	0.2 mL	.1 mL	50	1000
7078-1CN	1.0 mL	0.2 mL	.01 mL	50	1000
7078-2N	2.0 mL	0.2 mL	.01 mL	35	700
7078-5N	5.0 mL	1.0 mL	.1 mL	30	960
7078-10N	10.0 mL	2.0 mL	.1 mL	20	720

PYREX Disposable Glass Wide-Tip Serological Pipet, TD, Multi-Pack, Sterile, Plugged Ordering Information

Cat. No.	Capacity	Negative Grad.	Grad. Interval	Qty/ Pk	Approx. Qty/Cs
7078A-5	5.0 mL	1.0 mL	.1 mL	25	400
7078A-10	10.0 mL	2.0 mL	.1 mL	20	400

PYREX Disposable Glass Shorty Serological Pipets, TD, Individually Wrapped, Sterile, Plugged Ordering Information

Cat. No.	Capacity	Negative Grad.	Grad. Interval	Qty/ Pk	Approx. Qty/Cs
7077B-1	1.0 mL	0.5 mL	.01 mL	250	500
7077B-5	5.0 mL	3.0 mL	.1 mL	200	400
7077B-10	10.0 mL	4.0 mL	.2 mL	200	400
7077B-25	25.0 mL	5.0 mL	.2 mL	100	400
7077B-50	50.0 mL	3.0 mL	.5 mL	25	50

PYREX® Disposable Glass Serological Pipets (Continued)**PYREX Disposable Glass Shorty Serological Pipets, TD, Multi-Pack, Sterile, Plugged Ordering Information**

Cat. No.	Capacity	Negative Grad.	Grad. Interval	Qty/Pk	Approx. Qty/Cs
7078B-1	1.0 mL	0.5 mL	.01 mL	10	500
7078B-5	5.0 mL	3.0 mL	.1 mL	10	400
7078B-10	10.0 mL	4. mL	.2 mL	10	400
7078B-25	25.0 mL	5.0 mL	.2 mL	5	400
7078B-50	50.0 mL	3.0 mL	.5 mL	5	50

PYREX Disposable Glass Serological Pipets, Disposable, TD, Bulk Pack, Non-Sterile, Unplugged Ordering Information

Cat. No.	Capacity	Negative Grad.	Grad. Interval	Qty/Pk	Approx. Qty/Cs
7079-5X	0.5 mL	0.2 mL	.01 mL	125	500
7079-1N	1.0 mL	0.2 mL	.01 mL	50	1000
7079-2N	2.0 mL	0.2 mL	.01 mL	35	700
7079-5N	5.0 mL	1.0 mL	.1 mL	30	960
7079-10N	10.0 mL	2.0 mL	.1 mL	20	720


PYREX Disposable Bacteriological Pipets, Multi-Pack, Sterile, Plugged

These PYREX borosilicate glass pipets meet the requirements recommended by the American Public Health Association as shown in “Standard Methods for the Examination of Dairy Products.”

- Both sizes are designed for gravity feed (blow-out is not required)
- These cotton-plugged pipets are packaged in sterile bags
- New I.S.O. color-coding is printed directly on the tear strip for easy identification by size

PYREX Disposable Bacteriological Pipets, Multi-Pack, Sterile, Plugged Ordering Information

Cat. No.	Capacity	Grad. Interval	Qty/Bag	Qty/Cs
7058-1X	1.1 mL	0.5, 1, 1.1	25	500
7058-2X	2.2 mL	1, 2, 2.1, 2.2	25	500





Corning® and PYREX® Disposable Pasteur Pipets

These disposable glass Pasteur unplugged pipets are available in convenient bulk-packs.

- ▶ Corning Pasteur pipets are made from economical soda lime glass
- ▶ PYREX Pasteur pipets are made from heat resistant borosilicate Code 7800 glass
- ▶ One-time use eliminates the danger of cross-contamination of valuable specimens or laboratory reagents
- ▶ Pipets have a constriction one inch below the top to facilitate plugging
- ▶ Ideal for rapid nonvolumetric transfer work

Corning and PYREX Disposable Pasteur Pipets Ordering Information

Cat. No.	Glass Type	Length	Qty/Pk	Qty/Cs
7095D-5X	Borosilicate	14.6 cm (5.75 inches)	200	1000
7095D-9	Borosilicate	22.9 cm (9 inches)	200	1000
7095B-5X	Soda lime	14.6 cm (5.75 inches)	200	1000
7095B-9	Soda lime	22.9 cm (9 inches)	200	1000
7095B-NMR	Soda lime	33.0 cm (13 inches)	100	100

Corning Microscope Slides

Corning microscope slides are made from water-white glass to maximize clarity.

- ▶ Corning frosted slides are frosted at one end on either one or both sides.
- ▶ Frosted slides are highly legible and easy to write on because they are sandblasted instead of etched.
- ▶ To minimize waste, all slides are inspected for chips and rough edges before being packaged.
- ▶ Products are packaged in one-half gross boxes (approx. 72 slides each)
- ▶ 5 gross = 10 boxes (720 slides total)
- ▶ 10 gross = 20 boxes (1440 slides total)

Corning Plain Microscope Slides Ordering Information


Cat. No.	Approx. Size	Thickness	Qty/Pk	Qty/Cs
2947-3x1	3 x 1 inch	0.90 - 1.10 mm	72	1440
2947-75x25	75 x 25 mm	0.90 - 1.10 mm	72	1440
2947-75x38	75 x 38 mm	0.90 - 1.10 mm	72	720
2947-75x50	75 x 50 mm	0.90 - 1.10 mm	72	720

Corning Frosted One Side Microscope Slides Ordering Information

Cat. No.	Approx. Size	Thickness	Qty/Pk	Qty/Cs
2948-3x1	3 x 1 inch	0.90 - 1.10 mm	72	1440
2948-75x25	75 x 25 mm	0.90 - 1.10 mm	72	1440

Corning Frosted Both Sides Microscope Slides Ordering Information

Cat. No.	Approx. Size	Thickness	Qty/Pk	Qty/Cs
2949-75x25	75 x 25 mm	0.90 - 1.10 mm	72	1440



Corning® Cover Glass

Corning cover glasses are made from optically clear water-white No. 0211 zinc titania glass.



- ▶ Resistant to surface attack or weathering
- ▶ Available in both rectangular and square formats in several thicknesses
- ▶ Packed in 1 ounce plastic boxes, 10 boxes/case

Corning Rectangular Cover Glasses Ordering Information

Cat. No.	Approx. Size	Thickness	Approx. Qty/Box	Boxes/Cs
2935-223	22 x 30 mm	#1 (0.13 - 0.16 mm)	118	10
2940-223	22 x 30 mm	#11/2 (0.16 - 0.19 mm)	97	10
2935-224	22 x 40 mm	#1 (0.13 - 0.16 mm)	88	10
2940-224	22 x 40 mm	#11/2 (0.16 - 0.19 mm)	73	10
2935-225	22 x 50 mm	#11/2 (0.16 - 0.19 mm)	72	10
2940-225	22 x 50 mm	#11/2 (0.16 - 0.19 mm)	58	10
2935-243	24 x 30 mm	#1 (0.13 - 0.16 mm)	108	10
2940-243	24 x 30 mm	#11/2 (0.16 - 0.19 mm)	89	10
2935-244	24 x 40 mm	#1 (0.13 - 0.16 mm)	81	10
2940-244	24 x 40 mm	#11/2 (0.16 - 0.19 mm)	67	10
2935-245	24 x 50 mm	#1 (0.13 - 0.16 mm)	65	10
2940-245	24 x 50 mm	#11/2 (0.16 - 0.19 mm)	54	10
2935-246	24 x 60 mm	#1 (0.13 - 0.16 mm)	54	10
2940-246	24 x 60 mm	#11/2 (0.16 - 0.19 mm)	45	10

Corning Square Cover Glasses Ordering Information

Cat. No.	Approx. Size	Thickness	Approx. Qty/Box	Boxes/Cs
2865-18	18 x 18mm	#1 (0.13 - 0.16 mm)	240	10
2870-18	18 x 18mm	#11/2 (0.16 - 0.19 mm)	198	10
2875-18	18 x 18mm	#2 (0.19 - 0.25 mm)	159	10
2865-22	22 x 22mm	#1 (0.13 - 0.16 mm)	160	10
2870-22	22 x 22mm	#11/2 (0.16 - 0.19 mm)	135	10
2875-22	22 x 22mm	#2 (0.19 - 0.25 mm)	106	10
2865-25	25 x 25mm	#1 (0.13 - 0.16 mm)	124	10
2870-25	25 x 25mm	#11/2 (0.16 - 0.19 mm)	102	10
2875-25	25 x 25mm	#2 (0.19 - 0.25 mm)	82	10



PYREX® Disposable Culture Tubes

These PYREX culture tubes are made from borosilicate glass to reduce pH changes and contaminants potentially leached from soda-lime glass.

- ▶ Choice of round bottom or flat bottom
- ▶ Available with or without a white marking spot
- ▶ Designed for both tissue culture and general bacteriological work.
- ▶ Available with or without cap threads
- ▶ Caps are sold separately; see back page

PYREX Disposable Glass Round Bottom Rimless Culture Tubes, Bulk Pack Ordering Information

Cat. No.	Approx. Capacity	Approx. O.D. x Length	G.P.I. Thread Finish	Qty/ Pk	Qty/ Cs
99445-10	4.0 mL	10 x 75 mm	No threads	250	1000
99445-12	6.0 mL	12 x 75 mm	No threads	250	1000
99445-13	10.0 mL	13 x 100 mm	No threads	250	1000
99445-15	11.0 mL	15 x 85 mm	No threads	250	1000
99445-16	15.0 mL	16 x 100 mm	No threads	250	1000
99445-16X	19.0 mL	16 x 125 mm	No threads	250	1000
99445-16XX	23.0 mL	16 x 150 mm	No threads	250	1000
99445-18	28.5 mL	18 x 150 mm	No threads	250	500
99445-20	36.0 mL	20 x 150 mm	No threads	250	500

PYREX Disposable Glass Round Bottom Screw Cap* Culture Tubes, White Marking Spot, Bulk Pack Ordering Information

Cat. No.	Approx. Capacity	Approx. O.D. x Length	G.P.I. Thread Finish	Qty/ Pk	Qty/ Cs
99447-13	7.5 mL	13 x 100 mm	13-415	250	1000
99447-16	15.5 mL	16 x 125 mm	15-415	250	1000
99447-16X	19.0 mL	16 x 150 mm	15-415	250	1000
99447-161	11.5 mL	16 x 100 mm	15-415	250	1000
99447-20	24.0 mL	20 x 125 mm	18-415	250	500
99447-20X	30.0 mL	20 x 150 mm	18-415	250	500


PYREX Disposable Glass Round Bottom Screw Cap* Culture Tubes, Bulk Pack Ordering Information

Cat. No.	Approx. Capacity	Approx. O.D. x Length	G.P.I. Thread Finish	Qty/ Pk	Qty/ Cs
99449-13	7.5 mL	13 x 100 mm	13-415	250	1000
99449-16	11.5 mL	16 x 100 mm	15-415	250	1000
99449-16X	15.0 mL	16 x 125 mm	15-415	250	1000
99449-16XX	19.0 mL	16 x 150 mm	15-415	250	1000
99449-20	24.0 mL	20 x 125 mm	18-415	250	500
99449-20X	30.0 mL	20 x 150 mm	18-415	250	500

PYREX Disposable Glass Flat Bottom Screw Cap* Culture Tubes Ordering Information

Cat. No.	Approx. Capacity	Approx. O.D. x Length	G.P.I. Thread Finish	Qty/ Pk	Qty/ Cs
99448-16X	11.5 mL	16 x 100 mm	15-415	250	1000
99448-16	17.0 mL	16 x 125 mm	15-415	250	1000
99448-19	29.5 mL	19.5 x 145 mm	18-415	250	500

*Caps are sold separately; see Cat. No. 99999.



PYREX® Disposable Glass Centrifuge Tubes

PYREX disposable glass centrifuge tubes are manufactured from Type I borosilicate glass.

- ▶ Available in four sizes: 5, 10, 15 and 50 mL
- ▶ Choice of cap styles: either threaded tops for disposable phenolic screw caps or with a tooled rim for secure snap cap sealing
- ▶ Caps must be purchased separately; see back page

PYREX Disposable Glass Snap Cap Centrifuge Tubes Ordering Information

Cat. No.	Approx. Capacity	Snap Cap† Part. No.	Approx. O.D. x Length	Max. RCF	Qty/ Cs
99501-5	5 mL	99500-1	13 x 110 mm	3,000	125
99501-10	10 mL	99500-2	16 x 114 mm	3,000	125
99501-15	15 mL	99500-2	17 x 126 mm	3,000	125
99501-50	50 mL	Tooled Rim*	29 x 137 mm	3,000	72

*Snap cap not available in the 50 mL size.

†Snap caps are sold separately.

PYREX Disposable Glass Screw Cap Centrifuge Tubes Ordering Information

Cat. No.	Approx. Capacity	Snap Cap† Part. No.	Approx. O.D. x Length	Max. RCF	Qty/ Cs
99502-5	5 mL	9998-13 99999-13	13 x 110 mm	3,000	125
99502-10	10 mL	9998-15 99999-15	16 x 114 mm	3,000	125
99502-15	15 mL	9998-15 99999-15	17 x 126 mm	3,000	125
99502-50	50 mL	9999-24 9998-24	29 x 137 mm	3,000	72

†Screw caps are sold separately.

Caps for PYREX Disposable Glass Centrifuge Tubes

Corning offers the following caps that fit the disposable glass centrifuge tubes.

Caps for PYREX Disposable Glass Centrifuge Tubes Ordering Information

Cat. No.	Description	G.P.I. Thread Style	Qty/ Cs
99500-1	Polyethylene Snap cap - fits 5mL Tube	n/a	500
99500-2	Polyethylene Snap cap - fits 10 mL and 15 mL Tube	n/a	500
9998-13	Reusable Phenolic Caps with Teflon® Liners - fits 5 mL Tube	13-415	288
9999-132	Reusable Phenolic Caps with Rubber Liners - fits 5mL Tube	13-415	288
99999-13	Disposable Phenolic Caps with Rubber Liners - fits 5mL Tube	13-415	1000
9998-15	Reusable Phenolic Caps with Teflon Liners - fits 10 and 15 mL Tube	15-415	288
9999-152	Reusable Phenolic Caps with Rubber Liners - fits 10 and 15 mL Tube	15-415	288
99999-15	Disposable Phenolic Caps with Rubber Liners - fits 10 and 15 mL Tube	15-415	1000
9998-24	Reusable Phenolic Caps with Teflon Liners - fits 50 mL Tube	24-400	12
9999-24	Reusable Phenolic Caps with Rubber Liners - fits 50mL Tube	24-400	12



Digital Stirring Hot Plates



PC-400D



PC-410D



Model PC-420D



Hot Plate Accessories

From left to right:
 440129 - Support Rod
 440141 - Holding Rod
 440140 - Boss Head Clamp
 6970SR - Stir Bar Retriever
 6795PR - Temperature Controller

The new Corning® digital hot plates and stirrers are made with the same durability and quality Corning has been putting into Pyroceram® top products since 1964. Digital displays are only available on 5" x 7" and 10" x 10" models. All models include Pyroceram top and a two year warranty.

- ▶ **Twin display stirring hot plate** – Truly affordable unit shows both temperature and stirring speed on separate digital displays to ensure accuracy time and time again.
- ▶ **Precise liquid temperatures and stirring speeds** – Temperature accuracy within $\pm 2^{\circ}\text{C}$ and stirring speeds at $\pm 5\%$ of the speed setting you choose
- ▶ **Digital displays** – Indicate when set temperature and stirring speeds are reached. Facilitate repeatability of standard operating procedures.
- ▶ **Safety indicators** – Power button and hot top indicators warn if unit is plugged in or is too hot to touch. Products have been independently tested and evaluated for compliance to product safety standards.
- ▶ **Angled front panel and large knobs** – Minimizes damage from spills and puts controls right at your fingertips

Corning Digital Hot Plates, Stirrers, and Stirring Hotplates Ordering Information

All models include Pyroceram® top and a two year warranty.

Hot Plates

Cat. No.	Model Number	Voltage	Temp. Range	Top Plate Size
6796-400D	PC-400D	230V/50Hz	5-550°C	5" x 7" (12.7 x 17.8 cm)
6796-600D	PC-600D	230V/50Hz	5-550°C	10" x 10" (25.4 x 25.4 cm)

Stirrers

Cat. No.	Model Number	Voltage	Stir Range	Top Plate Size
6796-410D	PC-410D	230V/50Hz	60-1150 RPM	5" x 7" (12.7 x 17.8 cm)
6796-610D	PC-610D	230V/50Hz	60-1150 RPM	10" x 10" (25.4 x 25.4 cm)

Stirring Hot Plates

Cat. No.	Model Number	Voltage	Temp./Stir Range	Top Plate Size
6796-420D	PC-420D	230V/50Hz	5-550°C/60-1150 RPM	5" x 7" (12.7 x 17.8 cm)
6796-620D	PC-620D	230V/50Hz	5-550°C/60-1150 RPM	10" x 10" (25.4 x 25.4 cm)

Accessories

Cat. No.	Description
6795PR*	Temperature Controller for digital hot plates and stirring hot plates
6970SR	Stir Bar Retriever
6795KIT	Accessory Kit (includes 6795PR - External temperature controller, 440129 - Support rods, 18," 6970SR - Stir bar retriever)
440129	Stainless Steel Support Rod – two 9" (45.7 cm) rods, $\frac{1}{16}$ " diameter. Rods can be screwed together to 18"
400084	Stainless Steel Support Rod – 12" (30.5 cm), $\frac{1}{16}$ " diameter.
440140†	Boss Head Clamp for use with $\frac{1}{16}$ " support rod and holding rod
440141	Holding Rod – thermocouple, stainless steel
400430	Stir Bar, 3/8" x 2," for use with models PC-610/620
401435	Stir Bar, 3/8" x 1," for use with models PC-210/220/410/420
409000	Stir Bar, 3/4" x 3," for use with model PC-611 stirrers
409021	Stir Bar, 5/16" x 1," for use with model PC-171 stirrers

*Temperature controllers are not available for use with stirrer only units.

†Bosshead clamps to support rod, holding rod connects to bosshead. Temperature controller slides into holding rod.



Alphabetical Index

8 Well PCR Tube Strips	82
96 Well Collection Plates	76
96 Well Filter Plates	76
96 Well Microplates	50-54, 75, 76, 79, 81
96 Well Polypropylene Blocks	75
96 Well Polypropylene Storage Microplates	55
96 Well Polystyrene Stripwell™ Microplates	54
96 Well UV Microplates	79
245 mm Square BioAssay Dishes	75
384 Thermowell® Gold Microplates	81
384 Well Collection Plates	76
384 Well Filter Plates	76
384 Well Microarray Printing Plates	93
384 Well Microplates	57-58, 60, 76, 81, 93
384 Well Polypropylene Storage Microplates	60
384 Well Polystyrene Microplates	58
384 Well UV Microplates	79
1536 Well Microplates	61

A

Aspirator	105
Assay Microplates	50-54, 57, 80

B

Background Reduction Kit	96
Bar Code Customization	84
BioAssay Dishes	75
Bottles	109

C

Caps	27, 134
Cell Culture Dishes	13
Cell Culture Flasks	9-12
Cell Lifters	38
Cell Scrapers	38
CellCube® Systems	30
CellSTACK® Culture Chambers	28
Centrifuge Tubes	113-114, 134
Centrifuge Tubes Filters	78
Cluster Tubes	55
Coated Slides	88-90
Collection Plates	76
Containers	110
Corning® CellBIND® Surface	3
Corning Cover Glass	92, 132
Cover Glass	92, 132
Cryogenic Accessories	116
Cryogenic Vials	116
Culture Dishes	13
Culture Flasks	75
Culture Tubes	5, 133
Cylinders	111

D

Digital Hot Plates	135
Digital Stirrers	135
Digital Stirring Hot Plates	135
Dishes	158-160
DNA-BIND™ Assay Microplates	80

E

E-Cube® Culture System	30
Epoxide Coated Slides	88
Epoxide Spotting Solution	92
Equipment	135
Erlenmeyer Flasks	36, 111

F

Flasks	9-12
--------	------

G

GAPS™ II Coated Slides	90
------------------------	----

H

Hot Plates	135
Hot Plates Accessories	135
Hybridization Chambers	96

L

Lids	65
Liquid Transfer System	105

M

Microarray Printing	88
Microarray Processing	94
Microarray Slide Mailers/Storage Boxes	93
Microarray Storage Pouches	94
Microcentrifuge Tubes	115
Microplate Accessories	65
Microscope Slides	131
Microspatulas	112
Multiple Well Plates, Cell Culture	15

P

PCR Microplates	81
PCR Tube Strips	82
PCR Tubes	82
Pipet Tips	101
Pipets	99, 100, 129-131
Pipetting Aids	100
Pipettors	100
Polycarbonate Microplates	81
Polypropylene Blocks	75
Polystyrene Microplates	58, 60
Printing Plates	93
Pronto!™ Background Reduction Kit	96
Pronto! Epoxide Spotting Solution	92
Pronto! Universal Hybridization Kits	94
Pronto! Universal Spotting Solution	91
Protein Crystallization Microplates	63

Alphabetical Index (CONTINUED)

R

Reagent Reservoirs	104
Robolids	66
Roller Bottles	26

S

Sealing Mats	82
Sealing Tape	66, 82
Spatulas	112
Spinner Flasks	31
Spin-X® Centrifuge Tube Filters	78, 108
Spotting Solution	91, 92
Stirrers	135
Stirring Hot Plates	135
Storage Bottles	109
Storage Mats and Accessories	65
Storage Microplates	55, 60
Stripwell Microplates	54
Syringe Filters	108

T

Thermowell® Gold Microplates	81
Transtar-96® Well Liquid Transfer System	105
Transwell® Permeable Supports	20
Tubes, Cluster	55
Tubes, Culture	133

U

UltraGAPS™ Coated Slides	89
Ultra-Low Attachment CellSTACK® Chambers	6
Ultra-Low Attachment Dishes	6
Ultra-Low Attachment Flasks	6
Ultra-Low Attachment Plates	6
Ultra-Web™ Synthetic Surface	6
Universal Hybridization Kits	94
Universal Spotting Solution	91
UV Microplates	79

V

Vacuum Filters	105-107
Volume Adapter	76
Volume Applicator	76

Catalog Number Index

Cat. No.	Page No.	Cat. No.	Page No.	Cat. No.	Page No.	Cat. No.	Page No.	Cat. No.	Page No.
2865-2949									
2865	132	3165	30	3340	5, 17, 52	3470	22	3575	58
2870	92, 132	3166	30	3341	65	3471	6, 15	3575BC	58
2875	132	3200	30	3342	60, 75	3472	22	3582	162
2935	132	3201	30	3343	55	3473	6, 15	3583	24
2940	132	3206	115	3344	55	3474	6, 17, 51	3584	76
2947	131	3207	115	3345	17, 66	3477	25	3585	17, 51
2948	131	3208	115	3346	65	3478	25	3590	51
2949	131	3209	115	3347	60, 75	3479	25	3591	51
		3213	115	3355	55	3480	25	3595	17, 51
		3260	14	3356	55	3491	22	3596	17, 51
3000-3998		3261	6, 13	3357	55	3492	22	3598	17, 51
		3262	6, 13	3359	55	3493	22	3599	17, 51
3000	13	3264	30	3360	51	3494	22	3600	52
3001	13	3268	28	3361	51	3495	22	3601	52
3003	112	3269	28	3362	17, 52	3496	22	3603	52
3004	112	3270	28	3363	55	3504	76	3603	17, 52
3005	112	3271	28	3364	55	3505	76	3604	52
3006	112	3272	28	3365	55	3506	15	3605	52
3007	112	3275	10	3366	51	3507	76	3610	17, 52
3008	38	3276	10	3367	51	3510	76	3614	17, 52
3010	38	3281	29	3368	51	3511	76	3615	52
3011	38	3282	29	3369	51	3512	15	3620	115
3012	112	3283	29	3370	51	3513	15	3621	115
3013	112	3284	29	3371	55, 76	3514	76	3622	115
3020	30	3286	30	3372	52	3516	15	3628	51
3021	30	3287	30	3374	24	3517	25	3631	52
3022	30	3289	4, 10	3378	23	3519	25	3632	52
3040	30	3290	4, 10	3379	23	3520	25	3635	53, 79
3041	30	3291	4, 11	3380	24	3521	25	3640	58
3055	10	3292	4, 13	3381	24	3524	15	3640BC	58
3056	10	3293	4, 13	3382	24	3526	15	3641	51
3059	11	3294	5, 13	3383	24	3527	15	3642	52
3067	5, 11	3295	5, 13	3384	24	3531	76	3643	59
3068	5	3296	5, 13	3385	24	3533	76	3650	52
3068	11	3298	4, 13	3386	24	3540	59	3651	52
3069	11	3300	5, 17, 51	3387	24	3542	18, 59	3653	59
3070	11	3302	30	3388	24	3544	59	3655	59
3071	11	3303	6, 28	3391	24	3548	15	3656	60, 93
3073	4, 11	3304	30	3392	24	3550	63	3657	60, 76
3080	65	3310	5, 28	3395	23	3551	63	3658	60
3081	65	3311	5, 28	3396	23	3552	63	3662	18, 58
3083	65	3312	5, 28	3397	23	3553	63	3663	18, 59
3085	65, 93	3313	28	3398	23	3554	63	3664	18, 59
3087	82	3319	28	3399	23	3555	63	3665	17, 51
3089	66	3320	5, 28	3401	21	3556	63	3666	17, 52
3090	66	3321	5, 28	3402	21	3557	63	3667	17, 52
3095	66	3324	29	3407	23	3560	162	3672	60, 93
3098	65	3328	29	3412	21	3561	31	3673	58
3099	17, 65, 93	3329	29	3413	21	3562	31	3674	58
3100	157	3330	5, 28	3414	21	3563	31	3675	59
3101	30	3331	29	3419	21	3564	31	3676	58
3135	30	3332	29	3420	21	3565	31	3677	58
3136	30	3333	29	3421	21	3567	31	3678	58
3138	30	3334	29	3422	21	3569	31	3679	79
3139	30	3335	5	3428	21	3570	18, 58	3679	53, 79
3144	30	3335	15	3450	22	3571	18, 58, 82	3680	58
3150	11	3336	15	3452	22	3572	58	3682	53
3151	11	3337	5, 15	3458	22	3573	58	3683	5, 18, 59
3152	31	3338	5, 15	3460	22	3574	58	3686	52
3153	31	3339	29	3462	22	3574BC	58	3688	17

Catalog Number Index (CONTINUED)

Cat. No.	Page No.	Cat. No.	Page No.	Cat. No.	Page No.	Cat. No.	Page No.	Cat. No.	Page No.
3690	51	3814	6, 1	3957	55, 75	4713	102	4962	100
3693	52	3815	6, 1	3958	55, 75	4714	102	4963	100
3694	52	3816	11	3959	55, 75	4715	102	4964	100
3695	51	3820	58	3960	55, 75	4780	102	4975	100
3696	17, 51	3821	58	3961	55, 75	4781	102		
3697	17, 51	3821BC	58	3964	60, 75	4783	103	6509-6970	
3699	80, 82	3822	18, 58	3965	60, 75, 76	4785	103	6509	81
3700	58	3824	58	3968	29	4786	102	6511	80, 81
3701	18	3824BC	58	3969	29	4787	102	6524	66
3701	58	3826BC	58	3985	18, 59	4801	103	6530	82
3702	58, 76	3849	65	3990	52	4803	102	6531	82
3702BC	58	3850	62	3991	52	4804	102	6542	82
3706	59	3851	62	3992	52	4806	102	6547	82
3707	18, 59	3852	62	3993	52	4807	103	6551	81
3711	59	3853	19, 62	3994	53	4808	103	6555	82
3712	18, 59	3854	19, 62	3995	53	4809	103	6569 ...	66, 80, 82, 93
3721	53	3855	19, 62	3997	17, 51	4810	103	6570	66, 82
3723	58	3857	62	3998	53	4815	104	6575 ...	64, 66, 80, 82
3724	61	3858	62			4821	103	6795	135
3724BC	61	3870XX1	8, 13	4010-4975		4823	103	6796	135
3725BC	61	3871XX1	8, 13	4010	99	4826	104	6970SR	135
3726	19, 61	3872XX1	8, 17	4011	99	4830	104		
3727	19, 61	3873XX1	8, 17	4012	99	4834	104	7000-7610	
3727BC	61	3875	17, 52	4020	99	4840	104	7000	99
3728	61	3880	53	4021	99	4844	102	7007	6, 17
3728BC	61	3881	53	4050	99	4845	102	7015	99
3729	61	3882	17, 53	4051	99	4846	102	7016	99
3729BC	61	3883	53	4100	99	4853	103	7017	99
3730	62	3884	53	4101	99	4854	103	7041	99
3740	82	3885	17, 53	4250	99	4860	101	7042	99
3741	82	3886	53	4251	99	4862	102	7045	99
3742	82	3887	53	4401	56	4863	101	7058	130
3743	82	3891	61	4408	56	4864	101	7077	129
3744	82	3891BC	61	4410	56	4865	101	7078	129, 130
3745	82	3893	19, 61	4411	56	4866	102	7079	130
3746	82	3893BC	61	4412	56	4867	101	7078B	131
3747	82	3895	61	4413	56	4868	102	7095D	131
3748	82	3895BC	59	4418	56	4870	104	7605	105
3749	82	3896	51	4484	99	4871	104	7606	105
3750	82	3897	51, 76	4485	99	4872	104	7610	105
3751	81	3898	51	4486	99	4873	104		
3752	81	3903	17, 53	4487	99	4876	105	8160-8950	
3753	81	3904	17, 53	4488	99	4877	105	8160	78, 108
3754	81	3907	4, 26	4489	99	4878	105	8161	78, 108
3755	81	3912	52	4490	99	4880	101	8162	78, 108
3756	81	3913	54	4491	99	4884	103	8163	78, 108
3757	81	3914	54	4492	99	4888	101	8169	78, 108
3773	64	3915	52	4500	32, 99	4894	104	8170	78, 108
3775	64	3916	17, 52	4501	99	4901	104	8388	109
3785	64	3917	17, 52	4502	32	4910	100	8390	109
3788	51	3922	52	4504	32	4911	100	8393	109
3789	52	3923	54	4510	32	4914	100	8396	109
3790	17	3924	54	4512	32	4922	100		
3792	52	3930	17, 65	4514	32	4923	100	9017-9999	
3795	51, 76	3931	17, 65	4515	35	4930	105	9017	51
3797	51	3935	65	4519	33, 34	4931	105	9018	51
3798	51	3936	61	4558	113	4958	101	9032	101
3799	17, 51	3937	61	4710	102	4959	100	9099	100
3801	23	3950	167	4711	102	4960	100	9102	17, 54
3802	23	3956	55, 75	4712	102	4961	100		

Cat. No.	Page No.	Cat. No.	Page No.	Cat. No.	Page No.	Cat. No.	Page No.	Cat. No.	Page No.
9186	100	401435	135	430492	117	430897	114	431225	108
9998	134	401654	30	430499	118	430909	115	431227	108
9999	134	401661	34	430512	107	430915	115	431229	108
10010-10024		402576-409831		430513	107	430917	115	431231	108
10010	5, 13, 27	402576	35	430514	107	430921	114	431252	36, 112
10024	5, 13, 27	402577	35	430515	106	430944	105	431253	36, 112
40001-40090		402579	35	430516	106	430945	105	431255	36, 112
40001	96	402614	35	430517	106	431079-431449		431256	36, 112
40003	90	402645	35	430518	109	431079	13	431272	75
40004	90	402648	34	430521	107	431080	13	431301	75
40005	90	402649	34	430525	118	431081	13	431306	13
40006	90	402650	34	430526	118	431082	13	431318	26
40015	89	409000	135	430588	14	431085	13	431321	26
40016	89	409021	135	430589	14	431096	106	431328	4, 13
40017	89	430015-430945		430591	14	431097	106	431329	4, 26
40018	89	430015	107	430597	14	431098	106	431339	37, 112
40019	89, 91	430049	107	430599	13	431110	13	431340	37, 112
40024	89, 94	430052	113	430624	107	431111	14, 75	431344	4, 26
40025	89, 91	430053	113	430625	107	431117	107	431346	5, 13
40026	94	430055	113	430626	107	431118	107	431355	113
40027	91	430157	25	430627	107	431119	118	431363	37, 112
40028	94	430165	13	430639	10	431120	118	431364	37, 112
40029	96	430166	13	430641	10	431121	118	431372	37
40030	94	430167	13	430656	117	431123	114	431386	117
40040	88, 92	430168	10	430658	116	431124	114	431401	36
40041	88	430172	25	430659	116	431131	118	431403	36
40042	88	430179	110	430661	116	431132	27	431404	36
40043	88	430180	110	430662	116	431133	26	431405	36
40044	88	430181	110	430663	116	431134	4, 26	431406	36
40047	92	430182	111	430669	26	431135	26	431407	36
40048	94	430183	36, 75, 111	430698	27	431143	36, 75, 111	431408	36
40080	96	430186	106	430720	10	431144	36, 75, 111	431430	109
40081	93	430195	26	430725	10	431145	36, 75, 111	431431	109
40082	93	430196	13	430756	106	431146	36, 75, 111	431432	109
40085	94	430236	114	430758	106	431147	36, 75, 111	431433	109
40086	94	430281	109	430766	113	431153	106	431444	37
40090	94	430282	109	430767	106	431154	106	431445	37
99445-99502		430289	117	430768	106	431155	106	431446	37
99445	133	430289	117	430769	106	431160	107	431447	37
99447	133	430290	113	430770	106	431174	107	431448	37
99448	133	430291	113	430771	106	431175	109	431449	37
99449	133	430293	13	430773	106	431191	26	440129-440141	
99500	134	430304	113	430776	114	431198	26	440129	135
99501	134	430314	106	430790	113	431205	106	440140	135
99502	134	430320	106	430791	113	431206	106	440141	135
400430-401661		430372	10	430823	11	431212	108		
400430	135	430421	36, 75, 111	430824	11	431215	108		
400640	35	430422	36, 75, 111	430825	11	431218	108		
400649	35	430487	117	430828	113	431219	108		
401392	34	430488	117	430829	113	431220	108		
		430489	117	430849	26	431221	108		
		430490	117	430851	26	431222	108		
		430491	117	430852	26	431224	108		
				430853	26				

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.