

Thermo Scientific Nunc Immobilizer

Streptavidin and Amino

Excellent binding capacity

A unique surface chemistry has been created for the Thermo Scientific Nunc Immobilizer plates and strips which provides excellent covalent coupling for optimal target molecule orientation and stability.

The concentrations of target molecule required for coupling is usually very low compared with what is needed for passive binding. The plate surface provides a high affinity surface especially for the coupling of small ligands.




Fig. 1.

Immobilizer Streptavidin plates were compared with 4 other suppliers (A, B, C, and D) of streptavidin plates, in an ELISA. Low concentrations of biotinylated anti-human IgG, (in 6 concentrations) was coupled to the streptavidin surface. Following addition of human IgG, the amount of bound human IgG was detected by adding HRP conjugated anti-human IgG and visualized by incubation with OPD substrate and measuring at 492 nm. The result shows that the required amount of added human IgG appeared to be approx. half of the amount that is required for the other suppliers' plates when measured at an OD of 1.

High signal to noise ratio

Assays conducted on the Nunc™ Immobilizer™ plates and strips are characterized by high signal readouts (Fig. 1).

This, in combination with no or very low non-specific binding and thus negligible background noise, ensures superior signal to noise ratios (Fig. 2).

Highly reproducible results

Immobilizer plates and strips are manufactured for consistent and reproducible assay results.





Fig. 2.

The signal to noise ratio was obtained on streptavidin plates from 4 different suppliers (A, B, C, and D). They were compared to the signal to noise ratio obtained on Immobilizer Streptavidin plates in an assay, where a biotinylated 15-mer DNA oligonucleotide bound to the streptavidin, was hybridized to a 98 bp DIG-incorporated pUC19 amplicon (400 fmol/well). The amount of hybridized DNA was detected by adding an anti-DIG antibody conjugated to HRP and reading of the signal at 450 nm following addition of TMB. The noise signal was obtained from all plates by following an identical procedure, but replacing the pUC19 plasmid template with water.

The covalently bound monolayer coating of the plates is highly uniform. Leaching is minimal, even when using very stringent washing steps. These characteristics make the Immobilizer plates and strips highly suitable for many applications.

Mechanism of action Immobilizer Streptavidin



Coupling of A) biotinylated DNA oligonucleotide, B) biotinylated peptide, and C) biotinylated protein to the streptavidin conjugated photoprobes on the solid phase of the Streptavidin Immobilizer plates.

No blocking

In contrast to most other coated plates, the Immobilizer plates and strips do not require separate costly, labor intensive and time-consuming blocking steps to prevent non-specific binding.

This is also due in part to the placement and density of the active sites and their spacer arms. These create a “molecular shield” which prevents direct access to the polystyrene hydrophobic surface, where non-specific binding could otherwise take place.

Higher signal using low sample concentration

The Immobilizer products are characterized by the high signals obtained even with very low sample concentrations.

As illustrated in Fig. 3, in certain cases when samples at one tenth concentration are used, the signals obtained on the Immobilizer surface are greater than those obtained using normal high binding ELISA plates and passive absorption.

Short incubation time


Coupling of target molecules to the Immobilizer surface takes place very quickly, compared to the time required for stable coating by passive adsorption.

The Immobilizer surface can be saturated after approximately 4 hours incubation with gentle shaking.


The coupling efficiency is typically unaffected by pH. However, one must be aware that pH can conformationally change the biomolecules that are bound. The presence of non-ionic detergents suppresses the coupling reaction.

Long term stable at room temperature


The shelf life of the Immobilizer plates and strips exceeds one year, when they are stored in their unopened packaging at room temperature. The manufacturing date is printed on the package. The Immobilizer Streptavidin has a minimum shelf life post manufacturing of 18 months. For Immobilizer Amino surface the period is at least 2 years.



Mechanism of action Immobilizer Amino



Coupling of D, peptides or E, proteins to the Immobilizer Amino surface.



www.thermoscientific.com/oemdiagnostics

© 2010 Thermo Fisher Scientific Inc. All rights reserved. “Immobilizer” is a trademark of Exiqon A/S, Vedbaek, Denmark. The product is produced under license from Exiqon A/S - EP 0820483 and foreign applications and patents. All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Austria

+43 1 801 40 0

Belgium

+32 53 73 42 41

China

+86 21 68654588

Denmark

+45 4631 2000

France

+33 2 2803 2180

Germany

+49 6184 90 6940

India

+91 22 6716 2200

Italy

+39 02 02 95059 or
434-254-375

Japan

+81 3 3816 3355

Netherlands

+31 76 571 4440

Nordic/Baltic countries

+358 9 329 100

North America

+1 585-586-8800

Russia/CIS

+7 (812) 703 42 15

Spain/Portugal

+34 93 223 09 18

South America

+1 585 899 7298

Switzerland

+41 44 454 12 12

UK/Ireland

+44 870 609 9203

Other Asian countries

+852 2885 4613

Countries not listed

+49 6184 90 6940 or
+33 2 2803 2180

TILSPNUNCTN64 0510