

PRODUCT INFORMATION

Thermo Scientific Maxima Hot Start Taq DNA Polymerase

#EP0603 5 x 500 U

Lot: _ Expiry Date: _

Concentration: 5 U/ μ L

Supplied with: 10 x 1.25 mL of 10X Hot Start
PCR Buffer
10 x 1.25 mL of 25 mM MgCl₂

Store at -20°C

www.thermoscientific.com/onebio

Description

Thermo Scientific Maxima Hot Start *Taq* DNA Polymerase is designed to enhance the specificity, sensitivity and yield of DNA amplification (1-4). In addition, the enzyme provides the convenience of reaction set-up at room temperature. Maxima® Hot Start *Taq* DNA Polymerase is a recombinant *Taq* DNA polymerase which has been chemically modified by the addition of heat-labile blocking groups to its amino acid residues. The functional activity of the enzyme is restored during a short 4-minute incubation at 95°C. The activated enzyme maintains the same functionality as *Taq* DNA polymerase: catalyzes 5'→3' synthesis of DNA, has no detectable 3'→5' proofreading exonuclease activity, but possesses 5'→3' exonuclease activity. It exhibits deoxynucleotidyl transferase activity, which frequently results in the addition of extra adenines at the 3'-end of PCR products. Before activation, the two activities are not detectable.

Applications

- Hot start PCR.
- RT-PCR.
- Highly specific amplification of complex genomic and cDNA templates up to 3 kb.
- Amplification of low copy DNA targets.
- Real-time PCR.
- Multiplex PCR.
- Generation of PCR products for TA cloning.

Source

E. coli cells carrying a cloned *pol* gene from *Thermus aquaticus*.

Definition of Activity Unit

One unit of the enzyme catalyzes the incorporation of 10 nmol of deoxyribonucleotides into a polynucleotide fraction in 30 min at 74°C.

Storage Buffer

Enzyme is supplied in: 20 mM Tris-HCl (pH 9.0), 1 mM DTT, 0.1 mM EDTA, 100 mM KCl, 0.5% (v/v) Tween 20, 0.5% (v/v) Nonidet P40 and 50% glycerol.

10X Hot Start PCR Buffer

200 mM Tris-HCl (pH 8.3 at 25°C), 200 mM KCl, 50 mM (NH₄)₂SO₄.

Inhibition and Inactivation

- Inhibitors: ionic detergents (deoxycholate, sarkosyl and SDS) at concentrations higher than 0.06, 0.02 and 0.01%, respectively (5).
- Inactivated by phenol/chloroform extraction.

CERTIFICATE OF ANALYSIS

Endodeoxyribonuclease Assay

No conversion of covalently closed circular DNA to nicked DNA was detected after incubation of 25 units of Maxima Hot Start *Taq* DNA Polymerase with 1 µg of pUC19 DNA for 4 hours at 37°C.

Deoxyribonuclease Assay

No degradation of single-stranded and double stranded [³³P]-labeled oligonucleotide was observed after incubation with 10 units of Maxima Hot Start *Taq* DNA Polymerase for 4 hours at 37°C.

Ribonuclease Assay

No contaminating RNase activity was detected after incubation of 80 ng of 2 kb RNA transcript with 25 U of Maxima Hot Start *Taq* DNA Polymerase for 4 hours at 37°C.

Functional Assay (PCR)

Maxima Hot Start *Taq* DNA Polymerase was tested in PCR amplification of a 354 bp DNA fragment of single copy ApoE gene from 50 pg of human genomic DNA in a 20 µL PCR mixture containing 1X Hot Start PCR Buffer, 2 mM MgCl₂, 0.2 mM dNTP, 0.5 µM primers and 0.5 units of Maxima Hot Start *Taq* DNA Polymerase. The PCR products were analyzed on 1% agarose gel.

Quality authorized by:

 Jurgita Zilinskiene

(continued on reverse page)

PROTOCOL

To prepare several parallel reactions and to minimize the possibility of pipetting errors, prepare a PCR master mix by adding water, buffer, dNTPs, primers and *Taq* DNA Polymerase. Prepare enough master mix for the number of reactions plus one extra. Aliquot the master mix into individual PCR tubes and add template DNA. Keep all reaction components on ice. Reaction set up can be performed at room temperature.

1. Gently vortex and briefly centrifuge all solutions after thawing.
2. Add the following components for each 50 μ L reaction at room temperature:

10X Maxima Hot Start <i>Taq</i> buffer	5 μ L
dNTP Mix, 2 mM each (#R0241)	5 μ L (0.2 mM of each)
Forward primer	0.1-1 μ M
Reverse primer	0.1-1 μ M
25 mM MgCl₂	1-4 mM*
Template DNA	1 pg - 1 μ g
Maxima Hot Start <i>Taq</i> DNA Polymerase	1.25-2 U
Water, nuclease-free (#R0581)	to 50 μ L
Total volume	50 μL

*Volumes of 25 mM MgCl₂ solution to be used:

Final concentration of MgCl ₂ , mM	1	1.5	2	2.5	3	4
Volume of 25 mM MgCl ₂ for 50 μ L reaction, μ L	2	3	4	5	6	8

3. Gently vortex and spin down the samples.
4. When using thermal cyclers without a heated lid, overlay the reaction mixture with 25 μ L mineral oil.
5. Perform PCR using recommended thermal cycling conditions:

Step	Temperature, °C	Time	Number of cycles
Initial denaturation/ enzyme activation	95	4 min	1
Denaturation	95	0.5-1 min	25-40
Annealing	T _m -5	0.5-1 min	
Extension	72	1 min/kb	
Final extension	72	5-15 min	1

References

1. D'Aquila, R.T., et al., Maximizing sensitivity and specificity of PCR by preamplification heating, *Nucleic Acids Res.*, 19, 3749, 1991.
2. Kellog, D.E., et al., TaqStart antibody: "Hot start" PCR facilitated by a neutralizing monoclonal antibody directed against *Taq* DNA polymerase, *BioTechniques*, 16, 1134-1137, 1994.
3. Horton, R.M., Hoppe, B.L., and Conti-Tronconi, B.M., AmpliGrease: "Hot Start" PCR using petroleum jelly, *BioTechniques*, 16, 42-43, 1994.
4. Dang, C. and Jayasena S.D., Oligonucleotide inhibitors of *Taq* DNA polymerase facilitate detection of low copy number targets by PCR, *J. Mol. Biol.*, 264, 268-278, 1996.
5. Weyant, R.S., et al., Effect of ionic and nonionic detergents on the *Taq* polymerase, *Biotechniques*, 9, 309-308, 1990.

NOTICE TO PURCHASER:

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 6,127,155, 5,677,152 (claims 1 to 23 only) and 5,773,258 (claims 1 and 6 only). The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

LIMITED USE LABEL LICENSE: Internal Research and Development Use Only.

The purchase of this product conveys to the buyer the limited, non-exclusive, non-transferable right (without the right to resell, repackage, or further sublicense) to use this product for internal research and development purposes. No other license is granted to the buyer whether expressly, by implication, by estoppel or otherwise. In particular, the purchase of the product does not include or carry any right or license to use, develop, or otherwise exploit this product commercially and no rights are conveyed to the buyer to use the product or components of the product for purposes including but not limited to provision of services to a third party, generation of commercial databases or clinical diagnostics. This product is sold pursuant to authorization from Thermo Fisher Scientific and Thermo Fisher Scientific reserves all other rights. For information on purchasing a license for uses other than internal research and development purposes, please contact outlicensing@lifetech.com or Out Licensing, Life Technologies Inc., 5791 Van Allen Way, Carlsbad, California 92008.

PRODUCT USE LIMITATION

This product is developed, designed and sold exclusively *for research purposes and in vitro use only*. The product was not tested for use in diagnostics or for drug development, nor is it suitable for administration to humans or animals. Please refer to www.thermoscientific.com/onebio for Material Safety Data Sheet of the product.

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.