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
Measurement

Information: To measure the sample volume, usually overflow measuring beakers or measuring cylinders are used.

The volume is chosen according to the measured value expected. Choose the optimum volume for the value you expect to be measured. Too large measuring ranges will lead to inaccurate results. To estimate the measured value you can reckon with approx. 80% of the CBS value.

- **Sample preparation and filling of the measuring bottles** compare following instruction: DIN 38409 part 52: "Measurement of the oxygen consumption".
- **Rinse measuring bottle with sample.** Empty thoroughly.
- **Exactly measure the required oxygen-saturated (thoroughly homogenized) quantity of the sample according to information.**
- **Put the magnetic stirring rod into the bottle.**
- **Insert rubber quiver in the neck of the bottle.**
- **Put 2 sodium hydroxide tablets into the rubber quiver with tweezers.** (Caution: The tablets must never come into the sample!)
- **Tightly close the measuring bottle with the OxiTop.**
- **Start measurement:** Press the M and S keys simultaneously for 2 seconds until the display shows 00.
- **Keep the measuring bottle with the OxiTop put on for 5 days at 20°C** (e.g. in an incubator). After the meas. temperature has been reached (after 1 hour at the earliest, after 3 hours at the latest; **AutoTemp function**), the OxiTop automatically starts the measurement of the oxygen consumption.
- **During the 5 days the sample is continuously stirred.** The OxiTop automatically stores one value every 24 hours for 5 days. To have the current value shown press the M key.
- **Readout of the stored values after the 5 days have passed:** Press the S key.
- **Convert the displayed measured value (digits) into the BOD value with the following table (Digits x Factor = BOD5 in mg/l):**

Key	Display		BOD5 Meas. range mg/l	Sample- volume ml	Factor
●	Measuring day				
S	either stored value, e.g.:	or system message "memory empty" (= no measured value):	0-40	432	1
S	1 sec		0-80	365	2
S	1 sec		0-200	250	5
S	1 sec		0-400	164	10
S	1 sec		0-800	97	20
S	1 sec		0-2000	43,5	50
S	1 sec		0-4000	22,7	100

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- *The scale shows zero or a too low value.
The measuring equipment is not water-tight.
Check rubber quiver, screw top and bottle.*
- *Measured values below zero.
The temperature of the sample had not sufficiently been adjusted ($< 20^{\circ}\text{C}$).*
- *Measuring range exceeded.
The measuring range chosen is too small. With very high values ($> 4000 \text{ mg/l}$)
predilute the sample.*

Errors due to procedure have not been mentioned.

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*Value remains below
measuring range*

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*Change batteries
(approx. every 3 years)*

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
*Value exceeds
measuring range*

Batteriewechsel Battery replacement

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- 2** - Remove batteries.
- Put in new batteries,
pay attention to polarity!
(e.g. WTW-Modell Batt/Oxitop,
Order No. 209 012).
- Place upper lid with tap on locking
pin (lower lid). Attention! Do not
squeeze cable connection.
- Close upper lid
(make the hook click).

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Disposal note:

Please dispose of the lithium batteries at a battery collection point.

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Labelling:

To label the bottles, WTW offers marking rings with numbers, which can be put on the neck of the bottle.
Model **MARK-6** (6 pcs) Order No. 209 013.

Cleaning:

OxiTop: with soft cloth and aqueous soap solution
(no alcohol or acetone!)

Checking the system:

see application report 895233.

Further literature:

- Application report 895230: "BOD measurement in household wastewater"
- Application report 895231: "BOD measurement in organically heavily contaminated wastewater"
- Application report 895232: "BOD measurement with presence of impeding or toxic substances"

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