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## **Operating instructions for electrodes for measuring pH values and redox voltages (ORP)\***

### **Preparations and general information**

If present, remove the watering cap covering the membrane and diaphragm. This cap contains potassium chloride solution\*\* (water for pH glass electrodes). The electrode is ready for making measurements. Electrodes which have been stored in the dry state must be soaked for at least 24 hours in potassium chloride solution\*\* before use. The water in the solution keeps the diaphragm in a properly functioning state and produces the hydrated layer on the membrane glass surface which is essential for making pH measurements. If an air bubble is present inside the electrode it must be removed by shaking in a manner similar to resetting a medical thermometer. If there is insufficient potassium chloride solution\*\* in the electrolyte compartment of the reference system, top-up as necessary.

This refilling operation is unnecessary for low maintenance electrodes with gel filling or REFERID system. Watering with potassium chloride solution is particularly important for low maintenance electrodes to prevent drying out of the diaphragm.

The cap on the refilling port must be taken off for calibration and for making measurements. The membrane and the diaphragm must be immersed in the sample solution whereby a hydrostatic overpressure of at least 5 cm solution column should be present in the electrolyte compartment.

When using combinations with separate electrodes (separate pH glass electrode and reference electrode), make sure that the internal reference systems are compatible, e.g. it is essential to use an Ag/AgCl reference electrode if a glass electrode with Ag/AgCl internal reference element is being used.

### **Measuring pH values**

#### **Calibration\***

1. Set the temperature of the buffer solution on the pH meter. (Valid when using measuring chains with installed Pt 1000, e.g. N 1042 A.)
2. Take off the cap of the refilling port of the electrode. Rinse the electrode with distilled water.

\* Please pay attention too to the the parating instructions for the pH meters from SCHOTT GERÄTE.

\*\* 3 mol/l for Ag/AgCl, 3.5 mol/l for Ti/TiCl<sub>3</sub>, 4.2 mol/l for Hg/Hg<sub>2</sub>Cl<sub>2</sub>

3. Dip the electrode into the buffer solution which has a pH value in the vicinity of the zero point (pH = 6.87 or pH = 7.00). For indication of pH value of the buffer solution press "CAL" button or "zero point" control knob ( $\Delta$  pH).
4. Rinse the electrode with distilled water and then dip it into the second buffer solution, e.g. pH = 4.01 or pH = 4.00.
5. For indication of pH value of the buffer solution press "CAL" button or "slope" control knob (mV/pH). The pH meter is now matched to the characteristic of the electrode (calibrated). After rinsing with distilled water, the electrode is now ready for making measurements.
6. If a reproducibility of  $\Delta$  pH =  $\pm 0.01$  is required, buffer solutions according to DIN 19 266 and NIST must be used, e.g. pH = 6.87 and pH = 4.01. Furthermore, the buffer solution, sample solution and electrode must be temperature equilibrated to  $\Delta T = \pm 0.2$  K, e.g. at 25 °C.

#### Please note

pH electrodes from SCHOTT-GERÄTE have a linear characteristic over the range from pH = 0 to the onset of the alkali error which, depending on the concentration of  $\text{Na}^+$  ions, usually occurs above pH = 11.

Alkaline buffer solutions are very unstable due to absorption of carbon dioxide from the air. This can produce errors. Thus it is advisable to calibrate preferentially in the acidic range.

#### pH measurement\*

1. Set the temperature of the buffer solution on the pH meter. (Valid when using measuring chains with installed Pt 1000, e.g. N 1042 A.)
2. Rinse the electrode and dip then into the sample solution.
3. Take the pH reading on the pH meter after the reading has become stable.

#### Measuring redox voltages (ORP)

(mV measurements with platinum or gold combination electrodes)

An Ag/AgCl reference system is normally used in metal combination electrodes (e.g. platinum combination electrodes). No calibration is made. Redox test solutions are available for checking metal combination electrodes, e.g. +470 mV platinum against Ag/AgCl, 3 mol/l KCl or +430 mV for platinum against calomel. The voltage reading with the combination electrode should agree with the nominal value to within  $\pm 10$  mV in this solution. Otherwise the electrode metal and/or the diaphragm must be cleaned as described under "important instructions".

For making redox voltage measurements, switch the pH meter to the mV range. Dip the electrode into the sample solution and take the reading.

For correct recording of the results, the type of reference system (Ag/AgCl) and the temperature should be noted in addition to the mV reading.

\* Please also pay attention to the operating instructions for the pH meters from SCHOTT-GERÄTE.

## Important instructions

### Customer application

SCHOTT electrodes are designated with the type number and serial numbers. The type number is located on the glass shaft or on the foil of pH glass and metal electrodes. The serial number is on the plug head or on the cable close to the electrode head.

### Electrode connection

The SCHOTT plug head system is reliable and simple to use. The provided 10.5–1.5 O-ring is necessary to protect the plug contact against moisture.

When the connection to the measuring unit is to be made via terminals or if a plug is to be attached, make sure that the black semiconductor layer over the cable core insulation is removed for at least 3 mm back from the bare inner conductor. Use resin-cored solder for making soldered connections.

When using screw plug head electrodes, make sure that a 10–3.5 O-ring and a 18/13/4 washer (contained in the scope of delivery) are used as gasket material.

### Marking

Suitable electrodes are available from SCHOTT-GERÄTE for all customary applications (pH, temperature and pressure measurements). The right electrode for a particular purpose can easily be found in our current catalog. The service life of electrodes will be shortened if they are used at the limits or outside the range of normal operating conditions.

It is essential to use electrodes with glass shaft for making measurements in organic solvents. Special electrodes must be used for making measurements at low temperatures ( $< 0$  °C). These electrodes contain low temperature electrolyte (3 mol/l KCl solution with glycerol 1:1). pH glass electrodes for operation at low temperatures have low resistance membrane glass.

### Maintenance

The electrolyte inside combination electrodes and reference electrodes must be topped up or replaced from time to time. To ensure proper functioning of the electrode combination, observe the specified concentration of the potassium chloride solution exactly. If crystals have formed in the electrolyte compartments after long-term measurements or storage at elevated temperatures due to evaporation of water, carefully warm the electrode in a water bath to redissolve the crystals. Then draw off the electrolyte while it is still warm and replace it with fresh electrolyte. Potassium chloride may crystallize in the diaphragm of electrodes which are stored in the dry state. This clogs the diaphragm leading to hard sensitivity, unsteady readings, prolonged settle time and measuring errors. The crystallized potassium chloride can be redissolved by prolonged warming in diluted potassium chloride solution (1–2 mol/l). The provided filled watering cap should be used when storing the electrode.

## Cleaning

Contamination at the membrane and diaphragm leads to measurement errors and thus must be avoided or removed.

**pH glass membranes** can be cleaned by rubbing with a moistened paper tissue. Use a suitable solvent (e.g. acetone) to remove organic contamination. Treatment with dilute hydrochloric acid (1:1), pepsin solution acidified with hydrochloric acid (to remove protein contamination) or chrome sulfuric acid is possible too. SCHOTT electrodes with type S membrane glass should be cleaned in warm sodium hydroxide solution.

**Diaphragms** of combination electrodes and reference electrodes can be cleaned by treating them with dilute hydrochloric acid (1:1), pepsin solution acidified with hydrochloric acid (to remove protein contamination) or chrome sulfuric acid. Use a suitable solvent (e.g. acetone) to remove organic contamination. Proper functioning of externally clogged ceramic diaphragms can be restored by careful abrasion with fine grade sandpaper or with a fine grade diamond file. CAUTION! The pH glass membrane must not be scratched!

Never treat platinum diaphragms with mechanical methods! After chemical cleaning the diaphragm can be rinsed clear by sucking under vacuum or pushing on the watering cap so that air bubbles rise in the inner electrolyte solution.

Taper diaphragms should be made ready for operation before commencing measurements by slightly raising and then plugging the taper sleeve onto the drilled taper core. This produces increased outflow of electrolyte ensuring proper wetting of the taper surface. The refilling port should be open during this procedure.

**Metal electrodes** can be cleaned with scouring powder or with fine grade emery paper, exercising due caution during all such operations.

## The following points require your special attention

For making measurements in the medical and foodstuff chemistry fields and in all cases where physiological effects due to mercury or thallium compounds could be produced, use electrodes with Ag/AgCl reference system. This is particularly important too as a precaution against poisoning in the case of mechanical breakage of the electrodes.

When disposing of electrodes, observe the relevant legislation for prevention of environmental pollution.

All information contained in these operating instructions constitutes orientation data valid at the time of printing. SCHOTT GERÄTE may undertake changes and modifications at any time for technical or economic reasons as well as to comply with the legal stipulations of the various federal states.