Fisher
Bioblock Scientific

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Contents

Safety precautions	4-6
Warranty terms	7
Product	8
Operating panel	9
Rear panel	9
Mains voltage	10-11
Voltage setting	11
Changing the fuses	11
Starting the pump	12-13
Pumping against back-pressure	14
Overload protection	15
When the pump is not in use	16
Running-in period for tubing	17
Tubing life	17

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
Contents

Analog interface	18–20
Accessories	21–22
Pump heads	23–24
Flow-rates per channel for multi channel pump-heads	25–27
Flow-rates per channel for single-channel pump-heads	28–30
Maintenance	31
Repairs	31
Spare-parts	31
Technical specifications	32

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Please note

We recommend you to read this operating manual carefully.

When operating a pump, certain hazards cannot be excluded.

ISMATEC SA does not take liability for any damage resulting from the use of an ISMATEC® pump.

ISMATEC SA does not admit responsibility for the handling of chemicals.

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Safety precautions

ISMATEC® tubing pumps are designed for pumping applications in laboratories and industry. As such it is assumed that Good Laboratory Practice (GLP) and our following recommendations will be observed.

- The circuit between mains supply and pump has to be earthed.
- The pump must not be operated outside the destined operating and environmental conditions.
- Place the unit in a well ventilated position at least 10 cm away from walls, partitions etc. Ensure that curtains and similar materials do not cover the ventilation slits.
- The pump **must not** be used:
 - for medical applications on human beings
 - in explosion proof chambers or in the presence of flammable gases or fumes

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Safety precautions

- The pump must be switched off when cassettes or tubing are inserted or changed.
- The permeability of pump-tubing depends on the material used and pressure conditions. Tubing can also become electro-statically charged. Please be aware of possible hazards when laying tubing in explosion-proof chambers.
- Tubing can tear and burst during operation. If this could cause damage, the necessary safety measures based on the specific situation must be taken.
- Tubing pumps consist of revolving parts (rotor of the pump-head). Therefore, the pump must not be operated before the cassettes are fully snapped-in or the tube-bed completely shut.
- Do not manipulate the pump-head before the pump is switched off and disconnected from the mains.

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Safety precautions

- Be particularly cautious that no parts of your body such as fingers, long hair, etc. or jewellery, or loose objects such as cables or tubing, etc. can be trapped by the revolving rotor.
- Only new fuses, according to the specifications stated on page 10 in this manual, must be used. The fuse-holder must not be short-circuited.
- Do not open or remove the housing while the pump is operating.
- Repairs may only be carried out by a skilled person who is aware of the hazard involved.
- For service and repairs carried out by the customer or by third-party companies ISMATEC SA denies any responsibility.

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Warranty

For all parts manufactured by ISMATEC®
from date of delivery: 2 years
All other parts, excluding consumables: 1 year
We warrant the perfect functioning of our products, provided they have been installed and operated correctly according to our operating instructions.

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Warranty Terms

If production or material faults can be proved, the defective parts will be repaired or replaced free of charge at our discretion. A defective pump must be returned in the original ISMATEC® packing or in a packet of equal quality. For pump-heads from manufacturers other than ISMATEC SA the warranty terms of the specific manufacturer are valid. The duration of the warranty is not affected by making a claim for warranty service. Further claims are excluded. Shipping costs are charged to the customer.

Our warranty becomes invalid in case of

- improper operation by the user, or if the pump is diverted from its proper use,
- unauthorized modification or misuse by the user or by a third party,
- improper site preparation and maintenance,
- operation outside environmental and electrical specifications for the product,
- use of third-party software, hardware, interfacing or consumables purchased by the user and which do not comply with our specifications,
- damages caused by contamination or leaks due to torn or bursted tubing

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Product

Package contents

- BVP drive
- Pump head as ordered
- Power cord
- Operating instructions

➔ For operating the selected pump-head, please refer to the instruction manual separately supplied with the corresponding pump-head.

Please check the package and its contents for transport damage. If you find any signs of damage, please contact your local ISMATEC® representative immediately.

Complaints can only be accepted within 8 days from receipt of the goods.

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


Operating panel

- 1 Mains switch
- 2 Control light
- 3 Speed selector
999 = 240 rpm (max. speed)
adjustable in 0.1% steps
- 4 RUN/STOP
- 5 Rotation direction
- 6 MAX/SLOW key for drive speed
(ideal for fast filling or emptying
the system)

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⚠ Before starting-up

Check if the voltage setting visible
in the window of the fuse-holder
complies with your local mains
voltage.

If necessary, the voltage setting
must be changed and the 2 fuses
must be replaced.

⚠ Socket/Power cord

Use exclusively the original power
cord supplied with the pump. The
socket must be earthed (protective
conductor contact).

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

⚡ Mains voltage

Voltage setting 115V / 230V and changing the fuses

➡ Switch the pump off, pull out the mains plug.

1 Pull out the fuse-holder by opening it with a small screw-driver (size 0).

2 Take out the voltage selector plate. Turn it and re-insert it into the fuse-holder so that the required voltage rating is facing the window of the fuse-holder.

3 Insert 2 new fuses

➡ 230 V: 2 x 1,25 A (slow-blow)
➡ 115 V: 2 x 2,50 A (slow-blow)

⚠ Use always 2 slow-blow fuses of the same type complying with the local mains voltage.

4 Shut the fuse-holder. The voltage rating is visible in the window.

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Starting the pump

1 Press in the mains switch, the control light is on.

2 Set the required speed on the 3-digit speed selector.
(999 = 240 rpm)
→ The rotation speed can be adjusted while the pump is running.

3 Choose the rotation direction
→ Never change the rotation direction while the pump is running.

4 The MAX/SLOW- button allows rapid filling (or emptying) of the tubes while the pump is running.

5 Start the pump by switching the RUN/STOP button.

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
Pumping against pressure

The BVP can be used for continuous duty at a differential pressure of max. 1,5 bar. For short-time operation higher differential pressures may be managed depending on the mounted pump-head and tubing i.d. (smaller i.d. is preferable).

Depending on the application conditions the FMI pump-heads allow operations at differential pressures of up to 7 bar.


In case of any doubts please contact your ISMATEC® agent.

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Overcurrent protector

The tubing pump BVP *Standard* features an overload protector.

When the overcurrent protector is activated, the pump is stopped. In case of an overload, pin 12 of the analog interface (open-collector) is connected to the ground (pin 1) (see page 19).

Whenever this situation occurs, **the pump must be switched off immediately.**

➔ Let the pump cool down (it takes about 2 minutes until the pump is ready again)

Before the pump is re-started, it is most important to check the reason for the overload (eg. too high differential pressure).

➔ Only when the cause of the overload has been detected and the failure corrected accordingly may the pump be started again.

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When the pump is not in use

When the pump is idle, we recommend you to release the tubing from pressure.

This helps to protect the tubing from unnecessary strain and prolongs its service-life.

⚠ Syphoning effect

When the tubing is released from squeezing, the fluid can flow back to the reservoir.

a Tubing used with cassettes do not need to be removed. Releasing the cassette on one side is sufficient.

b The single-channel pump-heads 360, 380 and Easy Load allow opening the tube-bed in order to release the tubing pressure.

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
△ Please note

Please refer to our detailed tubing documentation.

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Running-in period for tubing

Every new tube requires a running-in period. If constant and reproducible flow rates are required, we recommend you to run new tubing in with water or the medium to be pumped for at least 1 to 3 hours before you start the application.

Tubing life

The service-life of the tubing depends on the application and the tubing material used. Example: chemicals, rotation speed, differential pressure, temperature, viscosity, pressure on tubing, etc.

General information on the service-life is stated in our ISMATEC® catalog (without obligation!).

In order to improve the lubrication and service-life of the tubing, we recommend users to spray both the tubing and the pump rollers with our silicone oil spray (Order No. SC0179).

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Analog interface

Pin 1, (GND) Ground

Reference potential for all other inputs.

Pin 3, start

In remote operation (pin 2 to GND) the pump starts when connected to pin 1 (GND)

Pin 4, direction

In the open position the pump turns counter-clockwise; when connected to pin 1 (GND) it turns clockwise.

Pin 4, direction

If this pin is left open, the pump runs in clockwise direction. Connecting pin 4 to the ground (pin1) makes the pump run in counter-clockwise direction.

Pin 5, speed IN

External speed control (0–5V, 0–10V, 0–20mA, 4–20mA) Input impedance and input range can be selected via a dip-switch inside the pump (see page 20).

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

Pin 7, +36 V_{dc}
About +36 V_{dc} (unstabilised) are available the pump motor (max. load 1 A).

Pin 9, speed OUT
The default setting is 0–10 V_{dc}, proportionally to the motor speed 2,4–240 rpm. Alternatively a frequency range from 0 to 12 kHz is available. The input range can be selected via a dip-switch inside the pump (see page 20).

Pin 10, +5V_{dc}
This pin is connected to the stabilised +5V_{dc} voltage supply of the pump.

Pin 12, overload
In case of an overload, this pin (open-collector) is connected to the ground.

➔ Switch off the mains voltage and wait 2 minutes for cooling down.

Pin 13, internal speed
If pin 13 is connected with pin 1 (ground), the pump speed can be adjusted via the speed selector on the operating panel instead of using the signal at pin 5.

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Settings of switch S1

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Pins		DIP-Switch 1	DIP-Switch 2	DIP-Switch 3	DIP-Switch 4	DIP-Switch 5	DIP-Switch 6
Pin 5 speed IN	0-5V	20k Ω	ON*	OFF*	OFF*	OFF*	ON*
	0-10V	38k Ω	OFF	ON	OFF	OFF	ON
	0-20mA	130 Ω	OFF	OFF	ON	OFF	ON
	4-20mA	130 Ω	OFF	OFF	OFF	ON	OFF

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
→ Make sure that the pump is disconnected from the mains supply.

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⚠ The instrument should only be opened by a qualified technician! Capacitors inside the pump may still be charged even though the mains plug has been disconnected some time ago.

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Accessories

Foot switch
Order No. IS 10039

This foot switch serves as a start/stop device. It is very useful when using the pump as a dispenser for filling tubes, bottles, etc.. Both hands are free for handling the bottles and tubing.

Spare-cassettes in PVDF
MS/CA pressure lever*
Order No. IS 3629

CA pressure lever*
Order No. IS 3820

* The cassettes with pressure lever are available on request. Depending on the application, this type of cassette may provide better results at elevated differential pressure conditions.

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21

Zubehör

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MS/CA Click'n'go Cassette
MS/CA Click'n'go cassette

MS/CA Click'n'go
Bestell-Nr. IS 3510
MS/CA Anpresshebel*
Bestell-Nr. IS 0649

Accessories

Spare-cassettes in POM

MS/CA Click'n'go
Order No. IS 3510
MS/CA pressure lever*
Order No. IS 0649

Accessoires

Cassettes de rechange en POM

MS/CA Click'n'go
No de commande IS 3510
MS/CA levier de pression*
No de commande IS 0649

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CA Click'n'go
Order No. IS 3710
CA pressure lever*
Order No. IS 0122

Adaptor for type CA
(order 2 adaptors per cassette)

* The cassettes with pressure lever
are available on request.
Depending on the application, this
type of cassette may provide
better results at elevated
differential pressure conditions.

CA Click'n'go cassette
No de commande IS 3710
CA levier de pression*
No de commande IS 0122

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Mounting

Each pump-head is supplied with an instruction manual.

For selecting the correct tubing, please refer to our detailed tubing documentation.

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Interchangeable Pump-heads

These interchangeable pump-heads can also be ordered separately and mounted on the BVP drive.

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

0.004–100 ml/min
MS 3
ISM 720

Interchangeable Pump heads

➔ for the piston pump heads the
suffix „-W“ resp. „-WT“ means: with
isolation (W) or temperature gland
(WT) for crystallizing media.

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


Flow rates

In the tables listed subsequently the values indicated are only approximate and determined as follows:

ml/min, per channel, with water and Tygon® tubing, without differential pressure

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Please note

We will be pleased to help the user to select the tubing material. However, the user himself has the final responsibility for the selection of the correct tubing material.

For reproducible results we recommend you to use tubing from ISMATEC®

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⚠ Piston pump-heads

After use, the piston pump-heads require thorough flushing
- especially after pumping media containing salt, protein or particles.
Do not pump media containing particles that exceed a diameter of 0,8 mm.

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Maintenance

Provided the BVP drive is operated properly and in compliance with this manual, the tubing is the only part that is subject to wear and tear.

Repairs

For repairs please send the defective BVP drive to your ISMATEC® representative.

Please give information on defect and date of purchase.

Spare-parts

After the warranty period your ISMATEC® representative will be pleased to send you on request:

- spare-parts
- parts lists
- wiring diagrams

Please give information on defect, date of purchase, serial-no., model.

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Please note

We also recommend you to observe our Warranty Terms as well as our Terms and Conditions of Sale.

In case of any queries, please contact your local ISMATEC® representative.

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Ask for the latest
ISMATEC catalogue.

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ISMATEC® - Your competent partner for demanding metering and
dispensing applications.

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Laboratoriumstechnik GmbH
A Unit of IDEX Corporation