
Revisions				SO 0046-2
Rev From	Rev To	ECO #	Date	Appr.
0599	0803	1942-03		

Notes

- BD Cat. No. 231746
- Blank (Sheet) Size : Length: 12" Width: 13"
 Number of Pages: 8 Number of Sheets: 1
 Page Size: Length 12" Width 3.25" Final Folded Size: n/a
- Style (see illustrations below): #4

- See Specification Control No. n/a for Material Information
- Ink Colors: Printed two sides ☒ Yes ☐ No
 No. of Colors: 1 PMS# Black
- Graphics are approved by Becton, Dickinson and Company. Supplier has the responsibility for using the most current approved revision level.

VS Controlled by BD Caribe, LTD

Label Design	Date	COMPANY CONFIDENTIAL. THIS DOCUMENT IS THE PROPERTY OF BECTON, DICKINSON AND COMPANY AND IS NOT TO BE USED OUTSIDE THE COMPANY WITHOUT WRITTEN PERMISSION	 Becton, Dickinson and Company 250 Schilling Circle Cockeysville, MD. 21030-0243 USA
Proofer	Date		
Checked By	Date		
Part Number: L000146		Category and Description Package Insert Dryslide Oxidase	Sheet: 1 of 9 Scale: 1:1
			A

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INTENDED USE
BBL™ DrySlide™ Oxidase is used for determining the oxidase reaction of bacteria.

SUMMARY AND EXPLANATION
The oxidase test is a qualitative procedure for determining the presence or absence of cytochrome c oxidase activity in bacteria. This activity is dependent on the presence of an intracellular cytochrome oxidase system that catalyzes the oxidation of cytochrome c by molecular oxygen, which then serves as the terminal electron acceptor in the organism's electron transport system.¹
Organisms containing cytochrome c as part of their respiratory chain are oxidase positive and turn the reagent purple; organisms lacking cytochrome c as part of their respiratory chain do not oxidize the reagent, leaving it colorless within the time limits of the test, and are oxidase negative.

PRINCIPLES OF THE PROCEDURE
N,N,N',N'-tetramethyl-*p*-phenylenediamine dihydrochloride is an aromatic amine that is dimethylated at its two amino groups. The reagent is colorless in the reduced state and dark purple in the oxidized state.
In the oxidase test, cytochrome oxidase produced by the microorganism does not directly oxidize the *p*-phenylenediamine reagent. It oxidizes cytochrome c, which in turn oxidizes the reagent to form a purple-colored compound.

REAGENTS
BBL DrySlide Oxidase is a disposable slide having four filter paper reaction areas containing N,N,N',N'-tetramethyl-*p*-phenylenediamine dihydrochloride and ascorbic acid as a reducing agent and stabilizer.

Warnings and Precautions:
For *in vitro* Diagnostic Use.
Handle **BBL DrySlide** Oxidase only by the slide mount. Do not touch the reaction area.
Follow established laboratory procedures in handling and disposing of infectious materials.
Storage Instructions: Store **BBL DrySlide** Oxidase at 15 – 30°C.
Store opened, resealed **BBL DrySlide** Oxidase pouches with unused slides at 15 – 30°C.
Use **BBL DrySlide** Oxidase from opened, resealed pouches within 1 week. Discard unused slides after 1 week.
The reaction area of the **BBL DrySlide** Oxidase is light sensitive. Protect from light.
The expiration date applies to the product in its intact pouch when stored as directed.
Product Deterioration: Do not use a product if it fails to meet specifications for identity and performance.

SPECIMEN COLLECTION AND PREPARATION
Collect specimens or samples in sterile containers or with sterile swabs and transport immediately to the laboratory according to recommended guidelines.²⁻⁸
Process each specimen using procedures appropriate for that sample.²⁻⁸
Obtain the test organism from a primary isolation culture streaked on an appropriate solid medium (i.e., chocolate agar or tryptic soy blood agar).
Select well-isolated colonies. Mixed cultures require subculturing for purification. Insufficient growth indicates the need for further incubation.
Pick isolated colonies or a sweep of confluent growth from the pure culture.

PROCEDURE
Material Provided: **BBL DrySlide** Oxidase.
Materials Required But Not Provided: Inoculation loop (platinum or plastic) or wooden applicator stick and quality control organisms.

- Test Procedure:**
- Open the **BBL DrySlide** Oxidase pouch and remove the required number of slides.
After opening, the remaining **BBL DrySlide** slides may be stored in the original pouch for up to 1 week at room temperature. Fold the top of the pouch over once and seal tightly with a self-adhesive sticker (provided). Discard unused slides 1 week after the package is opened.
Do not place used or partially used **BBL DrySlide** Oxidase (with organisms) back into the pouch.
 - Using an appropriate inoculation device, pick isolated colonies or a sweep of confluent growth from the culture to be tested.
 - Smear the organism directly onto the reaction area of the **BBL DrySlide** Oxidase. Each reaction area can accommodate up to 4 tests. To ensure a proper reaction, spread the inoculum on the **BBL DrySlide** reaction area to a 3 – 4 mm size.
 - Examine the reaction area for appearance of a dark purple color within 20 sec. Disregard color development after 20 sec.

User Quality Control:
Identity Specifications – 2" x 2" slides with four reaction areas. The slides are contained in foil pouches. Reaction area should be white paper with no blue tint. Each pouch contains 3 slides.
Performance Response – Perform the oxidase test following the above procedure with the test organisms listed below.

Organism	ATCC™	Oxidase Reaction	Color Response
<i>Escherichia coli</i>	25922	–	Colorless (no change) to gray
<i>Neisseria sicca</i>	9913	+	Dark purple
<i>Pseudomonas aeruginosa</i>	27853	+	Dark purple

Quality control requirements must be performed in accordance with applicable local, state and/or federal regulations or accreditation requirements and your laboratory's standard Quality Control procedures. It is recommended that the user refer to pertinent NCCLS guidance and CLIA regulations for appropriate Quality Control practices.

RESULTS

Oxidase-positive organisms produce a purple or dark color within 20 sec. Oxidase-negative organisms produce no color change or a change to light gray within the 20 sec test period. As noted, disregard color development after 20 sec.

LIMITATIONS OF THE PROCEDURE

The oxidase test is used in the presumptive identification of *Neisseria* and in the differentiation and identification of gram-negative bacilli. Examine all other oxidase-positive organisms by Gram stain to determine cellular morphology and Gram reaction. Additional biochemical testing using pure culture is recommended for complete identification. For further information, consult appropriate references.²⁻⁸

Weak oxidase producers, such as *Pasteurella*, may appear negative within the time limits of the test.⁸

Mucoid *Pseudomonas aeruginosa* may require longer than 10 – 20 sec to become positive.

Repeat oxidase reactions that are inconsistent with other biochemical reactions.

Determine oxidase reactions of gram-negative bacilli on colonies obtained from nonselective and nondifferential media to ensure valid results.⁹

Media containing high levels of glucose may inhibit oxidase activity, causing false-negative oxidase reactions.

PERFORMANCE CHARACTERISTICS

Prior to release, all lots of **BBL DrySlide** Oxidase slides are tested for specific product specifications. Samples of the lot are inoculated in the reaction area of the slide with a loopful of 18- to 24-h cultures of *Escherichia coli* (ATCC 25922), *Neisseria sicca* (ATCC 9913) and *Pseudomonas aeruginosa* (ATCC 27853). The slides are examined 20 sec after inoculation for a color change.

AVAILABILITY

Cat. No.	Description
231746	BBL™ DrySlide™ Oxidase, 25 x 3 slides.

REFERENCES

1. Steel, K.J. 1961. The oxidase reaction as a taxonomic tool. *J. Gen. Microbiol.* 25: 297-306.
2. Baron, E.J., L.R. Peterson, and S.M. Finegold. 1994. *Bailey & Scott's diagnostic microbiology*, 9th ed. Mosby-Year Book, Inc., St. Louis.
3. Pezzlo, M. (ed.) 1994. Aerobic bacteriology, p. 1.0.1-1.20.47. *In* H.D. Isenberg (ed.), *Clinical microbiology procedures handbook*, vol. 1 American Society for Microbiology, Washington, D.C.
4. Cunniff, P. (ed.). 1995. Official methods of analysis of AOAC International, 16th ed. AOAC International, Gaithersburg, Md.
5. FDA bacteriological analytical manual, 8th ed. AOAC International, Gaithersburg, Md.
6. Vanderzant, C., and D.F. Splittstoesser (ed.). 1992. *Compendium of methods for the microbiological examination of foods*, 3rd ed. American Public Health Association, Washington, D.C.
7. Flowers, R.S., W. Andrews, C.W. Donnelly, and E. Koenig. 1992. Pathogens in milk and milk products, p. 103 – 212. *In* R.T. Marshall (ed.), *Standard methods for the examination of dairy products*, 16th ed. American Public Health Association, Washington, D.C.
8. Murray, P.R., E.J. Baron, M.A. Pfaller, F.C. Tenover, and R.H. Tenover (ed.). 1995. *Manual of clinical microbiology*, 6th ed. American Society for Microbiology, Washington, D.C.
9. MacFaddin, J.F. 1980. *Biochemical tests for the identification of medical bacteria*, 2nd ed. Williams and Wilkins, Baltimore.

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© 2003 BD.
ATCC is a trademark of the American Type Culture Collection.
Dryslide is a trademark of Difco Laboratories, a subsidiary of Becton, Dickinson and Company.
BD, BD Logo and BBL are trademarks of Becton, Dickinson and Company.
BENEX Limited